


CytoFLEX mosaic

Spectral Detection Module for CytoFLEX S and CytoFLEX LX Flow Cytometers For Research Use Only. Not for use in diagnostic procedures.



D17052AA December 2024





CytoFLEX mosaic Spectral Detection Module for CytoFLEX S and CytoFLEX LX Flow Cytometers

PN D17052AA (December 2024)

© 2024 Beckman Coulter, Inc. All rights reserved.

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

All other trademarks are the property of their respective owners.

Contact Us

If you have any questions, contact our Customer Support Center.

- Worldwide, find us via our website at www.beckman.com/support/technical.
- In the USA and Canada, call us at 1-800-369-0333.
- In Austria, call us at 0810 300484.
- In Germany, call us at 02151 333999.
- In Sweden, call us at +46(0)8 564 859 14.
- In the Netherlands, call us at +31 348 799 815.
- In France, call us at 0825838306 6.
- In the UK, call us at +44 845 600 1345.
- In Ireland, call us at +353 (01) 4073082.
- In Italy, call us at +39 0295392 456.
- In other locales, contact your local Beckman Coulter Representative.

EC REP

Beckman Coulter Eurocenter S.A. 22, rue Juste-Olivier Case Postale 1044 CH - 1260 Nyon 1, Switzerland Tel: +41 (0) 22 365 36 11

Glossary of Symbols is available at beckman.com/techdocs (PN C24689).

May be covered by one or more pat. - see www.beckman.com/patents

Original Instructions

Revision History

Initial Issue AA, 12/2024 Software Version 1.0

This document applies to the latest software listed and higher versions. When a subsequent software version affects the information in this document, a new issue will be released to the Beckman Coulter Web site. For labeling updates, go to www.beckman.com and download the latest version of the manual or system help for your instrument.

D17052AA iii

iv D17052AA

Safety Notices

Read all product manuals and consult with Beckman Coulter-trained personnel before attempting to operate the instrument. Do not attempt to perform any procedure before carefully reading all instructions. Always follow product labeling and manufacturer's recommendations. If in doubt as to how to proceed in any situation, contact us.

Beckman Coulter, Inc. urges its customers to comply with all national health and safety standards such as the use of barrier protection. This may include, but it is not limited to, protective eyewear, gloves, and suitable laboratory attire when operating or maintaining this or any other automated laboratory analyzer.

This manual assumes that users have basic knowledge of the Windows operating system, as well as experience working with laboratory testing technology. Users are invited to consult the appropriate documentation for such information.

Alerts for Danger, Warning, and Caution



DANGER indicates an imminently hazardous situation which, if not avoided, will result in death or serious injury.



WARNING indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.



CAUTION indicates a potentially hazardous situation, which, if not avoided, may result in minor or moderate injury. It may also be used to alert against unsafe practices.

D17052AA V

Safety Precautions



Risk of operator injury if:

- All doors, covers and panels are not closed and secured in place prior to and during instrument operation.
- The integrity of safety interlocks and sensors is compromised.
- Instrument alarms and error messages are not acknowledged and acted upon.
- You contact moving parts.
- You mishandle broken parts.
- Doors, covers and panels are not opened, closed, removed and/or replaced with care.
- Improper tools are used for troubleshooting.

To avoid injury:

- Keep doors, covers and panels closed and secured in place while the instrument is in use.
- Take full advantage of the safety features of the instrument.
- Acknowledge and act upon instrument alarms and error messages.
- Keep away from moving parts.
- Report any broken parts to your Beckman Coulter Representative.
- Open/remove and close/replace doors, covers and panels with care.
- Use the proper tools when troubleshooting.



System integrity could be compromised and operational failures could occur if:

- This equipment is used in a manner other than specified. Operate the instrument as instructed in the Product Manuals.
- You introduce software that is not authorized by Beckman Coulter into your computer. Only operate your system's computer with software authorized by Beckman Coulter.
- You install software that is not an original copyrighted version. Only use software that is an original copyrighted version to prevent virus contamination.



If you purchased this product from anyone other than Beckman Coulter or an authorized Beckman Coulter distributor, and, it is not presently under a Beckman Coulter service maintenance agreement, Beckman Coulter cannot guarantee that the product is fitted with the most current mandatory engineering revisions or that you will receive the most current information bulletins concerning the product. If you purchased this product from a third party and would like further information concerning this topic, contact us.

Vİ D17052AA



Risk of instrument damage. This device is intended for indoor use only. To avoid device damage, do not install the instrument outdoors.



Risk of personal injury. Safety protection can be impaired if used in a manner not specified by the manufacturer. To avoid personal injury, use the instrument according to the manufacturer's instructions only.

Symbol Explanations

Symbol	Warning Condition	Action
	Biohazard Symbol	Consider all materials (specimens, reagents, controls, and monoclonal antibodies) and areas these materials come into contact with as being potentially infectious.
		Wear appropriate barrier protection and follow safe laboratory procedures when handling any material in the laboratory.
	Caution	To indicate that caution is necessary when operating the device or control close to where the symbol is placed, or to indicate that the current situation needs operator awareness or operator action in order to avoid undesirable consequences.
4	Caution, Risk of Electric Shock	To indicate hazards arising from dangerous voltages.
	Warning, Sharp elements	To warn of sharp element.
	Warning, Crushing of hands	To warn of a closing motion of mechanical parts of equipment.

D17052AA vii

Symbol	Warning Condition	Action
	Recycling Symbol WEEE Wheeled Bin Symbol	The symbol of a crossed-out wheeled bin on the product is required in accordance with the Waste Electrical and Electronic Equipment (WEEE) Directive of the European Union. The presence of this marking on the product indicates:
		1. that the device was put on the European Market after August 13, 2005 and
		2. that the device is not to be disposed via the municipal waste collection system of any member state of the European Union.
		For products under the requirement of WEEE directive, please contact your dealer or local Beckman Coulter office for the proper decontamination information and take-back program which will facilitate the proper collection, treatment, recovery, recycling, and safe disposal of device.
		For the Japan market:
		This system is considered an industrial waste, subject to special controls for infectious waste. Before disposal of the system, refer to the Waste Disposal and Public Cleaning Law for compliance procedures.
	Laser Radiation	The laser radiation symbol indicates that there can be laser light radiation in the area. Take precautions to prevent exposure.
50	RoHS Caution Symbol	This label indicates that the electronic information product contains certain toxic or hazardous substances. The center number is the Environmentally Friendly Use Period (EFUP) date, and indicates the number of calendar years the product can be in operation. Upon the expiration of the EFUP, the product must be immediately recycled. The circling arrows indicate the product is recyclable. The date code on the label or product indicates the date of manufacture.
CLASS 1 LASER PRODUCT COMPLIES WITH 21 OFR 1040, 10 AND 1040, 11 EXCEPT FOR CONFERMANCE WITH 166 06255-16 d. 3, AS DESCRIBED IN LASER NOTICE NO.56, DATED MY 8, 2019, MANUFACTURED 7	Laser Compliance	This label indicates that the product is a Class 1 Laser Product and is in compliance with international standard and US requirements.
CLASS 1 LASER PRODUCT PRODUIT LASER CLASSE 1 1 失激光产品	Laser Class I	This label indicates that this product is a Class I Laser product. Take precautions to prevent exposure.

VIII D17052AA

Symbol	Warning Condition	Action
	Consider all materials (specimens, controls, monoclonal antibodies, and so forth) as being potentially infectious.	Wear standard laboratory attire and follow safe laboratory procedures when handling any material in the laboratory.
EC REP	EC Representative	Indicates the Authorized representative in the European Community.
CE	European Conformity (CE Mark) Regulatory Mark	A "CE" mark indicates that a product has been assessed before being placed on the market, and has been found to meet European Union safety, health, and/or environmental protection requirements.
UK	UKCA Mark	A "UKCA" mark indicates that a product has been assessed before being placed in UK market, and has been found to meet UK safety, health, and/or environmental protection requirements.
TÜVRheinland c us	TUV symbol	This test mark is proof of compliance with Canadian National standards and U.S. National standards. It demonstrates an electrical product has been successfully tested and certified by a Certification Body accredited by the Standards Council of Canada (SCC) and by a Nationally Recognized Testing Laboratory (NRTL) accredited by the Occupational Safety and Health Administration (OSHA).
	Off (power)	To indicate disconnection from the mains, at least for mains switches or their positions, and all those cases where safety is involved.
	On (power)	To indicate connection to the mains, at least for mains switches or their positions, and all those cases where safety is involved.
	Fuse	To identify fuse boxes or their location
	Manufacturer	Indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.
<u>i</u>	Consult the Instructions for Use Symbol	Indicates the need for the user to consult the instructions for use.

D17052AA ix

Symbol	Warning Condition	Action
<u></u>	Date of Manufacture	Indicates the date when the medical device was manufactured.
SN	Serial Number	Indicates the manufacturer's serial number so that specific medical device can be identified.
#	Model Number	Indicates the model number or type number of a product.
	RCM Symbol	The "RCM" (Regulatory Compliance Mark) is depicted as a triangle with a partial circle and check. The mark is applied to products that comply with the EMC requirements of the Australian Communications Media Authority (ACMA) for use in Australia and New Zealand.
Info for USA only: California Proposition 65 WARNING Cancer & Reproductive Harm Sawww P65Warnings.ca.gov	California Proposition 65	This product can expose you to chemicals known to the State of California to cause Cancer and Reproductive Harm. For more information go to https://www.P65Warnings.ca.gov.

X D17052AA

Electromagnetic Compatibility (EMC)

This device complies with the emissions and immunity requirements as specified in the EN/IEC 61326 series of Product Family Standards for a "basic electromagnetic environment". Such equipment is supplied directly at low voltage from public mains network. This equipment is not intended for residential use.



This device generates, uses, and can radiate unintentional radio-frequency (RF) energy. If this device is not installed and operated correctly, this RF energy can cause interference with other equipment. It is the responsibility of the end user to be sure that a compatible electromagnetic environment for the device can be maintained so that the device operates as intended.

In addition, other equipment can radiate RF energy to which this device is sensitive. If one suspects interference between this device and other equipment, Beckman Coulter recommends the following actions to correct the interference:

- Evaluate the electromagnetic environment before installation and operation of this device.
- 2. Do not operate this device close to sources of strong electromagnetic radiation (for example: unshielded intentional RF sources), as these can interfere with proper operation. Examples of unshielded intentional radiators are handheld radio transmitters, cordless phones, and cellular phones.
- 3. Do not place this device near medical electrical equipment that can be susceptible to malfunctions caused by close-proximity to electromagnetic fields.
- 4. This device has been designed and tested to CISPR 11, Class A emission limits. In a domestic environment, this device may cause radio interference, in which case, you may need to take measures to mitigate the interference.

D17052AA xi

Safety Notices Electromagnetic Compatibility (EMC)

D17052AA χij

Contents

Revision History, iii Safety Notices, v Alerts for Danger, Warning, and Caution, v Safety Precautions, vi Symbol Explanations, vii Electromagnetic Compatibility (EMC), xi Introduction, xxvii Overview, xxvii How to Use Your Manual, xxvii About this Manual, xxvii Conventions, xxix Graphics, xxx Advancing Your Skills, xxx System Overview, 1-1 Overview, 1-1 Product Description, 1-1 Main Components, 1-2 Optical Components, 1-4 Wavelength Division Multiplexer (WDM), 1-4 Optical Fiber, 1-4 Fluidics System, 1-6 Fluid Containers/Cubitainers, 1-6 Fluidics Module, 1-6 Sample Station, 1-6 Sample Tube Holder Positions, 1-6 Plate Loader Components, 1-6 Plate Holder Components, 1-6 System Configuration, 1-7 System Configuration [CytoFLEX S], 1-7 System Configuration [CytoFLEX LX], 1-9

System Configuration [CytoFLEX mosaic Spectral Detection

CHAPTER 1:

CHAPTER 2:

```
Module], 1-11
Consumables and Supplies, 1-12
         Reagents, 1-12
         Ordering Information, 1-12
Instrument Specifications, 1-13
         Dimensions [CytoFLEX S], 1-13
         Dimensions [CytoFLEX LX], 1-13
         Dimensions [CytoFLEX mosaic Spectral Detection Module], 1-13
         Installation Category, 1-13
         Maximum Altitude, 1-13
         Pollution Degree, 1-13
         Acoustic Noise Level, 1-13
         Electrical Ratings, 1-14
         Cytometer, 1-14
         CytoFLEX mosaic Spectral Detection Module, 1-14
Performance Characteristics, 1-16
         Performance Characteristics [CytoFLEX S], 1-16
         Performance Characteristics [CytoFLEX LX], 1-16
         Performance Characteristics [CytoFLEX mosaic Spectral Detection
            Module], 1-17
Reagent Limitations, 1-18
Safety Data Sheets (SDS/MSDS), 1-18
Using the CytExpert for Spectral Software, 2-1
Overview, 2-1
Launching the Software, 2-1
Main Software Screen, 2-2
         Start Page, 2-2
         Acquisition Screen, 2-4
            Acquisition Screen Navigation, 2-5
            Collection, 2-6
            Collection [With Plate Loader], 2-7
            Test Tubes, 2-8
            Test Tube Management Module, 2-9
            Plot Area, 2-10
            Status Bar, 2-12
         Analysis Screen, 2-13
         Spectral Unmixing Experiment Screen, 2-14
             Spectral Unmixing Controls, 2-15
         QC Experiment Screen, 2-16
             QC Report Screen, 2-16
            QC Experiment Screen, 2-17
             QC Screen Navigation, 2-17
         Software Menu, 2-18
            Acquisition and Analysis Screen Menu, 2-20
```

```
User Management, 2-27
         Creating, Deleting, and Modifying Users in User Manager, 2-28
         Unlocking a User Account, 2-30
         Resetting a User Password, 2-31
         Changing a User Password, 2-32
         Forgot Username or Password, 2-33
         Exporting the Credential File, 2-35
Role Management, 2-36
         Creating, Deleting, and Modifying User Roles in Role Manager, 2-36
Account Policies, 2-36
User Management Operation Log, 2-39
Graphic and Gating Styles, 2-39
         Plots, 2-39
         Gates, 2-40
Plate Type Library, 2-40
         Adding a Plate Type, 2-40
         Editing a Plate Type, 2-40
         Duplicating a Plate Type, 2-40
         Deleting a Plate Type, 2-41
Fluorescent Tag Library, 2-41
         Adding a Fluorescent Tag, 2-42
         Deleting a Fluorescent Tag, 2-42
         Restoring the Default Settings, 2-43
         Exporting Fluorescent Tag Files, 2-43
         Importing Fluorescent Tag Files, 2-44
         Editing a Fluorescent Tag, 2-45
Software Settings, 2-45
Operation Principles, 3-1
Overview, 3-1
Sample Flow, 3-1
         Sample Loading, 3-1
         Hydrodynamic Focusing, 3-1
Laser Beam Shaping, 3-2
Cell Illumination, 3-2
         Forward Scatter, 3-2
         Side Scatter and Fluorescent Light, 3-2
Light Collection, Separation and Measurement, 3-2
         Forward Scatter Collection, 3-2
         Side Scatter and Fluorescent Light Collection, 3-2
             Side Scatter, 3-2
             Fluorescent Light, 3-3
Signal Processing, 3-3
```

CHAPTER 3:

```
Data Storage, 3-4
                     Automated Software Features, 3-4
                     Parameters, 3-4
                              TIME Parameter, 3-4
                     Plot Display, 3-4
                     Statistics, 3-5
CHAPTER 4:
                     Daily Startup, 4-1
                     Overview, 4-1
                     Pre-Startup Inspection, 4-1
                              Checking Waste and Reagent Levels [4 L Fluid Containers], 4-1
                              Checking Waste and Reagent Levels [10 L Fluid Cubitainers], 4-1
                              Power Source Inspection, 4-2
                              Workstation Connections Inspection, 4-2
                     Turning On the Instrument, 4-2
                     Logging into the Software, 4-3
                              Logging Out of the Software, 4-6
                              Locking the Account, 4-6
                              Selecting the Proper Sample Injection Mode, 4-6
                              Selecting the Plate Loader Sample Injection Mode [With Plate
                                 Loader, 4-7
                              Running the System Startup Program [with the Single Tube
                                 Loader, 4-7
                              Running the System Startup Program [With Plate Loader], 4-7
                              Selecting Experiments from the Start Page, 4-7
                     Initializing the Instrument, 4-7
CHAPTER 5:
                     Instrument Quality Control and Standardization, 5-1
                     Overview, 5-1
                     Quality Control, 5-2
                              Preparing the QC Sample, 5-3
                                 Required Materials, 5-3
                                 CytoFLEX Ready to Use Daily QC Fluorospheres Preparation
                                     Process, 5-3
                                 CytoFLEX Daily IR QC Fluorospheres Preparation Process, 5-3
                              Preparing the QC Sample [With Plate Loader], 5-3
                                  Required Materials, 5-3
                                 CytoFLEX Ready to Use Daily QC Fluorospheres Preparation
                                     Process, 5-4
                                 CytoFLEX Daily IR QC Fluorospheres Preparation Process, 5-4
                              Importing Lot-Specific Target Values, 5-4
                              Collecting QC Data, 5-7
                              Collecting QC Data [With Plate Loader], 5-11
                              Confirming Results, 5-15
```

Creating Levey-Jennings Charts, 5-18 QC Result Manager, 5-18
Standardization, 5-18 Preparing the Standardization Sample, 5-18 Required Materials, 5-18 Preparation Process, 5-18 Generating Target Median Values, 5-19 Creating a New Standardization Item, 5-23 Adding a New Standardization Item, 5-23 Creating a New Standardization Item from Acq. Setting Catalog, 5-27 Performing the Standardization, 5-32 Applying the Standardized Acquisition Settings, 5-35 Standardization Target Library, 5-37
Data Acquisition and Sample Analysis, 6-1
Overview, 6-1
Creating an Experiment, 6-1 Creating an Experiment [Without Plate Loader], 6-1 Creating an Experiment [With Plate Loader], 6-2 Setting Sample Wells, 6-2 Modifying Well Settings, 6-7 Setting the Label, 6-7 Creating Plots and Gates, 6-8 Creating and Adjusting Auto Gates, 6-17 Setting Customized Parameters, 6-17 Setting Custom Statistics, 6-19
Configuring Acquisition Settings, 6-19 Changing the Tube Name, 6-19 Laser Settings, 6-19 Setting Laser Target Power Settings [CytoFLEX LX Only], 6-19 Adjusting the Gain, 6-19 Link to Acq. Setting, 6-21 Adjusting the Threshold, 6-23 Setting Collection Conditions, 6-24 Setting Plot Display Conditions, 6-24 Viewing Saturation [Without Plate Loader], 6-24 Viewing Saturation [With Plate Loader], 6-28
Load Sample and Record Data, 6-33 Before Running Samples, 6-33 Sampling and Collecting Data [Without Plate Loader], 6-33 Sampling and Collecting Data [With Plate Loader], 6-33
Analyzing and Exporting Data, 6-33 Exporting FCS Files, 6-39 Exporting Single Tube Files, 6-39 Exporting Multiple FCS Files, 6-41 Exporting Plots or the Statistics Table of Multiple Tubes as Picture

CHAPTER 6:

CHAPTER 7:

CHAPTER 8:

```
Files, 6-43
         Importing and Exporting Instrument Settings, 6-44
            Importing Instrument Settings, 6-44
            Exporting Instrument Settings, 6-45
            Importing and Exporting Unmixing Settings, 6-46
         Printing Graphics, 6-46
Saving the Experiment, 6-46
         Concluding the Experiment, 6-46
Unmixing, 7-1
Overview, 7-1
Preparing the Spectral Unmixing Sample, 7-1
         Preparing the Spectral Unmixing Sample [Without Plate
             Loader, 7-2
         Preparing the Spectral Unmixing Sample [With Plate Loader], 7-2
Creating a Spectral Unmixing Experiment, 7-2
         Creating a Spectral Unmixing Experiment [Without Plate
            Loader, 7-2
         Creating a Spectral Unmixing Experiment [With Plate Loader], 7-5
Loading Sample and Recording Data, 7-9
Calculating Unmixing Values, 7-10
Confirming Unmixing Results, 7-12
Saving the Spectral Unmixing Experiment, 7-14
         Concluding the Experiment, 7-14
Adjusting Spectral Unmixing, 7-15
         Importing and Exporting Unmixing Matrix, 7-15
            Importing Unmixing Settings from Unmixing Matrix Files, 7-15
            Importing Unmixing Settings from the Spectral Library, 7-17
            Exporting Unmixing Settings to Unmixing Matrix Files, 7-20
            Exporting Unmixing Settings to Unmixing Matrix Files from the
                Spectral Library, 7-22
         Managing Spectral Library, 7-23
            Importing Spectral Library Files, 7-23
            Exporting Spectral Library Files, 7-25
            Deleting Spectral Library Files, 7-26
Fine Adjustment, 7-27
         Auto Fine Adjusting, 7-27
         Manual Fine Adjustment, 7-27
         Fine Adjustment Matrix, 7-28
Data Review, 8-1
Overview, 8-1
```

Copying Experiments and Importing Data, 8-1

Copying a Previously Acquired Experiment, 8-1

```
Importing Previously Acquired Data, 8-1
Setting the Plots and Statistics, 8-2
Calculating Sample Volume and Concentration, 8-2
Adjusting Unmixing Settings, 8-2
Exporting Results, 8-3
Daily Shutdown, 9-1
Overview, 9-1
Preparing the Cleaning Solution, 9-1
Running Daily Clean, 9-1
Shutting Down the Instrument, 9-1
Auto Shutdown [CytoFLEX LX Only], 9-2
Troubleshooting, 10-1
Overview, 10-1
Precautions/Hazards, 10-1
         Laser Related Hazards, 10-1
            Laser Beam Hazards, 10-2
            Laser Warning Labels, 10-2
Disposal Precaution, 10-6
Troubleshooting Table, 10-6
Backup and Restore, 10-8
         Backup, 10-8
         Restore, 10-8
         Log Cleanup, 10-8
Cleaning Procedures, 11-1
Overview, 11-1
Preparing the Instrument for Transport or Storage, 11-1
         Lifting and Carrying Instructions [CytoFLEX S and CytoFLEX
            LX], 11-2
         Lifting and Carrying Instructions [CytoFLEX mosaic Spectral
            Detection Module, 11-2
Replacement/Adjustment Procedures, 12-1
Overview, 12-1
Switching between Spectral Mode and Conventional Mode, 12-2
```

Switching to Conventional Mode, 12-3 Switching to Spectral Mode, 12-7

CHAPTER 9:

CHAPTER 10:

CHAPTER 11:

CHAPTER 12:

APPENDIX A:

APPENDIX B:

Instrument Installation, A-1 Overview, A-1 Instrument Transportation and Storage, A-1 Installation Environment Validation, A-2 Worktable, A-2 Ventilation and Cleaning, A-2 Power Source, A-3 Temperature and Humidity, A-4 Waste Disposal, A-4 Unpacking the Instrument and Inspecting the Materials for Defects or Omissions, A-5 Installing the Instrument and Connecting the Equipment, A-6 Installing the Instrument and Connecting the Equipment [CytoFLEX S], A-6 Installing the Instrument and Connecting the Equipment [CytoFLEX LX], A-6 Installing the CytoFLEX mosaic Spectral Detection Module, A-6 CytExpert for Spectral Software Installation Options, A-7 Installing the CytoFLEX for Spectral Software, A-7 Reinstalling the CytExpert for Spectral Software, A-8 Uninstalling the Software, A-8 Reinstalling the CytExpert for Spectral Software, A-14 Restoring the Cytometer and CytoFLEX mosaic Spectral Detection Module Configuration File, A-19 CytExpert for Spectral Electronic Record Management, B-1 Overview, B-1 Software Menu, B-1 Experiment Management, B-2 Closed File System, B-2 Experiment Directory Management, B-2 Folder Hierarchy Management, B-2 Experiment Related Operations, B-3 Importing an Experiment/Template, B-3 Exporting an Experiment/Template, B-3 Log, B-3 Electronic Signature, B-4 User Management, B-4 User Administration, B-4 Logging In and Out of the Software, B-4 Locking the Account, B-4 Role Management, B-4

User Management, B-4 Account Policies, B-4

APPENDIX C: Sample Injection Mode Control Kit, C-1

Overview, C-1

APPENDIX D: Deep Well Plate, D-1

Specimen Collection Plate Specifications, D-1

APPENDIX E: Cyber Security, E-1

Overview, E-1

Drive Encryption, E-2

BitLocker Keys, E-2

Changing Encryption Strength Level for BitLocker, E-2

Enabling BitLocker, E-6

Backing Up Your Recovery Key, E-13

BitLocker Decryption, E-17

Suspending BitLocker Drive Encryption, E-21

Recovering Your BitLocker Key, E-25

Protection from Malware Software, E-27

Installing Trellix Application Control, E-28

Enabling Trellix Application Control, E-35

Exporting the Certificate, E-39

Importing BEC Certificate, E-47

Switching Trellix Application Control to Update Mode, E-48

Disabling Trellix Application Control, E-50

Security Updates, E-52

Enabling System Protection, E-52

Creating a Restore Point, E-55

Downloading Security Updates, E-57

Notification Options, E-58

Installing the Security Updates, E-58

Recovering from Failed Security Updates, E-62

User Management, E-63

System Hardening, E-64

Enabling Automatic Log-Off, E-64

Setting Account Lockout Policy, E-74

Setting Password Policy, E-76

Disabling Remote Access, E-79

Disabling Auto Play, E-81

Disabling Unnecessary Services, E-83

Disabling Unauthorized Applications, E-88

Enabling SMB Digitally Sign Communications, E-91

Enabling Audit Object Access, E-94

Viewing Windows Security Logs, E-104

Enabling Installation Restriction, E-105

Enabling Microsoft Defender Firewall, E-108 Enabling Microsoft Defender, E-110 Enabling Network Time Protocol, E-111

Remote Access, E-116

APPENDIX F: Table of Hazardous Substances, F-1

Table of Hazardous Substances, F-1

Abbreviations

Beckman Coulter, Inc.

SOFTWARE END-USER LICENSE AGREEMENT

Related Documents

Illustrations

1.1	Main Components [CytorLEX 3 Without Plate Loader], 1-2
1.2	Main Components [CytoFLEX LX], 1-3
1.3	System Connections, 1-7
1.4	Power Entry Module, 1-8
1.5	Front of Cytometer [CytoFLEX S Without Plate Loader Shown], 1-8
1.6	System Connections, 1-9
1.7	Power Input Module, 1-10
1.8	Front of Cytometer, 1-10
1.9	Front and Left Sides of the CytoFLEX mosaic Spectral Detection Module, 1-11
1.10	Back Cover Connections of the CytoFLEX mosaic Spectral Detection Module, 1-12
2.1	Plot Area [Raw Worksheet Shown], 2-10
2.2	Plot Area [Unmixed Worksheet Shown], 2-11
2.3	Drawing Controls Toolbar (Top of Screen), 2-14
2.4	QC Report Screen, 2-16
2.5	Software Menu Tree*, 2-18
2.6	QC Software Menu Tree, 2-19
2.7	User Manager (Card View), 2-27
2.8	User Manager (Grid View), 2-27
2.9	Account Policies - Password Policy, 2-37
2.10	Account Policies - Account Lockout Policy, 2-38
2.11	Account Policies - Application Inactivity Policies, 2-38
2.12	Fluorescent Tag Library, 2-41
3.1	Light Path through the WDM with a Single Port, 3-3
6.1	Plate Window [CytoFLEX LX Shown], 6-2
6.2	All Gates - Example Experiment, 6-15
6.3	Circular Gating Logic - Example Experiment, 6-15
10.1	Laser Warning Label on the Laser Optical Bench [CytoFLEX LX], 10-3
10.2	Laser Warning Label within the Optical Bench (Located Inside the Cytometer) [CytoFLEX S], 10-3
10.3	Laser Warning Label within the Optical Bench (Located Inside the Cytometer) [CytoFLEX LX], 10-4

10.4	Laser Warning Label on the 355 nm Laser [CytoFLEX LX], 10-4
10.5	Laser Warning Labels on the Cytometer Back Cover [CytoFLEX S], 10-5
10.6	Laser Warning Labels on the Cytometer Back Cover [CytoFLEX LX], 10-5
10.7	Laser Warning Labels on the CytoFLEX mosaic Spectral Detection Module Back Cover, 10-6

Tables

10.1	Troubleshooting [With CytoFLEX mosaic Spectral Detection Module], 10-7
12.1	Optical Fibers with Indicators [CytoFLEX S], 12-3
12.2	Optical Fibers with Indicators [CytoFLEX LX], 12-3
E.1	Role-based Access Control Policies, E-64
E.2	List of Unnecessary Services, E-83
F.1	Table of Hazardous Substances Name and Concentration [CytoFLEX mosaic 63], F-2
F.2	Table of Hazardous Substances Name and Concentration [CytoFLEX mosaic 88], F-3

Introduction

Overview

This is an introductory chapter that gives a brief description of each document. It contains the following topics presented in this order:

- How to Use Your Manual
- About this Manual
- Conventions
- Graphics
- Advancing Your Skills

How to Use Your Manual

Your CytoFLEX mosaic Spectral Detection Module includes the manuals listed below:

- Use this **Instructions for Use** manual for information on the day-to-day operation of your CytoFLEX flow cytometer equipped with the CytoFLEX mosaic Spectral Detection Module. You can find detailed step-by-step procedures for Daily Startup and Quality Control, running samples, analyzing data, and performing Startup and Shutdown. This manual also contains physical and system specifications, safety and troubleshooting information, as well as information about what your CytoFLEX mosaic Spectral Detection Module does and the methods guiding operation. It also contains procedures for cleaning and maintenance.
- The CytoFLEX Platform Instructions for Use manual provides information on the day-to-day operation of your CytoFLEX flow cytometer not equipped with a CytoFLEX mosaic Spectral Detection Module.

NOTE To switch between Spectral Mode and Conventional Mode, refer to CHAPTER 12, Switching between Spectral Mode and Conventional Mode.

About this Manual

The information in the Instructions for Use manual is organized as follows:

CHAPTER 1, System Overview

Provides information regarding the individual components of the CytoFLEX flow cytometer and the CytoFLEX mosaic Spectral Detection Module, and the corresponding functions of these components.

D17052AA XXVII

CHAPTER 2, Using the CytExpert for Spectral Software

Provides an overview of each aspect of the software's functions.

CHAPTER 3, Operation Principles

Describes how the instrument measures scattered light and fluorescence as cells pass through the laser beam.

CHAPTER 4, Daily Startup

Provides instructions for starting your CytoFLEX flow cytometer equipped with a CytoFLEX mosaic Spectral Detection Module and navigating to the sample testing standby state.

CHAPTER 5, Instrument Quality Control and Standardization

Provides instructions for performing daily quality control (QC) on your CytoFLEX flow cytometer equipped with a CytoFLEX mosaic Spectral Detection Module to confirm the instrument is working correctly and to ensure accurate experimental data measurement. Quality control allows you to determine whether your instrument can provide adequate signal strength and precision.

CHAPTER 6, Data Acquisition and Sample Analysis

Provides instructions for operating the instrument, including data acquisition, analyzing and exporting results.

CHAPTER 7, Unmixing

Describes how to create an spectral unmixing experiment and automatically calculate unmixing values after acquiring the single color data. It also explains how to use these calculations for other experiments.

CHAPTER 8, Data Review

Describes how to use the Analysis screen to analyze data that has already been acquired.

CHAPTER 9, Daily Shutdown

Describes how to keep the instrument in optimal condition through daily cleaning during the shutdown procedure.

CHAPTER 10, Troubleshooting

Describes some common problems and their solutions in a basic troubleshooting matrix.

CHAPTER 11, Cleaning Procedures

Describes how to carry out certain routine and nonscheduled cleaning procedures.

CHAPTER 12, Replacement/Adjustment Procedures

Describes how to carry out certain routine and nonscheduled replacement and adjustment procedures.

APPENDIX A, Instrument Installation

Provides the instrument installation procedures for your CytoFLEX flow cytometer.

APPENDIX B, CytExpert for Spectral Electronic Record Management

Provides the instructions for using the CytExpert Electronic Record Management software option.

APPENDIX C, Sample Injection Mode Control Kit

Provides the instructions for using the CytoFLEX Sample Injection Mode Control Kit.

XXVIII D17052AA

APPENDIX D, Deep Well Plate

Provides a list of deep well plates suggested for use on the Plate Loader DW.

APPENDIX E, Cyber Security

Provides instructions for cyber security.

APPENDIX F, Table of Hazardous Substances

Provides the table of hazardous substances with the hazardous substance name and concentration.

Conventions

This document uses the following conventions:

- **Bold face** font indicates buttons or selections that appear on the workstation screen.
- The term "select" is used to indicate the following action:
 - To click with a mouse.

NOTE The verb "press" is reserved for mechanical buttons, such as keys on the keyboard.

- The software path to a specific function or screen appears with the greater than (>) symbol between screen options.
- Links to information in another part of the document for additional information are in blue and are underlined. To access the linked information, select the blue, underlined text.
- The information in your Instructions for Use manual applies to the CytoFLEX S instruments equipped with and without a plate loader, and the CytoFLEX LX, unless otherwise specified.

NOTE When information in this document only applies to the plate loader configuration, it is marked [With Plate Loader]. When information in this document only applies to the configuration not equipped with a plate loader, it is marked [Without Plate Loader].

When information in this document only applies to one instrument in the platform, it is marked either [CytoFLEX S] or [CytoFLEX LX].

When information in this document only applies to the CytoFLEX mosaic Spectral Detection Module, it is marked [CytoFLEX mosaic Spectral Detection Module]. When information in this document only applies to one type of the CytoFLEX mosaic Spectral Detection Module, it is marked either [CytoFLEX mosaic 63] or [CytoFLEX mosaic 88].

NOTE There are two kinds of plate loaders, the standard CytoFLEX Plate Loader and the CytoFLEX Plate Loader Deep Well (DW). The CytoFLEX Plate Loader DW (hereinafter called as Plate Loader DW) is compatible with CytoFLEX S and CytoFLEX LX instruments. When information in this document only applies to one plate loader, it is marked either [Plate Loader DW] or [Standard Plate Loader].

• The CytExpert for Spectral Electronic Record Management software installation screens are shown in all instances unless otherwise specified.

D17052AA XXIX

IMPORTANT IMPORTANT is used for comments that add value to the step or procedure being performed. Following the advice in the IMPORTANT adds benefit to the performance of a piece of equipment or to a process.

NOTE NOTE is used to call attention to notable information that should be followed during use, or maintenance of this equipment.

Graphics

All graphics, including screens and printouts, are for illustration purposes only and must not be used for any other purpose. For example, software screens that show the CytoFLEX for Spectral software in the background may not depict the latest production version of the system.

Advancing Your Skills

To learn more about using your CytoFLEX mosaic Spectral Detection Module system, e-Learnings are available through the Beckman Coulter Learning Center.

- Instructor-Led Coursed
 - CytoFLEX mosaic In-Person Training
 - CytoFLEX mosaic Remote Training
- E-Learning Modules
 - CytoFLEX mosaic Basic Operator Training

For further information about the available courses, visit the Beckman Coulter Learning Center at http://www.beckman.com/service/training, and enter CytoFLEX mosaic in the Search field or contact us.

XXX D17052AA

CHAPTER 1 System Overview

Overview

This chapter describes the individual components of the CytoFLEX flow cytometer, CytoFLEX mosaic Spectral Detection Module and the corresponding functions of these components.

This chapter contains information on:

- Product Description
- Main Components
- Optical Components
- Fluidics System
- Sample Station
- Plate Loader Components
- System Configuration
- Consumables and Supplies
- Instrument Specifications
- Performance Characteristics
- Reagent Limitations
- Safety Data Sheets (SDS/MSDS)

Product Description

The CytoFLEX mosaic Spectral Detection Module is an optional accessory which can be attached to the CytoFLEX S and CytoFLEX LX flow cytometers to switch the cytometers from Conventional Mode to Spectral Mode. Spectral Mode enables the spectral flow functionality and allows users to run more complex assays and use spectral unmixing technology for spectral overlap correction.

The CytoFLEX S and CytoFLEX LX flow cytometers equipped with a CytoFLEX mosaic Spectral Detection Module are used for the qualitative and quantitative measurement of biological and physical properties of cells and other particles. These properties are measured when the cells pass through one or multiple laser beams in a single-file.

For Research Use Only. Not for use in diagnostics procedures.

D17052AA 1-1


Main Components



Risk of instrument damage and/or instrument stability. Do not place any objects on top of the instrument, as this could cause warping of the top cover or affect the stability of the optical path.

The instrument consists of four main components: Fluid Containers/Cubitainers, Cytometer, the Workstation and CytoFLEX mosaic Spectral Detection Module.

Figure 1.1 Main Components [CytoFLEX S Without Plate Loader]



1. **CytoFLEX mosaic 63.** Provides signal collection with up to 63 channels using a CytoFLEX S configuration.


NOTE To switch between Spectral Mode and Conventional Mode, refer to CHAPTER 12, Switching between Spectral Mode and Conventional Mode.

NOTE Refer to the *CytoFLEX Platform Instructions for Use* manual if the instrument is being used without the CytoFLEX mosaic Spectral Detection Module connected.

- Fluid Containers. Accommodates sheath fluid and waste liquid as required for operation of the instrument.
- 3. Cytometer. Provides signal generation.
- 4. Workstation. Displays the content of the workstation and data generated by the Cytometer.

1-2 D17052AA

Figure 1.2 Main Components [CytoFLEX LX]



1. CytoFLEX mosaic 88. Provides signal collection with up to 88 channels using a CytoFLEX LX configuration.

NOTE To switch between Spectral Mode and Conventional Mode, refer to CHAPTER 12, Switching between Spectral Mode and Conventional Mode.

NOTE Refer to the *CytoFLEX Platform Instructions for Use* manual if the instrument is being used without the CytoFLEX mosaic Spectral Detection Module connected.

- 2. Cytometer. Provides signal generation.
- 3. Fluid Cubitainers. Accommodates sheath fluid and waste liquid as required for operation of the instrument.

NOTE The CytoFLEX LX does not have a fluid container holder.

4. Workstation. Displays and manipulates the contents of the Workstation and displays data generated by the Cytometer.

D17052AA 1-3

Optical Components

NOTE If you are using Conventional Mode, refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 1, Optical Components for detailed instructions.

NOTE To switch between Spectral Mode and Conventional Mode, refer to CHAPTER 12, Switching between Spectral Mode and Conventional Mode.

Wavelength Division Multiplexer (WDM)

Each WDM corresponds to a different laser, or in some cases two lasers. The WDMs are contained within the CytoFLEX mosaic Spectral Detection Module and cannot be altered to change filter configuration. Note that the CytoFLEX S equipped with a CytoFLEX mosaic 63 allows for up to 63 channel while the CytoFLEX LX equipped with a CytoFLEX mosaic 88 allows for 88 channels.

NOTE If you are using Conventional Mode, refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 1, Wavelength Division Multiplexer (WDM) for detailed instructions.

NOTE To switch between Spectral Mode and Conventional Mode, refer to CHAPTER 12, Switching between Spectral Mode and Conventional Mode.

Optical Fiber



Risk of data integrity damage.


- During use, verify that the optical fibers are securely connected between the fiber intake bracket at the back of the instrument and the back of the CytoFLEX mosaic Spectral Detection Module. A loose connection can alter the optical path and affect fluorescence detection.
- Take care when disconnecting the fibers as this could contaminate the tip and weaken the signal. Always ensure the fiber is reconnected right away.
- Do not kink the optical fibers.

NOTE If you are using Conventional Mode, refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 1, Optical Fiber for detailed instructions.

NOTE To switch between Spectral Mode and Conventional Mode, refer to CHAPTER 12, Switching between Spectral Mode and Conventional Mode.

1-4 D17052AA

Fluorescence emitted by laser-excited fluorochromes is picked up and delivered by each optical fiber to the corresponding detector module. CytoFLEX S optical fiber has a colored ring on the end that connects to the WDM, indicating the color of the corresponding laser. CytoFLEX LX optical fiber has a numbered label on the end that connects to the WDM, indicating the number of the corresponding laser. Ensure that the correct fiber is properly connected to the corresponding WDM.



- 1. Violet laser fiber
- 2. NUV or UV laser fiber
- 3. Red laser fiber

- 4. Yellow laser fiber
- 5. Blue laser fiber
- 6. Infrared laser fiber

D17052AA 1-5

Fluidics System

Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 1, Fluidics System for detailed instructions.

Fluid Containers/Cubitainers

Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 1, Fluid Containers/Cubitainers for detailed instructions.

Fluidics Module

Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 1, Fluidics Module for detailed instructions.

Sample Station

Refer to the CytoFLEX Platform Instructions for Use manual, CHAPTER 1, Sample Station for detailed instructions.

Sample Tube Holder Positions

Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 1, Sample Tube Holder Positions for detailed instructions.

Plate Loader Components

Refer to the CytoFLEX Platform Instructions for Use manual, CHAPTER 1, Plate Loader Components for detailed instructions.

Plate Holder Components

Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 1, Plate Holder Components for detailed instructions.

1-6 D17052AA


System Configuration

CAUTION

Risk of data loss and/or instrument damage. Never shut off the power or disconnect a data cable while the instrument is in the process of performing a task. This could cause data loss or damage to the system.

System Configuration [CytoFLEX S]

Figure 1.3 System Connections




- 1. Monitor
- 2. Mouse
- 3. Keyboard
- 4. Computer

- 5. Cytometer
- 6. Fluid Container holder
- 7. CytoFLEX mosaic 63


D17052AA

Figure 1.4 Power Entry Module



- 1. Power switch. Turns Cytometer on and off. An indicator light glows when the power is on.
- 2. Fuse. Protects the internal system from damage by high electrical current.
- **3. Power line socket.** Supplies the power to the Cytometer.

Figure 1.5 Front of Cytometer [CytoFLEX S Without Plate Loader Shown]



1. **Load button.** In addition to the software controls, this button can be used for automatic sample loading and data recording.


NOTE This function is not available in the Plate Loader sample injection mode.

1-8 D17052AA

System Configuration [CytoFLEX LX]

IMPORTANT The fluid cubitainers must be on the same level as the Cytometer. Do not place the fluid cubitainers on the floor.

Figure 1.6 System Connections




- 1. Monitor
- 2. Keyboard
- 3. Mouse
- 4. Computer

- 5. Cytometer
- 6. Fluid Cubitainers
- 7. CytoFLEX mosaic 88


D17052AA

Figure 1.7 Power Input Module



- 1. Power switch. Turns Cytometer on and off. An indicator light glows when the power is on.
- **2. Fuse.** Protects the internal system from damage by high electrical current.
- **3. Power line socket.** Supplies the power to the Cytometer.

Figure 1.8 Front of Cytometer




1. **Load button.** In addition to the software controls, this button can be used for automatic sample loading and data recording.

NOTE This function is not available in the Plate Loader sample injection mode.

1-10 D17052AA

System Configuration [CytoFLEX mosaic Spectral Detection Module]

Figure 1.9 Front and Left Sides of the CytoFLEX mosaic Spectral Detection Module



1. **Soft button.** In addition to the **Turn Off CytoFLEX mosaic** in the Cytometer menu on the software, this button can be used to turn the CytoFLEX mosaic Spectral Detection Module on/off. An indicator light glows when the power is on.

NOTE To fully remove the instrument from power, the power switch on the left side of the instrument must be in the off position.

- 2. Power switch. Turns the CytoFLEX mosaic Spectral Detection Module on and off. An indicator light glows when the power is on.
- **3. Power line socket.** Supplies the power to the Cytometer.




Figure 1.10 Back Cover Connections of the CytoFLEX mosaic Spectral Detection Module

- 1. Optical fiber connectors. Connects to the optical fibers.
- **2. Ethernet connector.** Connects to the Ethernet cable of the workstation.
- 3. FSC cable connector. Connects to the FSC cable.

Consumables and Supplies

Reagents

Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 1, Reagents for detailed instructions.

NOTE The CytoFLEX Daily QC Fluorospheres are not intended to be used with the CytoFLEX mosaic Spectral Detection Module.

Ordering Information

Your instrument may be upgraded to a more highly configured model. For information on specific upgrades, replacement parts, or supplies, visit:

• www.beckman.com/supplies/cytoflex-platform-upgrades

Otherwise, contact us.

1-12 D17052AA

Instrument Specifications

Dimensions [CytoFLEX S]

Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 1, Dimensions [CytoFLEX S] for detailed instructions.

Dimensions [CytoFLEX LX]

Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 1, Dimensions [CytoFLEX LX] for detailed instructions.

Dimensions [CytoFLEX mosaic Spectral Detection Module]

Dimensions						
Instrument	CytoFLEX mosaic 63	42.5 cm x 34.5 cm x 53.5 cm				
dimensions (Length x Width x Height)	CytoFLEX mosaic 88	42.5 cm x 34.5 cm x 53.5 cm				
Weight	CytoFLEX mosaic 63	30 kg				
	CytoFLEX mosaic 88	33 kg				

Installation Category

Installation Category 2

Maximum Altitude

Do not operate at an altitude greater than 2000 m (6561 ft). Beckman Coulter assumes no responsibility for any problems resulting from operating instrument at an altitude greater than 2000 m (6561 ft).

Pollution Degree

Pollution Degree 2

Acoustic Noise Level

Measure Level: <65 dBA

Electrical Ratings

Voltage: AC 100-240 V ± 10%

Frequency: 50/60 Hz

Power:

• CytoFLEX mosaic Spectral Detection Module: 450 VA

CytoFLEX S: 250 VACytoFLEX LX: 300 VA

Cytometer

Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 1, Cytometer for detailed instructions.

CytoFLEX mosaic Spectral Detection Module

Optics							
Excitation optics	The optical system is alignment free. The laser delays are automatically adjusted by the daily QC system, if required. No user intervention is required to ensure optimum system performance.						
	CytoFLEX S : The system can be configured with up to four spatially-separated lasers.						
	CytoFLEX LX : The system can be configured with up to six spatially-separated lasers.						
Emission optics	Alignment-free integrated optics quartz flow cell design with >1.3 NA.						
	Flow Cell dimensions: 430 μm x 180 μm internal diameter.						

1-14 D17052AA

Optics						
Laser devices	CytoFLEX S standard wavelengths	Violet laser • Wavelength: 405 nm, 80 mW • Beam spot size: 5 μm x 80 μm Blue laser • Wavelength: 488 nm, 50 mW • Beam spot size: 5 μm x 80 μm Yellow laser • Wavelength: 561 nm, 30 mW • Beam spot size: 5 μm x 80 μm				
		• Wavelength: 638 nm, 50 mW				
		Beam spot size: 5 μm x 80 μm				
	CytoFLEX LX standard	Ultraviolet (UV) laser				
	wavelengths	Wavelength: 355 nm, 20 mW				
		Beam spot size: 5 μm x 80 μm				
		Near Ultraviolet (NUV) laser				
		Wavelength: 375 nm, 60 mW				
		• Beam spot size: 5 μm x 80 μm				
		Violet laser				
		Wavelength: 405 nm, 80 mW				
		• Beam spot size: 5 μm x 80 μm				
		Blue laser				
		Wavelength: 488 nm, 50 mW				
		• Beam spot size: 5 μm x 80 μm				
		Yellow laser				
		Wavelength: 561 nm, 30 mW				
		• Beam spot size: 5 μm x 80 μm				
		Red laser				
		Wavelength: 638 nm, 50 mW				
		Beam spot size: 5 μm x 80 μm				
		Infrared (IR) laser				
		Wavelength: 808 nm, 60 mW				
_ , ,, ,, ,,	G	Beam spot size: 5 μm x 80 μm				
Forward scatter detection		uilt-in 488/8 band-pass filter.				
Fluorescence and side scatter detection	Fluorescence and side scatter light collected by the objective delivered by fiber optics to a patent-pending design with high performance, solid-state, high efficiency, low-noise detector a					
	Reflective optics with a sindetector.	ngle transmission band-pass filter in front of each				
Side scatter (SSC) configuration	The system offers the abit to better resolve particles	ility to collect side scatter lights from each laser s from noise.				

Electronics				
Signal processing 7 decade data display				
Digital sampling rate 6.25 MHz				
Signal	Pulse area and height for every channel, width for two selectable channels			

Data Management					
Software	CytExpert for Spectral	CytExpert for Spectral software			
Language	English				
FCS format	FCS 3.1 and FCS 3.0				
Operating system	Windows 10 Enterprise LTSC 2019 64 bit				
Minimum	Processor	Intel Core i7-11700K Processor, 3.6GHz			
workstation/ computer	Memory	64 GB RAM			
requirements	Storage	1 TB SSD			
[CytoFLEX S and	Ethernet	2 Ethernet interface			
CytoFLEX LX]	USB	4 USB 3.0 ports			
	Monitor resolution	tion 1920 x 1080 or 2560 x 1440			

Performance Characteristics

Performance Characteristics [CytoFLEX S]

Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 1, Performance Characteristics [CytoFLEX S] for detailed instructions.

Performance Characteristics [CytoFLEX LX]

Refer to the CytoFLEX Platform Instructions for Use manual, CHAPTER 1, Performance Characteristics [CytoFLEX LX] for detailed instructions.

1-16 D17052AA

Performance Characteristics [CytoFLEX mosaic Spectral Detection Module]

	P	erformance				
Sensitivity	MESF	FITC: <30 molecules of equivalent soluble fluorochrome (MESF-FITC)				
		PE: <10 molecules of equivalent soluble fluorochrome (MESF-PE)				
		PB450: <50 molecules of equivalent soluble fluorochrome (MESF-PB450)				
		APC: <25 molecules of equivalent soluble fluorochrome (MESF-APC)				
Fluorescence	Detection range	Ultraviolet (UV) laser:				
detection		Detectors: 20				
configuration		Wavelength: 365-950 nm				
		Near Ultraviolet (NUV) laser				
		Detectors: 18				
		Wavelength: 420-950 nm				
		Violet laser:				
		Detectors: 20				
		Wavelength: 420-950 nm				
		Blue laser:				
		Detectors: 16				
		Wavelength: 498-950 nm				
		Yellow laser:				
		Detectors: 12				
		Wavelength: 567-950 nm				
		Red laser:				
		Detectors: 10				
		Wavelength: 649-950 nm				
		Infrared (IR) laser:				
		Detectors: 3				
		Wavelength: 830-950 nm				
Fluorescence detection	Up to 79 fluorescent detectors are designed to cover a broad wavelength spectrum					
Fluorescence	rCV <3%	rCV <3%				
resolution	Daily QC Fluorospheres (The CytoFLEX Flow cytometer is capable of achieving <3% rCV. Using CytoFLEX Daily QC Fluorospheres, CytoFLEX Ready to Use Daily QC Fluorospheres or Daily IR QC Fluorospheres (for 808 nm Laser) for daily QC, the pass criteria is \leq 5% for the Violet, Blue, Yellow, and Red lasers while the pass criteria is \leq 7% for NUV, UV and IR lasers.				

Performance						
Blue/Yellow/Red/ Infrared/Ultraviolet/ Near Ultraviolet side scatter resolution	<300 nm					
Violet side scatter resolution	80 nm relative to polystyrene particles					
Forward and blue side scatter resolution	•	Scatter performance is optimized for resolving lymphocytes, monocytes, and granulocytes as well as nanoparticles.				
Carryover	Single Tube Loader format ≤1.0%					
Signal acquisition speed	30,000 events/second					

Reagent Limitations

Refer to the CytoFLEX Platform Instructions for Use manual, CHAPTER 1, Reagent Limitations for detailed instructions.

Safety Data Sheets (SDS/MSDS)

To obtain an SDS or MSDS for CytoFLEX Platform reagents used on the CytoFLEX Platform systems:

- 1. On the Internet, go to www.beckman.com:
 - a. Select Safety Data Sheets (SDS/MSDS) from the Support menu.
 - **b.** Follow the instructions on the screen.
 - **c.** Contact us if you have difficulty locating the information.
- **2.** If you do not have Internet access, contact us.

1-18 D17052AA

Using the CytExpert for Spectral Software

Overview

The CytExpert for Spectral software is a full-feature software package that controls the instrument operation, collection of experiment data, and analysis of the results. This chapter will explain the software functions and features.

NOTE The Cytexpert for Spectral software only applies to the CytoFLEX S and CytoFLEX LX instruments equipped with the CytoFLEX mosaic Spectral Detection Module. If your CytoFLEX instrument is not equipped with a CytoFLEX mosaic Spectral Detection Module or if you are using Conventional Mode, please use the CytExpert software. Refer to the CytoFLEX Platform Instructions for Use manual for detailed instructions on using the CytExpert software.

IMPORTANT The CytExpert software and the CytExpert for Spectral software can be installed on the same computer, but you cannot use the CytExpert for Spectral software and the CytExpert software together.

This chapter contains information on:

- Launching the Software
- Main Software Screen
- User Management
- Role Management
- **Account Policies**
- **User Management Operation Log**
- Graphic and Gating Styles
- Software Settings

Launching the Software



Select the desktop shortcut to launch the CytExpert for Spectral software.

If there is no desktop shortcut, run the "CytExpert.exe" software directly from the software installation directory. The default installation path is C:/Program Files/CytExpert for Spectral. Or,

> All Programs > CytExpert for Spectral.


Refer to CHAPTER 4, Logging into the Software, for detailed instructions on opening the software and confirming the connection status.

Main Software Screen

Hover your cursor over any button to display a text pop-up of the button's function.


Start Page

The login window automatically opens after the software has been launched. The start page automatically opens after logging into the software.



2-2 D17052AA

The following operations can be selected from the start page:




- **New Experiment.** For creating a new experiment. The process creates a file with the .xit extension and a folder with the same file name where the data (.fcs files) are kept.
- **New Experiment from Template.** For creating an experiment using a template saved from a previously saved experiment.
- **New Spectral Unmixing.** For creating a new spectral unmixing experiment. The process creates a file with the .xitsp extension and a folder with the same file name where the data (.fcs files) are kept.
- **Open Experiment.** For opening a previously created experiment.
- Open Spectral Unmixing. For opening a previously created spectral unmixing experiment.
- **Exit.** For exiting CytExpert for Spectral.

The Experiment, Template, and Spectral Unmixing tabs below give you the option of opening one of the 10 most recently opened experiments.

Acquisition Screen

Selecting New Experiment, New Experiment from Template, or Open Experiment automatically opens the Acquisition screen. The Acquisition screen can be accessed by selecting upon the left side of the page.




- 1. Navigation. Gives the option of accessing the Acquisition screen or Analysis screen.
- **2. Menu.** Allows you to configure settings for sample acquisition, instrument operation, and software options.
- **3. Instrument operation controls.** Controls sample loading/unloading and data acquisition and recording.
- **4. Collection.** Establishes control over data recording options, displays the acquisition status, and controls the sample flow rate.
- **5. Test tubes.** Allows you to configure and duplicate sample tubes, set display attributes, manage experimental data and unmixing.
 - **NOTE** The Tube section of the screen can be expanded or retracted by dragging the top border of the Tube section of the screen. Expanding this section covers other elements of the screen, including: Events to Display, Events/Sec, and the Acquisition buttons.
- **6. Plot area.** Includes plot and gating controls, as well as an area for creating plots and generating graphs. There are two worksheets display in this area, Raw worksheet and Unmixed worksheet.
- **7. Status bar.** Displays instrument connection status and system information.

2-4 D17052AA

Acquisition Screen Navigation

The Acquisition screens have two navigation icons, one for the Acquisition screen and the other for the Analysis screen.




- 1. Acquisition screen icon. Accesses the Acquisition screen.
- **2. Analysis screen icon.** Accesses the Analysis screen.

Collection

Standby state Acquisition Restart Initialize Record (Standby N Backflush Boost ♦ Next Tube * Acq. Setting... Events/Sec: 0.0 Abort(%): 0.00 Events: 00:00:00 Time: Events to Display: 1000 Events ▼ Events to Record: Events 10000 3 All Events Sec ▼ Time to Record: 600 ■ Volume to Record: 10 μL Sample Flow Rate: 10 µL/min Slow Medium Fast Custom

Initialized state



- 1. Acquisition control. Controls sample loading/unloading and data acquisition and recording.
- 2. Acquisition status. Displays such information as the acquisition rate (Events/Sec), event count, duration, and abort (%).
- **3.** Acquisition conditions. Sets the necessary conditions for recording data.
- 4. Sample flow rate. Sets the acquisition rate for data collection.

NOTE High acquisition rate may increase the abort rate and measurement CVs.


Custom: The flow rate can be adjusted in 1 μ L increments.

2-6 D17052AA


Collection [With Plate Loader]

Ct.





Initialized state



- 1. Acquisition control. Controls sample loading/unloading and data acquisition and recording.
- **2. Acquisition status.** Displays such information as the acquisition rate (Events/Sec), event count, duration, and abort (%).
- **3.** Acquisition conditions. Sets the necessary conditions for recording data.
- **4. Sample flow rate.** Sets the acquisition rate for data collection.

NOTE High acquisition rate may increase the abort rate and measurement CVs.

Test Tubes

Manual/Semi-Automatic Sample Injection Mode





Plate Loader Sample Injection Mode



Test Tube Status



- **1. Tube management controls.** Manages sample tubes. Used to add, copy, or delete attributes, and open the tube property.
- **2. Test tube status indication.** Displays a colored symbol in front of each tube indicating the status of the tube processing.
 - O indicates that the tube data was not acquired.
 - indicates that the tube data was acquired by selecting **Run**, and the tube data can be overwritten.
 - indicates that the tube data was acquired by selecting **Record** or **Auto Record**, and the tube data cannot be overwritten.
 - **(F)** indicates that the tube data was imported from FCS file.
 - **(F)** indicates that the tube data was not acquired, and the unmixed data was imported from FCS file.
 - 10 indicates that the tube data imported from FCS file is damaged.
 - ① indicates that the tube data was not acquired, and the unmixed data imported from FCS file is damaged.
 - 📵 indicates that tube data was not acquired, and was linked to the Acquisition Catalog.
 - indicates that the tube data was acquired by selecting **Run**, and was linked to the Acquisition Catalog.
 - indicates that the tube data was acquired by selecting **Record** or **Auto Record**, and was linked to the Acquisition Catalog.
 - **G** indicates that the tube data was acquired by selecting **Run**, and the tube data was saturated.

2-8 D17052AA

• **(a)** indicates that the tube data was acquired by selecting **Record** or **Auto Record**, and the tube data was saturated.

NOTE In spectral unmixing experiments, \triangle to the right of the test tube status indication symbol indicates that there is an unmixing gap.

NOTE Icons to the left of the test tube status indication symbol indicates the test tube's unmixing and fine adjustment status.

- Indicates that the tube data does not contain the unmixing matrix.
- indicates that the unmixing matrix was not applied to the tube.
- indicates that the unmixing matrix and fine adjustment matrix were not applied to the tube.
- indicates that the unmixing matrix was applied to the tube, but the unmixing matrix is missing.
- indicates that the unmixing matrix and fine adjustment matrix were applied to the tube, but the unmixing matrix is missing.
- indicates that the unmixing matrix was applied to the tube.
- // indicates that the unmixing matrix and fine adjustment matrix were applied to the tube.
- **3. Test tube list.** Displays the sample tubes used in the experiment. Right-click a tube in the list to perform additional operations.

NOTE In the Plate Loader Sample Injection mode the well number displays at the end of the tube name.

Test Tube Management Module

NOTE The is only available if the CytExpert for Spectral User Management software option is installed.

Select to open a previously created experiment.

New Experiment Test Tube Management Module [Manual/Semi-Automatic Sample Injection Mode]



New Experiment Test Tube Management Module [Plate Loader Sample Injection Mode]



New Spectral Unmixing Experiment Test Tube Management Module [Manual/Semi-Automatic Sample Injection Mode]




New Spectral Unmixing Experiment Test Tube Management Module [Plate Loader Sample Injection Mode]



Plot Area

The CytExpert for Spectral software supports two types of worksheets, Raw worksheet and Unmixed worksheet. Select the **Raw** or **Unmixed** tab in the bottom of the screen to switch between worksheets. Refer to Figure 2.1 and Figure 2.2.

Figure 2.1 Plot Area [Raw Worksheet Shown]



- 1. **Plot controls.** For creating single or multiple plots, such as dot plots, histograms, density plots, pseudo color plots, contour plots and spectrum plots.
- 2. Statistics and hierarchy controls. For creating statistical and hierarchical charts.
- **3. Graphical gating controls.** For creating graphical gates.
- **4. Zoom controls.** For zooming in and out within a plot.
- **5. Pan axis display controls.** For scaling axis ranges in the plots.
- **6. Gain adjustment control.** For increasing and lowering gain adjustments on the plots.


NOTE The gain adjustment control only works when a sample is running.

7. Threshold control. For setting the minimum particle size limit, scatter value, or fluorescence intensity that acquisition will allow.

2-10 D17052AA

- 8. Undo and redo controls. For undoing or redoing an action in the drawing area.
- **9. Display controls.** For controlling how plots and tables are aligned and arranged.
- **10.** Rearrange. For restoring the plots to the default positions.
- **11. Printing controls.** For printing and previewing the plot area.
- **12. Plot area.** For creating plots and displaying statistics and hierarchy tables.
- **13.** Raw worksheet. For editing and displaying raw data, and exporting or printing the results.
- **14. Unmixed worksheet.** For editing and displaying unmixed data, and exporting or printing the results.





- 1. **Plot controls.** For creating single or multiple plots, such as dot plots, histograms, density plots, pseudo color plots, and contour plots.
- 2. Statistics and hierarchy controls. For creating statistical and hierarchical charts.
- **3. Graphical gating controls.** For creating graphical gates.
- **4. Zoom controls.** For zooming in and out within a plot.
- **5. Pan axis display controls.** For scaling axis ranges in the plots.
- **6. Auto scaling control.** For scaling axis ranges in the plots automatically.
- **7. Fine Adjustment control.** For using fine adjustment functions.
- 8. Undo and redo controls. For undoing or redoing an action in the drawing area.
- **9. Display controls.** For controlling how plots and tables are aligned and arranged.
- **10.** Rearrange. For restoring the plots to the default positions.
- **11. Printing controls.** For printing and previewing the plot area.

- **12. Plot area.** For creating plots and displaying statistics and hierarchy tables.
- 13. Raw worksheet. For editing and displaying raw data, and exporting or printing the results.
- **14. Unmixed worksheet.** For editing and displaying unmixed data, and exporting or printing the results.

Status Bar



- **1. Communication connection status.** Displays whether the Cytometer and the Workstation are connected.
- **2. Instrument status information.** Displays the status of the Cytometer.
- **3. CytoFLEX mosaic status information.** Displays the status of the CytoFLEX mosaic Spectral Detection Module.
- **4. Sampler status.** Displays the sample injection mode state. There are two sample injection modes: Semi-automatic sample injection mode and manual sample injection mode.
 - **NOTE** CytoFLEX Cytometers equipped with a plate loader have three sample injection modes: Semi-automatic sample injection mode, manual sample injection mode, and plate loader sample injection mode.


5. Fluid status information. Displays the liquid level of the Fluid Containers/Cubitainers.

2-12 D17052AA

Analysis Screen

The Analysis screen is similar to the Acquisition screen, without the acquisition control modules.

NOTE Only the unmixed data can be analyzed in the analysis screen.



The Tube management module cannot add new sample tubes. Return to the Acquisition screen to add new sample tubes.


Tube Management



NOTE The is only available if you have the User Management software option installed. Select open a previously created experiment.


Drawing controls (see Figure 2.3) include the multi-data histograms and graphical display data controls.

Figure 2.3 Drawing Controls Toolbar (Top of Screen)



Spectral Unmixing Experiment Screen

The Spectral Unmixing Experiment screen appears when you open or create a new Spectral Unmixing experiment.




- 1. **Tube management.** Displays sample tubes required for the spectral unmixing experiment.
 - **NOTE** The Tube management section of the screen can import saved data (.fcs) files for computational purposes.
- 2. Plot area. Displays plots and spectrum plot used for preview, and unmixing plots and gating.

NOTE When the unmixing matrix is generated and applied, unmixing plots appear.

2-14 D17052AA

Spectral Unmixing Controls

The control area includes the spectral unmixing controls, coordinate pan axis display controls, gain adjustment controls, fine adjustment controls, and the undo and redo controls. The spectral unmixing controls gives you the option of calculating the spectral unmixing value, displaying the spectral unmixing matrix, or changing the spectral unmixing parameters.




QC Experiment Screen

The Quality Control (QC) Experiment screen appears when you access a QC experiment.

QC Report Screen

Before starting the QC routine, a Settings screen appears.

Figure 2.4 QC Report Screen




- 1. Menu. Allows you to configure settings related to QC experiments.
- 2. Acquisition control. Controls sample loading/unloading and data recording.
- 3. Lot selection. Allows you to select the lot number of the QC reagent.
- **4. QC results list area.** Displays the time and results of completed QC runs.
- 5. QC reports area. Displays detailed reports for the selected QC experiment.

2-16 D17052AA

QC Experiment Screen


When acquiring QC samples, the software opens the QC screen.



- 1. QC experiment progress indicator. Displays the QC stage.
- 2. Plot area. Displays the QC plots.

QC Screen Navigation

The Analysis screens have two navigation icons, one for the QC screen and the other for the Levey-Jennings (LJ) charts. Refer to CHAPTER 5, Creating Levey-Jennings Charts.



- 1. QC screen icon. Accesses the QC screen.
- **2. LJ screen icon.** Accesses the Levey-Jennings (LJ) screen.

Software Menu

The CytExpert for Spectral software contains the following selectable menu items:

Figure 2.5 Software Menu Tree*

File	Cytometer		Settings	QC/Standardization	Advanced	Account	Log	Signature ‡‡	Backup/Restore	
				Start			Experiment			View Help
New Experiment	Acq. Setting		Set Label ††	QC/Standardization	Delay Setting	User Manager	Operation Log ‡‡	Sign	Backup	File
New Experiment from Template			Set Customized Parameter ††		Laser Setting	Role Manager	System Operation Log ‡‡	Reject	Restore	Third Party Notice
New Spectral Unmixing	Boost †		Set Retention Period ‡‡		Maintenance	Account Policies	User Management Operation Log	Signature Records	Log Clean-up ##	About
Open Experiment	Initialize		Retention Period Options ‡‡		Event Rate Setting	Change Password		Signature Setting	2	
Open Spectral Unmixing	Standby		Spectral Unmixing Setup ***		Plate Type Library ‡	Export Credential File				
Save	Prime		Spectral Unmixing		Restore Configuration					
Save As	Deep Clean		Unmixing Matrix		Saturation Warning]				
Save As Template ††	Calibrate Sample Flow Rate		Spectral Library							
Import FCS File ††	System Startup Program		Fluorescent Tag Library							
Export FCS File (Raw) ††	Daily Clean		Events Display Setting							
Export FCS File (Unmixed) ††	Sample	Manual	Set Experiment Directory ‡‡							
Recent	Injection Mode	Semi Automatic	Options							
Recent Template		Plate Loader **								
Recent Spectral Unmixing	Sampler Reset									
Close Experiment	Turn On †††									
Experiment Explorer ‡‡	Turn Off †††									
Exit	Acq. Setting Cat	alog								
	Cytometer Confi	guration								
	Cytometer Information		1							

^{*} The menu options for **File**, **Cytometer**, **Settings**, and **QC/Standardization** change when you select Start QC/Standardization. Refer to Figure 2.6.

- **† Boost** is only active in the Manual Sample Injection mode.
- **‡ Plate Type Library** is only an option if the Plate Loader module is installed and the Plate Loader Sample Injection mode is selected.
- ** Plate Loader is only an option if the Plate Loader module is installed.
- *** Only available when the spectral unmixing experiment is open.
- †† Only available when a standard experiment is open.
- ‡‡ These options are only available if the CytExpert for Spectral Electronic Record Management software option is installed.
- ††† These options are only available on the CytoFLEX LX flow cytometer.

The Cytobank icon is hyper linked to the Cytobank spotlight page where you can login to the Cytobank Premium server, request a 30-day free trial and access additional information. The Cytobank platform allows you to analyze, manage, and securely share flow cytometry data. FCS files can be uploaded to the platform with the related attachments such as PDF files.

2-18 D17052AA

Figure 2.6 QC Software Menu Tree


File	Cytometer		Settings	Advanced	Account	Log	Signature ‡‡	Backup/Restore	Help
			Retention Period			Experiment Operation			View Help
New Experiment	Backflush		Options ##	Delay Setting	User Manager	Log ‡‡	Sign	Backup	File
									Third
New Experiment from			Fluorescent Tag		3	System Operation Log			Party
Template	Boost †		Library	Section 2 to the second section 2	Role Manager	##	Reject	Restore	Notice
			QC/Standardization			User Management			
New Spectral Unmixing	Initialize		Setting	Maintenance	Account Policies	Operation Log	Signature Records	Log Clean-up ##	About
20 120 10 100	10.00		28 6602	Event Rate	a a a		20 10 15 10 1		
Open Experiment	Standby		Target Library	Setting	Change Password		Signature Setting	l	
	**********		Standardization	Plate Type	Export Credential				
Open Spectral Unmixing	Prime		Target Library	Library ‡	File				
			Set Experiment	Restore					
Recent	Deep Clean		Directory ##	Configuration					
	Calibrate Sample Flow			Saturation					
Recent Template	Rate		Options	Warning]				
Recent Spectral									
Unmixing	System Start	tup Program							
Close									
QC/Standardization	Daily Clean								
Experiment Explorer ##	-	Manual							
	Sample	Semi							
Exit	Injection	Automatic							
	Mode	Plate							
		Loader **							
	Sampler Reset								
	Turn On †††								
	Turn Off †††		1						
	Acq. Setting Catalog		ľ						
	Cytometer Configuration		1						
	Cytometer Information		J						

- **† Boost** is only active in the Manual Sample Injection mode.
- **‡ Plate Type Library** is only an option if the Plate Loader module is installed and the Plate Loader Sample Injection mode is selected.
- ** Plate Loader is only an option if the Plate Loader module is installed.
- # These options are only available if the CytExpert for Spectral Electronic Record Management software option is installed.
- **†††** These options are only available on the CytoFLEX LX flow cytometer.

The Cytobank icon is hyper linked to the Cytobank spotlight page where you can login to the Cytobank Premium server, request a 30-day free trial and access additional information. The Cytobank platform allows you to analyze, manage, and securely share flow cytometry data. FCS files can be uploaded to the platform with the related attachments such as PDF files.

Acquisition and Analysis Screen Menu


CytExpert for Spectral - User Management




File Menu

For creating new experiments, opening existing experiments, saving new experiments and data, and importing/exporting FCS data files.

User Management



Electronic Record Management




2-20 D17052AA

Cytometer Menu

For configuring Cytometer settings and controlling Cytometer functions. Depending on the Cytometer state, certain functions may not be available.

NOTE The **Turn On, Turn Off, Turn On CytoFLEX mosaic** and **Turn Off CytoFLEX** selections are only available on the CytoFLEX LX.


Cytometer On and CytoFLEX mosaic Spectral Detection Module On




Cytometer On and CytoFLEX mosaic Spectral Detection Module Off




Cytometer On and CytoFLEX mosaic Spectral Detection Module Power Off




NOTE If the software restarts, the **Turn On CytoFLEX mosaic** selection is not available when the Cytometer is on and the CytoFLEX mosaic Spectral Detection Module is off.




Cytometer Off and CytoFLEX mosaic Spectral Detection Module On



Cytometer Off and CytoFLEX mosaic Spectral Detection Module Off




Cytometer Off and CytoFLEX mosaic Spectral Detection Module Power Off




2-22 D17052AA


Cytometer Power Off and CytoFLEX mosaic Spectral Detection Module On




Cytometer Power Off and CytoFLEX mosaic Spectral Detection Module Off



Cytometer Power Off and CytoFLEX mosaic Spectral Detection Module Power Off




NOTE If the software restarts, the **Turn On CytoFLEX mosaic** selection is not available when the Cytometer is power off and the CytoFLEX mosaic Spectral Detection Module is turned off.




Settings Menu

Used to select and/or change software options and settings.

User Management



Electronic Record Management



QC/Standardization Menu

Select **Start QC/Standardization** from the QC/Standardization menu to start the QC routine.




NOTE The QC/Standardization menu is the same for the CytExpert for Spectral User Management and the CytExpert for Spectral Electronic Record Management software options.

2-24 D17052AA


Advanced Menu

Used to access advanced settings for experienced users. Includes laser time delay settings.


User Management - Semi-Automatic/Manual Sample Injection Mode




Electronic Record Management -Semi-Automatic/Manual Sample Injection Mode



User Management - Plate Loader Sample Injection Mode




Electronic Record Management - Plate Loader Sample Injection Mode



NOTE To restore the Cytometer and CytoFLEX mosaic Spectral Detection Module configuration file, refer to APPENDIX A, Restoring the Cytometer and CytoFLEX mosaic Spectral Detection Module Configuration File.

Account Menu

Used to for user account management settings.




NOTE The Account menu is the same for the CytExpert for Spectral User Management and the CytExpert for Spectral Electronic Record Management software options.


Log Menu

Used to access logs including the Experiment Operation Log, the System Operation Log, and the User Management Operation Log.

User Management




Electronic Record Management



Signature Menu

Used to sign experiment and view signature details.



NOTE The Signature menu is only available in the CytExpert for Spectral Electronic Record Management software option.


Backup/Restore Menu

Used to backup/restore databases. Refer to CHAPTER 10, Backup and Restore.

User Management




Electronic Record Management



Help Menu

For displaying software version information and system Instructions for Use.



NOTE The Help menu is the same for the CytExpert for Spectral User Management and CytExpert for Spectral Electronic Record Management software options.

2-26 D17052AA

User Management

IMPORTANT Only the Initial Administrator or authorized users can manage users.

User Management is used to create and manage user accounts.

Select **Account > User Manager**. The User Manager window appears.

Figure 2.7 User Manager (Card View)





Figure 2.8 User Manager (Grid View)



1. Search text box. Filters users by username and display name.

- 2. View drop-down. Toggles between Card View (see Figure 2.7) and Grid View (see Figure 2.8).
- **3.** New. Used to create a new user profile.
- **4.** Modify. Used to modify an existing user profile.
- **5. Delete.** Used to delete an existing user profile.
- **6. Unlock.** Used to unlock an existing account that has been locked.

NOTE An account locks after 3 failed password attempts. The number of attempts can be changed by the administrator. Refer to Account Policies.

NOTE An account automatically unlocks after 30 minutes. The duration can be changed by the administrator. Refer to Account Policies.


- 7. Reset password. Used to reset an existing user password to a random password.
- **8.** Close. Closes the User Manager window.

Creating, Deleting, and Modifying Users in User Manager

The initial system administrator indicating with \boxed{b} in the username cannot be deleted or modified.

Creating a New User in User Manager

Select **New** in the User Manager window. The New window appears.




- **2** Fill in the new user information.
 - a. Enter the Username.
 - **b.** Enter the Full Name.
 - c. Select the User Role.
 - **d.** Select the Enabled checkbox to enable the user.

NOTE The Enabled checkbox can only be changed by an administrator.

2-28 D17052AA

3 Select **OK**. The Create User Result window appears.



NOTE Select Copy Username and Password to Clipboard to copy the username and password to a clipboard and inform the user of the change. The user is required to change the new password immediately upon the initial login.

- 4 Select ok.
- 5 Select Close.

Deleting Users in User Manager


IMPORTANT If an account has been used and log information has been generated related to it, the account cannot be deleted, but it can be disabled.

- 1 Select the user to be deleted in the User Manager window then select **Delete**.
- **2** The following message appears: *Are you sure to you want to delete the user?* Select **Yes** to confirm.
- 3 Select Close to exit the User Manager window.

Modifying Users in User Manager

IMPORTANT If an account has been used and log information has been generated related to it, the username cannot be modified.

1 Select **Modify** in the User Manager window. The Modify window appears.



NOTE The Initial System Administrator is a system default user and you can only modify its Full Name.

2 Modify the user information as necessary.

NOTE Uncheck the enabled box to disable a user.

- 3 Select **OK**.
- **4** Select **Close** to exit the User Manager window.

Unlocking a User Account

Select a Locked user in the User Manager window and select **Unlock**.


NOTE You cannot unlock an active user.

2-30 D17052AA

Resetting a User Password

IMPORTANT Only the administrator or authorized users can reset a password.


1 Select the user to be reset in the User Manager window and select Reset Password.



The confirm window appears.



2 Select **Yes**. The Reset Password Result window appears. A new random password is automatically generated.




NOTE Select Copy Username and Password to Clipboard to copy the username and password to a clipboard and inform the user of the change. The user is required to change the new password immediately upon the initial login.

- 3 Select **OK**.
- **4** Select **Close** to exit the User Manager window.

Changing a User Password

Beckman Coulter recommends changing your password on a regular basis.

1 Select **Account > Change Password**. The Change Password window appears.



2 Enter the current password, the new password, and confirm the new password.

2-32 D17052AA

3 Select **οκ**.

Forgot Username or Password

The Forgot username or password function allows only the initial system administrator to reset a password using the system administrator credential file. If you are a common user and forget your password, contact the administrator to reset a password for you. For instructions, refer to Resetting a User Password.

IMPORTANT The initial system administrator password cannot be reset without the system administrator credential file. Beckman Coulter is not responsible for and will not be able to recover your initial system administrator account if the system administrator password is forgotten and the credential file is lost. For instructions on exporting the credential file, refer to Exporting the Credential File.

1 Select Forgot username or password on the login window.




The Reset Password window appears.



IMPORTANT Find the system administrator credential file from the backup folder when installing the CytExpert or Spectral software. You can also export the credential file by selecting **Export Credential File** from the **Account** menu if you are logged in. Refer to Exporting the Credential File.

2 Navigate to the system administrator credential file and select **Open**.



3 Select **OK**. The Change Password window appears.



2-34 D17052AA

4 Enter the new password, and confirm the new password.

NOTE The new password must contain at least ten digits and at least three of the following character types by default:


- Uppercase characters of European languages (A through Z, with diacritic marks, Greek and Cyrillic characters);
- Lowercase characters of European languages (a through z, sharp-s, with diacritic marks, Greek and Cyrillic characters);
- Base 10 digits (0 through 9);
- Non-alphabetic characters (for example, !, \$, #, %);
- Any Unicode character recognized as an alphabetic character but not uppercase or lowercase.
- **5** Select **οκ**.

Exporting the Credential File

1 Select Account > Export Credential File. The Export Credential File window appears.



- **2** Select and browse to the desired backup directory to export the credential file.
- **3** Select **OK**. The following system prompt appears.



4

Select **OK**.

Role Management

Refer to the CytoFLEX Platform Instructions for Use manual, CHAPTER 2, Role Management.

Creating, Deleting, and Modifying User Roles in Role Manager

Refer to the CytoFLEX Platform Instructions for Use manual, CHAPTER 2, Creating, Deleting, and Modifying User Roles in Role Manager.

Account Policies

IMPORTANT Only an Administrator or an account has the Access account policies permission can manage account policies.

Account policies is used to define the default properties for the password policy, account lockout policy, and application inactivity policy.

Select **Account > Account Policies**. The Account Policies window appears.

2-36 D17052AA

Figure 2.9 Account Policies - Password Policy



NOTE The allowable range for each entry is as follows:

Password History: 12-24 times


Password Length: 10-12 characters

Minimum Age for Password: 0-89 days

Password Expiration: 1-90 days

Reminder for Expiration: 1-90 days

Figure 2.10 Account Policies - Account Lockout Policy




NOTE The allowable range for each entry is as follows:

• Invalid Login Attempts: 3-5 times

• Lockout Duration: 30-1440 minutes

Figure 2.11 Account Policies - Application Inactivity Policies



NOTE The allowable range for each entry is as follows:

Inactivity Duration: 1-15 minutes

2-38 D17052AA

User Management Operation Log

Refer to the CytoFLEX Platform Instructions for Use manual, CHAPTER 2, User Management Operation Log.

Graphic and Gating Styles


Plots

The CytExpert for Spectral software offers a variety of plot formats including:


- Single-parameter plots and histogram overlays
- Dual-parameter plots: dot plots, density plots, pseudo color plots, contour plots, and dot plot overlays
- Spectrum plots

NOTE Histogram Overlays and Dot Plot Overlays can only be created from multiple samples in the Analysis screen. A maximum of 10 samples can be overlaid.

NOTE Spectrum option is only an option in the Raw worksheet.



D17052AA



Gates

Refer to the CytoFLEX Platform Instructions for Use manual, CHAPTER 2, Gates.

NOTE Gating is not available in spectrum plots.

Plate Type Library

Refer to the CytoFLEX Platform Instructions for Use manual, CHAPTER 2, Plate Type Library.

Adding a Plate Type

Refer to the CytoFLEX Platform Instructions for Use manual, CHAPTER 2, Adding a Plate Type.

Editing a Plate Type

Refer to the CytoFLEX Platform Instructions for Use manual, CHAPTER 2, Editing a Plate Type.

Duplicating a Plate Type

Refer to the CytoFLEX Platform Instructions for Use manual, CHAPTER 2, Duplicating a Plate Type.

2-40 D17052AA

Deleting a Plate Type


Refer to the CytoFLEX Platform Instructions for Use manual, CHAPTER 2, Deleting a Plate Type.


Fluorescent Tag Library

IMPORTANT The maximum number of fluorescent tags allowed in the Fluorescent Tag Library is 500.

The Fluorescent Tag Library is used to manage fluorescent tags. Fluorescent tag can be added, deleted, and edited from the Fluorescent Tag Library.

Select **Settings** > **Fluorescent Tag Library** to access the Fluorescent Tag Library. Refer to Figure 2.12.





- **1. Search text box.** Filters fluorescent tag by tag name and peak channel.
- **2.** Add. Creates a new fluorescent tag.
- **3. Delete.** Deletes an existing fluorescent tag.
- **4. Default.** Applies the default settings.
- **5. Export.** Exports fluorescent tag files.

- **6. Import.** Imports fluorescent tag files.
- **7. OK.** Saves the fluorescent tag.
- **8.** Cancel. Cancels the settings.

Adding a Fluorescent Tag


1 Select **Settings** > **Fluorescent Tag Library**. The Fluorescent Tag Library window appears. Refer to Figure 2.12.

IMPORTANT The maximum number of fluorescent tags allowed in the Fluorescent Tag Library is 500.

- 2 Select Add. Enter the fluorescent tag name and select the peak channel from the Peak Channel drop-down menu.
- **3** Select **OK** to exit the Fluorescent Tag Library window.

Deleting a Fluorescent Tag

- 1 Select **Settings** > **Fluorescent Tag Library**. The Fluorescent Tag Library window appears. Refer to Figure 2.12.
- 2 Select the fluorescent tag to be deleted, and select **Delete**. The following message appears.




NOTE The default fluorescent tags, shown in bold, cannot be deleted from the Fluorescent Tag Library.

- 3 Select Yes.
- 4 Select **OK** to exit the Fluorescent Tag Library window.

2-42 D17052AA

Restoring the Default Settings

- 1 Select **Settings** > **Fluorescent Tag Library**. The Fluorescent Tag Library window appears. Refer to Figure 2.12.
- **2** Select **Default**. The following message appears.



- 3 Select Yes.
- 4 Select **OK** to exit the Fluorescent Tag Library window.

Exporting Fluorescent Tag Files

- 1 Select **Settings** > **Fluorescent Tag Library**. The Fluorescent Tag Library window appears. Refer to Figure 2.12.
- 2 Select **Export** to specify a path and filename for the fluorescent tag file you are saving.
- 3 Select Save.

NOTE The generated file ends in .ft.

4 Select Close.

Importing Fluorescent Tag Files


- 1 Select **Settings** > **Fluorescent Tag Library**. The Fluorescent Tag Library window appears. Refer to Figure 2.12.
- 2 Select **Import** and locate the path where the fluorescent tag file is saved. Select the corresponding fluorescent tag file (.ft) to import.

NOTE If the fluorescent tag exists in the fluorescent tag library, the importing will overwrite the existing fluorescent tag. The system will prompt to ask you to confirm.

IMPORTANT The import failures could occur if:

- The number of fluorescent tags exists in the fluorescent tag library and the imported fluorescent tag file exceeds 500.
- The integrity of the imported fluorescent tag file is compromised.
- 3 Select Open to import.

NOTE The fluorescent channels imported that are not supported by the software are shown in gray.



4 Select Close.

2-44 D17052AA

Editing a Fluorescent Tag

- 1 Select **Settings** > **Fluorescent Tag Library**. The Fluorescent Tag Library window appears. Refer to Figure 2.12.
- 2 Select the fluorescent tag to be edited and select the desired peak channel from the Peak Channel drop-down menu.

NOTE The default fluorescent tags' name, shown in bold, cannot be edited. Select **Default** to restore the default settings. Refer to Restoring the Default Settings.


3 Select Close.


Software Settings

Select **Settings** > **Options** to configure the software settings.

In the **Experiment** settings, you can set the experiment's default save path.


NOTE The Experiment setting is only available if the CytExpert User Management software option is installed.







In the **Tube** settings, you can select the columns that display in the tube section of the screen.

In the **Plot** settings, you can define the background of the graphics display area, configure the histograms, and set the default signal parameters to either the channel's area or the channel's height. The default is **Area**. You can also set the default axis display range, and set the spectrum plot resolution.




2-46 D17052AA






In the **Page Setup** settings, you can change the page size, orientation, margin size, and display options. Select **Show page breaks** to display page boundaries within the Acquisition or Analysis views for simplifying plot arrangement for printing.



In the **Plate Loader** settings, you can select the plate type, sampling sequence, mix, and backflush settings for the plate loader.

NOTE This setting is only available in the Plate Loader sample injection mode.



2-48 D17052AA

CHAPTER 3 Operation Principles

Overview

This chapter explains how the Cytometer measures scattered light and fluorescence as cells pass through the laser beam.

The illustrations in this chapter are not exact representations of the inside of the Cytometer. They are for explanatory purposes only.

This chapter contains information on:

- Sample Flow
- Laser Beam Shaping
- Cell Illumination
- Light Collection, Separation and Measurement
- Signal Processing
- Data Storage
- Automated Software Features
- Parameters
- Plot Display
- Statistics

Sample Flow

Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 3, Sample Flow for detailed instructions.

Sample Loading

Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 3, Sample Loading for detailed instructions.

Hydrodynamic Focusing

Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 3, Hydrodynamic Focusing for detailed instructions.

Laser Beam Shaping

Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 3, Laser Beam Shaping for detailed instructions.

Cell Illumination

Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 3, Cell Illumination for detailed instructions.

Forward Scatter

Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 3, Forward Scatter for detailed instructions.

Side Scatter and Fluorescent Light

Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 3, Side Scatter and Fluorescent Light for detailed instructions.

Light Collection, Separation and Measurement

Forward Scatter Collection

Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 3, Forward Scatter Collection for detailed instructions.

Side Scatter and Fluorescent Light Collection

Both side scatter and fluorescence are measured 90 degrees from the laser axis.

NOTE If you are using Conventional Mode, refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 3, Side Scatter and Fluorescent Light Collection for detailed instructions.

Side Scatter

The CytoFLEX mosaic Spectral Detection Module collects SS from each individual laser. It is much more intense than FL.


Side scatter light collected by the objective lens is delivered by fiber optics to a patent-pending design with high performance, solid-state, high efficiency, and low-noise detector arrays.

3-2 D17052AA

Fluorescent Light

Fluorescence and scattered light are transmitted by optical fibers to the Wavelength division multiplexer (WDM). Each WDM is a unique detector array that corresponds to a different laser. Refer to CHAPTER 1, Wavelength Division Multiplexer (WDM). Each WDM contains optical filters and detectors for detecting channel fluorescence or scatter from a particular laser. It is necessary to ensure that the filter and software settings match for each channel.

Figure 3.1 Light Path through the WDM with a Single Port



- 1. Fiber array photo detectors (FAPD)
- 2. SSC fiber array photo detector (SSC FAPD)
- 3. Lens
- 4. Lens
- 5. SSC filter
- **6.** ND filter
- 7. Filters

- 8. Filter
- 9. 30-degree reflector
- 10. Filter
- 11. Mirror
- 12. Mirror
- **13.** Lens
- 14. Light path

Signal Processing

Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 3, Signal Processing for detailed instructions.

Data Storage

Refer to the CytoFLEX Platform Instructions for Use manual, CHAPTER 3, Data Storage for detailed instructions.

NOTE Sample results can be printed out, saved to removable media, saved to a local hard drive or saved to a network drive. You can store sample results in Flow Cytometry Standard (FCS 3.1 or FCS 3.0) files.

Automated Software Features

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 3, Automated Software Features for detailed instructions.

Parameters

TIME Parameter

Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 3, TIME Parameter for detailed instructions.


Plot Display


The results of sample analysis appear on the Workstation screen as graphs called plots. Refer to CHAPTER 2, Graphic and Gating Styles.


3-4 D17052AA

Statistics

The Statistics Setting window allows you to change the display of the header, statistical elements and cell populations included.







Operation Principles Statistics

D17052AA 3-6


CHAPTER 4 Daily Startup

Overview

IMPORTANT Verify that the two correct USB configuration keys are securely connected to the computer USB ports. If the USB configuration keys are not connected, the following error message appears: CytExpert cannot find the license. Please check whether the correct USB configuration key has been plugged in.

This chapter describes the instrument startup procedure.

Workflow:



This chapter contains information on:

- Pre-Startup Inspection
- Turning On the Instrument
- Logging into the Software
- Initializing the Instrument

Pre-Startup Inspection

Before using the Cytometer with the CytoFLEX mosaic Spectral Detection Module, perform the following system checks.

Checking Waste and Reagent Levels [4 L Fluid Containers]

Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 4, Checking Waste and Reagent Levels [4 L Fluid Containers] for detailed instructions.

Checking Waste and Reagent Levels [10 L Fluid Cubitainers]

Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 4, Checking Waste and Reagent Levels [10 L Fluid Cubitainers] for detailed instructions.

Power Source Inspection

Check the power cable located below the power switch on the back of the Cytometer and the CytoFLEX mosaic Spectral Detection Module, and verify it is securely connected to the Cytometer, the CytoFLEX mosaic Spectral Detection Module and the appropriate power source.

Workstation Connections Inspection

Check that the monitor, mouse, keyboard, the Cytometer, and the CytoFLEX mosaic Spectral Detection Module are properly connected to the computer. Refer to Figure 1.3.

Turning On the Instrument



- 1. If the Cytometer, CytoFLEX mosaic Spectral Detection Module or Workstation fails to start properly, check first to see whether the power cable and connection cables are properly connected.
- 2. Never shut off the power or disconnect a data cable while the instrument is performing a task. Doing so can result in data loss or damage to the system.

NOTE If you are using Conventional Mode, refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 4, Turning On the Instrument for detailed instructions.

- 1 Turn on the main power switch located on the back of the Cytometer.
- Wait for the Cytometer to finish powering on, then turn on the CytoFLEX mosaic Spectral Detection Module using the power switch located on the left side of the CytoFLEX mosaic Spectral Detection Module followed by the power button on the front.
- **3** Wait for the CytoFLEX mosaic Spectral Detection Module to finish powering on, then turn on the Workstation.

NOTE The CytoFLEX mosaic Spectral Detection Module takes approximately 40 seconds to power on.

4-2 D17052AA

Logging into the Software

1 Log in to the Windows operating system and double-click the CytExpert for Spectral desktop icon to open the software, the login window appears.



NOTE The default software shortcut appears on the desktop. If you do not see the icon, the default installation path is under C:/Program Files/CytExpert for Spectral. Double-click CytExpert.exe to run the software.


- 2 Enter your username and password.
- 3 Select .

NOTE The display name of the user that is currently logged in displays in the top, right corner of the software screen.




D17052AA

- **4** Confirm that the software, the CytoFLEX mosaic Spectral Detection Module and the Cytometer are properly connected.
 - **a.** Open the software. The Startup screen appears.



b. Verify that the connection indicator light in the lower left corner of the software screen is green, and *Connected* is displayed. The left side shows the connection status, the middle shows the instrument status, and the right side shows the status details.



c. Verify that the *Sheath* and *Waste* flow indicators in the lower right corner of the software screen are green indicating that the fluidics system is normal.



4-4 D17052AA

d. Verify that the CytoFLEX mosaic Spectral Detection Module status indicator in the lower right corner of the software screen has a green checkmark indicating that the CytoFLEX mosaic Spectral Detection Module is ready.



NOTE SytoFLEX mosaic indicates that the CytoFLEX mosaic Spectral Detection Module is warming up. This process takes approximately 4 minutes.

NOTE

• A red connection indicator light indicates that there is a faulty connection. Ensure that the instrument is properly turned on and connected. If necessary, restart both the Cytometer and the Workstation.




 After the instrument initializes, a warning beep sounds if there is a problem with the fluidics system. If a flow indicator is red and blinking, it means that the fluidics system requires attention.



- When the waste fluid sensor is disconnected, the waste flow indicator shows that the waste container is full or nearly full.
- A gray CytoFLEX mosaic Spectral Detection Module status indicator indicates that there is a faulty connection. Ensure that the CytoFLEX mosaic Spectral Detection Module is properly turned on and connected. If necessary, restart the CytoFLEX mosaic Spectral Detection Module.



 Select the status information in the lower left to open the system log. Send a copy of the system log to your Beckman Coulter Representative for support if a service call is requested.



Logging Out of the Software

Select the username displayed in the top-right corner of the software screen and select Log out.



Locking the Account

Select the username displayed in the top-right corner of the software screen and select **Lock**.



The account locks automatically if it remains inactive for a specified duration. Refer to CHAPTER 2, Account Policies - Application Inactivity Policies.

Selecting the Proper Sample Injection Mode

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 4, Selecting the Proper Sample Injection Mode for detailed instructions on the following procedures:

- Using Semi-Automatic Injection Mode
- Using Manual Injection Mode

4-6 D17052AA

Selecting the Plate Loader Sample Injection Mode [With Plate Loader]

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the CytoFLEX Platform Instructions for Use manual, CHAPTER 4, Selecting the Plate Loader Sample Injection Mode [With Plate Loader] for detailed instructions on the following procedure:

Using Plate Loader Injection Mode

Running the System Startup Program [with the Single Tube Loader]

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 4, Running the System Startup Program [with the Single Tube Loader] for detailed instructions.

Running the System Startup Program [With Plate Loader]

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 4, Running the System Startup Program [With Plate Loader] for detailed instructions.

Selecting Experiments from the Start Page

Refer to CHAPTER 2, Start Page.

Initializing the Instrument

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 4, Initializing the Instrument for detailed instructions.

Daily StartupInitializing the Instrument

D17052AA 4-8

Instrument Quality Control and Standardization

Overview

This chapter provides information on performing daily Quality Control (QC) on the CytoFLEX flow cytometer equipped with the CytoFLEX mosaic Spectral Detection Module and how to confirm that the instrument is working properly within the specified parameters. Quality Control allows you to determine whether your instrument can provide adequate signal strength and precision.

This chapter also provides information on performing standardization. CytoFLEX Ready to Use Daily QC Fluorospheres, CytoFLEX Daily IR QC Fluorospheres or any other reference material that is relevant for your application(s) may be used as the standardization sample(s). The system can only recognize a single peak.

Standardization can be used to monitor the Median Fluorescent Intensities (MFI), or target values for scatter and fluorescent parameters that have been defined and optimized for a specific application. All channels in the current configuration, those with/without an assigned QC target, can be tracked as necessary via Standardization since this portion of the CytExpert for Spectral software is used to assess application specific settings. Standardization, however, does NOT replace QC as the Cytometer's optical alignment (rCV statistical analysis), Laser Power and Laser Delay outputs are not measured during the run.


NOTE Beckman Coulter recommends performing QC on a daily basis.

NOTE CytExpert for Spectral QC includes an automated QC routine with Levey-Jennings (LJ) charts tracking and logging.


NOTE CytExpert for Spectral Standardization allows for application-specific settings to be established and applied to future experiments.

NOTE Perform QC on the same day prior to performing the CytExpert for Spectral standardization.

QC Workflow:



Standardization Workflow:



This chapter contains information on:

- Quality Control
 - Preparing the QC Sample
 - Preparing the QC Sample [With Plate Loader]
 - Importing Lot-Specific Target Values
 - Collecting QC Data
 - Collecting QC Data [With Plate Loader]
 - Confirming Results
- Standardization
 - Preparing the Standardization Sample
 - Generating Target Median Values
 - Creating a New Standardization Item
 - Performing the Standardization
 - Applying the Standardized Acquisition Settings
 - Standardization Target Library

Quality Control

The QC process verifies important system functions. The system:

- 1. Measures the laser power of each individual laser and ensures that each laser meets the system specifications.
- **2.** Loads the QC sample and begins to acquire data.
- **3.** Verifies that the actual laser delays match those set in the software and will adjust the delay accordingly.
- **4.** Notifies you if laser delay is >2 µs from the previous setting. The software automatically changes the laser delay setting.

OR

Notifies you if laser delay is >5 μ s from the previous setting. Manual laser delay adjustments are required. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 12, Setting Laser Delay.

5. Verifies and calibrates the gain settings. If any of these parameters are outside of the operating limits, the system automatically adjusts these parameters. If the system is unable to adjust these parameters to fall within the operating limits, the system notifies you.

5-2 D17052AA

Preparing the QC Sample

Required Materials

The following materials are required to complete the QC process:

- CytoFLEX Ready to Use Daily QC Fluorospheres
- CytoFLEX Daily IR QC Fluorospheres (for systems configured with an IR laser)
- CytoFLEX Sheath Fluid or another nonionic antimicrobial sheath fluid
- Sample tubes (12 x 75 mm)
- Vortexer

CytoFLEX Ready to Use Daily QC Fluorospheres Preparation Process

Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 5, CytoFLEX Ready to Use Daily QC Fluorospheres Preparation Process for detailed instructions.

CytoFLEX Daily IR QC Fluorospheres Preparation Process

Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 5, CytoFLEX Daily IR QC Fluorospheres Preparation Process for detailed instructions.

Preparing the QC Sample [With Plate Loader]

Required Materials

The following materials are required to complete the QC process:

- CytoFLEX Ready to Use Daily QC Fluorospheres
- CytoFLEX Daily IR QC Fluorospheres (for systems configured with an IR laser)
- CytoFLEX Sheath Fluid or another nonionic sheath fluid
- Standard 96-well plate
 - 96-well flat-bottom
 - 96-well V-bottom
 - 96-well U-bottom
- 96-well deep well plate
 - 96-well V-bottom
 - 96-well U-bottom
- Vortexer

CytoFLEX Ready to Use Daily QC Fluorospheres Preparation Process

Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 5, Preparing the QC Sample [With Plate Loader], CytoFLEX Ready to Use Daily QC Fluorospheres Preparation Process for detailed instructions.

CytoFLEX Daily IR QC Fluorospheres Preparation Process

Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 5, Preparing the QC Sample [With Plate Loader], CytoFLEX Daily IR QC Fluorospheres Preparation Process for detailed instructions.

Importing Lot-Specific Target Values


Import lot-specific target values for each new lot of CytoFLEX Ready to Use Daily QC Fluorospheres, or CytoFLEX Daily IR QC Fluorospheres.

NOTE If you are using Conventional Mode, refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 5, Importing Lot-Specific Target Values for detailed instructions.



Risk of erroneous QC results. Different target value information correspond to different lot numbers. Selecting the wrong lot number will lead to erroneous QC results.


- 1 Select **Start QC/Standardization** from the QC/Standardization menu.
- **2** Select **Target Library** from the Settings menu. The Target Library window appears.



D17052AA


IMPORTANT The Beckman Coulter website may prompt you to select your Region and Country prior to the Beckman Coulter Technical Documents and Software page.

3 Select **Download Target File**. The Beckman Coulter Technical Documents and Software Downloads page appears.




NOTE If your CytoFLEX Workstation does not have access to the internet, navigate to https:// www.beckmancoulter.com/wsrportal/page/softwareDownloadSearch using a computer with access to the internet and save the file to a USB drive. If the website is not accessible, contact us.

- 4 If necessary, register and log in to the Beckman Coulter website.
- 5 In the Search By Product section of the screen, select the following:
 - a. Select Research & Discovery from the Market Segment drop-down menu.
 - **b.** Select **Flow Cytometry** from the Product Line drop-down menu.
 - **c.** Select **Instruments** from the Product Platform drop-down menu.
 - d. Select CytoFLEX mosaic from the Product drop-down menu.
 - e. Select CytoFLEX Ready to Use QC Fluorospheres Target Value, or CytoFLEX Daily IR QC Fluorospheres Target Values from the Software Name drop-down menu.
 - **f.** Select **English** from the Language drop-down menu.
- 6 Select Search.



7 The search results appear below the Search By Item/REF Number tab.

8 Select CytoFLEX Daily Ready to Use QC Fluorospheres Target Values, or CytoFLEX Daily IR QC Fluorospheres Target Values under the Software Name column. The CmpytoFLEX Daily Ready to Use QC Fluorospheres Target Values or CytoFLEX Daily IR QC Fluorospheres Target Values page appears.




5-6 D17052AA

- **9** Select **Download** next to the correct lot number from the CytoFLEX Daily Ready to Use QC Fluorospheres Target Values or CytoFLEX Daily IR QC Fluorospheres Target Values page.
- 10 If the File Download pop up window appears, select **Save** and browse to the desired file path.
- 11 Select Import from the Target Library window in the CytExpert for Spectral software.
- 12 Navigate to the file saved in Step 10 and select Open.
- 13 Select Close to exit the Target Library window.

Collecting QC Data

QC data and reports are saved by default. Select **QC/Standardization Setting** in the Settings menu to change the default save settings or modify the file path these files are saved to.



NOTE If you are using Conventional Mode, refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 5, Collecting QC Data for detailed instructions.

- 1 Double-click to start the CytExpert for Spectral software.
 - **a.** Ensure that the **Connected** icon on the Status Bar near the bottom-left side of the display is green.




NOTE If the icon is not green, ensure that the Cytometer USB is securely connected to the Workstation and restart the Workstation. Refer to Table 10.1, Troubleshooting [With CytoFLEX mosaic Spectral Detection Module].

b. Ensure that the **CytoFLEX mosaic** icon on the Status Bar near the bottom-right side of the display has a green checkmark.



- **NOTE** If the icon doesn't have a green checkmark, ensure that the network cable is securely connected to the Workstation and restart the Workstation. Refer to Table 10.1, Troubleshooting [With CytoFLEX mosaic Spectral Detection Module].
- **2** Verify the laser settings. Refer to CHAPTER 6, Laser Settings.
- 3 Select Start QC/Standardization in the QC/Standardization menu to access the QC experiment.

[CytoFLEX LX Shown]



Ensure that the QC bead lot number is selectable in the Lot No. drop down menu. If the lot number is not selectable, refer to Importing Lot-Specific Target Values, then select the proper lot number.


5-8 D17052AA

- 4 Select Initialize.
- Insert the prepared QC sample tube (see CytoFLEX Ready to Use Daily QC Fluorospheres Preparation Process) into the tube holder.
- 6 Select **Start** to load the sample and begin to run the QC procedure.

 Completed processes appear on the left. Plots appear on the right. The QC experiment sequentially detects the system configuration, laser power, laser delay, signal strength, and coefficient of variation.



During QC, the software automatically seeks the CytoFLEX QC Fluorospheres and computes the results. The software returns to the QC Report screen after the QC run is complete.



If the sampling rate is too low, the Cytometer stops the QC run and displays a prompt that the QC run fails to reach the required event flow rate. This is not considered a QC failure. If this situation occurs, increase the sample concentration by preparing a new tube of CytoFLEX Ready to Use Daily QC Fluorospheres (see CytoFLEX Ready to Use Daily QC Fluorospheres Preparation Process) and then perform the experiment.

NOTE Mix the fluorospheres dropper bottle thoroughly before dispensing it into a sample tube.

NOTE If the sample flow rate is not correct, calibrate the sample flow rate before running the QC. For instructions on calibrating the sample flow rate, refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 12, Calibrating the Sample Flow Rate or CHAPTER 12, Calibrating the Sample Flow Rate [With Plate Loader].

If the lot number of CytoFLEX QC Fluorospheres is new and QC fails, the following software message appears. Select **Yes**.



NOTE Target gain values must be established for each new lot number of CytoFLEX QC Fluorospheres. QC could fail up to 3 times upon running each new lot number for the first time until target gain values are established.

If the lot number of CytoFLEX QC Fluorospheres is NOT new and QC fails, refer to Step 2 of CHAPTER 5, Confirming Results, or CHAPTER 10, Troubleshooting.


If QC passes, proceed to Step 9.

9 Run Daily Clean to remove any residual fluorosphere particles if needed. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 11, Daily Clean.

5-10 D17052AA

Collecting QC Data [With Plate Loader]

QC data and reports are saved by default. Select **QC/Standardization Setting** in the Settings menu to change the default save settings or modify the file path these files are saved to.



NOTE If you are using Conventional Mode, refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 5, Collecting QC Data [With Plate Loader] for detailed instructions.

- 1 Double-click to start the CytExpert for Spectral software.
 - **a.** Ensure that the **Connected** icon on the Status Bar near the bottom-left side of the display is green.



NOTE If the icon is not green, ensure that the Cytometer USB is securely connected to the Workstation and restart the Workstation. Refer to Table 10.1, Troubleshooting [With CytoFLEX mosaic Spectral Detection Module].

b. Ensure that the **CytoFLEX mosaic** icon on the Status Bar near the bottom-right side of the display has a green checkmark.




NOTE If the icon doesn't have a green checkmark, ensure that the network cable is securely connected to the Workstation and restart the Workstation. Refer to Table 10.1, Troubleshooting [With CytoFLEX mosaic Spectral Detection Module].

2 Verify the laser settings. Refer to CHAPTER 6, Laser Settings.

IMPORTANT Ensure that the QC bead lot number is selectable in the Lot No. drop-down menu. If the lot number is not selectable, refer to CHAPTER 5, Importing Lot-Specific Target Values, then select the proper lot number.

3 Select Start QC/Standardization in the QC/Standardization menu to access the QC experiment.


- 4 Select Initialize.
- 5 Select **Eject**.
- **6** Insert the prepared QC well plate (see CytoFLEX Ready to Use Daily QC Fluorospheres Preparation Process or CytoFLEX Daily IR QC Fluorospheres Preparation Process) into the plate holder.



IMPORTANT Ensure the well position on the plate matches the well position selected in the software.

8 Select the appropriate QC well.

5-12 D17052AA




9 Select the desired plate type from the Plate Type dropdown menu.


NOTE The available plate types included in the dropdown menu depend on the settings selected in the Plate Library. Refer to CHAPTER 2, Plate Type Library.

- 10 Select the Mix and Backflush settings in the top of the Plate Settings window.
- 11 Select ok.
- 12 Select **Start** to load the sample and begin to run the QC procedure. The message *Please confirm* that the correct plate is placed properly and press OK appears. Select **OK**.

Completed processes appear on the left. Plots appear on the right. The QC experiment sequentially detects configuration, laser power, laser delay, signal strength, and coefficient of variation.



During QC, the software automatically seeks the CytoFLEX QC Fluorospheres, computes the results and determines whether the gain setting meets the assay target values. The software returns to the QC screen after the QC run is complete.




13 If the sampling rate is too low, the Cytometer stops the QC run and displays a prompt that the QC run fails to reach the required event flow rate. This is not considered a QC failure. If this situation occurs, increase the sample concentration by preparing a new tube of CytoFLEX Ready to Use Daily QC Fluorospheres (see CytoFLEX Ready to Use Daily QC Fluorospheres Preparation Process) and then perform the experiment.

NOTE Mix the fluorospheres dropper bottle thoroughly before dispensing it into a sample tube.

NOTE If the sample flow rate is not correct, calibrate the sample flow rate before running the QC. For instructions on calibrating the sample flow rate, refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 12, Calibrating the Sample Flow Rate or CHAPTER 12, Calibrating the Sample Flow Rate [With Plate Loader].

14 If the lot number of CytoFLEX QC Fluorospheres is new and QC fails, the following software message appears. Select **Yes**.



NOTE Target gain values must be established for each new lot number of CytoFLEX QC Fluorospheres. QC could fail up to 3 times upon running each new lot number for the first time until target gain values are established.

5-14 D17052AA

If the lot number of CytoFLEX QC Fluorospheres is NOT new and QC fails, refer to Step 2 of CHAPTER 5, Confirming Results, or CHAPTER 10, Troubleshooting.


If QC passes, proceed to Step 15.

15 Run Daily Clean to remove any residual fluorosphere particles. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 11, Daily Clean [With Plate Loader].

Confirming Results

Select **Start QC/Standardization** in the QC/Standardization menu to return to the QC Setting screen at any time to review completed experiment results.

1 Select a QC run from the QC Process list on the left and a QC report appears on the right.



NOTE The results column indicates a passing QC result with a and a failed QC result with a. QC results must meet the following criteria to pass:

- The gain differences must be ≤20% from the target gain.
- The median fluorescence intensity (MFI) differences must be ≤5% from the target MFI.
- The rCV must be within the target range.

The report area on the right displays detailed experiment results, including laser power, delay, testing conditions, and signal results. The same \bigcirc and \bigcirc symbols are used to indicate each result. For items that fail, values falling outside the prescribed range are displayed in red font. In the Comment area, an explanation appears for each failed item.

QC Report

		•	
Bead Lot No.:	RTU 15ARK Rev AA Spectral		
Bead Expires:	2025-10-04	QC Date:	2024-10-23 14:16
Cytometer Name:	CytoFLEX LX	Cytometer SN:	BG46032
CytoFLEX mosaic Name:	CytoFLEX mosaic 88	CytoFLEX mosaic SN:	MP8003
Loader Type:	Semi Automatic		

Threshold

Channel: FSC(Height) Mode: Manual Value: 150000

Laser

Laser	Delay(µs)	Default Delay(µs)	Difference Delay(µs)	Power(mW)	Target Power(mW)	Result
UV	80.80	83.04	-2.24	20	18-22	0
Violet	-91.20	-93.60	2.40	83	70-120	0
Blue	0.00	0.00	0.00	51	40-60	0
Yellow	-40.64	-41.60	0.96	30	20-40	0
Red	39.84	40.64	-0.80	53	40-60	0
IR	-123.52	-124.00	0.48	64	50-70	0

Signal Value

Parameter	Gain	Target Gain	%Difference Target Gain	Median	Target Median	%Difference Target Median	rCV(%)	Target rCV(%)	Width	Assay Gain	Result
Scatters Ø										0	
FSC	108	113	-4.42	643632.4	650536.1	-1.06	2.56	-	993.1	108	S
USSC	353	353	0.00	1625229.0	1605534.0	1.23	3.82	-	751.3	353	Ø
BSSC	46	44	4.55	2587672.0	2565757.0	0.85	5.41	-	969.9	46	Ø
YSSC	107	107	0.00	3453703.0	3577131.0	-3.45	2.78	-	1233.2	107	✓
RSSC	37	35	5.71	3198941.0	3120219.0	2.52	6.38	-	1189.8	37	✓
IRSSC	201	201	0.00	3828751.0	3877891.0	-1.27	4.22	-	1570.3	201	✓
UV Ø											
U1	428	428	0.00	205886.3	203763.4	1.04	1.79	-	760.5	428	Ø
U2	341	341	0.00	2840192.0	2813103.0	0.96	1.61	-	767.1	341	✓
U3	146	146	0.00	4448925.0	4398882.0	1.14	1.56	7.00	767.0	146	Ø
U4	208	208	0.00	3967010.0	3940995.0	0.66	1.38	7.00	765.0	208	Ø
U5	264	264	0.00	3511488.0	3495994.0	0.44	1.32	7.00	766.9	264	Ø
U6	357	357	0.00	3218966.0	3183130.0	1.13	1.44	7.00	767.7	357	Ø
U7	316	316	0.00	3722054.0	3690475.0	0.86	1.38	-	770.7	412	✓
U8	500	500	0.00	3555493.0	3504376.0	1.46	1.20	-	770.1	500	Ø
U9	444	444	0.00	3615184.0	3581448.0	0.94	1.29	-	772.6	661	✓
U10	587	587	0.00	3550787.0	3504628.0	1.32	1.49	-	773.8	587	Ø
U11	814	814	0.00	1961909.0	1949023.0	0.66	1.49	-	772.5	814	Ø
U12	661	661	0.00	1628309.0	1617955.0	0.64	1.16	-	775.2	661	✓
U13	528	528	0.00	2674078.0	2665883.0	0.31	1.21	-	781.0	528	Ø
U14	660	660	0.00	2893684.0	2875115.0	0.65	1.20	-	785.6	660	Ø

Specification

Delay: -5.00µs ≤ Difference Delay ≤ 5.00µs.
Gain: -20.00% ≤ %Difference Target Gain ≤ 20.00%.
Median: -5.00% ≤ %Difference Target Median ≤ 5.00%.
rCV: rCV(%) ≤ Target rCV(%).

Result

QC Passed.

Commen


The gain of VSSC is fixed to 1. No target value for IR1(IR,(830 - 850)),IR2(IR,(850 - 905)),IR3(IR,(905 - 950)).

5-16 D17052AA

2 If QC fails, follow the procedure below:

NOTE When the Confirm window appears, select **Yes** if you are sure that the QC failure is caused by the QC reagent batch difference or you need to use the adjusted gain.

Otherwise, select **No**.



- **a.** Verify whether the beads used were within their shelf life and stored in accordance with the appropriate instruction manual.
- **b.** Verify whether the allocated sample tube was prepared as required and correctly positioned.
- **c.** Prime the Flow Cell, and retest. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 12, Priming the Flow Cell.
- **d.** Run the Daily Clean, and retest. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 11, Daily Clean.
- **e.** Run the Deep Clean Procedure, and restart. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 11, Deep Clean Procedure.
- **f.** Repeat Steps c-d.

NOTE If QC fails two times in a row on the same day after repeating Steps a-f, contact us.

- If necessary, you can select (for CSV format) or (for PDF format) in the top left corner of the report area to export the QC results.
- 4 Select Close QC/Standardization in the File menu to exit the QC screen.

Creating Levey-Jennings Charts

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the CytoFLEX Platform Instructions for Use manual, CHAPTER 5, Creating Levey-Jennings Charts for detailed instructions.

QC Result Manager

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 5, QC Result Manager for detailed instructions.

Standardization

Ensure that the standardization sample has been run at optimized experiment settings to determine the standardization sample threshold setting as well as median values for all relevant channels.

Preparing the Standardization Sample

Use Beckman Coulter CytoFLEX Ready to Use Daily QC Fluorospheres or any other reference material that is relevant for your application.

Required Materials

The following materials are required to complete the QC process:

CytoFLEX Ready to Use Daily QC Fluorospheres, or other material applicable for your application

NOTE Note that the VSSC and NSSC channels are designed for small particle detection. It is not recommended to use the CytoFLEX Ready to Use Daily QC Fluorospheres when performing standardization for the VSSC and NSSC channels.

- CytoFLEX Daily IR QC Fluorospheres (for systems configured with an IR laser)
- CytoFLEX Sheath Fluid
- Sample tubes (12 x 75 mm)
- Vortexer

Preparation Process

For procedures, refer to CytoFLEX Ready to Use Daily QC Fluorospheres Preparation Process and/or CytoFLEX Daily IR QC Fluorospheres Preparation Process.

5-18 D17052AA

Generating Target Median Values

Generate lot-specific target values for each new lot of standardization beads.

NOTE If you want to generate target median values automatically for all channels from the Acq. Setting Catalog using existing settings in the Acq. Setting Catalog, refer to Creating a New Standardization Item from Acq. Setting Catalog.

NOTE If you are using Conventional Mode, refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 5, Generating Target Median Values for detailed instructions.

- 1 Double-click to start the CytExpert for Spectral software.
 - **a.** Ensure that the **Connected** icon on the Status Bar near the bottom-left side of the display is green.



NOTE If the icon is not green, ensure that the Cytometer USB is securely connected to the Workstation and restart the Workstation. Refer to Table 10.1, Troubleshooting [With CytoFLEX mosaic Spectral Detection Module].

b. Ensure that the **CytoFLEX mosaic** icon on the Status Bar near the bottom-right side of the display has a green checkmark.



NOTE If the icon don't have a green checkmark, ensure that the network cable is securely connected to the Workstation and restart the Workstation. Refer to Table 10.1, Troubleshooting [With CytoFLEX mosaic Spectral Detection Module].


- 2 Select **New Experiment** to create an experiment.
- Add a tube and edit the tube name. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 6, Changing the Tube Name for detailed instructions.
- Select **Acq. Setting** to input the optimized Gain and Threshold in the relevant channels in the Acq. Setting window. Refer to CHAPTER 6, Configuring Acquisition Settings.


NOTE Use the FSC channel as the trigger channel. The threshold may need to be adjusted to visualize the QC beads populations. If so, record this value for future reference.

5 Load the sample tube.

NOTE The sample tube holder accommodates 1.5 mL, 2.0 mL, and 12 x 75 mm sample tubes.


- **6** Create a FSC/SSC dot plot and gate on the singlet bead population. Refer to CHAPTER 6, Creating Plots and Gates.
- 7 Create histograms for each channel and apply the gate.
- 8 Select Run.
- **9** Add a statistics table and right-click the table and select **Statistics Setting**. The Statistics Setting window appears.







5-20 D17052AA

10 Select the **Population** tab and select the relevant population for the tube.



11 Select the **Statistics** tab then select the Median Fluorescence value for all parameters used and select **OK**.





NOTE The median values are the target settings that will be used for standardization.

12 Right-click the statistics table and select **Export to CSV File**.

If Excel is not available, manually record all the median values or take a screen shot.

5-22 D17052AA

13 Save the experiment.

NOTE Rerun the experiment and create a new standardization file if:

- · You change the standardization fluorosphere used.
- The Lot number for the standardization fluorosphere is changed.

Creating a New Standardization Item


The standardization item can be created in two ways:

- Add a new standardization item in the Standardization Target Library using the previously generated target median values. Refer to Generating Target Median Values and Adding a New Standardization Item.
- Create a new standardization item from the Acq. Setting Catalog using the previously saved acquisition settings. Refer to CHAPTER 6, Exporting Instrument Settings and Creating a New Standardization Item from Acq. Setting Catalog.


Adding a New Standardization Item

NOTE If you are using Conventional Mode, refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 5, Adding a New Standardization Item for detailed instructions.


Select **Start QC/Standardization** in the QC/Standardization menu to access the Standardization screen.



2 Select Settings > Standardization Target Library. The Standardization Target Library window appears.



3 Select **Add**. The Add Standardization Target Value window appears.



5-24 D17052AA

4 Enter the Item, Lot No. and Expiration date from the drop-down menu located at the top of the Add Standardization Target Value window.

NOTE A single Lot No. can include several Items, but you cannot add duplicate Items under the same Lot No.

NOTE If the Lot No. selected already exists, the Expiration date cannot be edited.

Select either Manual or Automatic threshold from the Standardization Test Setting or Acquisition Setting section of the window.

NOTE If you select **Manual** threshold, enter a value greater than 0, but less than 1,048,575.

NOTE Keep the threshold setting the same as previous Step 4 in Generating Target Median Values.

NOTE Unchecking the *Use the same threshold setting for acquisition setting* checkbox allows you to specify custom threshold settings and when to save the test item into the Acquisition Setting Catalog.

6 Enter the target median values saved in Step 11 in Generating Target Median Values into the corresponding channels.

Or copy the median values from the previously exported CSV file in Step 12 in Generating Target Median Values and paste into the corresponding median column,


NOTE Do not set the median tolerance range any lower than 5%.

NOTE FSC is a required channel.

NOTE Verify that the target values are entered correctly.

7 Select **OK** to save the target value.

The saved results display in the Standardization Target Library window. This item is ready to be run through the Standardization experiment. Refer to Performing the Standardization.




8 Select **Close** to exit the Standardization Target Library window,


5-26 D17052AA


Creating a New Standardization Item from Acq. Setting Catalog

1 Select **Start QC/Standardization** in the QC/Standardization menu to access the Standardization screen.



2 Select **Settings** > **Standardization Target Library**. The Standardization Target Library window appears.






3 Select Create from Acq. Setting Catalog.... The Create from Acq. Setting Catalog window appears.


5-28 D17052AA

4 Select Select... located at the top left of the Create from Acq. Setting Catalog window. The Acq. Setting Catalog window appears.



5 Select the desired Item and select **OK**.

NOTE If the following prompt appears, select **Yes** to proceed. If you do not want to proceed, select **No**.



6 Enter the Lot No. and Expiration date from the drop-down menu located at the top of the Create from Acq. Setting Catalog window.

NOTE A single Lot No. can include several Items, but you cannot add duplicate Items under the same Lot No.

NOTE If the Lot No. selected already exists, the Expiration date cannot be edited.

7 Select either Manual or Automatic threshold from the Standardization Test Setting or Acquisition Setting section of the window.


NOTE If you select **Manual** threshold, enter a value greater than 0, but less than 1,048,575.

NOTE Unchecking the *Use the same threshold setting for acquisition setting* checkbox allows you to specify custom threshold settings.

8 Select Initialize from the Collect Control Panel and select Record.

NOTE If you need to set the number of events to record in the Collect Control Panel, enter a value between 2,000 and 100,000.


9 Wait for the software to finish recording. The Save as FCS File window appears.



If you want to export the FCS file, navigate to the desired file path and select **Save**. If you do not want to export the FCS file, select **Cancel**.

5-30 D17052AA

10 The generated median values display in the Create from Acq. Setting Catalog window.




- 11 Select **OK** to save the target value and exit the Create from Acq. Setting Catalog window.
- 12 Select Close to exit the Standardization Target Library window.

Performing the Standardization

NOTE If you are using Conventional Mode, refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 5, Performing the Standardization for detailed instructions.

Select **Start QC/Standardization** in the QC/Standardization menu to access the Standardization screen.




2 Select the **Standardization** radio button.




5-32 D17052AA

3 Select the correct Lot No. from the Lot No. drop-down menu.



NOTE Ensure the Lot No. corresponds to the standardization sample that generated the target median values.


4 Select the Items to be standardized.




- 5 Select Initialize.
- **6** Select **Start**. The following software prompt appears. Select **Yes**.



The Process section of the screen displays the process details.



Once the process is complete, the standardization results displays in the Process list on the left and a Standardization Report displays on the right.




NOTE The result column indicates a passing standardization result with a and a failed standardization result with a.

NOTE The updated Standardization item is added to the Acquisition Catalog automatically and overwrites the previously existing standardized settings for this item.

5-34 D17052AA

7 Select **Acq. Setting Catalog** from the Cytometer menu to verify the gain settings. The Acq. Setting Catalog window appears.



NOTE S designates test items from Standardization.

- **8** Select **Close** to exit the Acq. Setting Catalog window.
- **9** Run Daily Clean. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 11, Daily Clean.
- 10 Select Close QC/Standardization.


Applying the Standardized Acquisition Settings

NOTE If you are using Conventional Mode, refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 5, Applying the Standardized Acquisition Settings for detailed instructions.


1 Open an experiment.

2 Select Acq. Setting from the Cytometer menu. The Acq. Setting window appears.

Acq. Setting - Tube1




3 Select Import from Catalog. The Acq. Setting Catalog window appears.



4 Browse for the item to import and select **Import**. The standardized settings are applied to the sample tube.

5-36 D17052AA

The Information window appears to notify of the corresponding channels with the changed gain as a result of the Standardization.



5 Select **OK**.

Standardization Target Library

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 5, Standardization Target Library for detailed instructions on the following procedures:

- Importing a Standardization Item
- Exporting a Standardization Item

NOTE The standardization items save as .mstgt file. This file can be used to standardize the settings between different instruments with the same laser configuration.

- Editing Standardization Item Parameters
- Duplicating Standardization Items
- Deleting Standardization Items

Instrument Quality Control and Standardization Standardization

5-38 D17052AA

Data Acquisition and Sample Analysis

Overview

This chapter contains information on how to use your Cytometer with CytoFLEX mosaic Spectral Detection Module, including data acquisition, analyzing and exporting results.

Workflow:



This chapter contains information on:

- Creating an Experiment
- Configuring Acquisition Settings
- Load Sample and Record Data
- Analyzing and Exporting Data
- Saving the Experiment

Creating an Experiment

Creating an Experiment [Without Plate Loader]

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 6, Creating an Experiment [Without Plate Loader] for detailed instructions.

NOTE CytExpert for Spectral software has two installation options, the CytExpert for Spectral User Management software option and the CytExpert for Spectral Electronic Record Management software option.

Creating an Experiment [With Plate Loader]

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 6, Creating an Experiment [With Plate Loader] for detailed instructions.

Setting Sample Wells

NOTE If you are using Conventional Mode, refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 6, Setting Sample Wells for detailed instructions.

Once the plate protocol is created, the plate window appears. Refer to Figure 6.1.






Figure 6.1 Plate Window [CytoFLEX LX Shown]

	Empty well.
	The well with color is set as sample well, but is not set for Auto Record.
3	The well is set as sample well, and is ready for Auto Record. The number labeled at the right bottom shows the order of auto record.
	Cleaning Agent well.
	Deionized water well.

6-2 D17052AA




- Left-click and drag your mouse to highlight the desired wells or hold the Control key and select each desired well.
- 2 Select or right-click the selected wells and select **Set As Sample Wells**. The Set As Sample Wells window appears.



NOTE Select or right-click and select **Set As Empty Wells** to reset selected wells as empty.

- **3** Enter the name in the Prefix box in the Naming Rules section of the window.
- 4 Select the desired Mix and Backflush duration from the Mix/Backflush section of the window.

- **5** Enter the Group Name in the Name box in the Group section of the window.
- **6** Select the sample well color using the color dropdown under the Group section of the window.
- **7** Select the desired acquisition settings under the Acq. Setting tab.




NOTE Select **Import from File** to import the settings from a FCS file.


NOTE If desired, import saved settings/standardization settings from the catalog.


6-4 D17052AA

8 Select the channels and adjust as needed.



9 Set unmixing under the Unmixing Details tab. Refer to CHAPTER 7, Unmixing for detailed instructions on setting unmixing.






10 Select Events to Record, Time to Record, or Volume to Record under the Stopping Rules tab.

NOTE Beckman Coulter recommends setting an acquisition time limit to stop the acquisition if the event limit cannot be reached.

11 Create the desired name and value under the Custom Metadata tab.



6-6 D17052AA

12 Select ox.

Modifying Well Settings

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 6, Modifying Well Settings for detailed instructions on the following procedures:

- Copying, Cutting, and Pasting Wells
- Applying Existing Well Settings to Additional Wells
- Moving the Location of a Well


Setting the Label

NOTE If you are using Conventional Mode, refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 6, Setting the Label for detailed instructions.


NOTE The Set Label function is not available in spectral unmixing experiments.

NOTE The Set Label function is only available when the tube contains the unmixing matrix. To import unmixing matrix, refer to CHAPTER 7, Importing and Exporting Unmixing Matrix.

- 1 Select **Set Label** in the Settings menu. The Set Label window appears.
- 2 In the Set Label window, modify the label as needed.



3 Select **Apply to.** The Apply Label Setting window appears.



4 Select the tubes to apply the label settings to and select **OK**.

NOTE The Select All checkbox allows you to select all of the tubes listed.

5 Select **Close** to exit the Set Label window.

Creating Plots and Gates

IMPORTANT The maximum number of elements allowed in an experiment is 200. Elements include plots, statistics tables, and gate hierarchy tables.

IMPORTANT The maximum number of gates allowed in an experiment is 200.


NOTE If you are using Conventional Mode, refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 6, Creating Plots and Gates for detailed instructions.

1 Use the plotting controls (see Figure 2.3) in the plot area to create plots and gates and to generate graphs as shown.

Use the icons to generate histograms, dot plots, density plots, pseudo color plots, contour plots, and spectrum plots.

6-8 D17052AA

The experiment uses scatter plots, histograms, polygon gating, four-quadrant gating, and line-segment gating.




- **a.** After selecting a plot, click and drag the mouse to adjust the position and select and drag the sizing handles at the edge of the graph to adjust the size of the graph.
- **b.** Select an axis name to change which channel is displayed. An "A" after the channel name indicates signal pulse area, while an "H" indicates height. The default setting is "A".

[Raw Worksheet Shown]

[Unmixed Worksheet Shown]


CD5 APC-A700-H CD56 PE-Cy5.5-H




NOTE To modify the default settings, select **Options** in the Settings menu. The Options window appears. Select **Plot** on the left side of the Options window. Under the Signal section of the window, change the Main Channel default by selecting the **Height** or **Area**.

NOTE When using both Height and Area signals, ensure the gain setting is set to where the Height signal does not reach its upper range.

c. Signal width can be used as a tool for doublet discrimination and to differentiate somatic cell adhesion. If necessary, select **\text{Acq. Setting...}* to open the Acq. Setting window.




d. Select the **Signal** tab, and select the channel with the required signal width.



6-10 D17052AA

- e. Plot properties can be configured to display axes in Log, Log-Linear, or Linear format.
 - 1) Double-click the plot or right-click the plot and select **Property** from the drop-down menu. The Plot Property screen appears.



- 2) Select whether to display the axes in logarithmic or linear format for both the X-axis and Y-axis. Enter a value for log-linear coefficient if the log-linear view is desired.
- 3) Select Close.


Or

Select the logarithmic axis on the plot. The slider appears. Drag the slider along the axis to change the log-linear coefficient and view events that are not shown, including events with negative values.


NOTE The log-linear slider is also available during data acquisition.

NOTE To reset the axis back to logarithmic, right-click on the axis and select **Property**. Select **X Axis Default** or **Y Axis Default** to reset the axis.

Histogram with logarithmic




Histogram with log-linear




D17052AA


Dot plot with logarithmic



Dot plot with log-linear



- **f.** You can adjust axis ranges using the pan axis display controls located at the top of the screen.
 - Select to zoom-in and define which area of a plot to enlarge. The selected area can be magnified to fill the entire graph. By selecting the zoom-out function, you can click on the graph and restore the plot to its original appearance before magnification.
 - Select to shift the axes. The mouse pointer appears as a hand. It allows you to drag the graph to reveal the axis segment you need.
 - Pan: Modifies the axis display range dimensions when panning both axes. When the pan control is selected, you can right-click the graph and select which axis you need to adjust when dragging. You can also pan directly to the default axis range.



— Single side pan: Modifies the axis display range dimensions when panning one axis.

NOTE Only the low end of the axis can be adjusted by the single side pan tool.

6-12 D17052AA


• Select to scale axis ranges in the plots automatically.


NOTE This feature is only available in the spectral unmixing experiment and in the unmixed worksheet opened in a standard experiment (see Figure 2.2).

- Double-click the border area of the plot to open the Plot Property window, or right-click the plot, then select **Property** to open the same Plot Property window.
- In the Plot Property window, manually enter the minimum and maximum display values for the X- and Y-axes. You can also select **Fit With Sample** to let the software automatically adjust the lower limit according to the signal and perform the corresponding log-linear transformation. The X- and Y-axes **Default** settings are the default parameters. The default parameters are 100-1,000,000.
 - **NOTE** Select **Fit With Sample** to identify the signal's lower limit, adjusting automatically as warranted. Selecting this item is recommended whenever the signal appears to be relatively low.
 - **NOTE** Select **Auto** to automatically set the upper and lower display limits of the axes based on the data already collected.

NOTE Select **Options** in the Settings menu, then select **Plot** to modify the default setting of the axis range under the Axis Default Setting section of the window.

2 To create gates, use the \Box \Box \Box \Box \Box \Box \Box \Box control buttons or right-click the plot and select the gate type required. Gates can be set according to different requirements to differentiate cell populations.






NOTE To add a vertex to a polygon gate:

- 1. Select the gate.
- 2. Hover your cursor over the perimeter of the gate until the cursor changes to the hand icon.
- **3.** Select the desired location for the new gate vertex.

NOTE A newly created gate becomes a subset of the plot where it appears. The relationship between parent and progeny/daughter gates can be changed when a displayed gate is subsequently modified.


The position of the same gate in different sample tubes may vary. To change the position of a gate and apply the change to all sample tubes accordingly, you can right-click the gate and select **Apply to All Tubes**.


You can also apply the change to select tubes by selecting **Apply to Tubes.**



6-14 D17052AA

- **3** Select the gates to display.
 - **a.** Select the heading area of the plot, select the parent population/gate to display in the plot from the drop-down menu. The selected parent gate appears in the heading area of the plot.





NOTE The CytExpert for Spectral software will not list gates which would create circular gating logic.

Figure 6.2 shows all gates defined in the example experiment below. Note that the only gate option in plot 1 of Figure 6.3 is P2 for the following reasons:

- Plot 1 cannot be gated on P1 because P1 is on that plot.
- Plot 1 cannot be gated on P2 because P2 is gated on P1.
- Plot 1 cannot be gated on the P2 OR P1 combo population because the gate logic contains P1.

Figure 6.2 All Gates - Example Experiment






Figure 6.3 Circular Gating Logic - Example Experiment



b. If necessary, you can select the **Combo Population** option from the drop-down menu to create a combination gate, using the Boolean relationships "AND", "OR", and "NOT" to produce a new gate. You can also select the population color or change the gate name.




- "AND" indicates that all selections must be satisfied. For example, "P1 AND P2" means that the data for the newly added gate represent the intersection of P1 and P2.
- "OR" indicates that only one of the selections needs to be satisfied. For example, "P1 OR P2" means that the data for the newly added gate represent the union of P1 and P2.
- "NOT" indicates exclusion from the selection. For example, "NOT P1" means that the data for the newly added gate represent the events that are not part of P1.

6-16 D17052AA

4 Select 🗀 to display the population hierarchy.

The Population Hierarchy function allows you to view how gates rank in relation to one another. To change the display color, double-click the default color and select the desired color from the drop-down color palette. To change the name of each gate, double-click the name of the desired gate. By hovering your mouse pointer over a combo population whose display name has just been changed, you can view its corresponding Boolean logical operation.



Creating and Adjusting Auto Gates

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 6, Creating and Adjusting Auto Gates for detailed instructions on the following procedures:


- Turning Auto Recalculate On/Off
- Adjusting Autogate Movement and Extent

Setting Customized Parameters


Set custom parameter to create fluorescence calculations.

NOTE If you are using Conventional Mode, refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 6, Setting Customized Parameters for detailed instructions.

1 Select **Set Customized Parameter** from the Settings menu. Or, right-click a test tube from the test tube menu and select **Set Customized Parameter**. The Set Customized Parameter window appears.



- 2 Select the Raw/Unmixed tab on the top of the Set Customized Parameter window as needed.
- **3** Enter a name for the parameter in the Name section.
- **4** Select the parameters for calculation from the Parameter dropdown menu.
- 5 Select the equation operations from the Open dropdown menu.
 The new parameter name is displayed in the list of parameters and statistic items.



6-18 D17052AA

Setting Custom Statistics

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the CytoFLEX Platform Instructions for Use manual, CHAPTER 6, Setting Custom Statistics for detailed instructions.

Configuring Acquisition Settings

NOTE Settings can be imported from the Acquisition Settings Catalog. Refer to CHAPTER 6, Importing Instrument Settings.

If unmixing settings are desired, import the unmixing matrix from the Spectral Library or import the unmixing file. Refer to CHAPTER 7, Importing and Exporting Unmixing Matrix. Unmixing matrix can be imported before or after performing the sample acquisition to generate unmixed data.

Changing the Tube Name

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the CytoFLEX Platform Instructions for Use manual, CHAPTER 6, Changing the Tube Name for detailed instructions.

Laser Settings

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the CytoFLEX Platform Instructions for Use manual, CHAPTER 6, Laser Settings for detailed instructions.

NOTE Lasers can only be enabled and disabled in Conventional Mode.

Setting Laser Target Power Settings [CytoFLEX LX Only]

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 6, Setting Laser Target Power Settings [CytoFLEX LX Only] for detailed instructions.

NOTE Lasers can only be enabled and disabled in Conventional Mode.

Adjusting the Gain

While the instrument is in use, the signal value can be increased or decreased by adjusting the instrument's gain configuration.

NOTE If you are using Conventional Mode, refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 6, Adjusting the Gain for detailed instructions.

1 Select ** Acq. Setting... on the left side of the screen. The Acq. Setting window appears.


2 Select the **Gain** tab in the Acq. Setting window.

Select or edit the instrument's default gain settings using one of the following methods:

- Edit the gain settings and select Set As Default to create a new default settings.
- Select **Default** to return to your saved default settings.
- Select **QC** Gain to use the instrument's QC settings.
- Select Assay Gain to use the instrument's assay settings.

NOTE In cases where you do not specify your own default parameters, the Assay settings and default settings are identical.

Adjust the gain setting of each channel under the **Gain** tab in the Acq. Setting window. Raising the gain increases the signal. Lowering the gain reduces the signal.



NOTE Gain adjustments have a predefined range between 1 and 3,000. For fine adjustments, use the text box and manually enter the gain values as needed.

NOTE The * symbol next to the channel indicates that:

- The gain of the channel has been modified.
- The gain of the channel is set to the minimum or maximum value.

NOTE If you adjust the gain of a laser, the gain of all corresponding channels changes simultaneously. If the gain of a channel has been modified, it cannot be modified by adjusting the corresponding laser's gain directly.

6-20 D17052AA

NOTE Optimize the gain settings according to your own experimental goals. The QC Gain values are only for reference.

NOTE Select **Link to Setting from Catalog** to apply the gain settings directly to skip adjusting the gain setting. Refer to Link to Acq. Setting.

Another option is to use the **Gain Control** button on the tool-bar in the graphic control area to adjust the gain values for cell population data to their desired levels, directly on the plots where the data appears during data collection.




4 If necessary, change the coordinate display range and the plot type.


Link to Acq. Setting

This operation will only affect gain settings, not threshold and width settings.

1 Select 🔀 Acq. Setting... on the left side of the screen. The Acq. Setting window appears.




2 Select **Link to Setting from Catalog** in the Acq. Setting window. The Acq. Setting Catalog window appears.



Or

Right-click the desired tube and select **Link to Setting from Catalog** from the drop-down menu. The Acq. Setting Catalog window appears.



6-22 D17052AA

3 Select the desired Item to link and select **OK**. The Information window appears.



NOTE Checking the **Convert Unmixing Matrix by the Gain** checkbox allows the unmixing matrix to automatically adjust by gain settings.

- 4 Select **OK**.
- **5** Select **Close** to exit the Acq. Setting window.


Adjusting the Threshold

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 6, Adjusting the Threshold for detailed instructions.

Setting Collection Conditions

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the CytoFLEX Platform Instructions for Use manual, CHAPTER 6, Setting Collection Conditions for detailed instructions.

NOTE In the CytExpert for Spectral software, the **Events to Record** also includes a drop down menu where you can select to record either raw data or unmixed data.



Setting Plot Display Conditions

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the CytoFLEX Platform Instructions for Use manual, CHAPTER 6, Setting Plot Display Conditions for detailed instructions.

Viewing Saturation [Without Plate Loader]

The saturation warning is an automatic iterative acquisition-and-check process. You can optimize the gain settings according to your own experiment goals to avoid signal saturation.


NOTE The Saturation Warning function is enabled by default. You cannot enable or disable the function during the data acquisition.


NOTE The Saturation Warning function only applies to the fluorescence channels.

NOTE Adjust the gain settings according to your own experimental goals as necessary to avoid signal saturation. Refer to Adjusting the Gain.


6-24 D17052AA

- 1 Confirm that the Saturation Warning function is enabled. If not, select **Saturation Warning** from the Advanced menu to enable it.
 - **NOTE** To disable the Saturation Warning function, select **Saturation Warning** from the Advanced menu The following system prompt appears, and select **Yes.** The saturation warning status displays on the top of the software screen.





- **2** Create a new experiment and new tube. Refer to Creating an Experiment.
- **3** Configure the acquisition settings. Refer to Configuring Acquisition Settings.
- 4 Load samples and record data as needed. Refer to Load Sample and Record Data.
- **5** When the acquired or recorded events reach 2,000, the confirmation window appears.



NOTE The default saturation threshold is 5%.

NOTE The acquisition continues when the confirmation window appears.

6 Select **Ignore** to close the confirmation window.

Or

Select **Stop** to stop the acquisition.


NOTE If you select **Ignore**, the confirmation window will no longer prompt.


7 When the acquisition completes, the Information window appears. Select **ok**.



6-26 D17052AA

f 8 Right-click a tube and select **View Saturation.** The View Saturation window appears.





9 Select Close.

Viewing Saturation [With Plate Loader]

The saturation warning is an automatic iterative acquisition-and-check process. You can optimize the gain settings according to your own experiment goals to avoid signal saturation.


NOTE The Saturation Warning function is enabled by default. You cannot enable or disable the function during the data acquisition.

NOTE The Saturation Warning function only applies to the fluorescence channels.

NOTE Adjust the gain settings according to your own experimental goals as necessary to avoid signal saturation. Refer to Adjusting the Gain.

Confirm that the Saturation Warning function is enabled. If not, select **Saturation Warning** from the Advanced menu to enable it.


NOTE To disable the Saturation Warning function, select **Saturation Warning** from the Advanced menu The following system prompt appears, and select **Yes.** The saturation warning status displays on the top of the software screen.



- **2** Create a new experiment and new tube. Refer to Creating an Experiment.
- **3** Configure the acquisition settings. Refer to Configuring Acquisition Settings.
- 4 Load samples and record data as needed. Refer to Load Sample and Record Data.
- **5** If you are running a single well, proceed to Step 6. If you are running a set of wells, skip to Step 9.

6-28 D17052AA

6 When the acquired or recorded events reach 2,000, the confirmation window appears.



NOTE The default saturation threshold is 5%.

NOTE The acquisition continues when the confirmation window appears.

7 Select **Ignore** to close the confirmation window.

Or

Select **Stop** to stop the acquisition.

NOTE If you select **Ignore**, the confirmation window will no longer prompt until the software completes the current sample tube acquisition.

f 8 When the acquisition completes, the saturation information window appears. Select f o K.



The saturated icon appears in the Tube List next to the tube if the tube data acquired is saturated. Skip to Step 11.



9 When the auto recorded events reach 2,000, the confirmation window prompts.



NOTE The default saturation threshold is 5%.

NOTE The acquisition continues when the confirmation window prompts.

6-30 D17052AA

10 Select Ignore current well or Ignore whole plate to close the confirmation window.


Or

Select **Stop** to stop the acquisition.


NOTE If you select **Ignore whole plate**, the confirmation window will no longer prompt.


NOTE If you select **Ignore current well**, the following confirmation window appears when the current sample tube acquisition completes. The acquisition stops.

- Select Ignore whole plate to continue the acquisition. The confirmation window will no longer prompt.
- Select Resume to continue the acquisition. The confirmation window will prompt when the next tube is saturated.
- Select Stop to stop the acquisition.



11 Right-click a tube and select **View Saturation**. The View Saturation window appears.





12 Select Close.

6-32 D17052AA

Load Sample and Record Data

Before Running Samples

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the CytoFLEX Platform Instructions for Use manual, CHAPTER 6, Before Running Samples for detailed instructions.

NOTE The Detector Configuration function is only available in Conventional Mode.

Sampling and Collecting Data [Without Plate Loader]

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 6, Sampling and Collecting Data [Without Plate Loader] for detailed instructions.

NOTE The Compensation function is only available in Conventional Mode.

Sampling and Collecting Data [With Plate Loader]

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 6, Sampling and Collecting Data [With Plate Loader] for detailed instructions on the following procedures:

- Creating a Heat Map [with Plate Loader]
- Refreshing a Heat Map
- Modifying Existing Heat Map Settings
- Deleting an Existing Heat Map
- Exporting a Heat Map


NOTE In the CytExpert for Spectral software, the Heat Map function is only available when the tube contains unmixed data.

Analyzing and Exporting Data

Select **Raw/Unmixed** located on the bottom of the plot area to open the raw/unmixed worksheet (see Figure 2.1 and Figure 2.2).



- **2** Select the sample tube to be analyzed.
- **3** Establish new gates or adjust the position of existing gates. Refer to Creating Plots and Gates.





NOTE Changing a gate's position does not affect the positions of other gates already established on a given sample tube. Each test tube individually records the positions of its associated gates. If you need to make a change that concerns all the tubes, you must select the gate, then right-click the correctly positioned gate and select **Apply to All Tubes**.

- 4 Select . The Gate Hierarchy screen appears.
- 5 Check the relationship between the parent and daughter gates in the Gate Hierarchy window.

NOTE Newly added gates become subsets of populations displayed in plots with existing gates. Name and display color can be modified. Right-click directly on a gate plot to change the name and color.


NOTE Select **No Color** to leave the gated events uncolored while retaining the color of the parent populations. By default, the populations defined by a vertical gate, hinged gate, or four-quadrant gate are uncolored.






6-34 D17052AA


6 Right-click the plot and select **Bring population to front** to make the display color of the specified gate appear in front of all other colors, or select **Send population to back** to hide the display color of the specified gate behind all other colors.





- 7 Select in the plot area to generate a statistical table.
- **8** Right-click the table and select **Statistics Setting** to modify the settings of the statistics display parameters. The Statistics Setting window appears.



The Statistics Setting window allows you to change the display of the header, statistical elements and cell populations included.









6-36 D17052AA

The final generated plots appear as below.



9 Right-click a plot and select **Export to Clipboard** or **Export to Graphic File** from the drop-down menu to select an image to export.




• **Export to Clipboard** copies the plot to the clipboard, allowing you to paste it directly into documents in common file formats.

NOTE Multiple plots can only be copied and pasted into Microsoft Word. If a single plot is copied, this can be pasted into both Microsoft Word or Microsoft PowerPoint.

• **Export to Graphic File** saves the plot as an image file.

NOTE Export to Graphic File can export plots in two selectable file formats. BMP bitmap format and EMF vector format.

10 To export statistics, right-click a statistical table to select any one of the available export options.



- **Export to CSV File** exports individual tube statistics as a single CSV file.
- Export All Samples to CSV File exports all tube statistics as a single CSV file.
- **Export to Clipboard** copies the statistics of an individual sample to the clipboard, allowing you to paste them directly into a Microsoft Excel file or other file formats.
- Export All Samples to Clipboard assembles the statistics for all the sample tubes of an experiment and copies them together to the clipboard. From there they can be pasted as a group into a Microsoft Excel file or other file formats.
- **Copy** converts a statistical table into an image format that can be pasted into documents.
- **11** Export the FCS file if necessary. Refer to Exporting FCS Files.

NOTE Ensure that any storage devices used with the instrument are free from viruses. To guard against data loss, Beckman Coulter recommends backing up data on a frequent and regular basis. Beckman Coulter is not liable for any loss of data resulting from computer viruses or damage to hardware.


6-38 D17052AA

Exporting FCS Files

Exporting Single Tube Files

NOTE If you are using Conventional Mode, refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 6, Exporting Single Tube Files for detailed instructions.


[Standard Experiment]: Right-click the desired tube from the test tube section of the screen and select Export FCS File (Raw) or Export FCS File (Unmixed). The Export FCS File window appears.



2 [Spectral Unmixing Experiment]: Right-click the desired tube and select Export FCS File from the drop-down menu. The Export FCS File window appears.



3 Select the population from the Population drop-down menu.



- 4 Select Area and/or Height.
- **5** Select the FCS format next to File Version.

NOTE The default setting is FCS 3.1.


6 Select **OK**. Specify the file path and file name, and select **Save**.

6-40 D17052AA

Exporting Multiple FCS Files

NOTE This function is not available in spectral unmixing experiments.

Select File > Export FCS File (Raw) or Export FCS File (Unmixed). The Export FCS File window appears.




2 Select the desired tubes to export.

NOTE The **Select All** checkbox allows you to select all of the tubes listed.

3 Select Next.

4 Select the population from the Population drop-down menu.




- 5 Select Area and/or Height.
- **6** Select the FCS format next to File Version.

NOTE The default setting is FCS 3.1.

6-42 D17052AA

7 Select ____ to specify the file path, and select **Next**. The following massage appears: *Export was successful*.



8 Select **Close** to exit the Export FCS File window.

Or

Select **Back** and repeat Steps 1-7 to continue exporting FCS files.

Exporting Plots or the Statistics Table of Multiple Tubes as Picture Files

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 6, Exporting Plots or the Statistics Table of Multiple Tubes as Picture Files for detailed instructions.


Importing and Exporting Instrument Settings

The CytExpert for Spectral software supports importing and exporting instrument settings to facilitate the experiment process. Only instrument settings identical to the current configuration can be imported with current detector settings.

Select Acq. Setting... to edit gain, threshold, and width. These can be imported from an experiment file or from a catalog of instrument settings.

NOTE If you are using Conventional Mode, refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 6, Importing and Exporting Instrument Settings for detailed instructions.

Importing Instrument Settings




NOTE Instrument settings can only be imported into tubes where data has not yet been recorded.

6-44 D17052AA


Select Import From File, locate the file with the required instrument settings, or select Import From Catalog to import the instrument settings.

[CytoFLEX LX Shown]



3 Select Close.

Exporting Instrument Settings



2 Select **Export to File** to export a current set of instrument settings. Specify the file path and file name and select **Save**.

NOTE The generated file ends in *.acq.

Or

Select **Export to Catalog**, give a name to the settings to be exported, and export the file to the software's Acquisition Setting Catalog, then select **OK**.

3 Select Close.

Importing and Exporting Unmixing Settings

The software supports unrestricted importing and exporting of unmixed data, regardless of whether the sample tube data has already been acquired. Imported unmixing values only cover channels identical with the current instrument configuration. The software automatically adjusts unmixing values according to differences in the gain level.

Printing Graphics

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the CytoFLEX Platform Instructions for Use manual, CHAPTER 6, Printing Graphics for detailed instructions.

Saving the Experiment

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the CytoFLEX Platform Instructions for Use manual, CHAPTER 6, Saving the Experiment for detailed instructions.

Concluding the Experiment

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 6, Saving the Experiment for detailed instructions.

6-46 D17052AA

Overview

Spectral unmixing is used to identify the fluorescence signal for each fluorospheres used in a given experiment. The spectral unmixing process is to disassemble sample raw data into single dye channels according to their spectral signatures. This chapter provides procedures for using the CytoFLEX mosaic Spectral Detection Module.

NOTE To switch between Spectral Mode and Conventional Mode, refer to CHAPTER 12, Switching between Spectral Mode and Conventional Mode.

Workflow:

Prepare spectral unmixing =	Create spectral unmixing experiment →	Acquire and record data	Calculate spectral unmixing	→	Apply spectral unmixing to sample for data analysis
-----------------------------	--	-------------------------	-----------------------------------	----------	---

This chapter contains information on:

- Preparing the Spectral Unmixing Sample
- Creating a Spectral Unmixing Experiment
- Loading Sample and Recording Data
- Calculating Unmixing Values
- Confirming Unmixing Results
- Saving the Spectral Unmixing Experiment
- Adjusting Spectral Unmixing
- Fine Adjustment

Preparing the Spectral Unmixing Sample

For the negative control sample and single positive control sample, you can use blood, cell lines, or dedicated compensation beads such as VersaComp Antibody Capture Beads. For details, refer to the appropriate reagent instructions for use. The negative control tube is used to determine the autofluorescence of the sample.

Preparing the Spectral Unmixing Sample [Without Plate Loader]

To perform a spectral unmixing experiment, prepare:

- A single positive control tube for each color present in the panel
- A negative control tube

NOTE Users can add more than one single positive control tube or negative control tube as necessary.

Preparing the Spectral Unmixing Sample [With Plate Loader]

To perform a spectral unmixing experiment, prepare:

- A single positive control well for each color
- A negative control well

NOTE Users can add more than one single positive control tube or negative control tube as necessary.


Creating a Spectral Unmixing Experiment

Creating a Spectral Unmixing Experiment [Without Plate Loader]

1 Select **File > New Spectral Unmixing** or **New Spectral Unmixing** on the start page to create a new spectral unmixing experiment.

NOTE The file name of the newly created spectral unmixing experiment has a ".xitsp" suffix.

7-2 D17052AA



2 Navigate to the desired file path and select **Save**. The Spectral Unmixing Setup window appears.

NOTE A **Control** is the reference material used for creating the unmixing matrix.



Risk of erroneous results. Select a negative control tube, according to which the fluorescence background will be set. If there is not a negative control tube available, then the single color tube must have a negative population.

- **3** Add a negative control tube.
 - a. Select Add.
 - **b.** Edit the Tube name and Remark as needed.
 - c. Select the desired control type from the **Set as Negative Control for** drop-down menu.

NOTE The default control type is **Sample and Single Stained**. If you select **Sample and Single Stained**, the tube will appear in the unmixing matrix and be used for unmixing calculation. If you select **Single Stained Only**, the tube will not appear in the unmixing matrix and not be used for unmixing calculation.

NOTE To add more negative control tubes, repeat Steps a.-c. A maximum of 10 negative control tubes can be added for the spectral unmixing experiments.

NOTE To delete a tube, select the desired tube and select **Delete**.

D17052AA

- **4** Add a single stained control tube.
 - a. Select Add.
 - **b.** Select the desired fluorescent tag from the **Tube** drop-down menu.

NOTE To manage the fluorescent tag, refer to CHAPTER 2, Fluorescent Tag Library.

c. Select the desired control type from the **Negative Control Type** drop-down menu.


NOTE The default control type is **Internal Negative Control**. The unstained control tube can be selected if needed.

- **d.** Edit the Label and Remark as needed.
 - **NOTE** Label information is retained in the unmixing matrix. Label information will be applied to the corresponding channels when the unmixing matrix is imported for the use in other experiments.

NOTE To add more single stained control tubes, repeat Steps a.-d. A maximum of 60 single stained control tubes can be added for the spectral unmixing experiments.

NOTE To delete a tube, select the checkbox to the left of the desired tube and select **Delete**.

NOTE To add multiple single stained control tubes, select **Add Multiple Controls**. The Add Single Stained Controls window appears.




Select the **Use** checkbox to add the desired tube. Select the desired control type from the Negative Control Type drop-down menu. Edit the Label and Remark as needed. Select **OK**.

NOTE The **Select All** checkbox allows you to select all of the tubes listed.

7-4 D17052AA

5 Select **OK**.

After confirmation, the software automatically generates the following unmixing spectral experiment.




6 Confirm that the gating and peak channel are appropriate. Use the pan tool to adjust the axis scale so that the sample signal appears in a suitable position. Adjust the gate so that it encloses the target cell population. Refer to CHAPTER 6, Creating Plots and Gates.

Creating a Spectral Unmixing Experiment [With Plate Loader]

1 Select **File** > **New Spectral Unmixing** or **New Spectral Unmixing** on the start page to create a new spectral unmixing experiment.

NOTE The file name of the newly created spectral unmixing experiment has a ".xitsp" suffix.



2 Navigate to the desired file path and select **Save**. The Spectral Unmixing Setup window appears.

NOTE A **Control** is the reference material used for creating the unmixing matrix.

∴ CAUTION

Risk of erroneous results. Select a negative control tube, according to which the fluorescence background will be set. If there is not a negative control tube available, then the single color tube must have a negative population.

- **3** Add a negative control tube.
 - a. Select Add.
 - **b.** Edit the Tube name and Remark as needed.
 - **c.** Select the desired control type from the **Set as Negative Control for** drop-down menu.
 - **NOTE** The default control type is **Sample and Single Stained**. If you select **Sample and Single Stained**, the tube will appear in the unmixing matrix and be used for unmixing calculation. If you select **Single Stained Only**, the tube will not appear in the unmixing matrix and not be used for unmixing calculation.
 - **d.** Select the plate type, sampling sequence, mix and backflush settings for the tube.
 - e. Select Set As Sample Well to assign well location.
 - **NOTE** To add more negative control tubes, repeat Steps a.-e. A maximum of 10 negative control tubes can be added for the spectral unmixing experiments.

7-6 D17052AA

- **4** Add a single stained control tube.
 - a. Select Add.
 - **b.** Select the desired fluorescent tag from the **Tube** drop-down menu.

NOTE To manage the fluorescent tag, refer to CHAPTER 2, Fluorescent Tag Library.

c. Select the desired control type from the **Negative Control Type** drop-down menu.

NOTE The default control type is **Internal Negative Control**. The unstained control tube can be selected if needed.

d. Edit the Label and Remark as needed.


NOTE Label information is retained in the unmixing matrix. Label information will be applied to the corresponding channels when the unmixing matrix is imported for the use in other experiments.

- **e.** Select the plate type, sampling sequence, mix and backflush settings for the tube.
- f. Select Set As Sample Well to assign well location.

NOTE To add more single stained control tubes, repeat Steps a.-f. A maximum of 60 single stained control tubes can be added for the spectral unmixing experiments.

NOTE To delete a tube, select the checkbox to the left of the desired tube and select **Delete**.


NOTE To add multiple single stained control tubes, select **Add Multiple Controls**. The Add Single Stained Controls window appears.



Select the **Use** checkbox to add the desired tube. Select the desired control type from the Negative Control Type drop-down menu. Edit the Label and Remark as needed. Select **OK**.

NOTE The **Select All** checkbox allows you to select all of the tubes listed.

5 If you need to add multiple single stained control tubes, select **Add Multiple Controls**. The Add Single Stained Controls window appears.



Select the **Use** checkbox to add the desired tube. Select the desired control type from the Negative Control Type drop-down menu. Edit the Label and Remark as needed. Select **OK**.


NOTE The Select All checkbox allows you to select all of the tubes listed

Select the plate type, sampling sequence, mix and backflush settings for all the tubes added.

7-8 D17052AA

6 Select ok.

After confirmation, the software automatically generates the following unmixing spectral experiment.



- Confirm that the gating and peak channel are appropriate. Use the pan tool to adjust the axis scale so that the sample signal appears in a suitable position. Adjust the gate so that it encloses the target cell population. Refer to CHAPTER 6, Creating Plots and Gates.
- **8** Before acquiring data, ensure the plate has been loaded properly. Data can be acquired as a single well or as a set of wells. Refer to CHAPTER 6, Sampling and Collecting Data [With Plate Loader].


Loading Sample and Recording Data


Refer to CHAPTER 6, Load Sample and Record Data for detailed instructions on loading sample and recording data.

NOTE If you change the acquisition settings of a negative control tube, the new acquisition settings will apply to all relevant single stained control tubes without data.


Calculating Unmixing Values

1 Select or select **Spectral Unmixing** from the Settings menu to calculate the spectral unmixing values. The following system prompt appears.





2 The Unmixing Result window appears.




7-10 D17052AA

3 Select Save As to export the unmixing matrix as a .unmix file and specify where to save it.

NOTE The unmixing matrix can also be imported for use in other experiments.

4 Select **Save to Library** to save the unmixing settings in the spectral library. The Save to Spectral Library window appears.



5 Specify the Keyword and select **OK**.


NOTE The settings stored in the spectral library are specific to the detector configuration. The spectral library can only be applied when the detector configuration are the same.

6 Select Close.

NOTE At any time, saved spectral unmixing experiments can be reopened and the unmixing values can be recalculated.

Confirming Unmixing Results

1 Select a sample tube from the tube list on the left side of the screen and the unmixing plots appear on the plot area.




NOTE The unmixing results appears on the right of the Name column. The A indicates that there is an unmixing gap. The indicates that the unmixing gap should be recalculated and refreshed. Refer to Fine Adjustment.

NOTE The software uses the positive-negative method to calculate the unmixing gap.

$$GAP_{(y \ to \ x)} = \frac{Y_{(Positive)} - Y_{(Negative)}}{X_{(Positive)} - X_{(Negative)}}$$

7-12 D17052AA

2 Select or select **Spectral Unmixing** from the Settings menu. The Unmixing Result window appears.



NOTE The **Change** button in the lower left of the Unmixing Result window allows you to change the unmixing algorithm. The two options currently offered are **LSM Unmixing** and **Poisson Unmixing**.

NOTE The **Normalization Setting** drop-down menu in the top right of the Unmixing Result window allows you to select the tube gain settings to apply to the unmixing matrix.

NOTE The **Signal** drop-down menu in the top right of the Unmixing Result window allows you to select which signal type to use, **Height** or **Area**.

NOTE In the Unmixing Details window:

- The spectral signature chart displays the unmixing matrix values based on the peak channel. Use the **Legends** checkbox to choose which unmixing control's spectral signature to show.
- The similarity index is a number between 0-1. The CytExpert for Spectral software displays the value in orange as a visual clue. The darker the shade, the higher the similarity.
- The complexity index is a measure of how distinguishable a collection of spectral signatures are from each other when unmixed together. The complexity index is used to assess the stability of a linear system.

NOTE In the Unmixing Matrix window, the CytExpert for Spectral software displays the unmixing values in either blue, or black, or orange, or red as a visual clue.

- Blue denotes the value falls below 0, which requires special attention.
- Black denotes the value falls within 0 100, which is acceptable.

NOTE In the Fine Adjustment Matrix window:

The Use checkbox applies the fine adjustment to the selected sample.

NOTE The Fine Adjustment range is from -500 to 500 depending on the spillover percentage calculated. The CytExpert for Spectral software displays the value in either black, or orange, or red as a visual clue.

- Black denotes the value falls within 0 99.99, which is acceptable.
- Orange denotes the value falls within 100 399.99, or -99.99 to -0.01, which requires special attention.
- Red denotes the value is within 400 500 or -500 to -100, which is unacceptable.
- **3** If there is an unmixing gap, follow the procedure below:
 - **a.** Verify whether the allocated single stained control tube was prepared as required and correctly positioned.
 - **b.** Verify that the gating and population hierarchy in each tube were correct.

NOTE The minimum recommended number of events in the specified gate is 50.

- **c.** Verify that the peak channel selected in each tube was appropriate.
- **d.** Repeat Calculating Unmixing Values to recalculate and export the unmixing results.
- **4** Optional: Use the Fine Adjustment function to adjust the unmixing gaps as needed. Refer to Fine Adjustment.

Saving the Spectral Unmixing Experiment

Selecting **Save** in the File menu allows you to save the experiment.

Selecting **Save As** and saving the experiment under a different name allows you to create a backup.

At any time, saved spectral unmixing experiments can be reopened and the unmixing values recalculated.

Concluding the Experiment

Conclude the experiment as follows:

- Select **Standby** to return the instrument to the standby state.
- Select **File > Close Experiment** to clear the experiment and return to the Start Page.

NOTE If changes were made to the experiment, the software prompts you to save the latest changes in the experiment before returning to the Start Page.

7-14 D17052AA


• Shut down the system. Refer to CHAPTER 9, Daily Shutdown.

Adjusting Spectral Unmixing

Importing and Exporting Unmixing Matrix

Importing Unmixing Settings from Unmixing Matrix Files

- 1 Create a new experiment or open a saved experiment. Refer to CHAPTER 6, Creating an Experiment.
- 2 Select the desired sample tube and select or select **Settings** > **Unmixing Matrix**. The Unmixing Matrix window appears.




- 3 Select **Import** and locate the path where unmixing matrix files are saved. Select the corresponding unmixing matrix file (.unmix) to import the unmixing values.
 - **NOTE** To import unmixing values from the spectral library, refer to Importing Unmixing Settings from the Spectral Library.

Both methods allow you to apply the gain independent unmixing capability. For example, importing the unmixing values with or without recalculating the unmixing matrix based on the gain settings.

NOTE Fine adjustment matrix is included in the unmixing matrix, and can be imported with the unmixing matrix files (.unmix).

4 After opening the desired unmixing file, the Import Unmixing Matrix window appears.



Select one of the following:


- Import unmixing matrix and convert based on the current gain.
- Import unmixing matrix only.
- Import unmixing matrix and gain.

NOTE

- If the tube does not have any data when importing unmixing values calculated from other
 instrument settings, the software prompts you to select whether the gain settings must be
 imported as well. Select Yes to import fluorescence channel gains settings along with the rest
 of the data. Select No to allow the CytExpert for Spectral software to adjust the unmixing matrix
 values based on the current gain settings.
- If the tube does have data when importing unmixing values from other instrument settings, the software prompts you to select whether the unmixing values are adjusted based on the current gain settings.
- It is important to note that automatic adjustments to unmixing values calculated from other
 instrument gain settings could result in incorrect unmixing. Always review the data after
 importing unmixing values to ensure the sample is unmixed properly.

7-16 D17052AA

5 Select **OK**. The imported unmixing matrix appears in the Unmixing Matrix window.





- **6** If necessary, select Apply to... to apply the unmixing values to the selected test tubes.
- 7 Select Close to exit the Unmixing Matrix window.

Importing Unmixing Settings from the Spectral Library

Create a new experiment or open a saved experiment. Refer to CHAPTER 6, Creating an Experiment.


2 Select the desired sample tube and select or select **Settings** > **Unmixing Matrix**. The Unmixing Matrix window appears.





7-18 D17052AA

- 4 In the Keyword/Controls column, the corresponding unmixing values can be selected for each channel. The unmixing values of the same keyword/controls can also be selected using the drop-down menus in the Keyword/Controls column.
- **5** Select **OK**. The imported unmixing matrix appears in the Unmixing Matrix window.




- **6** If necessary, select Apply to... to apply the unmixing values to the selected test tubes.
- 7 Select Close.

Exporting Unmixing Settings to Unmixing Matrix Files


Standard Experiment]: Select the desired tube that contains the unmixing matrix and select



or select **Settings > Unmixing Matrix**. The Unmixing Matrix window appears.




[Spectral Unmixing Experiment]: Select the desired tube and select or select Settings > Unmixing Matrix. The Unmixing Result window appears.



7-20 D17052AA

2 [Standard Experiment]: Select Export. Specify a path and filename for the unmixing file you are saving and select Save.

[Spectral Unmixing Experiment]: Select Save as. Specify a path and filename for the unmixing file you are saving and select Save.




NOTE The generated file ends in .unmix.


3 Select Close.

Exporting Unmixing Settings to Unmixing Matrix Files from the Spectral Library

1 Select **Settings > Spectral Library**. The Spectral Library window appears.



- 2 In the Keyword/Controls column, select the desired unmixing values to export.
- **3** Select **Save as**. The Select algorithm window appears.



Select one of the following:

- LSM Unmixing
- Poisson Unmixing
- 4 Select OK.

7-22 D17052AA

5 Specify a path and filename for the unmixing matrix file you are saving, and select **Save**.

NOTE The generated file ends in *.unmix.

6 Select Close.

Managing Spectral Library

Unmixing values can be managed in the Spectral Library.


Importing Spectral Library Files

1 Select **Settings** > **Spectral Library**. The Spectral Library window appears.



D17052AA 7-23

2 Select **Import** and locate the path where spectral unmixing settings are saved. Select the corresponding spectral library file to import.



- **3** Select **Open** to import.
- 4 Select Close.


7-24 D17052AA

Exporting Spectral Library Files

l Select **Settings > Spectral Library**. The Spectral Library window appears.



2 Select **Export** to specify a path and filename for the unmixing file you are saving.



D17052AA 7-25


3 Select Save.

NOTE The generated file ends in .cfg.

4 Select Close.

Deleting Spectral Library Files


1 Select **Settings** > **Spectral Library**. The Spectral Library window appears.



2 In the Keyword/Controls column, the corresponding unmixing values can be selected for each channel. The unmixing values of the same keyword/controls can also be selected using the drop-down menus in the Keyword/Controls column.

7-26 D17052AA

3 Select **Delete**. The following message appears:



Fine Adjustment

The fine adjustment can be adjusted in an experiment in three ways: auto fine adjustment, manual fine adjustment, and fine adjustment matrix.

Auto Fine Adjusting

Select the populations that need to be adjusted in the bivariate plot. Select from the graphic control area to use auto fine adjustment.

NOTE Select **Refresh All Unmixing Gaps** to recalculate and refresh all unmixing gaps after using the fine adjustment.

NOTE The auto fine adjustment function is only available in the spectral unmixing experiment.

Manual Fine Adjustment

Select the populations that need to be adjusted in the bivariate plot. Select from the graphic control area, then click and drag the mouse pointer inside the plot to adjust. In a standard experiment, you can drag the mouse pointer vertically and/or horizontally. In a spectral unmixing experiment, you can only drag the mouse pointer vertically.

NOTE Select **Refresh All Unmixing Gaps** to recalculate and refresh all unmixing gaps after using the fine adjustment.


D17052AA

Fine Adjustment Matrix


1 [Standard Experiment]: Select the desired tube that contains the unmixing matrix and select



or select **Settings** > **Unmixing Matrix**. The Unmixing Matrix window appears.




[Spectral Unmixing Experiment]: Select the desired tube and select or select Settings > Unmixing Matrix. The Unmixing Result window appears.



7-28 D17052AA

2 Select the **Fine Adjustment Matrix** tab, and adjust the fine adjustment value as needed.



- 3 Select Close.
- 4 Select Refresh All Unmixing Gaps to recalculate and refresh all unmixing gaps.

D17052AA 7-29

UnmixingFine Adjustment


7-30 D17052AA

CHAPTER 8 Data Review

Overview

This chapter discusses how to use the Analysis screen to analyze data. Data can be analyzed using any computer equipped with the CytExpert for Spectral software. No online connection is required.

Workflow:



This chapter contains information on:

- Copying Experiments and Importing Data
- Setting the Plots and Statistics
- Calculating Sample Volume and Concentration
- Adjusting Unmixing Settings
- Exporting Results

Copying Experiments and Importing Data

Copying a Previously Acquired Experiment

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 8, Copying a Previously Acquired Experiment.

Importing Previously Acquired Data

The CytExpert for Spectral software can import and analyze compatible FCS data files acquired by other CytoFLEX flow cytometers.

NOTE If you are using Conventional Mode, refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 8, Importing Previously Acquired Data for detailed instructions.

1 Create a new experiment or open a saved experiment. Refer to CHAPTER 6, Creating an Experiment.

D17052AA 8-1

2 [Standard Experiment]: Select File > Import FCS File, specify the file path and select Open.
Imported data files appear in the Tube screen.

The symbol in front of each data tube indicates that the data tube is an imported data file. Imported data files are copied and saved in the folder where the current experiment data files are saved.

3 [Spectral Unmixing Experiment]: Right-click the desired tube and select Import FCS Files from the drop-down menu. Specify the file path and select Open.

Imported data files appear in the Tube screen.

The **(F)** symbol in front of each data tube indicates that the data tube is an imported data file. Imported data files are copied and saved in the folder where the current experiment data files are saved.

Setting the Plots and Statistics

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 8, Setting the Plots and Statistics for detailed instructions on the following procedures:

- Opening the Analysis Screen
- Creating Histogram and Dot Plot Overlays

For configuring gates and generating statistics, refer to CHAPTER 6, Data Acquisition and Sample Analysis.

Calculating Sample Volume and Concentration

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the CytoFLEX Platform Instructions for Use manual, CHAPTER 8, Calculating Sample Volume and Concentration for detailed instructions.

Adjusting Unmixing Settings

Data unmixing can be carried out at any time. You can select the desired tube in the tube list on the

left side of the screen and select or select **Settings > Unmixing Matrix**. Refer to CHAPTER 7, Adjusting Spectral Unmixing for detailed instructions on adjusting unmixing settings.

8-2 D17052AA

Exporting Results

Refer to CHAPTER 6, Data Acquisition and Sample Analysis.

D17052AA 8-3

Data Review Exporting Results


8-4 D17052AA

CHAPTER 9 Daily Shutdown

Overview

This chapter provides procedures for shutting down the CytoFLEX instrument equipped with the CytoFLEX mosaic Spectral Detection Module.

Workflow:



This chapter contains information on:

- Preparing the Cleaning Solution
- Running Daily Clean
- Shutting Down the Instrument
- Auto Shutdown [CytoFLEX LX Only]

Preparing the Cleaning Solution

Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 9, Preparing the Cleaning Solution for detailed instructions.

Running Daily Clean

Refer to the CytoFLEX Platform Instructions for Use manual, CHAPTER 9, Running Daily Clean for detailed instructions.

Shutting Down the Instrument

- 1 Shut down the Cytometer. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 9, Shutting Down the Instrument for detailed instructions on shutting down the Cytometer.
- **2** Long press the soft button on the front of the Cytometer mosaic Spectral Detection Module to turn off the CytoFLEX mosaic Spectral Detection Module.

D17052AA 9-1

3 Optional: Turn off the main power switch on the left of the Cytometer mosaic Spectral Detection Module.

Auto Shutdown [CytoFLEX LX Only]

You can set up the system to automatically shutdown the Cytometer.

To schedule an auto shutdown after acquisition, refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 6, Creating an Experiment.

To schedule an auto shutdown during Daily Clean, refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 11, Daily Clean [With Plate Loader].

9-2 D17052AA

CHAPTER 10 Troubleshooting

Overview

IMPORTANT In addition to the information stated, never disassemble the instrument or have it repaired by unauthorized personnel. Beckman Coulter bears no responsibility for any problems arising from the unauthorized repair of the instrument.

This chapter introduces solutions to common problems. If there is a problem, follow the information in this chapter to carry out self-inspection. If the problem cannot be resolved, contact us.

This chapter contains information on:

- Precautions/Hazards
- Disposal Precaution
- Troubleshooting Table
- Backup and Restore

Precautions/Hazards

Laser Related Hazards

Beckman Coulter design and manufacture of the instrument complies with the requirements governing the use and application of a laser specified in regulatory documents issued by the:

- U.S. Department of Health and Human Services
- Center for Devices and Radiological Health (CDRH)
- International Electrotechnical Commission (IEC)

In compliance with these regulatory documents, every measure has been taken to ensure the health and safety of users and laboratory personnel from the possible dangers of laser use.

Use the instrument according to the information in the manuals.

Use controls or adjustments or performance of procedures other than those specified herein might result in hazardous radiation exposure.

To ensure your safety, the Cytometer lasers are covered with protective shields. Do not remove these shields.

No user-serviceable assemblies are accessible. Do not attempt to remove the laser or open it. The instrument has components that are dangerous to the operator. If any attempt has been made to

D17052AA 10-1

defeat a safety feature, or if the instrument fails to perform as described in its manuals, disconnect the power and contact us.

Laser Beam Hazards

The CytoFLEX Platform flow cytometer can be configured with up to 6 solid-state diode lasers that are capable of producing laser light at the following levels:

- 355 nm, 20 mW solid-state diode laser
- 375 nm, 60 mW solid-state diode laser
- 405 nm, 80 mW solid-state diode laser
- 488 nm, 50 mW solid-state diode laser
- 561 nm, 30 mW solid-state diode laser
- 638 nm, 50 mW solid-state diode laser
- 808 nm, 60 mW solid-state diode laser

A laser beam is a unique light source that shows characteristics different from conventional light sources. The safe use of the laser depends upon familiarity with the instrument and the properties of coherent, intense beams of light.



Risk of personal injury. The laser beam can cause eye damage if viewed either directly or indirectly from reflective surfaces (such as a mirror or shiny metallic surfaces). To prevent eye damage, avoid direct exposure to the laser beam. Do not view it directly or with optical instruments.

Indirect contact with the laser beam from reflective surfaces (such as jewelry or a screwdriver) is called specular reflection and might also cause damage.

For these reasons, it is important to:

- Limit access to the Cytometer to trained and experienced personnel.
- Never attempt to remove a shield housing a laser and defeat interlocks.
- Never remove a warning label.
- Contact us if a label is missing or unclear.

Laser Warning Labels



Risk of personal injury from radiation exposure. Never remove the shield surrounding a laser. Never remove covers and defeat interlocks.

CDRH-approved and IEC compliant labels are also placed near or on those covers that when removed might expose laser radiation. If necessary, a cover with a CDRH-approved or IEC compliant label must be removed by a qualified Beckman Coulter Representative only.

10-2 D17052AA

Refer to the following figures for the locations of the CDRH-approved and IEC compliant labels:

See Figure 10.1 for the Laser Warning Label on the Cytometer optical bench.

See Figure 10.2 and Figure 10.3 for the Laser Warning Label on the optical bench (located Inside the Cytometer).

See Figure 10.4 for the Laser Warning Label on the 355 nm Laser.

See Figure 10.5 and Figure 10.6 for the Laser Warning Labels on the Cytometer Back Cover.

See Figure 10.7 for the Laser Warning Label on the CytoFLEX mosaic Spectral Detection Module Back Cover.

The laser product is classified as CLASS 1 when all protective measures are in place. This product complies with 21 CFR Parts 1040.10 and 1040.11 as well as EN60825-1. See Figure 10.1.

Figure 10.1 Laser Warning Label on the Laser Optical Bench [CytoFLEX LX]





Figure 10.2 Laser Warning Label within the Optical Bench (Located Inside the Cytometer) [CytoFLEX S]




D17052AA 10-3

DANGER - LASER RADIATION
AVOID EYE OR SKIN EXPOSURE TO
DIRECT OR SCATTERED RADIATION
CLASS 4 LASER PRODUCT
危险 激光辐射
避免眼睛或皮肤暴露于
直接或散射的辐射中
4类激光产品

Figure 10.3 Laser Warning Label within the Optical Bench (Located Inside the Cytometer) [CytoFLEX LX]


Figure 10.4 Laser Warning Label on the 355 nm Laser [CytoFLEX LX]



10-4 D17052AA

Figure 10.5 Laser Warning Labels on the Cytometer Back Cover [CytoFLEX S]

Figure 10.6 Laser Warning Labels on the Cytometer Back Cover [CytoFLEX LX]



D17052AA 10-5

Figure 10.7 Laser Warning Labels on the CytoFLEX mosaic Spectral Detection Module Back Cover

Disposal Precaution





Risk of biohazardous contamination if you have skin contact with the waste container, its contents, and its associated tubing. The waste container and its associated tubing could contain residual biological material and must be handled with care. Clean up spills immediately. Dispose of the contents of the waste container in accordance with your local regulations and acceptable laboratory procedures.

Use universal precautions when working with pathogenic materials. Means must be available to decontaminate the instrument and to dispose of biohazardous waste.

Troubleshooting Table

Table 10.1 lists problems that you could encounter while running the CytoFLEX flow cytometer with the CytoFLEX mosaic Spectral Detection Module, the probable causes of each problem, and the corrective actions. These problems are listed alphabetically in the Index, under the primary entry "troubleshooting".

Refer to the CytoFLEX Platform Instructions for Use manual, CHAPTER 10, Troubleshooting Table for detailed instructions on the Cytometer troubleshooting.

10-6 D17052AA

 Table 10.1
 Troubleshooting [With CytoFLEX mosaic Spectral Detection Module]

Problem	Probable Cause	Corrective Action	
The CytoFLEX mosaic Spectral Detection Module cannot be turned on.	The CytoFLEX mosaic Spectral Detection Module is turned off in the Cytometer menu. [CytoFLEX mosaic 88] The power switch is in the off position and the Turn On selection will not function in the Cytometer menu. [CytoFLEX mosaic 88] The power cable is not securely connected. The fuse is blown.	 Ensure the power switch is in the on position on the back of the CytoFLEX mosaic Spectral Detection Module. Select Turn On in the Cytometer menu. [CytoFLEX mosaic 88] Ensure that the power cable is securely connected to the back of the CytExpert mosaic Spectral Detection Module. If the problem persists, contact us. 	
The CytExpert for Spectral software is not responsive after the software has been launched.	The digital signature validation is timing out and failing due to a malfunction in the DNS server or firewall.	 Type Registry Editor from the Windows search bar to locate the Register Editor. Select Register Editor to open the Register Editor window. Select HKEY_CURRENT_USER > SOFTWARE > Microsoft > Windows > CurrentVersion > WinTrust > Trust Providers > Software Publishing. Double click State to open the Edit DWORD Value window. Change the Value data to 23e00. Select OK. Close the Registry Editor window. NOTE Problems resulting from incorrect Registry Editor settings may require the 	
The connection indicator light in the lower right corner of the software screen is gray and displays <i>Disconnected</i> .	Data connection error FAPD data connection error The CytoFLEX mosaic Spectral Detection Module is not turned on. The CytoFLEX mosaic Spectral Detection Module's power cable is disconnected.	 operating system reinstallation. Ensure that the network cable is securely connected to the back of the Cytometer and the back of the Workstation. Refer to Figure 1.3 or Figure 1.6. Restart the software. Restart the Workstation. Refer to CHAPTER 4, Initializing the Instrument. Turn on the CytoFLEX mosaic Spectral Detection Module using the power switch on the back of the CytoFLEX mosaic Spectral Detection Module. Verify that the power cable is securely connected to the back of the CytoFLEX mosaic Spectral Detection Module. If the problem persists, contact us. 	

D17052AA 10-7

Backup and Restore

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 10, Backup and Restore.

Backup

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 10, Backup.

Restore

This software function is the same in both the CytExpert software and the CytExpert for Spectral software. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 10, Restore.

Log Cleanup

Use Log Cleanup to delete Experiment Operation Logs and System Operation Logs before a selected date.


NOTE These procedures are only available if you have the Electronic Record Management software option installed.

NOTE If you are using Conventional Mode, refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 10, Log Cleanup.

1 Ensure all experiment files are closed.

10-8 D17052AA


Select Backup/Restore > Log Clean-up. The Log Clean-up window appears.



Backup the data and try again if the operation system prompts the following message:




3 Select the desired date and time.



D17052AA


4 Select **Next**. The Confirm window appears.



Close the experiment and try again if the operation system prompts the following message:



5 Select **Yes**. The Clean-up succeeded window appears.



6 Select Finish.

10-10 D17052AA

Cleaning Procedures

Overview

This chapter describes how to carry out certain routine and nonscheduled cleaning procedures. Proper cleaning can help extend the service life of the instrument and ensure experimental accuracy. When conducting any cleaning, take all necessary biosafety precautions and use proper personal protective equipment.

Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 11, Cleaning Procedures for detailed instructions on the following procedures:

- Routine Cleaning
 - Daily Clean
 - Daily Clean [With Plate Loader]
 - Cleaning the Sample Station
 - Deep Clean Procedure
 - Cleaning the 4 L Sheath Fluid Container
 - Cleaning the 4 L Waste Container
- Nonscheduled Cleaning
 - Surface Cleaning and Disinfection
 - Preparing the Instrument for Transport or Storage

This chapter contains information on:

Preparing the Instrument for Transport or Storage

Preparing the Instrument for Transport or Storage





When the instrument is to be transported or is not to be used for 30 days or more, complete the emptying processes to prevent instrument damage and to reduce the possibility of biological contamination. Contact us if you have any questions.

1 Prepare the CytoFLEX for transportation or storage. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 11, Preparing the Instrument for Transport or Storage for detailed instructions.

D17052AA 11-1


- **2** Power down the CytoFLEX mosaic Spectral Detection Module and disconnect all the cables.
- If the CytoFLEX mosaic Spectral Detection Module is to be transported or stored, put the instrument into the Beckman Coulter supplied packing and comply with the requirements in APPENDIX A, Instrument Transportation and Storage, regarding correct placement during transportation and storage.

Lifting and Carrying Instructions [CytoFLEX S and CytoFLEX LX]


Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 11, Lifting and Carrying Instructions for detailed instructions.


Lifting and Carrying Instructions [CytoFLEX mosaic Spectral Detection Module]

1 Position a person on the left and right sides of the CytoFLEX mosaic Spectral Detection Module.




11-2 D17052AA





3 Gently lift the CytoFLEX mosaic Spectral Detection Module as shown in the figure below.



! WARNING

Risk of personal injury. Use caution when lowering the CytoFLEX mosaic Spectral Detection Module to avoid pinching fingers.

4 Lower the CytoFLEX mosaic Spectral Detection Module to its designated location.

D17052AA

Cleaning ProceduresPreparing the Instrument for Transport or Storage

11-4 D17052AA

Replacement/Adjustment Procedures

Overview

This chapter describes how to carry out certain routine and nonscheduled maintenance procedures. Proper maintenance can help extend the service life of the instrument and ensure experimental accuracy. When conducting any maintenance work, take all necessary biosafety precautions.

IMPORTANT In addition to parts specifically discussed, for all replacement parts, use only parts provided by Beckman Coulter to ensure proper functioning of the instrument. Never disassemble any part of the instrument without prior authorization. Beckman Coulter assumes no responsibility for any instrument problems resulting from the use of any part not authorized by Beckman Coulter for use with the instrument.

Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 12, Replacement/Adjustment Procedures for detailed instructions on the following procedures:

- Routine Replacement/Adjustment
 - Front Cover Removal and Reinstallation
 - Right-Side Cover Removal and Reinstallation
 - Filling the 4 L Sheath Fluid Container
 - Replacing the 10 L Sheath Fluid Cubitainer
 - Emptying the 4 L Waste Container
 - Emptying the 10 L Waste Cubitainer
 - Managing the Maintenance Reminder

NOTE If you have previously installed the CytExpert software, the Maintenance data will automatically synchronize between the CytExpert Software and the CytExpert for Spectral Software.

- Adding the Deep Clean Solution
- Replacing the Sheath Fluid Filter
- Unclogging the Sample Probe [With Sample Injection Mode Control Kit]
- Unclogging the Sample Probe [Without Sample Injection Mode Control Kit]
- Replacing the Sample Probe and/or the Sample Peristaltic Pump Tubing
- Replacing the Sample Probe Assembly [With Plate Loader]
- Changing the Sample Probe from the Single Tube Sample Station to the Plate Loader
 [CytoFLEX With Plate Loader]

D17052AA 12-1

- Changing the Sample Probe from the Plate Loader to the Single Tube Sample Station
 [CytoFLEX With Plate Loader]
- Inspecting the Liquid Flow Path for Leaks
- Priming the Flow Cell
- Replacing the Plate Holder [With Plate Loader]
- Plate Loader Module Removal and Reinstallation [With Plate Loader]
- Changing the Event Rate Setting
- Nonscheduled Replacement/Adjustment
 - Calibrating the Sample Flow Rate
 - Calibrating the Sample Flow Rate [With Plate Loader]
 - Setting Laser Delay
 - Replacing the Fuse
 - Replacing the Sheath Fluid Harness and/or Waste Harness
 - Changing Sample Mixing and Backflush Settings
 - Calibrating the Plate Position [With Plate Loader]

This chapter contains information on:

Switching between Spectral Mode and Conventional Mode

Switching between Spectral Mode and Conventional Mode

The CytoFLEX flow cytometer equipped with a CytoFLEX mosaic Spectral Detection Module supports two work modes: Conventional Mode and Spectral Mode.


- Spectral Mode: Attach the CytoFLEX mosaic Spectral Detection Module to the CytoFLEX flow cytometer to enable its spectral flow functionality and allow users to run more complex assays and use spectral unmixing technology for spectral overlap correction. Refer to Switching to Spectral Mode for detailed instructions on switching the instrument from Conventional Mode to Spectral Mode.
- Conventional Mode: Adjust how the optical fibers are connected to the Cytometer to use the
 Cytometer's optical components for collecting optical signals. Refer to Switching to
 Conventional Mode for detailed instructions on switching the instrument from Spectral Mode
 to Conventional Mode.

12-2 D17052AA

Table 12.1 Optical Fibers with Indicators [CytoFLEX S]

Optical Fibers	Indicator	Name
	Violet	Violet fiber
	Yellow	Yellow fiber
	Red	Red fiber
	Blue	Blue fiber

Table 12.2 Optical Fibers with Indicators [CytoFLEX LX]



Switching to Conventional Mode




Risk of data integrity damage.

- Verify that the optical fibers are securely connected to the WDM. A loose connection can alter the optical path and affect fluorescence.
- Do not kink the optical fibers.
- 1 Confirm that the instrument is in the standby state or that the instrument and CytoFLEX mosaic Spectral Detection Module are turned off.
- **2** Open the top cover.

D17052AA 12-3

Locate the Forward Scatter Switch and move the FSC connector from the Spectral Port (1) to the Conventional Port (2).


[CytoFLEX LX Shown]




D17052AA 12-4

4 Disconnect the rest of the optical fibers.

[CytoFLEX LX Shown]




[CytoFLEX S Shown]




D17052AA 12-5

Remove the connector protective caps from the WDM optical fiber connectors, and connect them to the adapter optical fiber connectors.

[CytoFLEX LX Shown]



6 Connect the correct optical fibers to the corresponding WDMs. Ensure the indicator number or color of the optical fiber matches the corresponding WDM. Refer to Table 12.1 and Table 12.2.



7 Close the top cover.

12-6 D17052AA


- **8** Log out of the CytExpert for Spectral software. Refer to CHAPTER 4, Logging Out of the Software.
- **9** Log into the CytExpert software. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 4, Logging into the Software for detailed instructions.
- **10** Initialize the instrument. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 4, Initializing the Instrument for detailed instructions.
 - **NOTE** If the instrument is turned off in Step 1, turn on the instrument. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 4, Turning On the Instrument for detailed instructions.
- **11** Run QC. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 5, Quality Control for detailed instructions.

Switching to Spectral Mode


- 1 Confirm that the instrument is in the standby state or that the instrument is turned off.
- **2** Open the top cover.

D17052AA 12-7

3 Locate the Forward Scatter Switch and move the FSC connector from the Conventional Port (2) to the Spectral Port (1).




4 Disconnect the rest of the optical fibers.




12-8 D17052AA

5 Remove the connector protective caps from the adapter optical fiber connectors, and connect them to the WDM optical fiber connectors.



6 Connect the correct optical fibers to the adapter optical fiber connectors. Ensure the indicator number or color of the optical fiber matches the corresponding adapter optical fiber connectors. Refer to Table 12.1 and Table 12.2



D17052AA 12-9

- **7** Close the top cover.
- **8** Turn on the instrument and the CytoFLEX mosaic Spectral Detection Module. Refer to CHAPTER 4, Turning On the Instrument.
- **9** Log out of the CytExpert software. Refer to the *CytoFLEX Platform Instructions for Use* manual, CHAPTER 4, Logging Out of the Software for detailed instructions.
- 10 Log into the CytExpert for Spectral software. Refer to CHAPTER 4, Logging into the Software.
- 11 Initialize the instrument. Refer to CHAPTER 4, Initializing the Instrument.
- 12 Run QC. Refer to CHAPTER 5, Quality Control.

12-10 D17052AA

Instrument Installation

Overview

[CytoFLEX S]: Your instrument may have been shipped directly to your laboratory, in which case you will need to set up and connect the Cytometer and the Workstation. Refer to this chapter for the instrument installation procedure.

Validate the environment → Open the package → Install the instrument → Prepare to boot

[CytoFLEX LX]: The CytoFLEX LX is installed by your Beckman Coulter Service Representative. Do not open the box or crate. Wait for a qualified Beckman Coulter Service Representative.

[CytoFLEX mosaic 63/88]: The CytoFLEX mosaic 63/88 is installed by your Beckman Coulter Service Representative. Do not open the box or crate. Wait for a qualified Beckman Coulter Service Representative.

This chapter contains information on:

- Instrument Transportation and Storage
- Installation Environment Validation
- Unpacking the Instrument and Inspecting the Materials for Defects or Omissions
- Installing the Instrument and Connecting the Equipment
- CytExpert for Spectral Software Installation Options
- Installing the CytoFLEX for Spectral Software
- Reinstalling the CytExpert for Spectral Software
- Restoring the Cytometer and CytoFLEX mosaic Spectral Detection Module Configuration File

Instrument Transportation and Storage

Refer to CHAPTER 11, Preparing the Instrument for Transport or Storage, prior to transportation or storage.

Attention to the following items is required when transporting or storing the instrument:

- Take caution to protect the instrument from exposure to rain or sunlight.
- Always place the instrument on a flat, stable surface, and take note of the symbol for this side up.
- Ambient temperature: -10 to 55 °C.
- Ambient humidity: < 80% RH.
- To prevent extrusion, the load on top cannot exceed 100 kg.

- CytoFLEX S Cytometer gross weight 27 kg; transport the instrument using only appropriate equipment so as to guard against personal injury.
- CytoFLEX LX Cytometer gross weight 103 kg; transport the instrument using only appropriate equipment so as to guard against personal injury.
- CytoFLEX mosaic 63 gross weight 68 kg; transport the instrument using only appropriate equipment so as to guard against personal injury.
- CytoFLEX mosaic 88 gross weight 71 kg; transport the instrument using only appropriate equipment so as to guard against personal injury.

Installation Environment Validation

IMPORTANT This instrument is intended for indoor use only.

Verify whether the installation environment satisfies the following requirements:

Worktable



Risk of instrument damage. Place the instrument on a level surface. Failing to do so places the system is in danger of toppling and can result in damage. Take all necessary precautions throughout the process of storing or transporting the instrument.

- The tabletop must be smooth and level.
- Minimum tabletop load bearing capacity [CytoFLEX S with CytoFLEX mosaic 63]: 90 kg.
- Minimum tabletop load bearing capacity [CytoFLEX LX with CytoFLEX mosaic 88]: 140 kg.
- The tabletop must not vibrate or shake.
- Minimum tabletop dimensions [CytoFLEX S with CytoFLEX mosaic 63]: 120 cm x 80 cm; minimum vertical space above tabletop: 100 cm.
- Minimum tabletop dimensions [CytoFLEX LX with CytoFLEX mosaic 88]: 200 cm x 80 cm; minimum vertical space above tabletop: 100 cm.
- Position the instrument in such a manner that it will facilitate disconnection of the power cable at the instrument end.

Ventilation and Cleaning

IMPORTANT If necessary, use ventilation equipment, but airflow must not be allowed to blow directly on the system, as it can affect the reliability of the data.

- Ensure that the working environment is well ventilated for proper heat dissipation.
- Maintain a clearance of at least 20 cm from the back of the instrument for heat dissipation.
- Keep the environment as dust free as possible.

A-2 D17052AA

A

- Avoid direct exposure to sunlight.
- Avoid placing near heat sources or exposing to drafts.
- Avoid corrosives or flammable gases.

Power Source



Risk of electric shock and/or instrument damage. Ensure that the power source is properly grounded. Improper grounding can cause electric shock and damage the system. Verify that the output voltage of the power outlet conforms to the system requirements and that a 5 A, time delay, T 5 AL 250 VAC fuse is installed. To prevent personal injury, Beckman Coulter recommends using a power source designed to protect against electrical shock.



Possible instrument damage could occur if you use an extension cord or a power strip to connect the Cytometer. Always plug the Cytometer into a dedicated outlet with an isolated ground.

The power source requirements are as follows:

- This instrument has been tested to and meets all applicable requirements for CE Marking.
- This instrument complies with the emission and immunity requirements described in IEC 61326-1.
- This equipment has been designated and tested to CISPR 11 Class A. In a domestic environment it may cause radio interference, in which case, you may need to take measures to mitigate the interference.
- It is advised that the electromagnetic environment should be evaluated prior to operation of the device.
- Do not use this device in close proximity to sources of strong electromagnetic radiation (unshielded intentional RF sources), as these may interfere with the proper operation.
- 100-240 volts, 50/60 Hz, 3-wire power cable, well grounded.
- Amperage not less than 10 A.
- The system requires a well-grounded power outlet (150 VA normal, 250 VA max) to provide the necessary power.
- Distance from system to socket less than 1.5 m.

Power consumption of the Plate Loader is <30 W.

Temperature and Humidity



Risk of instrument damage and/or erroneous results. To ensure reliability, the system must be operated in the specified environment, within the required temperature and humidity ranges. If the ambient temperature or humidity level falls outside the ranges mentioned above, use appropriate air conditioning.

- CytoFLEX S: Ambient temperature: 15-27°C with fluctuations of no more than <±2°C per hour.
- CytoFLEX LX: Ambient temperature: 15-30°C with fluctuations of no more than <±2°C per hour.
- **CytoFLEX mosaic Spectral Detection Module:** Ambient temperature: 15-30°C with fluctuations of no more than <±2°C per hour.
- Relative humidity: 15% RH-80% RH, non-condensing.

Waste Disposal





Risk of biohazardous contamination if you have skin contact with the waste container, its contents, and its associated tubing. The waste container and its associated tubing might contain residual biological material and must be handled with care. Clean up spills immediately. Dispose of the contents of the waste container in accordance with your local regulations and acceptable laboratory procedures.

The waste line from the Cytometer is connected to a waste container/cubitainer. Dispose of the system's waste in accordance with your local regulations and acceptable laboratory procedures.

The waste line supplied with the instrument can be connected to an open drain. If you use an open drain, mechanically secure the waste tube into the drain so the tube cannot accidentally come out of the drain. This prevents spillage.

A-4 D17052AA

Unpacking the Instrument and Inspecting the Materials for Defects or **Omissions**

Possible instrument damage could occur if you uncrate the instrument, install it, or set it up especially when you are in doubt of an instrument damage during shipment. Keep the instrument in its package until your Beckman Coulter Representative unpacks it for installation and setup.

Unpacking the Instrument and Inspecting the Materials for Defects or Omissions

When you receive your instrument, carefully inspect all cartons. If you see signs of mishandling or damage, file a claim with the carrier immediately. If separately insured, file the claim with the insurance company.

Take care to store the instrument in a suitable environment where it can remain in the proper position.

[CytoFLEX S and CytoFLEX LX]: Check that the following components on the packing list are present:

- Cytometer
- Cables
- Computer
- Mouse
- Keyboard
- Monitor
- Fluid Container holder
- Sheath fluid container [CytoFLEX S]
- Waste container
- Sheath fluid tubing
- Waste tubing
- USB configuration key
- Software USB

NOTE For the computer related components, refer to the manufacturer's manual.

[CytoFLEX mosaic Spectral Detection Module]: Check that the following components on the packing list are present:

- CytoFLEX mosaic
- Ethernet cable
- Power cable
- FSC cables
- USB configuration key
- Software USB

- Optical Fibers
- Top Cover and hinge board
- WDM Fixed block
- Cable ties
- FSC connect box
- FSC cable strain relief bush
- Connector protective covers
- Laser cover [CytoFLEX mosaic 88]
- Fiber digital wire clips [CytoFLEX mosaic 88]

Installing the Instrument and Connecting the Equipment



Risk of personal injury. Do not transfer the instrument from one location to another after the instrument installation is completed by service.

Installing the Instrument and Connecting the Equipment [CytoFLEX S]

Refer to the *CytoFLEX Platform Instructions for Use* manual, APPENDIX A, Installing the Instrument and Connecting the Equipment [CytoFLEX] for detailed instructions. For connecting the instrument, refer to Figure 1.3.

Installing the Instrument and Connecting the Equipment [CytoFLEX LX]

The CytoFLEX LX must be must be installed by a trained Beckman Coulter service engineer. For connecting the instrument, refer to Figure 1.6.


Installing the CytoFLEX mosaic Spectral Detection Module

The CytoFLEX mosaic Spectral Detection Module must be installed by a trained Beckman Coulter service engineer. For connecting the instrument, refer to Figure 1.3, Figure 1.6 and Figure 1.10.

A-6 D17052AA

CytExpert for Spectral Software Installation Options

CytExpert for Spectral software has two installation options upon install.



- CytExpert for Spectral User Management software option. User Login is required to run the system. Contains features and functionality that facilitates user and role management.
- CytExpert for Spectral Electronic Record Management software option. User Login is required to run the system. Contains features and functionality that facilitates compliance with 21 CFR Part 11 guidelines for Electronic Records and Signatures.

Installing the CytoFLEX for Spectral Software

The CytExpert for Spectral software is installed by Beckman Coulter service.

Reinstalling the CytExpert for Spectral Software

IMPORTANT The CytExpert software and the CytExpert for Spectral software can be installed on the same computer, but you cannot use the CytExpert for Spectral software and the CytExpert software together.

Use this procedure to:

- Change the software option installed. Refer to CytExpert for Spectral Software Installation Options for the differences between each software option available.
- Reinstall the same version of software.
- Upgrade your software to a higher version.

NOTE If you are using Conventional Mode, refer to the *CytoFLEX Platform Instructions for Use* manual, APPENDIX A, Reinstalling the CytExpert Software.



Risk of data loss. Reinstalling the CytExpert for Spectral software could overwrite your database. Ensure you backup your database prior to software reinstallation.

Uninstalling the Software

IMPORTANT Ensure you backup all of your data before uninstalling the software.

- Backup your CytExpert for Spectral data if you previously had the User Management or Electronic Record Management software option installed. Refer to CHAPTER 10, Backup and Restore.
- 2 Insert the software USB into the computer.

NOTE If the Autoplay window appears, select **Open folder to view files**.



A-8 D17052AA

3 Select **CytExpert for Spectral Setup_x.x.x.exe**. The Welcome to the CytExpert for Spectral Setup Wizard window appears.




A User Account Control message may display the message *Do you want to allow this app to make changes to your device?* Select **Yes** to continue.

4 Select **Next**. The Change your installation of CytExpert for Spectral screen appears.




5 Select **Remove**. The Begin remove of CytExpert for Spectral screen appears.




A-10 D17052AA

6 Select **Remove**. The Removing Progress window appears.



- **IMPORTANT** If you select **Yes** to remove the CytExpert for Spectral settings, CytExpert for Spectral database, cytometer configuration, and temporary files, your settings, database, configuration files, and temporary files will be overwritten and you will need to reinstall your cytometer configuration and restore any databases you might have.
- **7** The message Would you like to remove CytExpert for Spectral settings, CytExpert for Spectral database, cytometer configuration and temporary files? appears. Select **No**.



A-12 D17052AA

8 When software removal is complete, the software displays the message *CytExpert for Spectral has been removed successfully.*



9 Select Close.

Reinstalling the CytExpert for Spectral Software

1 Insert the software USB into the computer.

NOTE If the Autoplay window appears, select Open folder to view files.




2 Select **CytExpert for Spectral Setup_x.x.x.exe**. The Welcome to the CytExpert for Spectral Setup Wizard window appears.



A User Account Control message may display the message *Do you want to allow this app to make changes to your device?* Select **Yes** to continue.

A-14 D17052AA


3 Select **Next**. The Beckman Coulter License Terms window appears.




- 4 Read the Beckman Coulter Customer End User License Agreement.
- **5** Select the *I accept the terms of this agreement* checkbox.

NOTE The checkbox is not selectable until you scroll all the way to the end of the agreement.

6 Select **Next**. The Choose a file location window appears.




7 Select **Next**. The Installation Options window appears.




8 Select the desired installation option. Refer to CytExpert for Spectral Software Installation Options.

A-16 D17052AA

9 Select **Next**. The Begin installation of CytExpert for Spectral window appears.



 $\textbf{10} \ \ \textbf{Select Install} \ \textbf{to begin installing the software.} \ \textbf{The Installation Progress window appears.}$




11 If the following software prompt appears, select **OK**.



NOTE The term "device" in this message refers to the cytometer.

12 Wait for the software to finish installing. The install complete window appears.




13 Select Close to finish the CytExpert for Spectral software installation.

A-18 D17052AA


Restoring the Cytometer and CytoFLEX mosaic Spectral Detection **Module Configuration File**

Insert the software USB into the computer.

NOTE If the Autoplay window appears, select **Open folder to view files**.




- Open the CytExpert for Spectral software. Refer to CHAPTER 4, Logging into the Software.
- Select **Advanced** > **Restore Configuration**. The Restore Configuration window appears.




NOTE Ensure the experiment is closed.

- Select and browse to the desired configuration file folder from your Software USB.
- Open the Cytometer configuration file CytExpert for Spectral_x.x.xxx_Config_Backup_.bak.


6 Select **Execute**. The Restore Configuration window appears.



7 Select **Yes** to start restore the configuration.



8 Select Finish.



NOTE If the CytExpert software is installed, its configuration will automatically synchronized with the restored configuration. Restart the CytExpert software to apply the restored configurations.

9 Select OK.

A-20 D17052AA

APPENDIX B

CytExpert for Spectral Electronic Record Management

Overview

IMPORTANT You must have the CytExpert for Spectral Electronic Record Management software option installed to use the features listed below. Refer to APPENDIX A, CytExpert for Spectral Software Installation Options.

Beckman Coulter's CytoFLEX with CytExpert for Spectral software version 1.0 and higher contains features and functionality that facilitates compliance with 21 CFR Part 11 guidelines for Electronic Records and Signatures. This electronic record management includes controls for user identification, permissions, electronic signatures, data integrity, operation and experiment logs and audit trails. CytExpert for Spectral software version 1.0 and higher contains a database that uses checksum matching to prevent tampering of the records and files that are indexed in the Closed File System.

This chapter contains information on:

- Software Menu
- Experiment Management
- Log
- Electronic Signature
- User Management

Software Menu

The CytExpert for Spectral Electronic Record Management software option includes additional software menu items that are not available in the CytExpert for Spectral User Management software option. Refer to Figure 2.5 and Figure 2.6 for comprehensive software menu trees and details on which menu item applies to each software option.

NOTE If your CytoFLEX instrument is not equipped with a CytoFLEX mosaic Spectral Detection Module or if you are using Conventional Mode, refer to the *CytoFLEX Platform Instructions for Use* manual, APPENDIX B, Software Menu.

Experiment Management



Risk of file corruption. Do not add, delete, or modify data from the Windows Explorer directory. Manage all data changes using Experiment Explorer to ensure file indexing remains intact.

Closed File System

Refer to the *CytoFLEX Platform Instructions for Use* manual, APPENDIX A, Closed File System for detailed instructions.

Experiment Directory Management

Refer to the *CytoFLEX Platform Instructions for Use* manual, APPENDIX A, Experiment Directory Management for detailed instructions.


Folder Hierarchy Management

Refer to the *CytoFLEX Platform Instructions for Use* manual, APPENDIX A, Folder Hierarchy Management for detailed instructions.

B-2 D17052AA

Experiment Related Operations

The Experiment Explorer dialog appears in place of the Windows File Explorer in the following operations: New/Open Experiment, New/Open Spectral Unmixing, Save As, Save As Template, Recent Template, and New Experiment from Template.



Importing an Experiment/Template

Refer to the *CytoFLEX Platform Instructions for Use* manual, APPENDIX B, Importing an Experiment/Template for detailed instructions on importing experiment (.xit), spectral unmixing (.xitsp), or experiment from template (.xitm) files into the system.

Exporting an Experiment/Template

Refer to the *CytoFLEX Platform Instructions for Use* manual, APPENDIX B, Exporting an Experiment/ Template for detailed instructions on exporting experiment (.xit), spectral unmixing (.xitsp), or experiment from template (.xitm) files from the system.

Log

Refer to the *CytoFLEX Platform Instructions for Use* manual, APPENDIX B, Log for detailed instructions on the following:

- Experiment Operation Log
- System Operation Log
- User Management Operation Log

Electronic Signature

Refer to the *CytoFLEX Platform Instructions for Use* manual, APPENDIX B, Electronic Signature for detailed instructions on the following procedures:

- Signing Experiments
- · Rejecting Experiment
- Setting the Signature Retention Period
- Signature Setting
- Printing an Experiment Signature

User Management

User Administration

Logging In and Out of the Software

Refer to CHAPTER 4, Logging into the Software and CHAPTER 4, Logging Out of the Software.

Locking the Account

Refer to CHAPTER 4, Locking the Account.

Role Management

Refer to CHAPTER 2, Role Management.

User Management

Refer to CHAPTER 2, User Management.

Account Policies

Refer to CHAPTER 2, Account Policies.

B-4 D17052AA

Sample Injection Mode Control Kit

Overview

The Sample Injection Mode Control Kit is a mechanical knob installed by your service engineer that enables users to switch between the Plate Loader sample injection mode and the Semi-Automatic or manual sample injection mode. The switch eliminates the need to manually re-route the tubing.

Refer to the CytoFLEX Platform Instructions for Use manual, APPENDIX C, Sample Injection Mode Control Kit for detailed instructions on the following procedures:

- Performance Characteristics [With the Sample Injection Mode Control Knob]
- Sample Injection Mode Control Kit Components
- Switching the Sample Probe from the Single Tube Sample Station to the Plate Loader [With Sample Injection Mode Control Knob]
- Switching the Sample Probe from the Plate Loader to the Single Tube Sample Station [With Sample Injection Mode Control Knob]

Sample Injection Mode Control Kit Overview

C-2 D17052AA

APPENDIX D Deep Well Plate

Specimen Collection Plate Specifications

Refer to the *CytoFLEX Platform Instructions for Use* manual, APPENDIX D, Specimen Collection Plate Specifications for detailed instructions.

D-1

Deep Well PlateSpecimen Collection Plate Specifications

D-2 D17052AA

APPENDIX E Cyber Security

Overview

The following procedures are the only validated and approved recommendations for cyber privacy and security. Contact your IT professional for assistance.

- Changes below require Administrator level Windows access.
- Before connecting to an external device such as a hard drive, USB, or DVD/CD, verify it does not have a virus or malware.


NOTE Beckman Coulter recommends using multiple anti-virus software to ensure that all viruses or malware can be identified.

• Close all unnecessary applications while operating the CytExpert for Spectral software.



Risk of security vulnerabilities and privacy issues. Use the USB and Ethernet ports with caution. Unauthorized or improper use of these ports will lead to potential information security risks. Beckman Coulter recommends disabling all the unused USB and Ethernet ports on the instrument computer. Ensure that only trusted and authorized devices are connected to the active ports.

To prevent unauthorized access and data breaches, Beckman Coulter recommends checking the use of the USB and Ethernet ports on a regular basis.



System integrity could be compromised and operational failures could occur if:

- You introduce software that is not authorized by Beckman Coulter into your computer. Only operate your system's computer with software authorized by Beckman Coulter.
- You apply domain policies which alter the default configuration.
- You alter the system in a manner other than the approved changes outlined below.

IMPORTANT Beckman Coulter Life Sciences has assembled a global product security team to assess vulnerabilities and determine responses in a coordinated Vulnerability Disclosure (CVD) process. If you encounter a potential threat to cyber security, please contact us via

https://www.beckman.com/about-us/compliance/coordinated-vulnerability-disclosure.

This chapter contains information on:

- Drive Encryption
- Protection from Malware Software
- Security Updates
- User Management
- System Hardening
- Remote Access

Drive Encryption

Your Windows 10 system is equipped with BitLocker. Use BitLocker disk encryption software to prevent unauthorized access to your hard drive, refer to BitLocker Keys. For instructions on BitLocker Encryption or BitLocker Decryption, see Enabling BitLocker or BitLocker Decryption.

IMPORTANT BitLocker is the only approved encryption method.

BitLocker Keys



Risk of data loss. Store your BitLocker key in a safe place. For instructions on saving your BitLocker key, refer to Backing Up Your Recovery Key. If you lose your BitLocker Key, you will lose access to your entire hard drive including all of your data. Beckman Coulter is not responsible for security keys and will not be able to recover your data off your hard drive. Work with your IT department in creating and storing the keys in a safe place.

Changing Encryption Strength Level for BitLocker


BitLocker defaults to XTS-AES 128-bit encryption method. To alter the encryption, follow the steps below before enabling BitLocker. Changes to the encryption level will only impact the encryption of new volumes. Work with your IT department to change the Encryption Strength Level.

E-2 D17052AA


E-3

Changing the encryption strength for BitLocker

Type **gpedit.msc** from the Windows search bar to locate the gpedit Control Panel.




2 Click **gpedit.msc**. The Local Group Policy Editor window appears.




D17052AA

3 Navigate to Computer Configuration > Administrative Templates > Windows Components > BitLocker Drive Encryption.




4 Double-click "Choose drive encryption method and cipher strength (Windows 10.....)".



E-4 D17052AA

Select **Enabled** and use the dropdown arrows to select the encryption method for the operating system drive, fixed data drive, and for the removable data drives.




NOTE The default is set to XTS-AES 128-bit.

- 6 Select Apply.
- 7 Select **OK** and close all windows.
- **8** Restart the computer and proceed to Enabling BitLocker.


Enabling BitLocker

Type Manage BitLocker from the Windows search bar to locate the Manage BitLocker Control Panel.




E-6 D17052AA


2 Double-click Manage BitLocker to open the BitLocker Drive Encryption window.




3 Select Turn on BitLocker from the Control Panel screen.



4 BitLocker verifies that your PC meets the requirements.




5 Select **Print the recovery key** to back up the recovery key. For instructions, refer to Backing Up Your Recovery Key.




E-8 D17052AA

6 Select Next.




7 Select "Encrypt used disk space only (faster and best for new PCs and drives)", then select **Next**.



IMPORTANT The system check ensures that BitLocker can read the recovery and encryption keys correctly before encrypting the drive. This check might take some time, but it is recommended to ensure that your selected unlock method works without requiring the recovery key.

8 Select "Run BitLocker system check", then select **Start encrypting**.



NOTE BitLocker restarts your computer before encrypting.

9 Select **Restart now** to start the BitLocker Encryption.




10 Enter your password in the Password screen.


NOTE After logging in you will see a Window that says **Encrypting in Progress**.

E-10 D17052AA


11 Open the Manage BitLocker Control Panel from the Windows search bar.



12 Under Operating System Drive, you will see the BitLocker is Encrypting the drive.



13 Once the BitLocker has encrypted the drive, the BitLocker will be on.




E-12 D17052AA

Backing Up Your Recovery Key


IMPORTANT The Back Up Your Recovery Key to your Microsoft Account is not supported at this time.

IMPORTANT The Recovery Key can only be saved to a USB or external drive. You cannot save to a location on your PC.


Type Manage BitLocker from the Windows search bar to locate the Manage BitLocker Control Panel.



2 Double-click **Manage BitLocker** to open the BitLocker Drive Encryption window.




3 Select "Backup your recovery key". The following window appears.




E-14 D17052AA

IMPORTANT The Recovery Key can only be saved to a USB or external drive. You cannot save to a location on your PC.


- **4** Select a destination to save the Recovery Key.
 - **Save to a file**: Select "Save to a file". The following screen appears.





• **Print the recovery key**: Select "*Print the recovery key*". The following screen appears.

1. Select "Microsoft Print to PDF", and then select Print.



2. Select the location that you want the file saved. In the File name, type **BitLocker Key** and in the Save as type, ensure that it is listed as .pdf.

5 Select Save then select Finish.


E-16 D17052AA

BitLocker Decryption


IMPORTANT Turning off the BitLocker feature decrypts all the data on the drive. Decrypting your operating system means that you have removed BitLocker from your system. The BitLocker key will change each time BitLocker is re-enabled. Only the latest BitLocker key will work. Previous BitLocker keys will be invalid and will not unlock the system. Please note this is a time-consuming process. However, you can also suspend BitLocker instead of decrypting the drive, see Suspending BitLocker Drive Encryption.

Turning off BitLocker


Type Manage BitLocker from the Windows search bar to locate the Manage BitLocker Control Panel.



2 Double-click Manage BitLocker to open the BitLocker Drive Encryption.



3 Select Turn off BitLocker.




E-18 D17052AA

The following message appears.




IMPORTANT Drive decryption might take a long time, but you can keep using your PC during the decryption process.

4 Select **Turn Off BitLocker**. The BitLocker Drive Decryption screen displays.



After the decryption has completed, the following screen displays.




E-20 D17052AA

Suspending BitLocker Drive Encryption


IMPORTANT Suspending BitLocker Drive Encryption is a temporary method for removing BitLocker protection without decrypting the Windows drive. Use Suspend BitLocker only if you need to update the BIOS firmware or Windows boot loader/startup files that requires access to the drive in a pre-boot environment. This will help prevent BitLocker from locking the drive and can avoid a lengthy decryption process.

Suspending BitLocker Protection


Type Manage BitLocker from the Windows search bar to locate the Manage BitLocker Control Panel.



2 Double-click **Manage BitLocker** to open the BitLocker Drive Encryption window.



3 Select Suspend Protection.



E-22 D17052AA

A confirmation message "Do you want to suspend BitLocker protection?" displays.



IMPORTANT If you choose to suspend BitLocker, your data will not be protected.

4 If you want to proceed, select **Yes** and the following popup information will display.




If you do not want to proceed, select No.


NOTE There are certain situations where you may want to suspend BitLocker such as updating your PC's firmware, hardware, or operating system. If you forget to resume BitLocker after these updates are completed, BitLocker will automatically resume the next time you restart your PC.

Restarting BitLocker Protection

Type Manage BitLocker from the Windows search bar to locate the Manage BitLocker Control Panel.




2 Double-click Manage BitLocker to open the BitLocker Drive Encryption window.



E-24 D17052AA

3 Select Resume Protection.



Recovering Your BitLocker Key

There might be situations where the BitLocker Key is requested from your computer, such as:

- Your drive is placed onto another computer
- Computer requests your recovery key; it may be due to one of the following but is not limited to:
 - Hardware changes
 - Firmware changes
 - Other Windows updates


BitLocker Recovery

1 BitLocker Recovery window will be displayed upon powering on your system.



E-26 D17052AA

2 Enter your key from your saved text file or printed copy, then select **Enter**.



3 Windows will reboot your computer.

Protection from Malware Software



The instrument workstation will be vulnerable to malware, viruses, and unauthorized access if the instrument PC is connected directly to the laboratory network without installation of all approved OS and malware patches.

IMPORTANT Only the version of the Trellix® Application Control (Version 8.4.1) has been validated to work with the CytoFLEX instrument.


For more information on Configuration of Trellix Application Control, refer to Trellix Application Control Guide.

Installing Trellix Application Control

IMPORTANT The Trellix Application Control Software (Version 8.4.1), which has been proved compatible with the software, should be purchased separately and downloaded from Trellix.

The following procedure can assist you to install the Trellix Application Control software. For more instructions, access the Trellix for help.

- 1 Unzip the Trellix Application Control file **SOLIDCORxxx-xxx_WIN.zip**
- **2** Right-click **Setup-win-8-10-2012-amd64-xxx.exe**, and select **Run as administrator**. The Trellix Installer opens.




E-28 D17052AA

3 Select **Next**. The License Agreement displays.




Read the license agreement and select I accept the terms in the license agreement, then select Next.

4 Enter the username and the serial number. Select **Next**.



If you select **Install without license key**, follow the procedures below to add the license key after the installation.

a. Right-click solidifier on your desktop, and select Run as administrator.



E-30 D17052AA

The following window displays.

```
Administrator: Trellix Solidifier Command Line


Microsoft Windows [Version 10.0.17763.6293]
(c) 2018 Microsoft Corporation. All rights reserved.

C:\Windows\system32>_
```


b. Enter cd C:\Program Files\McAfee\Solidcore.



c. Enter the license key in the console.




The license key is added.





E-32 D17052AA

5 Select the destination folder and select **Next**.




6 Select **Install**. The system starts to install Trellix Solidifier





E-34 D17052AA

Wait for the Trellix Solidifier to finish installing. Select Finish.



Enabling Trellix Application Control

The enabled Trellix Application Control Software prevents your system from executing unauthorized applications on your instrument workstation.

IMPORTANT Ensure all the necessary applications have been installed prior to enabling the Trellix Application Control.

Right-click solidifier. on your desktop, and select **Run as administrator**. The following window displays.





2 Enter cd C:\Program Files\McAfee\Solidcore.




E-36 D17052AA

3 Enter **sadmin so** in the console to scan the system and create an allowlist. The solidification process takes about 30 minutes depending on your system size and CPU performance.






4 Enter **sadmin enable** to enable the Trellix Application Control.



IMPORTANT Beckman Coulter recommends setting a password to manage the authority of enabling/disabling the Trellix Application Control.

5 Enter **sadmin passwd** to set a password to lock the Trellix Application Control authority.




6 Restart the system.

E-38 D17052AA


Exporting the Certificate

Follow procedures below to export the certificate of CytExpert for Spectral software into your PC.


1 Right-click on the application **CytExpert.exe**, and select **Properties**. The following window displays.



2 Select the **Digital Signatures** tab. Then select the signature and select **Details** to view the signature information.




The Digital Signature Details window displays.




E-40 D17052AA


3 Select View Certificate.



4 Select the **Details** tab and then select **Copy to File**.




The Certificate Export Wizard window appears.




E-42 D17052AA

5 Select **Next** to export the certificate.




6 Select Base-64-encoded X.509 and then select Next.




E-44 D17052AA

7 Select a destination to save the certificate and then select **Save**.



8 Select Next.



E-46 D17052AA

9 Select **Finish** when the success message displays.



Importing BEC Certificate

Follow the procedures below to add the certificate of CytExpert for Spectral software into the Trellix Certificate list. All the certified applications will be regarded as trusted and authorized by the Trellix Application Control.

IMPORTANT Trellix Application Control recognize only X.509 certificates.

Copy the certificate *BEC.cer* from the CytExpert for Spectral software disk into your instrument PC.



Right-click solidifier... on your desktop, and select **Run as administrator**. The following window displays.

```
Administrator: Trellix Solidifier Command Line

Microsoft Windows [Version 10.0.17763.6293]
(c) 2018 Microsoft Corporation. All rights reserved.

C:\Windows\system32>_
```

3 Enter sadmin cert add C:\BEC.cer to add the certificate.

Switching Trellix Application Control to Update Mode

The update mode allows the Trellix Application Control to scan the instrument workstation and add the newly-installed application into the allowlist automatically. Switch the Trellix Application Control to the update mode when you execute scheduled or emergency changes on your instrument workstation.

E-48 D17052AA


1 Right-click solidifier... on your desktop, and select **Run as administrator**. The following window displays.

```
Administrator: Trellix Solidifier Command Line

Microsoft Windows [Version 10.0.17763.6293]
(c) 2018 Microsoft Corporation. All rights reserved.


C:\Windows\system32>_
```

2 Enter **sadmin bu** to switch the Trellix Application Control to Update Mode.



3 Execute scheduled or emergency changes, patch installations, or software updates for your instrument workstation.

4 Enter sadmin eu to exit the Update Mode.




Disabling Trellix Application Control

IMPORTANT Disabling Trellix Application Control is a temporary method in exceptional cases, for example, Windows or firmware changes. This will help stop Trellix Application Control from blocking unauthorized applications. Beckman Coulter recommends enabling the Trellix Application Control in most cases.


All the applications installed in this mode will not be included in the Trellix allowlist.

1 Right-click solidifier... on your desktop, and select **Run as administrator**. The following window displays.




E-50 D17052AA

2 Enter **sadmin disable** in the console.



- 3 Select Enter.
- **4** Restart the System.
- **5** Enter **sadmin status** to view the status of Trellix Application Control.



Security Updates



Risk of software incompatibility. Only security updates posted to the Beckman Coulter website should be installed on the instrument controller PC. Do not install patches from any other source.


Install the update on one instrument and verify functionality before installing on other instruments.

IMPORTANT For information on the currently supported Windows build and security updates, consult the Release Notes of the latest validated update posted to the website. Refer to Downloading Security Updates.

IMPORTANT Before installing the security updates, ensure you have a recent backup of the instrument database and experiment files (see CHAPTER 10, Backup and Restore) and the current BitLocker key (see Backing Up Your Recovery Key). Ensure the backups are saved to an encrypted external media or secured network drive.


Enabling System Protection

1 In the Start Menu, type **Restore** to locate the **Create a Restore Point**.




E-52 D17052AA


 ${\bf 2} \quad \text{Double-click $\sf Create \ a \ Restore \ Point$ to open the System Properties window.}$



3 Select Configure.



The following window appears.




E-54 D17052AA

- 4 Select Turn on system protection, then select Apply.
- **5** Close the System Properties window by selecting **OK**.


Creating a Restore Point

A restore point is used to remove all software and update installations that occurred after it was created. If installing an update on your PC causes any loss of functionality, the restore point can be used to place the system back into the operable state present before the update installation. In all cases, a system restore with any restore point would not cause loss of any data (i.e., LMD files) on the system.




Risk of unexpected software modification. To ensure a system recovery point is available that is indicative of the most recent software configuration, the user should create a system restore point immediately prior to running the security update installer. Otherwise, a user would need to utilize the most recent automatically created restore point which might be up to 72 hours old and not include all recent software modifications.

1 In the Start Menu, type **Restore** to locate the **Create a Restore Point**.




D17052AA

2 Double-click Create a Restore Point to open the System Properties window.



3 Select Create.



E-56 D17052AA

- **4** Type a name for the restore point and select **Create**.
- **5** Wait for the restore point creation to complete and then select **Close**.
- **6** Close the System Properties window by selecting **OK**.

Downloading Security Updates

In order to download and install security updates, you will need to access the Beckman Coulter website at www.beckmancoulter.com.

IMPORTANT Account registration is required for new users.

- 1 Go to www.beckmancoulter.com.
- 2 Using the menu bar in the in the upper right hand corner of the screen, hover over **Login**, then scroll down to select **Software Downloads**.
- **3** To search by product:
 - a. Select Research & Discovery from the Market Segment dropdown menu.
 - **b.** Select **Flow Cytometry** from the Product Line. dropdown menu.
 - **c.** Select **Instruments** from the Product Series dropdown menu.
 - **d.** Select **CytoFLEX mosaic** from the Product dropdown menu.
- 4 Select Search.

IMPORTANT Print out the Release Notes for each update as they will not be available when the next update is released. The Component Name section of the release notes identifies the Microsoft updates that are included in the download. This information may be requested by your IT team and will be needed to recover from a failed update.

5 Each security update is cumulative and only the most recent validated set of updates will be available for download. Select the first file in the results to view the release notes and then select **Download**.

If more than one update is listed in the search results, select seach update.

If the security updates cannot be not downloaded directly to the instrument controller system, transfer the downloaded files via an encrypted external media or secured network drive to the instrument controller PC.

Notification Options

To receive notification when a new security update is validated and available for download, you must sign up for an account on www.beckmancoulter.com and follow the procedure for setting up subscriptions.

Sign up for subscriptions as outlined below:

- Select the **Document Language** consistent with your IFU language.
- Select Research & Discovery from Market Segment.
- Select Flow Cytometry from Product Line.
- Select Instruments from Product Series.
- Select CytoFLEX mosaic from Product.

If you do not sign up for product notifications, you will have to check periodically for newly released security updates, see Downloading Security Updates.

NOTE Beckman Coulter will only validate and post updates for Windows 10 LTSC 2019.

Installing the Security Updates



Risk of software incompatibility. Beckman Coulter can only support security updates on unmodified systems. Application of domain policies or other configuration changes may impact compatibility.


IMPORTANT If more than one update is available for your product, install the updates in the order indicated by the sequence number as outlined in the file naming convention below.

RSU-<Operating System> <Product Name>-<Year><Month>_<Sequence #>_<Sequence Revision#> **Example:** RSU-W10LTSC-CytoFLEX for Spectral-2025Feb_1_0

- 1 Ensure you have a recent copy of your BitLocker key, refer to Backing Up Your Recovery Key.
- 2 Login to Windows as a user with Administrator access.
- **3** Close all program and save open files.

E-58 D17052AA

- **4** Power down the instrument.
- **5** Open the first security update file (For example: RSU-W10LTSC-CytoFLEX for Spectral-2025Feb_1_0) on the Instrument Controller PC by double-clicking on the file. The BOSS SUPP Package Setup window appears.



If the Microsoft Windows Defender Smart Screen prompt appears, select **More Info** and **Run Anyway** to launch the installer.


A User Account Control message may display the message *Do you want to allow this app to make changes to your device?* Select **Yes** to continue.

NOTE The security update installer will automatically detect and install only the necessary updates.

6 Select **Next**. The Begin installation of BOSS SUPP Package window appears.




7 Select **Install** to begin installing the security updates. The Installing Progress window appears. This process may take several minutes.




E-60 D17052AA

f 8 Wait for the security updates to finish installing. The install complete window appears.




9 Select **Finish**. When prompted to restart the system, select **Yes**.





10 Wait while the system completes the installation of the security updates. This process may take several minutes and will return you to the Windows login screen.



- **11** Install any additional downloaded security updates (sequence numbers above 1) by repeating Steps 5 to 9.
- **12** Verify functionality of instrument software before installing security updates on any other system. If functionality is impaired, see Recovering from Failed Security Updates.

Recovering from Failed Security Updates


These steps below should only be followed if the instrument software or operating system functionality is impaired after a security update was applied.

- 1 Power down the instrument and close all programs. Do not turn your workstation off.
- 2 In the start menu, type view update and select View your Update History.
- 3 Select Uninstall Updates.
- 4 Highlight the update listed under Microsoft Windows that matches the cumulative security update KB number in the release notes. Refer to Component Name section of the Release Notes, see Step 5, Downloading Security Updates.
- **5** Select **Uninstall** and select **Yes** when prompted to verify the uninstall. This process may take several minutes to complete.
- **6** When prompted to restart, select **Restart Now**.

E-62 D17052AA

- 7 Once the PC has been rebooted, login and verify Microsoft update has been removed by repeating Steps 2 and 3. Verify the KB selected in Step 4 is no longer listed.
- **8** Power on the instrument and verify restored software functionality.
 - If software functionality has been restored, continue using your instrument without the update installed. Contact us for troubleshooting the update incompatibility.
 - If software functionality has not been restored and the PC was not provided by Beckman Coulter, contact laboratory IT for system restore options.
 - If software functionality has not been restored and the PC was provided by Beckman Coulter, contact us for system restore.

User Management



Risk of unauthorized access and modification. CytoFLEX workstation transmits unencrypted data to any network printer. To avoid unauthorized access to transmitted data, ensure adequate network protection mechanism are in place.



Risk of unauthorized access and modification. To protect your workstation from unauthorized access and modification, please refer to the following recommendations.

Beckman Coulter recommends that you use role-based access control policies to manage user accounts on workstations. Assigning appropriate roles to users accessing workstations is the responsibility and obligation of workstation administrators. Table E.1 is a good practice recommended by Beckman Coulter to manage user accounts according to the use requirements of the workstation and combined with Windows built-in user groups.

Table E.1 Role-based Access Control Policies

Requirements for Using a Workstation	Windows Built-in User Groups	Role Example
Responsible for the daily maintenance of the workstation and have full control authority over the workstation. Can install, uninstall software, install security updates, restore the system, and manage user accounts on workstation.	Administrators	 IT Administrator Device Administrator
The workstation is used daily to control the cytometer to complete collection, QC, analysis and other work.	User	 Operator Experimenter


System Hardening

Enabling Automatic Log-Off

IMPORTANT Automatic log-off can protect your instrument PC from unauthorized operation when you leave the PC temporarily. Follow the steps below to enable the Screen Saver protection.


Enabling Screen Saver Protection

1 Type **Edit group policy** from the Windows search bar to locate the Edit group policy editor.




E-64 D17052AA

 $\begin{tabular}{ll} \bf 2 & {\tt Select} \begin{tabular}{ll} \bf Edit \begin{tabular}{ll} {\tt group} \begin{tabular}{ll} {\tt policy} \end{tabular}. \begin{tabular}{ll} {\tt Select} \begin{tabular}{ll} \bf Edit \begin{tabular}{ll} {\tt group} \begin{tabular}{ll} {\tt policy} \end{tabular}. \begin{tabular}{ll} {\tt select} \be$




3 Navigate to User Configuration > Administrative Templates > Control Panel > Personalization.




E-66 D17052AA

4 Select **Enable screen saver.** The following window appears.




- **5** Select **Enabled** and **Apply**.
- **6** Select **OK** to close the Enable screen saver window.

7 Select Password protect the screen saver from the Local Group Policy Editor window.




E-68 D17052AA

8 Select Enabled and Apply.




9 Select **OK** to close the window.

10 Select Screen saver timeout from the Local Group Policy Editor window.




E-70 D17052AA

The Screen saver timeout window appears.




11 Select Enabled and Apply.



12 Select **OK** to close the window.


Setting Lock Screen

1 Type Lock screen settings from the Windows search bar to locate Screen saver settings.




E-72 D17052AA

The Settings window appears.




2 Scroll down to the bottom and select **Screen saver settings**. The Screen Saver Settings window appears.



- 3 Select Apply.
- 4 Select **OK** to close the window.


Setting Account Lockout Policy

1 Type Local Security Policy from the Windows search bar to locate the Local Security Policy.




E-74 D17052AA

 ${\bf 2} \quad \hbox{Double-click ${\bf Local Security Policy}$ to open the Local Security Policy window.}$




3 Navigate to Security Setting > Account Policies > Account Lockout Policy. Configure the Security Settings as below.



NOTE You must configure the **Account lockout threshold** first.


Setting Password Policy


1 Type Local Security Policy from the Windows search bar to locate the Local Security Policy.




E-76 D17052AA

 ${\bf 2} \quad \hbox{Double-click ${\bf Local Security Policy}$ to open the Local Security Policy window.}$







- **4** Customize the password policy according to the following.
 - Enforce password history: 10 passwords remembered
 - Maximum password age: 90 days
 - Minimum password age: optional
 - Minimum password length: 10 characters
 - Password must meet complexity requirements: enabled


The password must contain all four of the following character types:

- letter in upper case (A -Z)
- letter in lower case (a -z)
- numbers (0-9)
- distinctive character (i.e.,,!, @, #, \$,%, ^, &).
- Store passwords using reversible encryption: optional


E-78 D17052AA

Disabling Remote Access


1 Type Control Panel from the Windows search bar to locate the Control Panel app.




2 Select Control Panel to open the Control Panel window.



3 Type System from the search bar to locate Allow remote access to your computer.




4 Select Allow remote access to your computer. The System Properties window appears.



E-80 D17052AA


5 Select **Don't allow remote connections to this computer** from the Remote tab.




6 Select **OK** to close the window.

Disabling Auto Play

1 Type Auto Play settings from the Windows search bar to locate the Auto Play settings.




2 Select **Auto Play settings**. The Settings window appears.



E-82 D17052AA

3 Select **off** to turn off Auto Play.



4 Close the window.


Disabling Unnecessary Services

It is recommended to disable all the unnecessary services while operating the CytExpert for Spectral software.


Table E.2 List of Unnecessary Services

Bluetooth Support Service	Microsoft FTP Service (If available)	Windows Media Player Network Sharing Service
Geo-location Service	Microsoft iSCSI Initiator Service	Windows Mobile Hotspot Service
Infrared Monitor Service	Remote Procedure Call (RPC) Locator	WinHTTP Web Proxy Auto-Discovery Service
Internet Connection Sharing	Routing and Remote Access	Xbox Live Game Save
Link-Layer Topology Discovery Mapper	SSDP Discovery	Xbox Live Networking Service
LxssManager (If available)	UPnP Device Host	

1 Type Computer Management from the Windows search bar to locate the Computer Management App.




2 Select Computer Management.




E-84 D17052AA


3 Navigate to Services and Applications > Services.



4 Double-click the service you want to disable.




The following window appears.



E-86 D17052AA

5 Select **Disabled** from the Startup type dropdown.




6 Select Apply.


7 Select **OK** to close the window.

Disabling Unauthorized Applications

1 Type **Edit group policy** from the Windows search bar to locate the **Local group policy editor** control panel.




The Local Group Policy Editor window displays.




E-88 D17052AA


2 Navigate to Computer Configuration > Administrative Templates > Windows Components > Store.



3 Select Disable all apps from Windows Store.



The following window appears.




- 4 Select Enabled and Apply.
- **5** Select **OK** to close the window.
- **6** Disable the following applications. Refer to Steps 3-5.
 - Turn off Automatic Download and Install of updates.
 - Turn off the offer to update to the latest version of Windows.
 - Turn off the Store application.

E-90 D17052AA


Enabling SMB Digitally Sign Communications

The Server Message Block (SMB) protocol provides the basis for file and print sharing and many other networking operations, such as remote Windows administration. To prevent man-in-the-middle attacks that modify SMB packets in transit, follow the instructions below to enable the SMB Digitally Sign Communications.


1 Type **Edit group policy** from the Windows search bar to locate the **Local group policy editor** control panel.




The Local Group Policy Editor window displays.



Navigate to Computer Configuration > Windows Settings > Security Settings > Local Policies > Security Options.




3 Select Microsoft network server: Digitally sign communications (always).




E-92 D17052AA

The following window appears.



4 Select **Enabled** and **Apply**. The Confirm window appears.




5 Select **Yes.** Then select **OK** to close the window.


Enabling Audit Object Access

This audit object access policy allows the Windows system to record an event for the CytoFLEX mosaic database if the security auditing setting for both the CytExpert for Spectral database and CytExpert for Spectral service database are enabled. When the CytExpert for Spectral database files and/or the CytExpert for Spectral service database are deleted, an event will be automatically created in the Windows security logs.

Type Local Security Policy from the Windows search bar to locate the Local Security Policy control panel.




2 Double-click Local Security Policy to open the Local Security Policy window.




E-94 D17052AA


3 Navigate to Security Settings > Local Policies > Audit Policy.




4 Double-click Audit object access to open the Audit object access Properties window.



5 Select Success and Failure, and select OK to confirm.




The security setting for audit object access is enabled.




6 Close the Local Security Policy window.


E-96 D17052AA

Navigate to Windows (C:) > ProgramData > Beckman Coulter > CytExpert for Spectral > Database.




8 Right-click the **CytExpert.db** database file, and select **Properties** from the drop-down menu to open the CytExpert.db Properties window.






E-98 D17052AA


9 Select the **Security** tab, and select **Advanced**.



The following window appears.




10 Select the Auditing tab and select Add.




NOTE You must be an administrator or have been given the appropriate privileges to view the auditing properties of this object.

The following window appears.




E-100 D17052AA


11 Configure the security settings as below.




a. Click **Select a principal**. The Select User or Group window appears.




b. Enter **Everyone** in the object name field and select **OK**.



c. Select the Type **All**.




d. Select **Show advanced permissions**. The following window appears.




E-102 D17052AA

e. Set the advanced permissions and Select **OK**.




The security settings for CytoFLEX for Spectral database is enabled.




12 Select ok.

Viewing Windows Security Logs

 $\label{eq:total_problem} \textbf{1} \quad \text{Type } \textbf{Event Viewer} \text{ from the Windows search bar to open the Event Viewer window.}$




2 Navigate to Windows Logs > Security.



E-104 D17052AA


3 Select **Find** from the right Actions area and type **CytExpert.db** to quickly find the security logs for CytExpert for Spectral database.




Type **CytoCtrlService.db** and select **Find Next** to quickly find the security logs for the CytoCtrlService database.

Enabling Installation Restriction


1 Type Edit group policy from the Windows search bar to locate the Edit group policy editor control panel.



2 Select Edit group policy editor. The Local Group Policy Editor window displays.




Navigate to Computer Configuration > Administrative Templates > Windows Components > Windows Installer.




E-106 D17052AA

4 Select Allow user control over installs.



The following window appears.



- 5 Select Disabled and Apply.
- **6** Select **OK** to close the window.

7 Disable Always install with elevated privileges. Refer to Steps 3-5.

Enabling Microsoft Defender Firewall

IMPORTANT The Instrument cannot work properly if the necessary network ports are disabled. The communication between the CytExpert for Spectral software and the Instrument requires the following network ports to be enabled.

Outbound:

- TCP/192.168.76.18:1521
- TCP/192.168.76.18:1522
- TCP/192.168.76.18:1523
- TCP/192.168.76.18:1524
- TCP/192.168.76.18:1525
- TCP/192.168.76.18:60298
- Type Windows Defender Firewall from the Windows search bar to locate the Windows Defender Firewall app.




E-108 D17052AA

2 Select Turn Windows Defender Firewall on or off.




3 Select Turn on Windows Defender Firewall for both Private network and Public network settings.



4 Select OK.


Enabling Microsoft Defender

- 1 Type Windows Defender settings from the Windows search bar.
- **2** Double-click **Windows Defender settings** to open the Windows Security window.
- 3 Select **Virus & threat protection**, and select **Manage settings**. The Virus & threat protection settings window appears.



E-110 D17052AA


4 Ensure the Real-tine protection is on.




Enabling Network Time Protocol

To ensure the accuracy of date and time information included in logs, follow the instructions below to enable the Network Time Protocol.

1 Type Change the date and time from the Windows search bar.




2 Double-click Change the date and time.




E-112 D17052AA

The Date and Time Settings window appears.




NOTE Ensure that the **Set time zone automatically** and **Set time zone automatically** are on.






The Date and Time window displays.




E-114 D17052AA

4 Select the Internet Time tab from the Date and Time window, and then select **Change settings** from the Internet Time tab.



The Internet Time Settings window displays.



D17052AA E-115

5 Select **Synchronize with an internet time server**, and select **OK**. The PC is set to automatically synchronize with the internet time server.

Remote Access

IMPORTANT Use BeckmanConnect for remote control over the CytoFLEX instrument. For instructions on using BeckmanConnect, refer to the BeckmanConnect website.

E-116 D17052AA

APPENDIX F

Table of Hazardous Substances

Table of Hazardous Substances

The Hazardous Substances Names and Concentration is shown in Table F.1 and Table F.2.

NOTE Refer to the *CytoFLEX Platform Instructions for Use* manual, APPENDIX H, Table of Hazardous Substances for the tables of hazardous substances of the CytoFLEX S and CytoFLEX LX.

D17052AA F-1

Table F.1 有毒有害物质名称及含量的标识格式 Table of Hazardous Substances Name and Concentration [CytoFLEX mosaic 63]

电子电气产品号码 EEP Part Number: D10641		产品名称 Product Name: CytoFLEX mosaic Spectral Detection Module 63 产品型号 Product Model Number: CytoFLEX mosaic								
部件名称				有毒 ^注	有害物质或元	素 Hazardous S	ubstances Nan	ne		
Component Name	铅 (Pb)	汞 (Hg)	镉 (Cd)	六价铬 (Cr ⁶⁺)	多溴联苯 (PBB)	多溴二苯醚 (PBDE)	邻苯二甲 二乙酯 (DEHP)	乙血压 (BBP)	舒张压 (DBP)	二异丁基苯酚 (DIBP)
金属件 Metal Part	0	0	0	0	0	0	0	0	0	0
塑料件 Plastic Part	0	0	0	0	0	0	0	0	0	0
玻璃制品 Glass Part	0	0	0	0	0	0	0	0	0	0
橡胶制品 Rubber Part	0	0	0	0	0	0	0	0	0	0
印刷电路板 Printed Circuit Assembly	0	0	0	0	0	0	0	0	0	0
金属(底盘/螺钉/组件/护罩) Metal (Chassis/Screw/ Subassemblies/Shields)	0	0	0	0	0	0	0	0	0	0
线缆及连接器 Wire/Cable/ Connectors	0	0	0	0	0	0	0	0	0	0
其他电子件 Other Electrical	0	0	0	0	0	0	0	0	0	0
标签 Labeling	0	0	0	0	0	0	0	0	0	0
包装材料 Package	0	0	0	0	0	0	0	0	0	0

This table is prepared in accordance with the provisions of SJ/T 11364

- 〇:表示该有毒有害物质在该部件所有均质材料中的含量均在GB/T 26572标准规定的限量要求以下
- x: 表示该有毒有害物质至少在该部件的某一均质材料中的含量超出GB/T 26572标准规定的限量要求
- (企业可在此处,根据实际情况对上表中打"×"的技术原因进行进一步说明)
- O: Indicates that the toxic or hazardous substances contained in all of the homogenous materials for this part is below the limit requirements in GB/T 26572.
- X: Indicates that the toxic or hazardous substance contained in at least one of the homogenous materials used for this part in above the limit requirement in GB/T 26572.
- (Enterprise may further provide in this box technical explanation for marking "X" based on their actual conditions)

Table F.2 有毒有害物质名称及含量的标识格式 Table of Hazardous Substances Name and Concentration [CytoFLEX mosaic 88]

电子电气产品号码 EEP Part Number: D10638		产品名称 Product Name: CytoFLEX mosaic Spectral Detection Module 88 产品型号 Product Model Number: CytoFLEX mosaic								
部件名称		有毒有害物质或元素 Hazardous Substances Name								
Component Name	铅 (Pb)	汞 (Hg)	镉 (Cd)	六价铬 (Cr ⁶⁺)	多溴联苯 (PBB)	多溴二苯醚 (PBDE)	邻苯二甲 二乙酯 (DEHP)	乙血压 (BBP)	舒张压 (DBP)	二异丁基苯酚 (DIBP)
金属件 Metal Part	0	0	0	0	0	0	0	0	0	0
塑料件 Plastic Part	0	0	0	0	0	0	0	0	0	0
玻璃制品 Glass Part	0	0	0	0	0	0	0	0	0	0
橡胶制品 Rubber Part	0	0	0	0	0	0	0	0	0	0
印刷电路板 Printed Circuit Assembly	0	0	0	0	0	0	0	0	0	0
金属(底盘/螺钉/组件/护罩) Metal (Chassis/Screw/ Subassemblies/Shields)	0	0	0	0	0	0	0	0	0	0
线缆及连接器 Wire/Cable/ Connectors	0	0	0	0	0	0	0	0	0	0
其他电子件 Other Electrical	0	0	0	0	0	0	0	0	0	0
标签 Labeling	0	0	0	0	0	0	0	0	0	0
包装材料 Package	0	0	0	0	0	0	0	0	0	0

This table is prepared in accordance with the provisions of SJ/T 11364

- 〇:表示该有毒有害物质在该部件所有均质材料中的含量均在GB/T 26572标准规定的限量要求以下
- x: 表示该有毒有害物质至少在该部件的某一均质材料中的含量超出GB/T 26572标准规定的限量要求
- (企业可在此处,根据实际情况对上表中打"×"的技术原因进行进一步说明)
- O: Indicates that the toxic or hazardous substances contained in all of the homogenous materials for this part is below the limit requirements in GB/T 26572.
- X: Indicates that the toxic or hazardous substance contained in at least one of the homogenous materials used for this part in above the limit requirement in GB/T 26572.
- (Enterprise may further provide in this box technical explanation for marking "X" based on their actual conditions)

Table of Hazardous Substances

Table of Hazardous Substances

F-4 D17052AA

Abbreviations

The following list is a composite of the symbols, abbreviations, acronyms, and reference designators either used in this manual or related to the information in it. When the same abbreviation (or reference designator) is used for more than one word (or type of component), all meanings relevant to this manual are included, separated by semicolons.

' — foot

" — inch

% — percent

°C — degrees Celsius

°F — degrees Fahrenheit

± — plus or minus

< — less than

> — greater than

 \leq — less than or equal to

µ — micron

μL — microliters

µm — micrometer

A — ampere

AC — alternating current

APC — Allophycocyanin

APC-A700 — Allophycocyanin-Alexa Fluor 700 tandem dye

APC-A750 — Allophycocyanin-Alexa Fluor 750

APC-Cy7 — Allophycocyanin-Cyanin 7

API — Application Programming Interface

Acq. — Acquisition

BCI — Beckman Coulter Incorporated

BMP — bitmap

BP — band-pass filter

CDRH — Center for Devices and Radiological Health

CFSE — carboxyfluorescein succinmidyl ester

cm — centimeters

CSV — comma separated value

CV — coefficient of variation

DNA — deoxyribonucleic acid

DW — deep well

ECD — Energy Coupled Dye

EFUP — Environmentally friendly Use Period

EMC — Electromagnetic Compatibility

EMF — enhanced metafile format

EMR — electromagnetic radiation

FAPD — Fiber Array Photo Detector

FCS — flow cytometry standard

FITC — Fluorescein isothiocyanate

FLn — Fluorescence Number

FSC — forward scatter

GB — gigabyte

GHz — gigahertz

Gr Wt — gross weight

H — humidity

Hz — hertz

IEC — International Electrotechnical Commission

IR — infrared

kg — kilograms

KO — Krome Orange

LED — light emitting diode

L — liter

LJ — Levey-Jennings

LWH — length, width, height

m — meter

MB — megabyte

MFI — median fluorescence intensity

MHz — megahertz

min — minute

mL — milliliter

mm — millimeter

mW — milliwatt

NA — numerical aperture

NaCIO — sodium hypochlorite solution

NaN₃ — sodium azide

nm — nanometer

Nt Wt - net weight

PB — Pacific Blue dye

PC5 — Phycoerythrin-Cy5 tandem dye

PC5.5 — Phycoerythrin-Cy5.5 tandem dye

PC7 — Phycoerythrin-Cy7 tandem dye

PE — Phycoerythrin

PEEK — polyether ether ketone

PerCP — Peridinin-Chlorophyll

PI — Propidium Iodide

PN — part number

QC — quality control

RAM — random access memory

rCV — robust coefficient of variation

RH — relative humidity

RoHS — Restriction of Hazardous Substances Directive

RPTM — real-time messaging protocol

S/N — serial number

SNR — signal to noise ratio

SSC — side scatter

USB — universal serial bus

UV — ultraviolet

V — volts

VA — volt-ampere

VAC — voltage alternating current

VSSC — violet side scatter

WDM — wavelength division multiplexer

W — watts

Abbreviations-2 D17052AA

Index

	routine maintenance, 11-1
create, 7-2, 7-5	4 L waste container (CytoFLEX)
data, 6-33	emptying, 12-1
routine maintenance, 12-1	
	Α
Symbols	Α
,	define, Abbreviations-1
define, Abbreviations-1	about this manual, xxvii
II	AC
define, Abbreviations-1	define, Abbreviations-1
°C	account
define, Abbreviations-1	locking, 4-6
°F	unlock, 2-30
define, Abbreviations-1	account menu, 2-25
μ	account policies, 2-36
define, Abbreviations-1	acoustic noise level, 1-13
μL	Acq.
define, Abbreviations-1	define, Abbreviations-1
μm	acquisition screen, 2-4
define, Abbreviations-1	collection, 2-6
<u>+</u>	navigation, 2-5
define, Abbreviations-1	plot area, 2-10
<	status bar, 2-12
define, Abbreviations-1	test tubes, 2-8
- laCon Alaboratetta on 1	acquisition settings
define, Abbreviations-1 ≤	to configure, 6-19 add
define, Abbreviations-1	
define, Abbreviations-1	plate type plate type
Numerics	add, 2-40 add a plate type, 2-40
10 L sheath fluid cubitainer	adding
replacing, 12-1	Deep Clean solution, 12-1
10 L waste cubitainer (CytoFLEX LX)	standardization item, 5-23
emptying, 12-1	adjust
4 L sheath fluid container	spectral unmixing, 7-15
to clean, 11-1	adjust compensation
4 L sheath fluid container (cleaning)	manually, 7-27
routine maintenance, 11-1	adjust spectral unmixing, 7-15
4 L waste container	adjusting
to clean, 11-1	autogate movement, 6-17
4 L waste container (cleaning)	gain, 6-19

spectral unmixing, 7-15, 7-27	beam shaping
threshold, 6-23	laser, 3-2
unmixing settings, 8-2	Beckman Coulter Customer Support Center
adjusting auto gates, 6-17	contacting, ii
adjustment	BMP
nonscheduled, 12-2	define, Abbreviations-1
routine, 12-1	bold face font
advanced menu, 2-25	conventions, xxix
alerts	BP
caution define, v	define, Abbreviations-1
important define, v	
note define, v	C
warning define, v	C
analysis	calculate
screen, 8-1	unmixing values, 7-10
analysis screen, 2-13	calculate unmixing values, 7-10
to open, 8-2	calculating
analyzing	concentration, 8-2
data, 6-33	sample volume, 8-2
APC	unmixing values, 7-10
define, Abbreviations-1	calibrating
APC-A750	plate position [with plate loader], 12-2
define, Abbreviations-1	sample flow rate, 12-2
APC-Cy7	sample flow rate with plate loader, 12-2
define, Abbreviations-1	call center, contact information, ii
API	caution
define, Abbreviations-1	define, v
applying	CDRH
standardization in QC, 5-32, 5-35	define, Abbreviations-1
auto fine adjustment, 7-27	CFSE
auto gates	define, Abbreviations-1
adjusting, 6-17	change
creating, 6-17	event rate setting, 12-2
auto recalculate	password, 2-32
turn off, 6-17	change the event rate setting, 12-2
turn on, 6-17	changing
auto shutdown, 9-2	mixer and backflush settings, 12-2
auto startup, 9-2	tube name, 6-19
autogate	changing sample probe
adjusting movement, 6-17	from plate loader to single tube sample
automated software feature, 3-4	station, 12-2
	from single tube sample station to plate loader, 12-1
В	channel
backflush settings	
	to set, 6-7 characteristics
to change, 12-2 backup, 10-8	
<u> </u>	performance, 1-16 check
band-pass filter, Abbreviations-1 BCI	waste and reagent levels (10 L fluid
define, Abbreviations-1	cubitainers), 4-1
MOTITIO, FIGURE VIGILIAND I	Cucitanici 0/1 1 1

waste and reagent levels (4 L fluid	acquisition settings, 6-19
containers), 4-1	confirm
cleaning	unmixing results, 7-12
4 L sheath fluid container, 11-1	confirm unmixing results, 7-12
4 L waste container, 11-1	confirming results, 5-15
nonscheduled, 11-1	connecting the instrument, A-6
routine, 11-1	contact information, Beckman Coulter
sample station, 11-1	Customer Support Center, ii
surfaces, 11-1	Conventional Mode, 12-2
cleaning solution	conventions
to prepare, 9-1	bold face font, xxix
cleaning the 4 L sheath fluid container	the term press, xxix
routine maintenance, 11-1	the term select, xxix
cleaning the 4 L waste container	conventions used, 1-xxix
routine maintenance, 11-1	copying
cleaning the sample station	experiments, 8-1
routine maintenance, 11-1	previously acquired experiment, 8-1
closed file system	create
electronic record management, B-2	spectral unmixing experiment, 7-2
cm	create a new experiment from a template
define, Abbreviations-1	start page operations, 2-3
collecting	create a new screen
data, 6-33	start page operations, 2-3
QC data, 5-7	create a spectral unmixing experiment
QC data with plate loader, 5-11	start page operations, 2-3
collection	create spectral unmixing experiment, 7-2
acquisition screen, 2-6	creating
collection conditions	dot plot overlays, 8-2
to set, 6-24	heat map (plate loader), 6-33
compensation	histogram overlays, 8-2
to manually adjust, 7-27	Levey-Jennings charts, 5-18
components	new experiment, 6-1
Cytometer, 1-2	new experiment with plate loader, 6-1, 6-2
Fluid Containers/Cubitainers, 1-2	plots and gates, 6-8
fluidics system, 1-6	spectral unmixing experiment, 7-2
instrument, 1-1	spectral unmixing experiment with plate
main, 1-2	loader, 7-5
optical, 1-4	users, 2-28
plate holder, 1-6	creating auto gates, 6-17
plate loader, 1-6	credential file
sample injection mode control kit, C-1	export, 2-35
Workstation, 1-2	CSV
concentration	define, Abbreviations-1
calculation, 8-2	custom statistics
sample volume, 8-2	setting, 6-19
concluding	customer-replaceable parts
experiment, 6-46, 7-14	ordering, F-1
configuration	customized parameters
system, 1-7	setting, 6-17
configuring	CV
CV1111EM1111E	~ 1

define, Abbreviations-1	routine maintenance, 11-1
cyber security, E-52	daily start up with plate loader
backing up your recovery key, E-13	routine maintenance, 11-1
BitLocker decryption, E-17	data
BitLocker encryption, E-6	analyzing, 6-33
BitLocker keys, E-2	collecting, 6-33
changing encryption strength, E-2	exporting, 6-33
disabling Trellix application control, E-50	sampling, 6-33
drive encryption, E-2	to import, 8-1
enabling Trellix application control, E-35	data sheets, material safety
exporting certificate, E-39	how to order, 1-18
good practices in cyber security, E-1	data storage
Importing certificate, E-47	description, 3-4
installing Trellix Application Control, E-28	Deep Clean
malware detection software, E-27	procedure, 11-1
operating system patching, E-57, E-62	routine maintenance, 11-1
recovering your BitLocker key, E-25	Deep Clean solution
remote access, E-116	to add, 12-1
suspending BitLocker drive encryption, E-21	to prepare, 12-1
switching Trellix Application control to	default password, 2-28
update mode, E-48	define
user management, E-63	%, Abbreviations-1
CytExpert	caution, v
installation options, A-7	important, v
CytExpert for Spectral software	note, v
reinstall, A-8	quality control (QC), 5-1
CytExpert nano certificate, E-39	symbols, -vii
Cytobank, 2-18, 2-19	warning, v
Dimensions, 1-13	delete
instrument specifications	plate type
dimensions, 1-13	plate type
CytoFLEX mosaic 63, 1-2, 1-3	delete, 2-41
CytoFLEX mosaic 88, 1-3	spectral library files, 7-26
Dimensions, 1-13	delete a plate type, 2-41
instrument specifications	delete spectral library files, 7-26
dimensions, 1-13	deleting
Dimensions, 1-13	heat map (plate loader), 6-33
instrument specifications	standardization items, 5-37
dimensions, 1-13	users, 2-29
Cytometer, 1-2, 1-3	detector (WDM)
Cytometer cannot be turned on, 10-7	optical components, 1-4
cytometer menu, 2-21	detector unit
	wavelength division multiplexor (WDM), 1-4
D	disinfection
	surfaces, 11-1
daily shutdown	disposal
routine maintenance, 11-1	precaution, 10-6
daily shutdown with plate loader	waste, A-4
routine maintenance, 11-1	DNA
daily start up	

define, Abbreviations-1	importing an experiment/template, B-3
document overview, xxvii	log, B-3
dot plot overlays	software menu, B-1
to create, 8-2	system operation log, B-3
downloading security updates, E-57	user management, B-4
drive decryption	user management operation log, B-3
enabling BitLocker, E-17	electronic signature
drive encryption, E-2	electronic record management, B-4
BitLocker keys, E-2	retract, B-4
BitLocker recovery, E-25	EMC
changing encryption strength level for	define, Abbreviations-1
BitLocker, E-2	EMF
enabling BitLocker, E-6	define, Abbreviations-1
dual-parameter plots, 2-39	emptying
duplicate	10 L waste cubitainer (CytoFLEX LX), 12-1
plate type	4 L waste container (CytoFLEX), 12-1
plate type	EMR
duplicate, 2-40	define, Abbreviations-1
duplicate a plate type, 2-40	environment validation
duplicating	installation, A-2
standardization item parameters, 5-37	environment validation (installation)
DW	power source, A-3
define, Abbreviations-1	temperature and humidity, A-4
	worktable, A-2
_	event rate setting
E	change, 12-2
ECD	exit the software
define, Abbreviations-1	start page operations, 2-3
edit	experiment
plate type	concluding, 6-46, 7-14
plate type	sign, B-4
edit,2-40	experiment (new)
edit a plate type, 2-40	creating, 6-1
editing	experiment directory management
standardization item parameters, 5-37	electronic record management, B-2
EFUP	experiment management
define, Abbreviations-1	electronic record management, B-2
electronic devices	experiment operation log
fuse, 1-8, 1-10	electronic record management, B-3
load button, 1-8, 1-10	experiment related operations
power switch, 1-8, 1-10	electronic record management, B-3
electronic record management	experiment saving, 6-46
closed file system, B-2	experiment selection
electronic signature, B-4	start page, 4-7
experiment directory management, B-2	experiment settings, 2-45
experiment management, B-2	experiment signature
experiment operation log, B-3	print, B-4
experiment related operations, B-3	experiments
exporting an experiment/template, B-3	to copy, 8-1
folder hierarchy management, B-2	export

credential file, 2-35	fluidics module
spectral library files, 7-25	fluidics system, 1-6
unmixing settings to unmixing matrix	fluidics system
files, 7-20	components, 1-6
unmixing settings to unmixing matrix files	Fluid Containers/Cubitainers, 1-6
from the spectral library, 7-22	fluidics module, 1-6
export spectral library files, 7-25	fluorescent light (FL)
exporting	collection, 3-2
data, 6-33	folder hierarchy management
FCS files, 6-39	electronic record management, B-2
heat map (plate loader), 6-33	four-quadrant gates
instrument settings, 6-44, 6-45	for dual-parameter plots, 2-40
plots of multiple tubes as picture files, 6-43	front cover
results, 8-3	removal, 12-1
standardization item, 5-37	FSC
statistics table of multiple tubes as a picture	define, Abbreviations-1
file, 6-43	functions
unmixing settings, 7-20	software, 2-1
exporting and experiment/template	fuse
electronic record management, B-3	electronic devices, 1-8, 1-10
exporting unmixing settings to unmixing	to replace, 12-2
matrix files, 7-20	-
exporting unmixing settings to unmixing	•
matrix files from the spectral	G
library, 7-22	gain
,,,	to adjust, 6-19
	gate settings, 2-47
F	gates
FAPD	to create, 6-8
define, Abbreviations-1	gates for dual-parameter plots
FCS	four-quadrant gates, 2-40
define, Abbreviations-1	lasso gates, 2-40
FCS files	polygon gates, 2-40
exporting, 6-39	rectangle gates, 2-40
files menu, 2-20	gates for single-parameter plots
fine adjustment, 7-27	line-segment gates, 2-40
auto fine adjustment, 7-27	vertical gates, 2-40
fine adjustment matrix, 7-28	
manual ting adjustment 7-27	GB
manual fine adjustment, 7-27	define, Abbreviations-1
fine adjustment matrix, 7-28	define, Abbreviations-1 Ghz
fine adjustment matrix, 7-28 FITC	define, Abbreviations-1
fine adjustment matrix, 7-28	define, Abbreviations-1 Ghz define, Abbreviations-1 Gr. Wt.
fine adjustment matrix, 7-28 FITC	define, Abbreviations-1 Ghz define, Abbreviations-1
fine adjustment matrix, 7-28 FITC define, Abbreviations-1	define, Abbreviations-1 Ghz define, Abbreviations-1 Gr. Wt.
fine adjustment matrix, 7-28 FITC define, Abbreviations-1 FLn	define, Abbreviations-1 Ghz define, Abbreviations-1 Gr. Wt. define, Abbreviations-1
fine adjustment matrix, 7-28 FITC define, Abbreviations-1 FLn define, Abbreviations-1	define, Abbreviations-1 Ghz define, Abbreviations-1 Gr. Wt. define, Abbreviations-1 graphic and gating styles, 2-39
fine adjustment matrix, 7-28 FITC define, Abbreviations-1 FLn define, Abbreviations-1 flow cell to prime, 12-2	define, Abbreviations-1 Ghz define, Abbreviations-1 Gr. Wt. define, Abbreviations-1 graphic and gating styles, 2-39 graphics
fine adjustment matrix, 7-28 FITC define, Abbreviations-1 FLn define, Abbreviations-1 flow cell to prime, 12-2 Fluid Containers, 1-2, 1-3	define, Abbreviations-1 Ghz define, Abbreviations-1 Gr. Wt. define, Abbreviations-1 graphic and gating styles, 2-39 graphics printing, 6-46
fine adjustment matrix, 7-28 FITC define, Abbreviations-1 FLn define, Abbreviations-1 flow cell to prime, 12-2 Fluid Containers, 1-2, 1-3 Fluid Containers/Cubitainers	define, Abbreviations-1 Ghz define, Abbreviations-1 Gr. Wt. define, Abbreviations-1 graphic and gating styles, 2-39 graphics
fine adjustment matrix, 7-28 FITC define, Abbreviations-1 FLn define, Abbreviations-1 flow cell to prime, 12-2 Fluid Containers, 1-2, 1-3	define, Abbreviations-1 Ghz define, Abbreviations-1 Gr. Wt. define, Abbreviations-1 graphic and gating styles, 2-39 graphics printing, 6-46

define, Abbreviations-1	standardization item, 5-37
hazard	importing an experiment/template
laser, 10-1	electronic record management, B-3
hazards	information
laser beam, 10-2	ordering, 1-12
radiation, 10-2	initializing the instrument, 4-7
hazards/precautions, 10-1	inspecting materials
heat map (plate loader)	unpacking, A-5
creating, 6-33	inspection
deleting, 6-33	liquid flow tubing, 12-2
exporting, 6-33	pre-boot, 4-1
modifying, 6-33	inspection when unpacking, A-5
refreshing, 6-33	installation
help menu, 2-26	environment validation, A-2
help, Beckman Coulter Customer Support	instrument, A-1
Center, ii	software, A-7
histogram	installation category, 1-13
types of display, 3-4	installation environment validation
histogram overlays	ventilation and cleaning, A-2
to create, 8-2	installation options
how to	CytExpert, A-7
use your manual, xxvii	installing
humidity and temperature	Plate Loader module, 12-2
installation environment validation, A-4	installing security updates, E-62
Hz	installing the instrument, A-6
define, Abbreviations-1	installing the instrument and connecting the
	equipment, A-6
I	unpacking inspection, A-6
	instrument
IEC	components, 1-1
define, Abbreviations-2	initializing, 4-7
import	inspecting materials, A-5
spectral library files, 7-23	installation, A-1
unmixing settings from the spectral	installing, A-6
library, 7-17	maintenance, 12-1
unmixing settings from unmixing matrix	shut off, 9-1
files, 7-15	transportation and storage, A-1
import and export unmixing matrix, 7-15	unpacking, A-5, A-6
import spectral library files, 7-23	instrument characteristics, 1-13
import unmixing settings from the spectral	instrument components
library, 7-17	Cytometer, 1-2
import unmixing settings from unmixing	Fluid Containers/Cubitainers, 1-2
matrix files, 7-15	Workstation, 1-2
important	instrument configuration
define, v	unpacking inspection, A-5
importing	instrument settings
data, 8-1	exporting, 6-45
instrument settings, 6-44	importing, 6-44
lot-specific target values, 5-4	importing and exporting, 6-44
previously acquired data, 8-1	instrument specifications, 1-13

instrument transport or storage	routine maintenance, 12-2
to prepare, 11-1	LJ
IR	define, Abbreviations-2
define, Abbreviations-2	load button
items	electronic devices, 1-8, 1-10
adding for standardization, 5-23	load sample and record data, 7-9
deleting for standardization, 5-37	locking
duplicating for standardization, 5-37	account, 4-6
editing for standardization, 5-37	log
exporting for standardization, 5-37	electronic record management, B-3
importing for standardization, 5-37	user management operation, 2-39
	log in
• •	software, 4-3
K	log menu, 2-26
kg	lot-specific target values
define, Abbreviations-2	to import, 5-4
KO	LWH
define, Abbreviations-2	define, Abbreviations-2
	define, nooreviations 2
L	NA.
	M
L 1.C All tree o	m
define, Abbreviations-2	define, Abbreviations-2
laser	main components, 1-2
hazards, 10-1	maintenance
settings, 6-19	instrument, 12-1
target power settings (CytoFLEX LX	See also routine maintenance; nonscheduled
only), 6-19	maintenance
laser beam	maintenance reminder
hazards, 10-2	to manage, 12-1
laser delay	malware detection software, E-27
to set, 12-2	manage
laser lines, 1-14	spectral library, 7-23
laser target power settings (CytoFLEX LX	manage spectral library, 7-23
only), 6-19	delete spectral library files, 7-26
laser wavelength, 1-14	export spectral library files, 7-25
lasso gates	import spectral library files, 7-23
for dual-parameter plots, 2-40	management
launching the software, 2-1	user, E-63
LED	managing
define, Abbreviations-2	maintenance reminder, 12-1
Levey-Jennings charts	manual
to create, 5-18	about, xxvii
light collection, separation and measurement	adjust compensation, 7-27
description, 3-2	manual fine adjustment, 7-27
line segment gates	manuals
for single-parameter plots, 2-40	updating, iii
liquid flow tubing	material safety data sheets, 1-18
inspection, 12-2	material safety data sheets (MSDS/SDS)
liquid flow tubing inspection	how to order, 1-18

maximum altitude, 1-13	N
MB	NA
define, Abbreviations-2	define, Abbreviations-2
menu	NaClO
account, 2-25	define, Abbreviations-2
acquisition and analysis screen, 2-20	NaN3
advanced, 2-25	define, Abbreviations-2
cytometer, 2-21	navigation
files, 2-20	acquisition screen, 2-5
help, 2-26	QC screen, 2-17
log, 2-26	new experiment
QC/Standardization, 2-24	creating, 6-1
settings, 2-24	new experiment with plate loader
signature, 2-26	
software, 2-18	creating, 6-1, 6-2
menu tree	nm define Abbreviations 2
software, 2-18	define, Abbreviations-2
MFI	nonscheduled cleaning, 11-1 nonscheduled maintenance
define, Abbreviations-2	
MHz	adding the Deep Clean solution, 12-1
define, Abbreviations-2	calibrating the sample flow rate, 12-2
min	calibrating the sample flow rate with plate
define, Abbreviations-2	loader, 12-2
mixer settings	preparing the instrument for storage or
to change, 12-2	transport, 11-1
mL	replacing the fuse, 12-2
define, Abbreviations-2	nonscheduled replacement/adjustment, 12-2
mm	note
define, Abbreviations-2	define, v
modifying	Nt. Wt.
heat map (plate loader), 6-33	define, Abbreviations-2
users, 2-30	
well settings with plate loader, 6-7	0
modifying well settings, 6-7	
applying existing well settings to additional	open a spectral unmixing experiment
wells, 6-7	start page operations, 2-3 open an experiment
moving well location [with plate loader], 6-7	1 1
movement	start page operations, 2-3 open software, 4-3
adjusting for autogates, 6-17	<u>.</u>
moving	opening
well location [with plate loader], 6-7	analysis screen, 8-2
MSDS (material safety data sheets)	operating system patching, E-52
how to order, 1-18	optical
multiple tube plots	components, 1-4
	optical bench
export as picture files, 6-43	optical component, 1-4
multiple tube statistics table	optical components
export as picture files, 6-43	detector (WDM), 1-4
mW define Abbrevietions ?	optical bench, 1-4
define, Abbreviations-2	optical fiber, 1-4

optical fiber, 1-4	Plate Loader module
ordering	installing, 12-2
customer-replaceable parts, F-1	removing, 12-2
ordering information, 1-12	plate loader PEEK tubing
overlay dot plots, 2-39	replacing, 12-1
overlay histograms, 2-39	plate loader PEEK tubing (replacing)
overview	routine maintenance, 12-1
document, xxvii	plate position [with plate loader]
	calibrating, 12-2
n	plate type, 2-40
P	add, 2-40
page setup settings, 2-47	plate type library, 2-40
parameters	add a plate type, 2-40
description, 3-4	delete a plate type, 2-41
duplicating standardization items, 5-37	duplicate a plate type, 2-40
editing standardization items, 5-37	edit a plate type, 2-40
TIME, 3-4	plot area
password	acquisition screen, 2-10
change, 2-32	plot display conditions
default, 2-28	to set, 6-24
forgot password, 2-33	plot settings, 2-46
reset, 2-31	plots
PB	dual-parameter, 2-39
define, Abbreviations-2	overlay dot plots, 2-39
PC5	overlay histograms, 2-39
define, Abbreviations-2	single-parameter, 2-39
PC5.5	to create, 6-8
define, Abbreviations-2	to set, 8-2
PC7	types of display, 3-4
define, Abbreviations-2	PN
PE	define, Abbreviations-2
define, Abbreviations-2	policies
PEEK	account, 2-36
define, Abbreviations-2	pollution degree, 1-13
PerCP	polygon gates
define, Abbreviations-2	for dual-parameter plots, 2-40
performance characteristics, 1-16	position
sample injection mode control knob, C-1	sample tube holder, 1-6
performing	power source
daily quality control, 5-1	installation environment validation, A-3
PI	power source inspection
define, Abbreviations-2	pre-boot, 4-2
picture files	power switch
to export, 6-43	electronic devices, 1-8, 1-10
plate holder	pre-boot inspection, 4-1
components, 1-6	check waste and reagent levels (10 L fluid
plate holder with plate loader	cubitainers), 4-1
replacing, 12-2	check waste and reagent levels (4 L fluid
plate loader	containers), 4-1
components, 1-6	

power source inspection, 4-2	define, Abbreviations-2
workstation connections inspection, 4-2	result manager, 5-18
precaution	standardization, 5-19
disposal, 10-6	QC data
precautions	collecting, 5-7
safety, -vi	QC data with plate loader
precautions/hazards, 10-1	collecting, 5-11
preparation process	QC experiment screen, 2-16, 2-17
QC sample (CytoFLEX Daily IR QC	QC menu tree
Fluorospheres), 5-3	software, 2-19
QC sample (CytoFLEX Daily IR QC	QC report screen, 2-16
Fluorospheres) with plate loader, 5-4	QC sample
QC sample (CytoFLEX Daily QC	preparation, 5-3
Fluorospheres), 5-18	QC sample (CytoFLEX Daily IR QC
QC sample (CytoFLEX Ready to Use Daily QC	Fluorospheres)
Fluorospheres), 5-3, 5-4	preparation process, 5-3
prepare	QC sample (CytoFLEX Daily IR QC
unmixing sample, 7-1	Fluorospheres) with plate loader
preparing	preparation process, 5-4
Deep Clean solution, 12-1	QC sample (CytoFLEX Daily QC Fluorospheres)
instrument for storage or transport, 11-1	preparation process, 5-18
QC sample, 5-3	QC sample (CytoFLEX Ready to Use Daily QC
QC sample with plate loader, 5-3	Fluorospheres)
spectral unmixing sample, 7-2	preparation process, 5-3, 5-4
Standardization sample, 5-18	QC sample with plate loader
the cleaning solution, 9-1	preparation, 5-3
preparing QC sample	QC screen
required materials, 5-3, 5-18	navigation, 2-17
preparing QC sample with plate loader	QC software menu tree, 2-19
required materials, 5-3	QC/Standardization menu, 2-24
press	Quality Control, 5-2
conventions, xxix	quality control, 5-2
previously acquired data	perform daily, 5-1
to import, 8-1	quality control (QC)
previously acquired experiment	define, 5-1
to copy, 8-1	
priming	n
flow cell, 12-2	R
print	radiation
experiment signature, B-4	hazards, 10-2
printing	RAM
graphics, 6-46	define, Abbreviations-2
procedure	rCV
Deep Clean, 11-1	define, Abbreviations-2
product description, 1-1	reagent
- · · · · · · · · · · · · · · · · · · ·	limitations, 1-18
0	reagent levels (10 L fluid cubitainers)
Q	check, 4-1
QC	reagent levels (4 L fluid containers)
applying standardization, 5-32, 5-35	check, 4-1

reagent limitations, 1-18	define, Abbreviations-2
recovery key, E-13	right-side cover
rectangle gates	removal, 12-1
for dual-parameter plots, 2-40	RoHS
refreshing	define, Abbreviations-2
heat map (plate loader), 6-33	role management, 2-36
reinstall	routine cleaning, 11-1
CytExpert for Spectral software, A-8	routine maintenance
removal procedures	cleaning the 4 L sheath fluid container, 11-1
front cover, 12-1	cleaning the 4 L waste container, 11-1
right-side cover, 12-1	cleaning the sample station, 11-1
removing	daily startup and shutdown, 11-1
Plate Loader module, 12-2	daily startup and shutdown with plate
replacement	loader, 11-1
nonscheduled, 12-2	Deep Clean, 11-1
routine, 12-1	liquid flow tubing inspection, 12-2
replacing	replacing the plate loader PEEK tubing, 12-1
10 L sheath fluid cubitainer, 12-1	replacing the sample peristaltic pump
fuse, 12-2	tubing for sample loading, 12-1
plate holder with plate loader, 12-2	replacing the sheath fluid filter, 12-1
plate loader PEEK tubing, 12-1	routine replacement/adjustment, 12-1
sample peristaltic pump tubing for sample	RPTM
loading, 12-1	define, Abbreviations-2
sheath fluid filter, 12-1	running
sheath fluid harness and/or waste	system startup program, 4-7
harness, 12-2	system startup program with plate
replacing the peristaltic pump for sample	loader, 4-7
loading	
routine maintenance, 12-1	C
replacing the sample probe for plate loader	S
routine maintenance, 12-1	S/N
replacing the sheath fluid filter	define, Abbreviations-2
routine maintenance, 12-1	safety
required materials	precautions, -vi
QC sample preparation, 5-3, 5-18	sample
QC sample preparation with plate loader, 5-3	checks before running, 6-33
reset	flow, description, 3-1
password, 2-31	sample acquisition position
restore, 10-8	sample tube holder position, 1-6
result manager	sample flow
QC, 5-18	description, 3-1
results	sample flow rate
confirming, 5-15	calibrating, 12-2
to export, 8-3	sample flow rate with plate loader
retention period	calibrating, 12-2
setting, B-4	sample injection control knob
retract	performance characteristics, C-1
electronic signature, B-4	sample injection mode
RH	selecting, 4-6
	to use manual injection mode, 4-6

to use plate loader injection mode, 4-7	acquisition, 2-4
to use semi-automatic injection mode, 4-6	analysis, 2-13, 8-1, 8-2
sample injection mode control kit	QC experiment, 2-16, 2-17
components, C-1	QC report, 2-16
sample loading position	software, 2-2
sample tube holder position, 1-6	Spectral Unmixing experiment, 2-14
sample peristaltic pump tubing	SDS (safety data sheets)
to replace, 12-1	how to order, 1-18
sample peristaltic pump tubing for sample	SDS/MSDS, 1-18
loading	select
to replace, 12-1	conventions, xxix
sample peristaltic pump tubing for sample	selecting
loading (replacing)	proper sample injection mode, 4-6
routine maintenance, 12-1	service, contact information, ii
sample probe	setting
to unclog, 12-1	channel, 6-7
sample probe change	collection conditions, 6-24
from plate loader to single tube sample	custom statistics, 6-19
station, 12-2	customized parameters, 6-17
from single tube sample station to plate	laser delay, 12-2
loader, 12-1	laser target power settings (CytoFLEX LX
sample probe for plate loader	only), 6-19
replacing, 12-1	plot display conditions, 6-24
sample probe for sample loading	retention period, B-4
to replace, 12-1	sample wells with sample wells, 6-2
sample probe for sample loading (unclogging)	signature, B-4
routine maintenance, 12-1	the plots and statistics, 8-2
sample probe switch	settings
from plate loader to single tube sample	acquisition, 6-19
station, C-1	experiment, 2-45
from single tube sample station to plate	gate, 2-47
loader, C-1	laser, 6-19
sample station, 1-6	page setup, 2-47
to clean, 11-1	plot, 2-46
sample station (cleaning)	software, 2-45
routine maintenance, 11-1	tube, 2-46
sample tube holder position	settings menu, 2-24
sample acquisition position, 1-6	sheath fluid cubitainer 10 L
sample loading position, 1-6	replacing, 12-1
standby position, 1-6	sheath fluid filter
sample wells with plate loader	to replace, 12-1
setting, 6-2	sheath fluid filter (replacing)
sampling	routine maintenance, 12-1
data, 6-33	shut off
save	the instrument, 9-1
spectral unmixing experiment, 7-14	shutdown
save spectral unmixing experiment, 7-14	automatic, 9-2
saving	shutdown daily
experiment, 6-46	routine maintenance, 11-1
screen	shutdown daily with plate loader

routine maintenance, 11-1	spectral unmixing experiment
side scatter (SS)	create, 7-2
collection, 3-2	save, 7-14
sign	to create, 7-2
experiment, B-4	spectral unmixing experiment screen, 2-14
signals	spectral unmixing controls, 2-15
processing of, operation principles, 3-3	spectral unmixing experiment with plate loader
signature	creating, 7-5
setting, B-4	spectral unmixing sample
signature menu, 2-26	prepare, 7-2
single-parameter plots, 2-39	SSC
SNR	define, Abbreviations-2
define, Abbreviations-2	Standardization
software	loading standardization sample, 5-19
functions, 2-1	Target value, 5-19
launching, 2-1	standardization, 5-18
log in, 4-3	applying in QC, 5-32, 5-35
main screen, 2-2	standardization item
open, 4-3	adding, 5-23
to start, A-19	deleting, 5-37
software installation, A-7	exporting, 5-37
starting the software, A-19	importing, 5-37
software menu, 2-18	standardization item parameters
acquisition and analysis screen menu, 2-20	duplicating, 5-37
electronic record management, B-1	editing, 5-37
software menu tree, 2-18	Standardization sample
software operations	preparation, 5-18
create a new experiment, 2-3	standardization target library
software screen, 2-2	deleting standardization target
start page, 2-2	settings, 5-37
start page operations, 2-3	duplicate standardization target
software settings, 2-45	settings, 5-37
specifications	editing standardization target settings, 5-37
characteristics, 1-13	exporting standardization target
instrument, 1-13	settings, 5-37
specimen	importing standardization target
See sample	settings, 5-37
spectral library	standby position
manage, 7-23	sample tube holder position, 1-6
spectral library files	start page
delete, 7-26	main software screen, 2-2
export, 7-25	selecting experiments, 4-7
import, 7-23	start page operations
spectral mode, 12-2	create a new experiment from a
spectral unmixing	template, 2-3
adjust, 7-15	create a spectral unmixing experiment, 2-3
to adjust, 7-15, 7-27	exit the software, 2-3
spectral unmixing controls	main software screen, 2-3
spectral unmixing experiment screen, 2-15	open a spectral unmixing experiment, 2-3

open an experiment, 2-3	target power settings (CytoFLEX LX only)
starting	laser, 6-19
the software, A-19	target value
startup	Standardization, 5-19
automatic, 9-2	temperature and humidity
startup daily	installation environment validation, A-4
routine maintenance, 11-1	storage, A-1
startup daily with plate loader	test tube
routine maintenance, 11-1	management module, 2-9
station	test tube management module, 2-9
sample, 1-6	test tubes
statistics	acquisition screen, 2-8
setting, 3-5	the term press
to set, 8-2	conventions, xxix
status bar	the term select
acquisition screen, 2-12	conventions, xxix
storage	threshold
temperature and humidity, A-1	to adjust, 6-23
storing the instrument, A-1	TIME parameter
support, Beckman Coulter customer, ii	description, 3-4
surface	to add
cleaning and disinfection, 11-1	plate type, 2-40
suspending BitLocker drive encryption, E-21	to use manual injection mode, 4-6
switch	to use plate loader injection mode, 4-7
Conventional Mode, 12-2, 12-3	to use semi-automatic injection mode, 4-6
Spectral Mode, 12-2	transporting the instrument, A-1
spectral mode, 12-7	Trellix application certificate, E-47
switch to Conventional Mode, 12-3	Trellix application control
switch to spectral mode, 12-7	disabling, E-50
Switching between Spectral Mode and	enabling, E-35
Conventional Mode, 12-2	update mode, E-48
switching sample probe	troubleshooting, 10-1
from plate loader to single tube sample	Cytometer cannot be turned on, 10-7
station, C-1	table, 10-6
from single tube sample station to plate	tube
loader, C-1	settings, 2-46
symbols	tube name
define, -vii	to change, 6-19
system configuration, 1-7	turn off
system operation log	auto recalculate, 6-17
electronic record management, B-3	turn on
system startup program	auto recalculate, 6-17
to run, 4-7	turning the power on, 4-2
system startup program with plate loader	
running, 4-7	U
Т	unclogging
	sample probe for sample loading, 12-1
table	unlock
troubleshooting, 10-6	user account, 2-30

unmixing	VAC		
adjust spectral unmixing, 7-15	define, Abbreviations-2		
calculate unmixing values, 7-10	ventilation and cleaning, A-2		
confirm unmixing results, 7-12	vertical gates		
create spectral unmixing experiment, 7-2	for single-parameter plots, 2-40		
fine adjustment, 7-27	VSSC		
load sample and record data, 7-9	define, Abbreviations-2		
prepare unmixing sample, 7-1	define, more viations 2		
save spectral unmixing experiment, 7-14			
unmixing results	W		
confirm, 7-12	W		
	define, Abbreviations-2		
unmixing sample	warning		
prepare, 7-1	define, v		
unmixing settings	waste container 4 L (CytoFLEX)		
adjust, 8-2	emptying, 12-1		
exporting, 7-20	waste cubitainer 10 L (CytoFLEX LX)		
import from the spectral library, 7-17			
import from unmixing matrix, 7-15	emptying, 12-1		
unmixing values	waste disposal, A-4		
to calculate, 7-10	waste harness		
unpacking inspection, A-5	to replace, 12-2		
check instrument configuration, A-5	waste levels (10 L fluid cubitainers)		
installing the instrument and connecting	check, 4-1		
the equipment, A-6	waste levels (4 L fluid containers)		
unpacking the instrument, A-5	check, 4-1		
USB	wavelength		
define, Abbreviations-2	laser, 1-14		
use your manual	wavelength division multiplexer (WDM), 1-4		
how to, xxvii	WDM		
user account	define, Abbreviations-2		
unlock, 2-30	well settings with plate loader		
user management, 2-27, E-63	modifying, 6-7		
electronic record management, B-4	collecting, 6-33		
user management operation	copying		
log, 2-39	wells, 6-7		
user management operation log	create		
electronic record management, B-3	spectral unmixing experiment, 7-5		
users	create spectral unmixing experiment, 7-5		
creating, 2-28	cutting		
deleting, 2-29	wells, 6-7		
modifying, 2-30	data		
UV	collecting, 6-33		
define, Abbreviations-2	sampling, 6-33		
define, neore viacions 2	modifying well settings		
	copying wells, 6-7		
V	cutting wells, 6-7		
V	pasting wells, 6-7		
define, Abbreviations-2	pasting		
VA	wells, 6-7		
define, Abbreviations-2	prepare		

unmixing sample, 7-2	prepare
sample injection mode	unmixing sample, 7-2
selecting, 4-7	spectral unmixing experiment, 7-2
sampling, 6-33	unmixing
selecting	create spectral unmixing experiment, 7-2
proper sample injection mode, 4-7	prepare unmixing sample, 7-2
spectral unmixing experiment, 7-5	routine maintenance
unmixing	unclogging the sample probe, 12-1
create spectral unmixing experiment, 7-5	unclogging the sample probe, 12-1
prepare unmixing sample, 7-2	Workstation, 1-2, 1-3
routine maintenance	component, 1-2
unclogging the sample probe, 12-1	workstation connections inspection
unclogging the sample probe, 12-1	pre-boot, 4-2
create	worktable
spectral unmixing experiment, 7-2	installation environment validation, A-2
create spectral unmixing experiment, 7-2	

Beckman Coulter, Inc. SOFTWARE END-USER LICENSE AGREEMENT

IMPORTANT: READ THIS END-USER LICENSE AGREEMENT:

THIS END-USER LICENSE AGREEMENT ("Agreement") GOVERNS THE USE OF CERTAIN COMPUTER PROGRAMS ASSOCIATED WITH THE CORRESPONDING PRODUCT ("Product") OF BECKMAN COULTER, INC. OR ITS AFFILIATES ("BCI"). YOUR USE OF SUCH COMPUTER PROGRAMS, YOUR PURCHASE OF THE PRODUCT, AND/OR YOUR USE OF THE PRODUCT CONSTITUTES ASSENT AND ACCEPTANCE OF THIS AGREEMENT BY YOU (collectively "You" or "Your") AND SIGNIFIES YOUR AGREEMENT TO BE BOUND BY THE TERMS AND CONDITIONS OF THIS AGREEMENT. IF YOU DO NOT AGREE WITH ALL OF THE TERMS AND CONDITIONS OF THIS AGREEMENT, YOU DO NOT HAVE THE RIGHT TO AND MAY NOT USE THE COMPUTER PROGRAMS PROVIDED.

- 1. GRANT OF LICENSE. Subject to all of the terms and conditions of this Agreement, BCI grants to You a non-exclusive, non-sub licensable license to use the computer programs provided by BCI ("Programs") in connection with and/or installed on the Product, and the associated user documentation ("Documentation") (together with the Programs referred to herein, collectively, as the ("Software")), only in accordance with the Documentation ("License"). You have no right to receive, use or examine any source code or design documentation relating to the Programs. The Software is licensed and not sold. As between the parties, BCI and its licensor retains all right, title and interest in and to the Software and any and all derivative works, except as expressly and unambiguously licensed herein, and BCI reserves all rights in the Software not granted to You. Except as otherwise expressly provided under the Agreement, You shall only use the Software for Your internal research and development purposes.
- **2. USE RESTRICTIONS.** You shall not (and shall not allow others to), directly or indirectly:
 - **a.** use (including make any copies of) the Programs beyond the scope of the license granted under Section 1.;
 - **b.** modify, translate, adapt, or otherwise create derivative works or improvements, whether or not patentable, of the Programs or any part thereof;
 - **c.** combine the Programs or any part thereof with, or incorporate the Programs or any part thereof in, any other programs;
 - **d.** reverse engineer, disassemble, decompile, decode, or otherwise attempt to derive or gain access to the source code of the Programs or any part thereof;
 - **e.** reverse engineer, disassemble, decompile, decode, or otherwise attempt to derive or gain access to the object code of the Programs or any part thereof;
 - **f.** remove, delete, alter, or obscure any trademarks or any copyright, trademark, patent, or other intellectual property or proprietary rights notices provided on or with the Software, including any copy thereof;
 - g. rent, lease, lend, sell, sub license, assign, distribute, publish, transfer, or otherwise make available the Programs, or any features or functionality of the Programs, to any third party for any reason, whether or not over a network or on a hosted basis, including in connection with the Internet or any web hosting, wide area network (WAN), virtual private network (VPN), virtualization, time-sharing, service bureau, software as a service, cloud, or other technology or service;
 - **h.** use the Software in, or in association with, the design, construction, maintenance, or operation of any hazardous environments or systems;
 - i. use the Software in violation of any law, regulation, or rule; or
 - **j.** use the Software for purposes of competitive analysis of the Programs, the development of a competing software product or service, or any other purpose that is to BCI's commercial disadvantage.

D17052AA Warranty-1

- 3. TERMINATION. The License is effective until terminated. You may terminate the Agreement and the License at any time by destroying all copies of the Programs. The License will terminate automatically without notice if You fail to comply with any provision of this Agreement. Within fourteen (14) days following termination, You shall cease all use of, and, destroy all copies of the Software in Your possession or control and so certify to BCI in writing. Except for the License, the terms of this Agreement shall survive termination. Termination is not an exclusive remedy and all other remedies will be available to BCI whether or not the License is terminated.
- **4.** SOFTWARE SUPPORT. Please refer to the BCI support number provided in the Documentation for the Product for information regarding support of the Product.
- 5. DISCLAIMER OF WARRANTY. THE SOFTWARE IS PROVIDED TO YOU "AS IS" AND WITH ALL FAULTS AND DEFECTS WITHOUT WARRANTY OF ANY KIND. TO THE MAXIMUM EXTENT PERMITTED UNDER APPLICABLE LAW, BCI, ON ITS OWN BEHALF AND ON BEHALF OF ITS AFFILIATES AND ITS AND THEIR RESPECTIVE LICENSORS AND SERVICE PROVIDERS, EXPRESSLY DISCLAIMS ALL WARRANTIES, WHETHER EXPRESS, IMPLIED, STATUTORY, OR OTHERWISE, WITH RESPECT TO THE SOFTWARE, INCLUDING ALL IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, TITLE, AND NON-INFRINGEMENT, AND WARRANTIES THAT MAY ARISE OUT OF COURSE OF DEALING, COURSE OF PERFORMANCE, USAGE, OR TRADE PRACTICE. WITHOUT LIMITATION TO THE FOREGOING, BCI PROVIDES NO WARRANTY OR UNDERTAKING, AND MAKES NO REPRESENTATION OF ANY KIND THAT THE LICENSED SOFTWARE WILL MEET YOUR REQUIREMENTS, ACHIEVE ANY INTENDED RESULTS, BE COMPATIBLE, OR WORK WITH ANY OTHER SOFTWARE, APPLICATIONS, SYSTEMS, OR SERVICES, OPERATE WITHOUT INTERRUPTION, MEET ANY PERFORMANCE OR RELIABILITY STANDARDS OR BE ERROR FREE, OR THAT ANY ERRORS OR DEFECTS CAN OR WILL BE CORRECTED.
- LIMITATION OF LIABILITY. TO THE FULLEST EXTENT PERMITTED UNDER APPLICABLE LAW:
 - a. IN NO EVENT WILL BCI OR ITS AFFILIATES, OR ANY OF ITS OR THEIR RESPECTIVE LICENSORS OR SERVICE PROVIDERS, BE LIABLE TO YOU OR ANY THIRD PARTY FOR ANY USE, INTERRUPTION, DELAY, OR INABILITY TO USE THE SOFTWARE; LOST REVENUES OR PROFITS; DELAYS, INTERRUPTION, OR LOSS OF SERVICES, BUSINESS, OR GOODWILL; LOSS OR CORRUPTION OF DATA; LOSS RESULTING FROM SYSTEM OR SYSTEM SERVICE FAILURE, MALFUNCTION, OR SHUTDOWN; FAILURE TO ACCURATELY TRANSFER, READ, OR TRANSMIT INFORMATION; FAILURE TO UPDATE OR PROVIDE CORRECT INFORMATION; SYSTEM INCOMPATIBILITY OR PROVISION OF INCORRECT COMPATIBILITY INFORMATION; OR BREACHES IN SYSTEM SECURITY; OR FOR ANY CONSEQUENTIAL, INCIDENTAL, INDIRECT, EXEMPLARY, SPECIAL, OR PUNITIVE DAMAGES, WHETHER ARISING OUT OF OR IN CONNECTION WITH THIS AGREEMENT, BREACH OF CONTRACT, TORT (INCLUDING NEGLIGENCE), OR OTHERWISE, REGARDLESS OF WHETHER SUCH DAMAGES WERE FORESEEABLE AND WHETHER OR YOU WERE ADVISED OF THE POSSIBILITY OF SUCH DAMAGES;
 - b. IN NO EVENT WILL THE COLLECTIVE AGGREGATE LIABILITY OF BCI AND ITS AFFILIATES, INCLUDING ANY OF THEIR RESPECTIVE LICENSORS AND SERVICE PROVIDERS, UNDER OR IN CONNECTION WITH THIS AGREEMENTOR ITS SUBJECT MATTER, UNDER ANY LEGAL OR EQUITABLE THEORY, INCLUDING BREACH OF CONTRACT, TORT (INCLUDING NEGLIGENCE), STRICT LIABILITY, AND OTHERWISE, EXCEED TWO HUNDRED AND FIFTY DOLLARS (US\$250); AND
 - **c.** THE LIMITATIONS SET FORTH IN THIS SECTION 6. SHALL APPLY EVEN IF YOUR REMEDIES UNDER THIS AGREEMENT FAIL OF THEIR ESSENTIAL PURPOSE.
- **7.** U.S. GOVERNMENT RIGHTS. The Software is commercial computer software, as such term is defined in 48 C.F.R. §2.101. Accordingly, if You are the US Government or any contractor therefor, You shall receive only those rights with respect to the Software as are granted to all other end users under license, in accordance with
 - **a.** 48 C.F.R. §227.7201 through 48 C.F.R. §227.7204, with respect to the Department of Defense and their contractors, or
 - **b.** 48 C.F.R. §12.212, with respect to all other US Government licensees and their contractors.

Warranty-2 D17052AA

- 8. EXPORT REGULATION. The Software may be subject to US export control laws, including the US Export Administration Act and its associated regulations. You shall not, directly or indirectly, export, re-export, or release the Software to, or make the Software accessible from, any jurisdiction or country to which export, re-export, or release is prohibited by law, rule, or regulation. You shall comply with all applicable federal laws, regulations, and rules, and complete all required undertakings (including obtaining any necessary export license or other governmental approval), prior to exporting, re-exporting, releasing, or otherwise making the Software available outside the US.
- 9. INTELLECTUAL PROPERTY RIGHTS. You acknowledge and agree that the Software is provided under license, and not sold, to You. You do not acquire any ownership interest in the Software under this Agreement, or any other rights thereto, other than to use the same in accordance with the license granted and subject to all terms, conditions, and restrictions under this Agreement. BCI shall retain its entire right, title, and interest in and to the Software and all intellectual property rights arising out of or relating to the Software, except as expressly granted to You in this Agreement. You shall use commercially reasonable efforts to safeguard all Software (including all copies thereof) from infringement, misappropriation, theft, misuse, or unauthorized access.
- 10. THIRD PARTY COMPONENTS; ADDITIONAL TERMS. The Software may contain or be delivered with one or more components, which may include third-party components, identified by BCI in the Documentation, readme.txt file, third-party click-accept or elsewhere (e.g., on www.beckman.com or attached to this Agreement as an appendix) (the "Identified Component(s)") as being subject to different license agreement terms, disclaimers of warranties, limited warranties or other terms and conditions (collectively, "Additional Terms") than those set forth herein. You agree to the applicable Additional Terms for any such Identified Component(s).
- 11. MISCELLANEOUS. You may not assign this Agreement or any of Your rights hereunder without BCI's prior written consent and any attempt to do so without such consent shall cause this Agreement and any of Your associated rights hereunder to be null and void. No failure to exercise any right hereunder will operate as a waiver thereof. If any provision of this Agreement shall be adjudged by any court of competent jurisdiction to be unenforceable or invalid, that provision shall be limited or eliminated to the minimum extent necessary so that this Agreement shall otherwise remain in full force and effect and enforceable. This Agreement shall be construed in accordance with the laws of the State of California and the United States without regard to conflicts of laws provisions thereof and without regard to the United Nations Convention on Contracts for the International Sale of Goods. The sole and exclusive jurisdiction and venue for any actions related to the subject matter hereof shall be the state and U.S. federal courts located in the County of Orange, California. You irrevocably submit to the jurisdiction of such courts and consent to venue in such forum with respect to any action or proceeding that relates to this Agreement. The prevailing party in any action to enforce this Agreement shall be entitled to recover its reasonable costs and expenses including reasonable attorneys' fees. No amendment to or modification of this Agreement will be binding unless in writing and signed by a duly authorized officer of BCI. This Agreement is in the English language only, which language shall be controlling and any revision of this Agreement in any other language shall not be binding. Both parties agree that this Agreement is the complete and exclusive statement of the mutual understanding of the parties relating to the subject matter of this Agreement.

D17052AA Warranty-3

Beckman Coulter, Inc. SOFTWARE END-USER LICENSE AGREEMENT

Warranty-4 D17052AA

Related Documents

CytoFLEX mosaic Instructions for Use PN D17052

CytoFLEX Platform Instructions for Use PN B49006

www.beckman.com

