User Manual

Cubic Vacuum Desiccators

Models: VDC-11, 21, 31, 41, 41D
11U, 21U, 31U, 41U
Manual No. D4001L002 Version 4.4

⚠ WARNING

Before using this product, read this entire Operator's Manual carefully. Users should follow all of the Operational Guidelines contained in this Manual and take all necessary safety precautions while using this product. Failure to follow these guidelines could result in potentially irreparable bodily harm and/or property damage.

Thank you for purchasing Jeio Tech's products.

Jeio Tech Co., Ltd. is committed to customer service both during and after the sale. If you have questions concerning the operation of your unit or the information in this manual, contact our Sales Department. If your unit fails to operate properly, or if you have questions concerning spare parts or Service Contracts, contact our Service Department.

Quality Management System

Jeio Tech Co, Ltd. is dedicated to providing world-best product quality and customer satisfaction. To ensure we maintain this commitment we have developed and implemented a total quality program, which conforms to the requirements according to DIN EN ISO 9001:2015 for the design, development, production, sales and servicing of biotechnology, environmental chemical engineering related products, and reliable measuring equipment for electric and electronics (ovens, incubators, constant temperature and Humidity chambers, constant temperature baths, refrigerating bath circulators, heat exchangers and shakers).

Visit our Web site at www.jeiotech.com to view a copy of our certificate.

Disclaimer

Jeio Tech Co., Ltd. is committed to a continuing program of product development and improvement, and reserves the right to change information, such as specifications, appearance, and dimensions, described in this document without notice.

Copyright

No part of this manual may be reproduced or transmitted in any form or by any means, including photocopying, recording, or using information storage and retrieval systems, for any purpose other than the purchaser's own use, without the express written permission of Jeio Tech Co., Ltd.

©2026. All Rights Reserved. Jeio Tech Co., Ltd.

Any other product names and services identified in this manual are trademarks or registered trademarks of their respective owners. No such use, or the use of any trade name, is intended to convey endorsement or other affiliation with Jeio Tech Co., Ltd.

Table of Contents

1	Safety Advice	1
1.1	User Manual	1
1.2	DANGER/CAUTION/NOTICE Alerts.....	1
1.3	Danger Alerts	2
1.4	Caution Alerts	2
2	Product Description	3
2.1	Introduction	3
2.2	Characteristics	4
2.3	Layout.....	6
3	Installation.....	8
3.1	Checking Shipment Damage and Unpacking.....	8
3.2	Checking the Contents of the Package	8
3.3	Installation Environment	1 0
3.4	Installation location	1 0
3.5	The labels attached to the instruments (VDC-41D)	1 1
3.6	VDC-41D Installation method (when used as Vacuum Defoaming Stirrer)	1 2
3.6.1	Vacuum Defoaming with Jeiotech Stand, MSA-0420	1 2
3.6.2	Vacuum Defoaming with another company product	1 6
3.6.3	How to Attach the Lab Sticker	1 6
4	Operation.....	1 7
4.1	Precautions	1 7
4.2	Operation of the 3-Way Valve	1 7
4.3	Vacuum formation	1 8
4.4	Multi-stacking of the Units	2 0
4.5	Storage	2 2
5	Maintenance	2 3
6	Troubleshooting.....	2 4
7	Accessories	2 5
8	Appendix.....	2 6
8.1	Technical Specifications	2 6
8.2	Disposal of the Unit	2 7
8.3	Warranty.....	2 7
8.3.1	Warranty period	2 7
8.3.2	Technical Service Contact Points	2 8
8.3.3	Certificate of Product Warranty	2 8

1 Safety Advice

1.1 User Manual

This manual contains important safety and operation information. You must carefully read, understand, and follow all the instructions in this manual prior to operating this instrument. Keep this manual in a safe place nearby for reference and make it easily available to all users.

1.2 DANGER/CAUTION/NOTICE Alerts

This manual highlights DANGER/CAUTION/NOTICE alerts to prevent injury or property damage and also to achieve optimum performance of your instrument.

These alerts are classified into three types in this manual depending on the importance and the risk levels as described below:

Symbols	Meaning
	Ignoring this alert could cause serious or even fatal injury.
	Ignoring this alert could cause serious injury or property damage.
	Ignoring this alert could cause operational problems.

Ensure to read the danger and caution alerts in 1.3 and 1.4 carefully before use.

1.3 Danger Alerts

 DANGER
Ignoring the following warnings could cause serious or even fatal injury.

- ◆ Never install or use this unit in explosive atmospheres.
- ◆ Never install this unit near to hazardous or flammable substances.
- ◆ Never disassemble, repair, or modify this unit on your own. Doing so will void your warranty and may result in injuries or product damages.

1.4 Caution Alerts

 CAUTION
Ignoring the following cautions could cause serious injury or property damage.

- ◆ Carefully read all the warning labels before using this unit. And do not remove or damage the warning labels.
- ◆ Do not move this unit during operation. When moving it, use the molded-in handles.
- ◆ Do not use this unit in wet or moist atmosphere or where water leakage is expected.
- ◆ Do not use this unit in contaminated atmosphere or where metallic dust exists.
- ◆ Do not expose this unit to any heat sources including direct sunlight.
- ◆ Do not allow moisture, organic solvent, dust, or corrosive gases to get into this unit.
- ◆ If you observe a strange smoke, odor or noise from this unit, stop using it. After the smoke or odor disappears, contact your dealer or Jeio Tech if any repair is required.
- ◆ Do not use chlorine bleach, ammonia-based cleaners, abrasives, ammonia, or metal scouring pads. Wipe with a soft damp cloth or a sponge soaked in water or diluted neutral detergent or ethanol diluted at 50% concentration.

2 Product Description

2.1 Introduction

Congratulations of your purchase of our high-quality cubic vacuum desiccator which is specially designed to provide maximum benefit for your investment with respect to performance, safety, ease of use, and durability.

This unit is a sturdily constructed 360° transparent vacuum desiccator ideal for storing, drying, cooling, and degassing various materials. Vacuum environment created by this unit providing a very low humidity, contamination-free, corrosion-free, and dust-free environment is ideal for:

- storing moisture-sensitive materials such as anhydrides and other hygroscopic compounds as well as materials which are highly reactive with atmosphere,
- removing residual solvents or air bubbles,
- cultivating anaerobic microbes,
- minimizing the reabsorption of humidity and the exposure to air during the drying or the cooling stages of biological/medical/food/materials experiments,
- conducting experiments at low pressure,
- drying at low pressure to minimize the heat effects.

Currently eight models of cubic vacuum desiccators are available from Jeio Tech (VDC-11/11U, 21/21U, 31/31U, and 41/41U) with the internal volume ranging from 11 liter to 45 liter. Note that UV blocking models (VDC-11U, 21U, 31U, and 41U) are offered to minimize damages or discoloration of light-sensitive samples by using UV-blocking resin which completely blocks all ranges of UV and the blue spectrum of the visible light.

2.2 Characteristics

◆ Key Vacuuming Features

(1) Outstanding Vacuum Sustainment

This unit maintains a 1 Torr (133 Pa) vacuum for over 72 hours at room temperature and stays gastight without a vacuum allowing reliable experiments over extended periods.

(2) Integral Molding of the Polycarbonate Body

Unlike conventional vacuum desiccators, the polycarbonate cubic body (except the door part) of this unit is integrally molded as a single piece for maximum sealing performance. This unique design without requiring any adhesive for body construction allows virtually no leakage.

Shock-resistant, high-tensile strength, and transparent polycarbonate body provides extreme durability as well as 360° visual observation from outside.

(3) Viton® 3-Way Valve

This unit comes standard with a Viton® 3-Way Valve offering consistent and uniform vacuum draw, vacuum release, or gas exchange without the inconvenience of connecting and disconnecting pump hoses to the unit.

If needed, up to 4 valves can be installed using the built-in ports at the base of the unit and the positioning of the valve(s) can be also easily changed.

When releasing the vacuum, If using the Pore plug on the 3-Way Valve, it allows air to slowly enter the small diameter inlet, which protects dust or powder samples from flying out.

(4) Airtight Door Sealing

The door is located in the front for convenient access. Equipped with a special door latch which is easily closed at normal pressure and a greaseless vacuum seal using high-quality silicon gasket, this unit enables airtight sealing between the door and the body and quick vacuum formation.

(5) Vacuum Defoaming Stirrer (VDC-41D)

VDC-41D model effectively removes air bubbles from viscous samples or mixtures through a combination of a mixing process and a vacuum process.

◆ Key Design Features

(1) Multi-stackable Design

The VDC models come in a variety of capacities to suit the user's needs with the internal volume ranging from 11 liter to 45 liter.

With the same foot print, different models can be stacked together in any convenient fashion allowing for lab space saving.

(2) Efficient Use of the Internal Space

Compared to round desiccators, cubic desiccators provide larger usable volumes with the same footprint. In addition, the internal space can be more efficiently utilized by using the standard perforated shelves designed for efficient airflow.

Two or three perforated shelves are provided and, if needed, up to four or six shelves can be used for efficient storage of the samples. (Refer to Product Checklist or Specifications for details.)


(3) Easy Vacuum Level Checking

With an built-in vacuum gauge installed on the front door, the vacuum level checking can be done conveniently.

(4) Desiccant Tray

This unit comes standard with a desiccant tray.


2.3 Layout


1. Vacuum Gauge
2. Door
3. Door Latch
4. 3-Way Valve
5. Non Slip Feet
6. Groove Handle
7. Valve Port Plug

8. Vacuum Seal
9. Shelf Rail
10. Perforated Shelves
11. Desiccant Tray
12. Stacking Groove
13. Air Inlet/Outlet Nozzle

[VDC-41D when impeller and vacuum mixing seal is connected]

- 1. Stand
- 2. Shaft (8 mm)
- 3. Vacuum mixing seal
- 4. Vacuum Gauge
- 5. Door
- 6. Body

- 7. Door Latch
- 8. 3-Way Valve
- 9. Gasket
- 10. Perforated Shelves
- 11. Desiccant Tray



3 Installation


3.1 Checking Shipment Damage and Unpacking

Upon receiving the instrument, check to ensure that no damage has occurred during shipment. It is important that any damage that occurred during shipment is detected before unpacking. If you find such damage, the carrier must be notified immediately.

3.2 Checking the Contents of the Package

After unpacking, check to ensure that all the parts and accessories described below are included in the package. If not, contact your dealer immediately.

Item		Figure	Quantity
Main Body			1
Vacuum Gauge			1
3-Way Valve			1
Pore plug			1
Port Plug	VCD-11/11U		1
	VDC-21/21U, 31/31U, 41/41D/41U		3
Vacuum Seal			1
Desiccant Tray			1

Item		Figure	Quantity
Perforated Shelves	VDC-11/11U, 21/21U		2
	VDC-31/31U, 41/41U/41D		3
User Manual			1
Stand	Only VDC-41D		1
Anchor impeller			1
Vacuum mixing seal			1
Lab sticker			1

3.3 Installation Environment

NOTICE
Avoid exposure to any heat sources including direct sunlight.

This unit is designed for indoor use in laboratory environments. The conditions for proper installation are as follows:

Permissible Ambient temperature	5 to 40°C (41 to 104° F)
Permissible Relative humidity	below 80%
Permissible Altitude	0 to 2,000 m (6,562 ft)

3.4 Installation location

NOTICE
Install the unit on a level surface.


3.5 The labels attached to the instruments (VDC-41D)

1) ID label

An ID label containing the basic information is attached to the bottom right panel of the instrument.

2) Safety label


Safety labels are attached to the instrument to provide important information about potential risks or how to avoid them.

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Due to the nature of the manufacturing process, small bubbles, flow marks or black spots may appear on the surface of the unit. But such micro defects do not affect on performance or quality of the unit. Please use with confidence.

Code No.
Description
Model
Voltage
Serial No.

JEIOtech Made in Korea
www.jeio.tech.com


3.6 VDC-41D Installation method

(when used as Vacuum Defoaming Stirrer)


3.6.1 Vacuum Defoaming with Jeiotech Stand, MSA-0420

- (1) Check the level of the floor where the product will be installed using a level.
- (2) Place the VDC-41D unit on the floor.
- (3) Attach the M14 nut to the provided support (see Fig. 1).
- (4) Assemble the Support to the bracket of the unit (see Fig. 1).
- (5) Attach the Clamp holder to the bracket (see Fig. 1, 2).

* Clamp holder is not included as a standard component.

(Fig 1. Stand Assembly)


(Fig 2. Completed Stand Assembly)


- (6) Install the O-ring onto the screw of the vacuum mixing seal.
- (7) Open the door of the VDC-41D, remove the shelf, and insert the impeller into the chamber (see Fig. 3).
- (8) Loosen the fixing nut of the vacuum mixing seal and pass the impeller through from

bottom to top (see Fig. 4).


(9) Fasten the vacuum mixing seal to the body using a spanner (see Fig. 5).

(Fig 3. Impeller Assembly 1)


(Fig 4. Impeller Assembly 2)


(Fig 5. Impeller Assembly 3)

(10) Mount the MSA-0420 onto the clamp holder, insert the impeller into the chuck of the MSA, and tighten using the chuck handle (see Fig. 6, 7).


(Fig 6. Securing the Impeller 1)

(Fig 7. Securing the Impeller 2)

- (11) When viewed from the side, ensure that the chuck, impeller, and vacuum mixing seal are aligned vertically, then tighten the clamp holder (see Fig. 8).

(Fig 8. Securing Chuck, Impeller, vacuum mixing seal)

- (12) Connect the power cable of the MSA-0420 and switch it on.
- (13) Run the impeller at 300 RPM and check for visible wobbling. If excessive wobbling is observed, recheck the vertical alignment of the chuck, impeller, and vacuum mixing seal.
- (14) Close the door, connect a vacuum pump to the valve, and check if vacuum is properly maintained. (Normal condition: After reaching a vacuum level of 1 Torr, the vacuum pump is stopped and mixing is performed at 300 RPM. Under this condition, the


vacuum level should be maintained below 10 mbar for 10 minutes; otherwise, it is considered abnormal.)

3.6.2 Vacuum Defoaming with another company product

- (1) Follow the installation procedure described in 3.5.1, but replace the MSA-0420 with a third-party overhead stirrer.
- (2) Ensure vertical alignment of the chuck, impeller, and vacuum mixing seal.
Only impellers with a shaft diameter of 8 mm can be installed.

3.6.3 How to Attach the Lab Sticker

- (1) Clean the surface of the shelf where the sticker will be applied and make sure it is completely dry.
- (2) Remove the protective film from the flat side of the Lab Sticker.
- (3) Place the flat side facing downward and firmly attach the sticker onto the shelf.
- (4) Make sure no air bubbles remain between the Lab Sticker and the shelf. Any trapped air may expand under vacuum conditions, so press out the air thoroughly during application.

(Fig. 9. Example of Lab Sticker attachment)

NOTICE

- The adhesive strength of the Lab Sticker may vary depending on the laboratory environment and maintenance conditions.
- Adhesion may be affected by factors such as temperature, humidity, and the surface condition of the beaker.
- If adhesion weakens, clean the sticker gently with a soft sponge and neutral detergent to remove dust and contaminants. Rinse thoroughly with water and allow it to dry completely before reuse.
- Do not use alcohol, strong acids, or strong alkaline cleaners, as they may reduce adhesion.
- Under normal usage conditions, the Lab Sticker is recommended to be replaced every 2–3 years, or immediately if adhesion decreases or wear is observed.


4 Operation




4.1 Precautions

⚠ CAUTION

- Do not form a vacuum when airtight containers are inside the unit.
- Note that this unit cannot be autoclaved.
- When attaching other equipments or additional 3-Way Valves after removing the port plugs, use proper O-rings and nuts for tight sealing.
- The vacuum desiccator is made of polycarbonate, so the thermal deformation temperature is 120°C. When putting high-temperature samples inside the unit, cool them thoroughly in advance or make sure that they are not in direct contact with the interior of the unit to avoid interior damage.

4.2 Operation of the 3-Way Valve

		Front Nozzle	Rear Nozzle
	Rear	Open	Closed
	Front		

	Closed	Open
 	Closed	Closed

4.3 Vacuum formation


The procedures for vacuum formation are as follows:

- (1) When using the unit for drying purposes, fill the tray with desiccants and put it inside prior to operation.
- (2) After putting the samples inside the unit, firmly close the door by locking the two latches properly. You will hear a click sound when the latch is properly engaged.
- (3) Connect the vacuum pump hose to one of the two nozzles of the 3-Way Valve.
- (4) After activating the vacuum pump, turn the valve handle towards the nozzle to open it and read the vacuum gauge to check whether the vacuum is forming properly.
- (5) When the desired vacuum level is achieved, turn the valve handle until it is perpendicular to the nozzle to close it and check the vacuum gauge whether there is any leakage.
- (6) When releasing the vacuum, slowly turn the valve handle towards the opposite nozzle to

return to the normal pressure.

→ At this time, using the Pore plug on the 3-Way Valve allows air to slowly enter the small diameter (Ø 1 mm) inlet, which protects dust or powder samples from flying out.

Here's how to use the Pore plug on the 3-Way Valve.

The cap of Pore plug

- ① Insert the loop, the end part of the Pore plug into the nozzle(Air inlet) of 3-Way Valve.
- ② Cap the Pore plug into the nozzle(Air inlet) of 3-Way Valve.
- ③ Turn the valve towards the capped air inlet to release the vacuum.

However, if you use the Pore plug, it takes 3 to 4 times longer than using 3-Way Valve only to recover to atmospheric pressure. The following table compares the recovery time to atmospheric pressure when using the Pore plug on 3-Way Valve for each model. Please, use it with reference.

Model		VDC-11	VDC-21	VDC-31	VDC-41
Internal volume		11 L	23 L	35 L	45 L
Time to recover to atmosphere pressure	3-Way Valve Air inlet (Ø, 4.5mm)	About 30 sec	About 1 min	About 1~2 min	About 1~2 min
	Pore plug (Ø, 1mm)	About 2~3 min	About 2~3 min	About 7~8 min	About 9~10 min


CAUTION

- Do not impact or drop the unit when a vacuum is formed inside.
- When forming or releasing a vacuum, turn the valve handle slowly to avoid rapid pressure change. If not, the stored samples can be damaged or dispersed.
- When forming a vacuum, install a cold trap in front of the vacuum pump. Moisture and chemicals from the samples may result in damages to the pump.
- Do not let the body of the unit into direct contact with acetone, benzene, toluene, chloroform, cresol, sodium hydroxide, highly concentrated nitric or sulfuric acid, acetic acids, or strong chlorine-based solvents.
- Ensure that the door seal is not damaged by strong acid or base. Damaged seal can cause malfunction of the unit.
- When a vacuum is formed inside the unit, do not apply excessive force to the 3-Way Valves, the vacuum gauge, or any of their vicinities. Any damages caused by external forces can cause malfunction of the unit.
- When the unit is not forming a vacuum, keep the valve handle perpendicular to the nozzles.
- When relocating the unit, use the built-in groove handles instead of holding on to other protruded parts such as the vacuum gauge or 3-Way Valves.

4.4 Multi-stacking of the Units

All VDC models are designed to be easily multi-stacked together by locking the feet of the unit into the stacking grooves located on the top surface of another unit.

VDC-21(U), 31(U), 41(U) models have the same width, so you can use a combination of several products as needed : as shown below

⚠ CAUTION

- Desiccator multi-stacking height can be inverted when it is more than 1m, causing damage to users. Therefore, it is recommended to be less than 1m.
- Desiccator load weight should not exceed maximum weight by referring to Tech data.
- Keep an appropriate distance from other products and avoid heat sources or direct sunlight and install them on a flat place.

4.5 Storage

If this unit is not used for an extended period of time, clean the instrument with soft cloth and store it in dry place after packing it properly.

4.6 Vacuum Defoaming Stirrer

- (1) Place the sample to be mixed and defoamed inside the desiccator, and attach the impeller at the proper height relative to the sample level.
- (2) Ensure that the chuck, shaft, and vacuum mixing seal are correctly aligned on a vertical axis.
- (3) Close the chamber door and set the desired RPM to mix the sample for a specified duration. (Adjust the mixing time according to the sample's viscosity and application, and verify that excessive shaking does not occur.)
- (4) Stop mixing once the sample is uniformly blended.
- (5) Verify the vacuum pump connection, then apply vacuum to remove air bubbles from the sample.
(Monitor for sample expansion and prevent overflow from the beaker.)
- (6) When bubble formation ceases, gradually introduce external air to return the chamber to atmospheric pressure.
- (7) If required, repeat the sequence of mixing → vacuum → venting.

CAUTION

- Exceeding this speed may cause Vacuum mixing seal deformation and malfunction.
- Ensure that the Chuck, Shaft, and Vacuum mixing seal are aligned along the vertical axis. Misalignment may lead to Vacuum mixing seal deformation and malfunction.
- Only 8mm diameter Shafts can be used.

5 Maintenance

CAUTION

- If the unit is contaminated, wear chemicals-resistant gloves before cleaning.
- Regularly check the 3-Way Valves, the vacuum gauge, or any of their vicinities for any damage.
- Regularly check the door seal as well as the surface where the door and the body are in contact and keep them always clean and undamaged.
- Make sure to clean the door seal regularly using diluted neutral detergent. For thorough cleaning, detach the seal from the door.
- Do not use chlorine bleach, ammonia-based cleaners, abrasives, organic solvents, or metal scouring pads when cleaning. Use a soft cloth all the time.
- Desiccator is made of polycarbonate and has good chemical resistance. However, it is recommended to dilute the ethanol concentration to 50% during washing because high purity ethanol may have some effects when used for a long time.
- If the 3-Way Valve does not move smoothly, apply Vacuum grease.
- The 3-Way Valve is a consumable part. If vacuum cannot be maintained due to damage to the 3-Way Valve, replace it.

6 Troubleshooting

Problem	Cause	Corrective Action
Vacuum Malfunction	Contaminated surface between the door and the body	Thoroughly clean the surface with a soft cloth using diluted neutral detergent.
	Damaged door seal	Replace the seal.
	Unsuitable room temperature	Check the temperature whether it is within the operating range.
	Misplaced or damaged 3-Way Valve	Check whether the valve is properly set. If damaged, replace it with a new one.
	Damaged vacuum gauge	Replace the gauge.
Vacuum or mixing malfunction (VDC-41D)	Vacuum Mixing Seal issues	Check whether the vacuum mixing seal is properly installed. Verify that the vacuum mixing seal is securely attached to the impeller.

7 Accessories


Component	Cat. No.	Description	
Vacuum Gauge	EDA4301	0 to -0.1 (Mpa)	
Vacuum Seal	VDC0032	VDC-11(U)	Silicone, 220 x 250 x 6.5
	VDC0033	VDC-21(U)	Silicone, 310 x 260 x 6.5
	VDC0015	VDC-31(U)	Silicone, 310 x 360 x 6.5
	VDC0034	VDC-41(U)(D)	Silicone, 310 x 460 x 6.5
Desiccant Tray	VDC0037	VDC-11(U)	PP, 158 x 24 x 194
	VDC0005	VDC-21(U)/31(U)/41(D)(U)	PP, 268 x 218 x 26
Perforated Shelves	VDC0035	VDC-11(U)	PC, 202 x 6.5 x 204
	VDC0003	VDC-21(U)/31(U)/41(D)(U)	PC, 289 x 11 x 304
3-Way Valve	AAAD4501	PP Ø9.5 x 25	
Silicagel (1kg)	BEA00078603	For all VDC models	
Beaker Dock Shelf Set (500ml)*	AAAD1507	VDC-41D	
MSA-0420*	AAH373325	VDC-41D	
Anchor impeller	AAA37572	VDC-41D	
Vacuum Mixing Seal	AAAD1506	VDC-41D	
Stand	AAA37508	VDC-41D	
Lab Sticker	AAA30551	VDC-41D	

※ When using the Beaker Dock Shelf Set, liquids with a viscosity of up to 30,000 cP can be mixed.

※ Third-party overhead stirrers may also be used. For details on the accessories included with the overhead stirrer, please refer to the user manual of the respective product.

8 Appendix

8.1 Technical Specifications

Model		VDC-11(U)	VDC-21(U)	VDC-31(U)	VDC-41(U)	VDC-41D
Figure						
Internal volume (L)		11	23	35	45	
Vacuum	Gauge range (Analog, MPa)	0 ~ -0.1				
	Max. Permissible	1.33 x 10 ⁻⁴ (1Torr)				
	Nozzle diameter (mm)	Ø 9.5				
Mixing	Max RPM	-				300 ¹⁾
	Max Viscosity (cPs)	-				6000 ¹⁾ (with Lab Sticker)
Material	Body, Shelves	Polycarbonate, transparent				
	3-Way Valve, Latch, Desiccant tray	Polypropylene				
	Vacuum Seal, Foot	Silicone Rubber				
Dimensions (mm)	Overall (W x D x H)	322 x 285 x 271	420 x 397 x 281	420 x 397 x381	420 x 397 x 491	420 x 397 x 501
	Door opening (W x H)	208x238	295x245	295x345	295x445	295x445
	Shelf (W x D x H)	202x204x6.5			289x304x11	
	Desiccant tray (W x D x H)	158x194x24			218x268x26	
Net weight (kg)		4.2	8.2	10.8	12.2	14.8

1) The specified performance may vary depending on the surface cleanliness of the Lab Sticker and the surrounding operating conditions.

Model		VDC-11(U)	VDC-21(U)	VDC-31(U)	VDC-41(U)	VDC-41D
Shelves	Number of standard/Max.	2 / 4	2 / 4	3 / 5	3 / 6	3 / 6
	Max. load/shelf(Kg)	3	5	5	5	5
Cat. No.		AAAD4001	AAAD4021	AAAD4011	AAAD4031	AAAD4431

Model	VDC-11U	VDC-21U	VDC-31U	VDC-41U
Color	Amber			
UV Protection (range)	UV A,B,C 100% (200~450nm)			
Cat. No.	AAAD4101	AAAD4121	AAAD4111	AAAD4131

8.2 Disposal of the Unit

Disposing of the unit must be done in an environmentally responsible way if it has been potentially exposed to bio-agents or radioactive samples. Failure to follow stringent requirements for appropriate disposal may lead to actions against you and your organization.

- (1) First, check with your laboratory or organization to ensure that you are following all the policies and procedures for disposal of laboratory equipments and wastes.
- (2) If not possible, contact your local governing body for regulations regarding disposal of laboratory equipments. Jeio Tech highly recommends you to find a local service provider that can properly dispose of your unit.

8.3 Warranty

8.3.1 Warranty period

- (1) If the product fails during normal and proper use, the warranty period for manufacturing liability shall be two years from the date of delivery.
- (2) If you request a repair, please check the items below and let us know so that you can get a quicker and more accurate repair.

- | |
|---|
| <ul style="list-style-type: none"> • Date of purchase: • Customer name / address / phone number / E-mail: • Fault status • Model • Serial number |
|---|

8.3.2 Technical Service Contact Points

Technical Services

Address : 153 (Youngsan-dong), Techno 2-ro Yuseong-gu, Daejeon-city, 34025, Republic of Korea

Tel : 1522-4298

Website : www.jeiotech.com

E-mail : help@jeiotech.com

8.3.3 Certificate of Product Warranty

This is a guide to the product quality warranty and service terms provided by JEIOTECH Co., Ltd

All products are guaranteed to operate normally when installed under proper conditions, as specified in the user manual for each product, and used according to the intended purposes and permissible operating conditions.

The warranty period for product quality starts from the delivery date and is indicated in the catalog provided for each product or in the purchase order agreed upon with the customer. If a product malfunctions within the warranty period, JEIOTECH Co., Ltd. provides free repair services. However, among the transportation costs incurred during the repair process, any special costs resulting from the customer's site conditions (e.g., door or wall modifications, use of special handling equipment) are not covered under the warranty. Furthermore, JEIOTECH's liability under this warranty is limited to replacing the product or refunding the original purchase price at JEIOTECH's sole discretion if reasonable attempts to repair or replace the product fail.

Parts replaced during free repairs within the warranty period become the property of JEIOTECH Co., Ltd. Additionally, the warranty period is not extended due to free repairs. JEIOTECH is not responsible for any losses incurred by the customer due to delays in repairs caused by issues such as part shortages or shipping delays.

The following cases will be handled as paid repairs, even within the warranty period:

1. Non-defective cases

- * Requests to resolve issues arising from improper operating environments or abnormal installations.
- * Reinstallation requests or connections with external equipment other than the initial installation.
- * Requests for usage instructions or simple adjustments that do not require disassembly.
- * Inspections of non-defective products.
- * Cleaning or debris removal service requests.

2. Failures caused by user errors:

- * Malfunctions due to usage not in accordance with the user manual.
- * Issues resulting from improper operating environments or abnormal installations.
- * Malfunctions caused by incompatible chemical usage.
- * Malfunctions caused by irregular or improper electrical or utility supplies.
- * Damage due to product movement, storage, external impact, or dropping by the customer.
- * Malfunctions resulting from unauthorized modifications or repairs without prior written approval from JEIOTECH Co., Ltd.
- * Malfunctions caused by using consumables or optional components not specified by JEIOTECH Co., Ltd.
- * Repairs performed by personnel other than JEIOTECH or its authorized service agents.

3. Other cases

- * When the product model or serial number has been damaged, removed, manipulated, or altered.
 - * Malfunctions caused by fire, flooding, or other natural disasters.
 - * Failure due to consumable parts reaching the end of their lifespan (e.g., lamps, filters, fuses, gaskets, packings, hoses, heat transfer fluids).
 - * Natural wear and tear (e.g., discoloration, screen printing degradation).
 - * Damage or malfunction of product accessories (e.g., power cables, data cables, glass components) or options.
 - * Malfunctions caused by errors in computers, computer components, or connected network systems.
-
- Before JEIOTECH performs free or paid repairs or accepts the return of a product, the customer must ensure the product has been decontaminated of chemical and biological hazards.
 - Replacement parts used during repairs may differ from the original but will function correctly and be equivalent in performance.
 - JEIOTECH does not bear any responsibility for damage to stored samples, data loss, or any other direct or indirect damages caused by product malfunctions, regardless of the warranty period.
 - JEIOTECH does not guarantee issues related to computers, computer components, or network systems connected to the product, regardless of the warranty period.
 - This warranty is the sole and exclusive warranty provided. Any implied warranties imposed by applicable laws are limited to the warranty period stated herein. This warranty cannot be modified by any verbal advice, statements, or documents on behalf of JEIOTECH unless explicitly labeled as a "Warranty Amendment" and signed by an authorized representative of JEIOTECH.
 - This warranty applies only to products purchased for direct use by the customer. It does not cover products purchased for resale purposes without prior written approval from JEIOTECH.
 - This warranty does not apply to consumables and is non-transferable.

- ✧ The contents of this manual can be changed or upgraded without prior notice.
- ✧ The copyright of this manual is reserved by Jeio Tech.

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

1. Stand
2. Shaft (8 mm)
3. Vacuum mixing seal
4. Vacuum Gauge
5. Door
6. Body

7. Door Latch
8. 3-Way Valve
9. Gasket
10. Perforated Shelves
11. Desiccant Tray

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





Item		Figure	Quantity
Main Body			1
Vacuum Gauge			1
3-Way Valve			1
Pore plug			1
Port Plug	VDC-11/11U		1
	VDC-21/21U, VDC-31/31U, VDC-41/41U/41D		3
Vacuum Seal			1

Item		Figure	Quantity
Desiccant Tray			1
Perforated Shelves	VDC-11/11U, VDC-21/21U		2
	VDC-31/31U, VDC-41/41D/41U		3
사용자 설명서			1
Stand			1
Anchor impeller			1
Vacuum mixing seal			1
Lab sticker			1

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