



Enrichment

Mouse Loft and Mouse Balcony

- Repeatable results from the experiments
- Easy logistics
- Good visibility



Mouse Loft:
environmental enrichment for GM500 IVC
equipped with autowatering system



Mouse Balcony:
environmental enrichment both
for GM500 IVC with water bottle
and with autowatering system



Mouse Loft and Mouse Balcony

The Mouse Loft and the Mouse Balcony are suspended environmental enrichments that create a more complex 3D environment in the cage in order to allow the mice to display normal behaviour patterns like jumping and climbing.

1/Repeatable results from the experiments

Mouse Loft and Mouse Balcony are durable environmental enrichments: they can be washed and autoclaved and this contributes to grant repeatable results from the experiments.

2/Easy logistics

The Mouse Loft and the Mouse Balcony are stackable, which makes handling, washing and storage operations much easier.

3/Good visibility

Both environmental enrichments are made of transparent polysulfone plastic for easy daily inspection and since they are placed in the front of the cage, they grant a very good visibility of the animals.

Repeatable results from the experiments

1



Easy logistics and good visibility

2&3



Mouse Loft



Mouse Balcony



HEADQUARTERS
ITALY | Tel. +39 0332 809711 • www.tecniplast.it • E-mail: tecnicom@tecniplast.it

SISTER COMPANIES
AUSTRALIA/NEW ZEALAND | Tel. + 61 2 8845 6500 • www.tecniplast.it • E-mail: info@tecniplast.com.au
FRANCE | Tel. 04 72 52 94 41 • www.tecniplast.fr • E-mail: info@tecniplast.fr
GERMANY | Tel. 08805 921320 • www.tecniplast.de • E-mail: info@tecniplast.de
UNITED KINGDOM | Tel. 0845 0504556 • www.tecniplast.it • E-mail: info@tecniplastuk.com
JAPAN | Tel. + 03 5770 5375 • www.tecniplast.it • E-mail: info@tecniplastjapan.co.jp
USA/CANADA | Toll Free: 877.669.2243 • www.tecniplastusa.com • E-mail: info@tecniplastusa.com

To find your local distributor, please visit www.tecniplast.it



ISO-9001 - Cert. n° 0875



ISO-14001 - Cert. n° 0256A/0

