

# NU-301 Horizontal Airflow Workstation

The AireGard NU-301 Horizontal Laminar Airflow Console Workstation provides a HEPA filtered clean work area that is ideal for operations requiring a particle-free, bacteria-free, clean air environment.

The NU-301 features true laminar airflow, a large 99.99% HEPA filter, disposable fiberglass pre-filters, high-capacity motor/blower, adjustable motor speed controller, Minihelic® pressure gauge, fluorescent lighting, and the HEPEX™ Zero Leak Airflow system.

## Applications

### Life Science Research

The clean work area provides an excellent work space for small laboratory appliances, microscopes, pipetting, or similar applications.

### Pharmacy

For use in compounding sterile drug preparations or IV admixture preparation.

### Animal or Plant

For procedures requiring work-in-progress protection from outside potentially harmful contaminants.

### Electronic Testing, Manufacturing, or Inspections



## Standard Sizes

	Cabinet Width, (ft)	Filter Height, (in)	Weight (lbs.)
NU-301-324	3	24	290
NU-301-330	3	30	300
NU-301-336	3	36	310
NU-301-424	4	24	330
NU-301-430	4	30	340
NU-301-436	4	36	350
NU-301-524	5	24	375
NU-301-530	5	30	385
NU-301-536	5	36	395
NU-301-624	6	24	420
NU-301-630	6	30	430
NU-301-636	6	36	440

