


Abeyance A440

Vial Inventory:

2 ml Slim Racks: 88
Boxes per Rack: 5
Total Boxes: 440
5.28 x 2.25 in (134 x 58 mm)

Square Racks: 32
Boxes per Rack: 13
Total Boxes: 416
5.50 x 5.50 in (140 x 140 mm)




Square Mini Racks: 92
Boxes per Rack: 13
Total Boxes: 1,196
3.18 x 3.34 in (81 x 85 mm)

SBS Racks: 48
5.50 x 3.63 in (140 x 92 mm)

Blood Bag Inventory:

25 ml Frames: 463
Bags per Frame: 7
Total Bags: 3,241
3.94 x 0.59 in (100 x 15 mm)

50 ml Frames: 276
Bags per Frame: 7
Total Bags: 1,932
6.30 x 0.59 in (160 x 15 mm)



250 ml Frames: 212
Bags per Frame: 4
Total Bags: 848
7.88 x 0.59 in (200 x 15 mm)

500 ml Frames: 172
Bags per Frame: 4
Total Bags: 688
9.30 x 0.59 in (236 x 15 mm)