



Cryo Store Pico Automated Sample Storage System



Keep samples (and staff) safe.
 Access needed samples without exposing innocent samples to uncontrolled transient warming or staff to cryogenic dangers.

Introducing the new Cryo Store Pico™ – an automated cryogenic storage system for high value samples that fits all spaces from discovery through to clinical use.

You can now get all the benefits of an automated cryogenic storage system that keeps samples and staff safe while making them accessible in a controlled way.

Designed to fit spaces with a standard 8-foot ceiling height and 32-inch doorways, the new Cryo Store Pico gives you the sample security and inventory management you need within a smaller footprint.

- 8 feet**
 Cryo that fits in your lab/office
- Stores up to **8,800**
 2mL vials
- Temperature **Stability**
 during retrieval
- 10 Day**
Emergency
 Hold Time

| | |
|--|---|
| Space Efficiency | Fits 8,800 2mL vials within small footprint |
| | User-friendly design improves ergonomics over manual freezers |
| Process control and risk mitigation | Precise cold chain handling with cryo-critical features and documentation |
| | Repeatable processes, secure login access, and onboard monitoring with alarms |
| Data collection and audit support | Real-time inventory control and reporting (supports 21 CFR, Part 11 compliance) |
| | Traceability and monitoring for sample safety and future regulatory compliance |