



# Barkey **vialguard pro**

Thawing your vials. Vertically.

► Upright dry thawing of vials

► For thawing of starting material, for cellbanks and for QC routine

► Full flexibility: vial types, volumes, quantity per run

► Automated and robust process

► GMP and GLP compliant

► Wide range of services such as IQ/OQ/PQ/PM



## Key features

- Barkey vialholder ensures proper thawing process (customizable for multiple types of vials)
- Ensures maximal product recovery for small volume thawing
- Processing multiple or single volumes at a time
- Vertical thawing as preferred thawing method
- GMP and GLP compliant
- Optimal for QC needs
- Tightly conforms to most commercially available cryovials and vials
- Versatile and low maintenance device
- For reproducible standardized process, thaw after thaw

## Reproducibility

- Automated and robust process
- Customized settings
- Reduced risk of human error

## Services

- Qualification package IQ, OQ, PQ
- Maintenance plan PM
- Expert training for users

## Accessories (optional)

- Documentation software Barkey TCP Logging Tool (via Ethernet connection)
- Bar code reader



|                        |   |
|------------------------|---|
| Product                | Barkey vialguard pro                            |
| Power supply           | 220-240 VAC 50/60 Hz or<br>100-120 VAC 50/60 Hz |
| Set temperature        | +25°C to +45°C                                  |
| Dimensions             | 340 x 600 x 320 mm (WxDxH)                      |
| Weight without filling | 21 kg   |
| Filling quantity       | 9 liters  |
| Safety class           | I   |
| Certification          | CE  |
| Barkey is certified    | DIN EN ISO 13485                                |

Technical amendments and errors expected.



Barkey GmbH & Co. KG  
Gewerbestrasse 8  
33818 Leopoldshoehe  
Germany

Tel +49 (0) 5202 / 9801-0  
Fax +49 (0) 5202 / 9801-99  
info@barkey.de  
www.barkey.de