



2D PLATFORM ROCKERS

2D Platform Rockers are ideal for gentle mixing of solutions, reagents and blotting and staining applications. Space saving design and robust construction provide durability, reliability and quiet operation.



NS-04A
Fixed Speed 2D Rocker



NS-05UA
Adjustable 2D Rocker



NS-14UA
Large Capacity Adjustable 2D Rocker

PRODUCT FEATURES:

- All models include a dimpled non-skid mat
- Large platforms with hooks for rubber strips to easily accommodate many tubes, flasks and bottles
- NS-05UA and NS-14UA models feature adjustable speed and timer
- NS-14UA model features adjustable speed, timer and tilt angle
- All models are suitable for use in cold room and standard incubator environments

Model	NS-04A	NS-05UA	NS-14UA
Product Name	24 RPM 2D Rocker	2D Rocker with Adjustable Speed and Time	Large Capacity 2D Rocker with Adjustable Speed, Tilt Angle, and Time
Speed	24 RPM	2~80 RPM	1~50 RPM
Speed Accuracy	NA	±1 RPM	
Tilt Angle	±10°		1~15° (Adjustable in 1° increments)
Timer	NA	0~99hr 59min	0~9999min
Platform Dimensions	307 x 267 mm (12.1 x 10.5 in)		
Max Load	2 kg (4.4 lb)		4.5 kg (9.9 lb)
Display	NA	LCD	LED
Accessories Included in Standard Packaging	NSA02 (Dimpled Rubber Mat) and 6 PCS of NSA01 (Rubber Strips)		
Operating Environment	Ambient Temperature: 5°C~40°C Ambient Relative Humidity: ≤80%		
Power Supply	AC110V±11V, 60Hz		
Max. Power	4W	10W	13.5W
Dimensions W×D×H	290 x 380 x 150 mm (15.0 x 11.4 x 5.9 in)		296 x 440 x 209 mm (17.3 x 11.7 x 8.2 in)
Net Weight	4.5 kg (9.9 lb)	6 kg (13.2 lb)	14 kg (30.9 lb)
IP Protection Rating	IP21		
Certifications	CE		cCSAus, CE, RoHS, WEEE

Specifications





2D Platform Rockers

2D Platform Rockers are ideal for gentle mixing of solutions, reagents and blotting and staining applications. Space saving design and robust construction provide durability, reliability and quiet operation.


Accessories and Specifications

For over 25 years Crystal Technology and Industries has been driven to deliver exceptional value to the scientific marketplace by developing innovative equipment and Cryo Storage Solutions for today's scientist.

Name	Part Number	Specifications	PCS/Pack
Rubber Strips for Securing Tubes, Flasks, Bottles to Platform (For NS-04A, NS-05UA, NS-14UA)	NSA01	WxD: 180 x 7 mm (7.1 x 0.3 in)	6
Dimpled Rubber Mat (For NS-04A, NS-05UA, NS-14UA)	NSA02	WxD: 298 x 258 mm (11.7 x 10.2 in)	1
Stacking Platform Assembly(For NS-04A, NS-05UA, NS14UA)	NSA04	WxDxH: 295 x 322 x 109 mm (11.6 x 12.7 x 4.3 in)	1



NSA01
Rubber Strips for Securing Tubes,
Flasks, Bottles to Platform



NSA02
Dimpled Rubber Mat



NSA04
Stacking Platform
Assembly