

Where Convenience Meets Compliance

EZ-DOCK SMARTBOX[®] Euthanasia Station



The EZ-Dock combines the best features of all our systems. A built-in flow circuit and times are already preset for the specific cage type to ensure AVMA compliancy. Plug in, connect the provided gas line and the system is ready for use! Once the cage is placed in dock properly, the cycles begin and shuts off when time has elapsed. The unit is designed to be as simple to use as possible.

Designed for specific cage types. More manufacturer's models coming soon.

Features:

- Compliant with AVMA Guidelines
- Cycle starts automatically when cage is docked
- Choose between adult and neonates
- LED indicator lights show when cycle is running and when complete
- Automatically shuts off if cage is removed for safety
- Uses house or bottled CO₂ supply at 15-100 psi
- Super E-Z



Shown with Tecniplast GM500 Cage.
Other Cage Types Available.



Scan to watch video

