



HEMOROAMI00

The HemoRoam line of mobile, battery-powered refrigerators enables healthcare providers to store blood between 2°C to 8°C and platelets at ambient temperatures. By reducing medical waste and saving time, the HemoRoam allows healthcare providers to transition from outdated disposable coolers to innovative, cost-effective solutions for the modern day.



THE FEATURES YOU NEED

- Outer case constructed with durable military-grade Aluminum (MIL-STD-810, MIL-STD-461, ATA 300 (Cat. 1), DIN EN 60 529)
- Ergonomic spring handles automatically stow for safe shipping
- Stackable outer case allows for efficient out-of-use storage.
- Smart change technology automatically switches charge between battery and electrical power plugs without the need for user intervention
- Digital microprocessor temperature controller with built-in low battery alarm
- Watertight, sealed, and anodized payload compartment
- Complies with cold-chain protocols of VFC, FDA, CDC, AABB & AATB
- Pad lockable

INCLUDED ACCESSORIES

- Castor wheels
- Handle bar
- File folder

OPTIONAL ACCESSORIES

- Cigarette lighter adapter
- Temperature data logger

SPECIFICATIONS

| |
|---------------------------------|
| OVERALL DIMENSIONS |
| WEIGHT |
| REFRIGERATED PAYLOAD DIMENSIONS |
| REFRIGERATED PAYLOAD CAPACITY |
| PROGRAMMABLE TEMPERATURE RANGE |
| TEMPERATURE UNIFORMITY |
| POWER INPUTS |

| |
|--|
| 31.75" X 23.25" X 33.25" |
| 135 LBS |
| 22" X 14" X 18" |
| 130 UNITS OF 300 ML RBC / 300 ML PLASMA |
| REFRIGERATOR: 1°C - 6°C |
| +/- 1.0°C AT SETPOINT 4.0°C |
| 12-24 VDC / 110-240 VAC |

BATTERY CAPACITY

| | NEW TOTAL WEIGHT | FRIDGE BATTERY LIFE | FREEZER BATTERY LIFE |
|--------|------------------|---------------------|----------------------|
| BBS50 | 180 LBS | 24-30 HRS | NA |
| BBS100 | 210 LBS | 48-60 HRS | NA |
| BBS200 | 280 LBS | 96-120 HRS | NA |