

Ultra & Precision Series Laboratory Refrigerators



	Ultra (LHUV)	Ultra Elite (LHEV)	Ultra Touch (LHTV)	Precision (LP)	Precision Touch (LPVT)
Controller					
Variable Speed Refrigeration	Bottom Mounted	Bottom Mounted	Bottom Mounted	Top Mounted	Top Mounted
Door Type Available	Glass or Solid	Glass or Solid	Glass or Solid	Glass or Solid	Glass or Solid
Thermal Performance	+/-2.5C	+/-1.5C	+/-1.0C	+/-1.0C	+/-1.0C
Advanced Airflow Technology	✓	✓	✓	✓	✓
High/Low Level Temp Alarms	✓	✓	✓	✓	✓
Remote Alarm Contacts	✓	✓	✓	✓	✓
12VDC Controller Battery Backup	X	✓	✓	✓	✓
Door Ajar Alarm	X	✓	✓	✓	✓
Touchscreen Display	X	X	✓	X	✓
Data & Event Logging and USB Download	X	X	✓	✓	✓
4-20mA	X	X	✓	X	✓
Stainless Steel Interior	X	X	X	✓	✓
Optional roll out stainless drawers	X	X	X	✓	✓
Warranty	2yr Parts & Labor, 5yr Compressor	2yr Parts & Labor, 7yr Compressor	3yr Parts & Labor, 7yr Compressor	2yrs Parts and Labor, 7yr Compressor	3yrs Parts and Labor, 7yr Compressor