

# TW 10LD

## LN<sub>2</sub> DEWAR



For safe storage and dispensing of small amounts of liquid nitrogen.



SPECIFICATIONS	10LD
LN <sub>2</sub> Capacity (L)	10
Static Hold Time (Days)	45
Evap. Rate (Liters/Day)	0.22
DIMENSIONS	
Necktube (in)	2.0
Overall Dia. (in)	11.4
Height (in)	23.5
WEIGHT	
Empty (Lbs)	14.5
Full (Lbs)	32.3

## FEATURES

- 10 Liter LN<sub>2</sub> capacity
- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Snap-on cap and necktube assures a tight closure
- Advanced insulation materials assure high thermal efficiency



855.750.8191 | [sales@ICBiomedical.com](mailto:sales@ICBiomedical.com)  
[www.ICBiomedical.com](http://www.ICBiomedical.com)  
Made In USA