

# New Product Announcement

Effective Nov 1, 2019. Subject to change without notice.

## LinkLabel™ BlueTooth Enabled Label Printer

- Compatible with BlueTooth devices
- Free app allows for limitless customization
- Hundreds of preloaded scientific symbols
- Labels in white, clear, colored and Ultra Low Temperature
- Lightweight and compact  
5.25 x 2.75 x 6.75in. (LxWxH)
- Prints barcodes, QR codes, serials and batches



The new LinkLabel™ Lab Printer brings a new level of freedom to laboratory labeling. MTC's free, downloadable LabelPro™ App allows users to endlessly customize and specialize their labels. Use the app to easily change text, font, positioning, size, and label length. The LinkLabel also has the ability to print batches, serial numbers, and a variety of barcodes.

With 12 available label cassettes, the customization doesn't end on your cell phone. Choose from multiple colors of tape, including transparent. Three label heights are available for different sized vessels. 6mm tapes are great for tissue culture plates, 12mm tapes are ideal from microtubes, and the extra tall 24mm tapes are ideal for large tubes and flasks. Special "ultra-low" label tapes are also available for storage in -80°C environments or for liquid nitrogen applications.

Cat. No.	Description	Quantity	List Price
L3000	LinkLabel™ Bluetooth Label Printer, w/ 115-230V power supply	1ea	\$299.00
<b>Label Tape Refill Cassettes</b>			
L9010-6WK	6mm lab tape, roll of 26 feet, white with black print	1ea	\$23.00
L9010-6CK	6mm lab tape, roll of 26 feet, clear with black print	1ea	\$28.00
L9010-12WK	12mm lab tape, roll of 26 feet, white with black print	1ea	\$25.00
L9010-12CK	12mm lab tape, roll of 26 feet, clear with black print	1ea	\$32.00
L9010-12BW	12mm lab tape, roll of 26 feet, blue with white print	1ea	\$25.00
L9010-12GW	12mm lab tape, roll of 26 feet, greb with white print	1ea	\$25.00
L9010-12RW	12mm lab tape, roll of 26 feet, red with white print	1ea	\$25.00
L9010-12YK	12mm lab tape, roll of 26 feet, yellow with black print	1ea	\$25.00
L9010-24WK	24mm lab tape, roll of 26 feet, white with black print	1ea	\$32.00
L9010-24CK	24mm lab tape, roll of 26 feet, clear with black print	1ea	\$41.00
L9010-12ULT	12mm lab tape, roll of 20 feet, Ultra-Low Temp, white with black print	1ea	\$28.00
L9010-24ULT	24mm lab tape, roll of 20 feet, Ultra-Low Temp, white with black print	1ea	\$41.00
<b>Also Available</b>			
L9010	LABeler™ handheld laboratory labeler, w/ 115-230V power supply	1ea	\$169.50
L9010-A	LABeler™ handheld laboratory labeler, w/out 115-230V power supply	1ea	\$129.50

Dealer Confidential:



www.mtcbio.com



(908) 279-0280



info@mtcbio.com