



The only vortexer available with Multi-Head ™
 Powerful Motor for INSTANT Vortexing
 Q-Drive™ Dynamic Balance System
 Touch or continuous operation





Unique Multi-Head allows for standard vortexing PLUS holds up to 8 microtubes

The $Mortexer^{TM}$ is MORE than your standard vortex mixer. The unique design of the Multi-Head TM (included) features a traditional cup head for general purpose mixing, plus holds up to eight microcentrifuge tubes.

The drive system is balanced with Q-Drive™ technology, providing a smooth vortexing action, while minimizing noise and excessive vibration. This results in quieter, more efficient and longer lasting operation. Increased speed, up to 3400 rpm produces an instant vortex of even the largest sample sizes, including full 50ml tubes.

Other features include: "no-walk" suction feet, cold room/incubator compatibility and spill-proof electronic controls. Two year warranty.

Technical Data:

Speed Range: 200 to 3400 rpm
Operating Modes: Touch or continuous
Operating Temp. Range: +4°C to +65°C
Dimensions: 13 x 16 x 17 cm

 5 x 6.3 x 6.7 in.

 Weight:
 3.8 kg/8.4 lb

 Electrical:
 115V AC.60 Hz

lectrical: 115V AC, 60 Hz, 1.5A 230V AC, 50 Hz, 0.75A

230V AC, 60 Hz, 0.75A

altituti

/arranty: **2 years**











Ordering Information:

BV1005* Mortexer™ Vortex Mixer (with Multi-Head)

Optional Accessories:

BV1000-FLAT Optional 3 inch flat head

BV1000-COMBO Optional head for tubes and microplate.

66 tubes (38 x 1.5ml and 28 x 0.5ml) or 1 microplate

BV1000-H15 Optional head for 12 x 1.5/2.0ml, horizontal BV1000-H150 Optional head for 4 x 15ml, horizontal

BV1000-H500 Optional head for 2 x 50ml, horizontal

To order a product in 230V / 50Hz, please add -E to item number.
 To order a product in 230V / 60Hz, please add -E2 to item number.



Benchmark Scientific Inc. PO Box 709 Edison, NJ 08817 Phone: 908-769-5555

Fax: 908-222-1864

Web: www.BenchmarkScientific.com Email: Info@BenchmarkScientific.com