



Neon Stool 9571, Bench-to-Counter Height, Freewheeling Casters, Black Vinyl, Foot Step, With 4D Arms

Article #:20025409-A*

**9-11 Week Factory Order*

Features

- The award-winning Neon Chair is known as the best lab chair globally for its design, superior quality, ergonomic features and overall comfort.
- Innovative Seat Cushion System allows for easy removal and replacements.
- All-inclusive ergonomic functionality with premium mechanism, back height lumbar adjustment, back lock, seat height adjustment, swivel-tilt, forward tilt, tilt tension and 4-D Arms.

Specifications

Height Adjustment:	Bench-to-Counter Height Range 25-1/2" - 36-1/4"
Seat Depth Range:	16-5/8" - 18-7/8"
Accent Trim Color:	In-Stock (Blue); Factory Orders (Orange, Green and Grey)
Frame Color:	Grey Back (Non-ESD Version)
Seat Material:	Black Vinyl and Enhanced Comfort Seat
Castors:	Freewheeling
Base:	Black Aluminum Base
Cleanroom Rating:	ISO Level 7
Bio-Safety BSL:	Level 2
BIFMA Rating:	Complies with Level e3 Platinum
Lab Grade Rating:	Complies with Scientific Equipment & Furniture Assn.- SEFA12
Weight Capacity:	300 lbs.
Product Weight:	40 lbs.
Warranty:	10-Year Warranty
Original Design:	German Engineering and Production
ESD Capable:	Available by Request