RBI-SHAKER

MAXIMUM CAPACITY, BENCHTOP MODEL



- + High Capacity, benchtop model
- → Circular 19mm orbit for aeration and mixing
- Exchange flask clamps in one second with MAGic Clamp™ platform.

Technical Data

Speed Range: 30 to 300 rpm Speed Increment: 1 rpm Orbit: 19mm (3/4") Max. Capacity: 5 x 2L or 9 x 1L Flasks Max. Load: 65 lbs (30kg)

Platform Dimensions: 18 x 18 in./45 x 45 cm (Included) 1 min. to 99 hours / continuous Timer:

65 lbs / 30 kg Weight:

18.5 x 21.25 x 5 in. / 47 x 54 x 12.5 cm

120V or 230V, 50-60Hz

The Orbi-Shaker XL offers flask capacity that far exceeds most benchtop models. With a large 18 x 18 inch platform and a powerful, beltless motor drive, the Orbi-Shaker XL is capable of mixing up to five 2L flasks or nine 1L flasks simultaneously. The included, non-slip rubber mat platform is ideal for low speed applications, while the convenient MAGic Clamp™ platform is also available for high speed mixing of Erlenmeyer flasks, tube racks or microplates. The MAGic Clamp™ allows for instant flask clamp exchange without the use of tools, screws or any other hardware.

H1000-MR-CMB

H1000-MR-T600 H1010-P-MP H1010-P-SP H1000-MR-T1550

Ordering Information

Orbi-Shaker XL, High capacity, benchtop shaker BT1011* Includes 18x18" platform with non-slip mat, 115V

Optional Platforms and Accessories:

H1010-MR MAGic Clamp universal platform, (18x18") Flask Clamp, 50ml (max. 49)
Flask Clamp, 125ml (max. 36)
Flask Clamp, 250ml (max. 25)
Flask Clamp, 500ml (max. 16)
Flask Clamp, 1000ml (max. 9) H1000-MR-50 H1000-MR-125 H1000-MR-250 H1000-MR-500 H1000-MR-1000 Flask Clamp, 2000ml (max. 5) Microplate Clamp (max. 8) H1000-MR-2000 H1000-MR-MP Tube Rack, 20x50 & 30x15 ml tubes (max. 4)
Tube Rack, 32x15 tubes (max. 8)
Tube Rack, 15x50 tubes (max. 8) H1000-MR-1550

H1000-MR-T15 H1000-MR-T50

Bottle Rack, 4x500/600ml conical bottles (max. 4) H1000-MR-T600

Tube Rack, 12 x 1.5/2.0ml horizontal (max. 6) Tube Rack, 4 x 15ml horizontal (max.8) H1000-MR-T15H H1000-MR-T150H

H1000-MR-T500H Tube Rack, 2 x 50ml horizontal (max.8)

Adjustable flask clamp for MAGic clamp platform, accepts 125ml (max. 36), 250ml (max. 25) or 500ml (max. 16)

Spring loaded platform, for uncommon vessels H1010-P-SP H1010-P-MP 6 x microplate platform (standard or deepwell)



PO Box 709 Edison, NJ 08818 Phone: 908-769-5555 Fax: 732-313-7007 Web: www.BenchmarkScientific.com Email: Info@BenchmarkScientific.com