

PHARMA ROAM25

The PharmaRoam line of mobile, battery-powered refrigerators and/or freezers enables healthcare providers to store vaccines and biologic products. By reducing medical waste and saving time, the PharmaRoam allows healthcare providers to transition from outdated disposable coolers to innovative, cost-effective solutions for the modern day.



THE FEATURES YOU NEED

- Outer case constructed with durable military-grade Aluminum (MIL-SPEC)
- Ergonomic spring handles automatically stow for safe shipping
- Stackable outer case allows for efficient out-of-use storage.
- Smart change technology automatically switches charge between battery and electrical power plugs without the need for user intervention
- Digital microprocessor temperature controller with built-in low battery alarm
- Watertight, sealed, an anodized payload compartment
- Complies with cold-chain protocols of VFC, FDA, CDC, AABB & AATB
- Pad lockable

SPECIFICATIONS

| | |
|---------------------------------|---|
| OVERALL DIMENSIONS | 23" X 15.5" X 23" |
| WEIGHT | 55 LBS |
| REFRIGERATED PAYLOAD DIMENSIONS | 15.5" X 7" X 9" |
| REFRIGERATED PAYLOAD CAPACITY | 4450 UNITS OF VACCINE DOSES OR 550 SINGLE USE SYRINGES |
| PROGRAMMABLE TEMPERATURE RANGE | REFRIGERATOR: 2°C - 8°C FREEZER: -20°C - 0°C |
| TEMPERATURE UNIFORMITY | +/- 1.0°C AT SETPOINT 4.0°C |
| POWER INPUTS | 12-24 VDC 110-240 VAC |

OPTIONAL ACCESSORIES

- Cigarette lighter adapter
- Mobility options
- Temperature data logger

BATTERY CAPACITY

| | NEW TOTAL WEIGHT | FRIDGE BATTERY LIFE | FREEZER BATTERY LIFE |
|--------|------------------|---------------------|----------------------|
| BBS22 | 67 LBS | 20 HRS | 4-8 HRS |
| BBS50 | 100 LBS | 40-48 HRS | 8-12 HRS |
| BBS100 | 130 LBS | 80-96 HRS | 20-24 HRS |