



ROEMER

PHARMAROAM3

The PharmaRoam line of mobile, battery-powered refrigerators and/or freezers enables healthcare providers to store vaccines and biologic products. By reducing medical waste and saving time, the PharmaRoam allows healthcare providers to transition from outdated disposable coolers to innovative, cost-effective solutions for the modern day.



THE FEATURES YOU NEED

- Outer case constructed with durable military-grade Aluminum (MIL-SPEC)
- Ergonomic spring handles automatically stow for safe shipping
- Stackable outer case allows for efficient out-of-use storage.
- Smart change technology automatically switches charge between battery and electrical power plugs without the need for user intervention
- Digital microprocessor temperature controller with built-in low battery alarm
- Watertight, sealed, an anodized payload compartment
- Complies with cold-chain protocols of VFC, FDA, CDC, AABB & AATB
- Pad lockable

SPECIFICATIONS

OVERALL DIMENSIONS	15.08" X 11.34" X 13.38"
WEIGHT	29 LBS
REFRIGERATED PAYLOAD DIMENSIONS	4.5" X 2.5" X 8"
REFRIGERATED PAYLOAD CAPACITY	480 UNITS OF VACCINE DOSES OR 60 SINGLE USE SYRINGES
PROGRAMMABLE TEMPERATURE RANGE	REFRIGERATOR: 2°C - 8°C FREEZER: -20°C - 0°C
TEMPERATURE UNIFORMITY	+/- 0.5°C AT SETPOINT 4.0°C
POWER INPUTS	12-24 VDC 110-240 VAC

BATTERY CAPACITY

*EXTERNAL BATTERY STORED BELOW UNIT, ADDS 7" IN HEIGHT

OPTIONAL ACCESSORIES

- Cigarette lighter adapter
- Mobility options
- Temperature data logger

	NEW TOTAL WEIGHT	FRIDGE BATTERY LIFE	FREEZER BATTERY LIFE
BBS22	43 LBS	24-28 HRS	4-6 HRS
BBS50	74 LBS	48-52 HRS	10-12 HRS