

PHARMAROAM60

The PharmaRoam line of mobile, battery-powered refrigerators and/or freezers enables healthcare providers to store vaccines and biologic products. By reducing medical waste and saving time, the PharmaRoam allows healthcare providers to transition from outdated disposable coolers to innovative, cost-effective solutions for the modern day.



THE FEATURES YOU NEED

- Outer case constructed with durable military-grade Aluminum (MIL-SPEC)
- Ergonomic spring handles automatically stow for safe shipping
- Stackable outer case allows for efficient out-of-use storage.
- Smart change technology automatically switches charge between battery and electrical power plugs without the need for user intervention
- Digital microprocessor temperature controller with built-in low battery alarm
- Watertight, sealed, an anodized payload compartment
- Complies with cold-chain protocols of VFC, FDA, CDC, AABB & AATB
- Pad lockable

SPECIFICATIONS

OVERALL DIMENSIONS

31.75" X 23.25" X 24.75"

WEIGHT

100 LBS

REFRIGERATED PAYLOAD DIMENSIONS

22" X 14" X 10"

REFRIGERATED PAYLOAD CAPACITY

**9800 UNITS OF VACCINE DOSES
OR 1230 SINGLE USE SYRINGES**

PROGRAMMABLE TEMPERATURE RANGE

**REFRIGERATOR: 2°C - 8°C
FREEZER: -20°C - 0°C**

TEMPERATURE UNIFORMITY

+/- 1.0°C AT SETPOINT 4.0°C

POWER INPUTS

**12-24 VDC
110-240 VAC**

BATTERY CAPACITY

	NEW TOTAL WEIGHT	FRIDGE BATTERY LIFE	FREEZER BATTERY LIFE
BBS50	145 LBS	36-40 HRS	8-10 HRS
BBS100	175 LBS	80-96 HRS	20-24 HRS
BBS200	245 LBS	160-196 HRS	36-48 HRS

OPTIONAL ACCESSORIES

- Cigarette lighter adapter
- Mobility options
- Temperature data logger