


MVE Research Dewars

The MVE Research Dewars are developed specifically for easy and safe transport of liquid nitrogen and samples within facilities. These units are vacuum insulated for superior thermal performance and long hold times.



Features include:

- All stainless steel construction
- Insulated lid
- Wide mouth for easy access
- Convenient carrying handle (except .5L)
- Meets laboratory safety requirements that prohibit glass lined dewars



	RD-6	RD-3	RD-2	RD-1	RD-1W	RD-0.5
Performance						
LN2 Capacity <i>L</i>	6	3	2	1	1	0.5
Handle	yes	yes	yes	yes	yes	no
Unit Dimensions						
Neck Opening <i>in. (mm)</i>	7.3 (185)	7.3 (185)	3.9 (99)	3.3 (84)	3.9 (99)	2.6 (66)
Inner Diameter <i>in. (mm)</i>	7.3 (185)	7.3 (185)	3.9 (99)	3.3 (84)	3.9 (99)	2.6 (66)
Outer Diameter <i>in. (mm)</i>	7.8 (198)	7.8 (198)	4.8 (122)	4.3 (109)	4.8 (122)	3.4 (86)
Usable Internal Height <i>in. (mm)</i>	10.6 (269)	6.3 (160)	7.1 (180)	6.2 (157)	6.1 (155)	7.1 (180)
Overall Height <i>in. (mm)</i>	11.8 (300)	7.5 (191)	12.3 (312)	9.1 (231)	7.0 (178)	8.0 (203)
NER without cover <i>L/hr</i>	0.4	0.2	0.1	0.1	0.2	0.1

90 Day Warranty