## **Rugged and Versatile Lab Rotator!**



## Another revolutionary product from the makers of the Vortex-Genie<sup>®</sup>!

- O Holds a wide variety of tubes, flasks, microplates, bags, and odd-shaped vessels at any angle.
- O Available in two INTERCHANGEABLE platform configurations.
- O Standard platform is designed to turn 90 degrees allowing different mixing motions.
- O MAGNETIC platform features unique spring clip plate system.
- O Designed for cold room and incubator use. Will easily fit in most incubators.
- O Uniquely designed for a wide selection of motions adjustable from gentle rolling to forceful end over end mixing.
- Ocontinuous-duty quality motor provides smooth, quiet operation.

## **Ordering Information**

Description
-------------

Catalog No.

Catalog No.	Description
SI-2100	Rotator Genie, Standard Platform, 120V
SI-2200	Rotator Genie, Magnetic Platform, 120V
SI-2101	Rotator Genie, Standard Platform, 230V, No Plug
SI-2201	Rotator Genie, Magnetic Platform, 230V, No Plug
SI-2102	Rotator Genie, Standard Platform, 230V, Euro Plug
SI-2202	Rotator Genie, Magnetic Platform, 230V, Euro Plug
SI-2103	Rotator Genie, Standard Platform, 230V, British Plug
SI-2203	Rotator Genie, Magnetic Platform, 230V, British Plug
SI-2104	Rotator Genie, Standard Platform, 230V, Swiss Plug
SI-2204	Rotator Genie, Magnetic Platform, 230V, Swiss Plug
SI-2106	Rotator Genie, Standard Platform, 240V, Aust. Plug
SI-2206	Rotator Genie, Magnetic Platform, 240V, Aust. Plug
SI-2105	Rotator Genie, Standard Platform, 100V
SI-2205	Rotator Genie, Magnetic Platform, 100V

# 



© 2014 Scientific Industries LTP0074 Rev. B

#### Home of the Multi-Task Genie<sup>™</sup> Family 80 Orville Drive, Suite 102 Bohemia, NY 11716 U.S.A. Phone: (631) 567-4700 • Fax: (631) 567-5896 Toll Free: 888-850-6208 Customer Service: custsvc@scientificindustries.com scientificindustries.com

### Supplied

Speed Overall Dimensions (W x D x H)

Platform Dimensions Standard Platform Magnetic Platform Weight Environmental

### **Specifications**

Standard Platform with Retaining Bands or Magnetic Platform 5-35 RPM 400 x 175 (290 platform horizontal) x 400mm (15.8 x 7.0 (11.5 platform horizontal) x 15.8 in)

260 x 260mm (10.25 x 10.25 in) 290 x 315mm (11.5 x 12.5 in) 5.2 kgs (11.5 lbs) Installation Category 2, 0° - 38°C (32° - 100°F), 95% max. relative humidity

## Rotator Genie<sup>™</sup> Accessories



Fig. No.	Catalog No.	Description
1	SI-2110	Standard Platform with Bands
2	SI-2210	Magnetic Platform
3	SI-1123(*)	1 Clip Plate for 12 each 10-13mm tubes (magnetic platform capacity: 12 clip plates)
4	SI-1124(*)	1 Clip Plate for 6 each 15-17mm tubes (magnetic platform capacity: 12 clip plates)
5	SI-1125(*)	1 Clip Plate for 3 each 28-30mm tubes (magnetic platform capacity: 12 clip plates)
6	SI-1120(*)	6 Clip Plates for 12 each 10-13mm tubes (magnetic platform capacity: 12 clip plates)
7	SI-1121(*)	6 Clip Plates for 6 each 15-17mm tubes (magnetic platform capacity: 12 clip plates)
8	SI-1122(*)	6 Clip Plates for 3 each 28-30mm tubes (magnetic platform capacity: 12 clip plates)
9	SI-1126(*)	Bag Mounting Strips (pack of 16)
10	SI-1130(*)	2 Clip Plates for Hybridization Tubes
11	SI-1134(*)	2 Universal Clip Plates
12	SI-1135(*)	1 Magnetic Covered 100 Microtube Box
13	SI-1136(*)	1 Clip Plate for 100ml Volumetric Flasks
14	SI-1137(*)	Dual Port Mixing Tube (use with SI-1130)

(\*) - Requires magnetic platform, (Cat. No. SI-2210), which is supplied with all SI-22XX models or can be purchased separately.



© 2014 Scientific Industries LTP0074 Rev. B

Home of the Multi-Task Genie<sup>™</sup> Family 80 Orville Drive, Suite 102 Bohemia, NY 11716 U.S.A. Phone: (631) 567-4700 • Fax: (631) 567-5896 Toll Free: 888-850-6208 Customer Service: custsvc@scientificindustries.com scientificindustries.com