



TUBE REVOLVERS

Powerful Tube Revolver with Five Unique Interchangeable Rotisseries



TR-01A



TR-02UA

PRODUCT FEATURES:

- Ergonomically designed for simple operation and minimal space usage
- Includes 5 unique interchangeable rotisseries
- Rotisseries can be mixed to accommodate a variety of tube sizes simultaneously
- Maintenance-free motor is powerful and quiet
- TR-02UA can be programmed to oscillate at every 40° to simulate hand mixing with rotation and shaking

For over 25 years Crystal Technology and Industries has been driven to deliver exceptional value to the scientific marketplace by developing innovative equipment and Cryo Storage Solutions for today's scientist.

Specifications

Model	TR-01A	TR-02UA
Product Name	Fixed Speed Tube Revolver	Variable Speed Tube Revolver
Speed	25 RPM	10~40 RPM
Speed Accuracy	NA	±1 RPM
Capacity	84x0.5/0.8mL / 60x1.5/2.0mL / 28x5/7mL / 28x10/15mL / 6x50mL Tubes	
Display	NA	LED
Accessories Included in Standard Packaging	2xTRA01 (10/15mL Paddle) + 2xTRA02 (5/7mL Paddle) + 2xTRA03 (1.5/2mL Paddle) + 2xTRA04 (0.5/0.8mL Paddle) + 1xTRA05 (50mL Tube Paddle)	
Operating Environment	Ambient Temperature: 4°C~60°C Relative Humidity: ≤ 80%	
Power Supply	AC110V, 60Hz	
Max. Power	4W	10W
Dimensions WxDxH	148 x 225 x 195 mm (5.8 x 8.9 x 7.7 in)	148 x 260 x 195 mm (5.8 x 10.2 x 7.7 in)
Net Weight	0.8 kg (1.8 lb)	1 kg (2.2 lb)
Protection Level	IP21	
Certifications	CE, RoHS, WEEE	cCSAus, CE, RoHS, WEEE



Available magnetic platform for mixing liquids in bags





TUBE REVOLVERS

Accessories



2x

TRA01
10ml/15 ml Tubes



2x

TRA02
5ml/7ml Tubes



2x

TRA03
1.5ml/2ml Tubes

Both models include the following 5 rotisserie paddles



2x

TRA04
0.5ml/0.5ml Tubes



1x

TRA05
50ml Tubes

Specifications

For over 25 years Crystal Technology and Industries has been driven to deliver exceptional value to the scientific marketplace by developing innovative equipment and Cryo Storage Solutions for today's scientist.

Name	Part Number	Capacity	PCS/Pack
Rotisserie Paddle for 10/15mL Tubes	TRA01	7 x 10/15mL Tubes	2
Rotisserie Paddle for 5/7mL Tubes	TRA02	7 x 5/7mL Tubes	2
Rotisserie Paddle for 1.5/2mL Tubes	TRA03	15 x 1.5/2mL Tubes	2
Rotisserie Paddle for 0.5/0.8mL Tubes	TRA04	21 x 0.5/0.8mL Tubes	2
Rotisserie Paddle for 50mL Tubes	TRA05	6 x 50mL Tubes	1
Magnetic Platform for Mixing Bags	TRA06	Platform Dimensions: 156x226 mm (6.1x8.9 in)	1



Available magnetic platform for mixing liquids in bags

