

TW 35HC

LN₂ HIGH-CAPACITY SMALL FREEZER

For long-term storage of large quantities at cryogenic temperatures.



FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Numbered index ring and tapered internal canister locator
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid



SPECIFICATIONS	35HC
LN ₂ Capacity (L)	35
Static Hold Time (Days)	130
Evap. Rate (Liters/Day)	0.27
DIMENSIONS	
Necktube (in)	4.7
Overall Dia. (in)	18.8
Height (in)	26.8
WEIGHT	
Empty (Lbs)	39.0
Full Weight (Lbs)	101.3
CAPACITY	
Number of Canisters	10
Canister Dimensions (in)	2.64 x 11
1.2ml & 2.0 ml vials (5/cane)	850
1.2ml & 2.0 ml vials (6/cane)	1020
.5 cc Straws (10/cane)	2800
.5 cc Straws- Bulk (2 Levels)	7300
.25 cc Straws (28/cane)	7280
.25 cc Straws- Bulk (2 Levels) in Goblets	16400



800.886.2796 | Fax 317. 297.7988
sales@ICBiomedical.com
www.ICBiomedical.com