

# TW 4LD

## LN<sub>2</sub> DEWAR



For safe storage and dispensing of small amounts of liquid nitrogen.



SPECIFICATIONS	4LD
LN <sub>2</sub> Capacity (L)	4
Static Hold Time (Days)	10
Evap. Rate (Liters/Day)	0.40
DIMENSIONS	
Necktube (in)	1.2
Outside Dia. (in)	7.6
Height (in)	17
WEIGHT	
Empty (Lbs)	6.6
Full (Lbs)	13.7

## FEATURES

- 4 Liter LN<sub>2</sub> capacity
- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Snap-on cap and necktube assures a tight closure
- Advanced insulation materials assure high thermal efficiency



855.750.8191 | [sales@ICBiomedical.com](mailto:sales@ICBiomedical.com)  
[www.ICBiomedical.com](http://www.ICBiomedical.com)  
Made In USA