

LABeler™ Lab Printer



- + Variety of label tapes for all common tube sizes
- + Prints a variety of 3D barcodes
- + Lightweight, for one handed operation
- + Built in cutter eliminates the need for scissors
- + "Ultra-Low" tape for liquid nitrogen storage



Perfect for labelling microtubes and cryovials



Eliminates the smearing/bleeding associated with markers & pens after autoclaving

Proper identification of scientific samples is often a key objective of the research laboratory. The process of writing on sample tubes/containers can be messy and unreliable. It can also be difficult to find the time to access a computer and design a layout for expensive ink jet labels. The new LABeler™ Lab Printer provides the convenience of a lab pen with the benefits of a large desktop printer. Simply type in the desired copy using the 50+ button key pad and press the print key... within seconds your label can be cut by pressing the cut lever and is ready for application.



With over 15 available sizes, colors & material combinations, the LABeler is ideal for label identification of any laboratory vessel for nearly any laboratory application. Standard label tapes can be autoclaved, boiled or even frozen to temperatures as low as -50°C. Special Ultra-Low label tapes are also available for storage in -80°C environments or for liquid nitrogen applications.



• Replaceable tape cassettes (26' length)



• Built in cutter, eliminates the need for scissors



• Available with 9V power supply (also runs off 6 AA batteries)

Ordering Information:

Item No.	Description:
L9010*	LABeler Lab Printer with 115V power supply
L9010-A*	LABeler Lab Printer without power supply (battery use)
<u>Standard Lab Tape (withstands temperatures of -50°C to 130°C)</u>	
L9010-6WK	26' cassette of 6mm lab tape, white with black print
L9010-6CK	26' cassette of 6mm lab tape, clear with black print
L9010-12WK	26' cassette of 12mm lab tape, white with black print
L9010-12CK	26' cassette of 12mm lab tape, clear with black print
L9010-12YK	26' cassette of 12mm lab tape, yellow with black print
L9010-12RW	26' cassette of 12mm lab tape, red with white print
L9010-12BW	26' cassette of 12mm lab tape, blue with white print
L9010-12GW	26' cassette of 12mm lab tape, green with white print
L9010-24WK	26' cassette of 24mm lab tape, white with black print
L9010-24CK130	26' cassette of 24mm lab tape, clear with black print
<u>Ultra-Low Temperature Tape, liquid nitrogen applications & temp. as low as -100°C</u>	
L9010-12ULT	20' cassette of 12mm ULT tape, white with black print
L9010-24ULT	20' cassette of 24mm ULT tape, white with black print