

D. For Autoclave (Mobile) Applications (continued)

Autoclavable Cage Racks

- Preconfigured, autoclavable cage racks available in two footprints – 18x60 (457x1524mm) and 24x60 (610x1524mm).
- Units are 69" (1753mm) high, including casters. **Casters must be ordered separately.**
- Carts come standard with swaged posts and aluminum split sleeves.

Width/Length/Height (in.)	Width/Length/Height (mm)	Approx. Pkd. Wt.		Description	Cat. No.
		(lbs.)	(kg)		
18x60x69	457x1524x1753	136	61.7	5-shelf, 18" (457mm) wide cage rack	LAB0681S5CL
24x60x69	610x1524x1753	166	75.3	5-shelf, 24" (610mm) wide cage rack	LAB0642S5CL

Note: To build a cage rack according to specific requirements, follow instructions (D) on previous page. Refer to Autoclave guidelines, below, for important checklist to follow.

See page 195 for Caster Selection.



LAB0642S5CL
(casters sold separately)

Autoclave Cage Rack Covers

Width/Length/Height (in.)	Width/Length/Height (mm)	Approx. Pkd. Wt.		Description	Cat. No.
		(lbs.)	(kg)		
20x62x61.5	508x1575x1562	6.0	2.7	Autoclavable cart cover for 18" (457mm) wide Cage Rack	LABCA0681CL
26x62x61.5	660x1575x1562	7.0	3.2	Autoclavable cart cover for 24" (610mm) wide Cage Rack	LABCA0642CL

Note: Autoclave cart covers are rated for approximately 70 cycles per cover. Color may vary from image shown.



LAB0642S5CL with
autoclave cart cover
(casters sold separately)

Autoclavable Lab Carts

- Preconfigured, autoclavable Lab carts are available in two footprints – 18x24 (457x610mm) and 24x36 (610x914mm).
- Choose between 2- and 3-shelf models
- Units are 39" (991mm) high, including casters. **Casters must be ordered separately.**
- Carts come standard with swaged posts and aluminum split sleeves.

Width/Length/Height (in.)	Width/Length/Height (mm)	Approx. Pkd. Wt.		Description	Cat. No.
		(lbs.)	(kg)		
18x24x39	457x610x991	25.0	11.3	2-shelf, 18" (457mm) wide Lab cart	LAB4281S2CL
18x24x39	457x610x991	37.5	17.0	3-shelf, 18" (457mm) wide Lab cart	LAB4281S3CL
24x36x39	610x914x991	46.5	21.1	2-shelf, 24" (610mm) wide Lab cart	LAB6342S2CL
24x36x39	610x914x991	65.5	29.7	3-shelf, 24" (610mm) wide Lab cart	LAB6342S3CL

Note: To build a lab cart according to specific requirements, follow instructions (D) on previous page. Refer to Autoclave guidelines, below, for important checklist to follow.

See page 195 for Caster Selection.



LAB4281S2CL
(casters sold separately)

Autoclave Lab Cart Covers

Width/Length/Height (in.)	Width/Length/Height (mm)	Approx. Pkd. Wt.		Description	Cat. No.
		(lbs.)	(kg)		
20x26x33.5	508x660x851	3.0	1.4	Autoclavable cart cover for 18" (457mm) wide Cage Rack	LABCA4281CL
26x38x33.5	660x965x851	3.5	1.6	Autoclavable cart cover for 24" (610mm) wide Cage Rack	LABCA6342CL

Note: Autoclave cart covers are rated for approximately 70 cycles per cover.

Autoclave Application Guidelines

When stainless solid shelving, as part of a mobile shelving unit, is used in an autoclave, the following Metro product guidelines are recommended to ensure optimal product performance and longevity:

Stainless, swaged posts with aluminum post caps must be utilized [part numbers 33/54/63UPS-SW]

Aluminum split sleeves [part number 9986S] with stainless rings should be used to mount solid stainless shelves.

Select either hi-temp phenolic casters or hi-temp nylon casters for application. These casters are designed specifically for autoclave environments. Bumpers should not be used as they are not designed to withstand autoclave temperature ranges.