

FREE
Equipment
with purchase of item below*

(Free with items: C2417, BT4001, B3D5000, or H5000-H/HC)

Your choice of:



Value: \$390

COMBI-Rotor™
Microcentrifuge

MyFuge™ 12

Item: C1012

Value: \$199.90



Vortex Speed:
2,800rpm (fixed)

Vornado™

Item: BV101-B



Value:
\$375

Tube Rotator
Accepts tubes 1.5 to 50ml

Roto-Mini™

Item: R2020



Touch
Control

- Full color touch screen
- mySpin™ 1 button program recall
- COMBI-Rotor™ accepts tubes and strips
- Compact, only 9 x 12 x 7.75 in.
- Optional rotors for 5ml tubes and 2ml Cryovials

MC-24™ Touch
High Speed Microcentrifuge

Max. Speed: 13,500rpm / 16,800xg
Max Capacity: 24 x 1.5/2.0ml
2 x PCR Strips (16 x 0.2ml)

C2417	MC-24 Touch Microcentrifuge with 24 place COMBI-Rotor	\$2179.00	\$1852.15
C2417-R5	Optional 8 x 5ml rotor for MC-24 Touch	\$922.00	\$783.70
C2417-R5-ADP2	Cryovial adapter for C2417-R5 rotor, 8/pk	\$115.00	\$97.75



- Designed to fit and maximize capacity in most CO₂ incubators and chambers
- All delicate electronics located in the remote controller, protected from harsh environments
- Adjust shaker settings without disturbing the CO₂ environment
- Instantly exchange flask clamps with MAGic Clamp Accessories

ORBI-SHAKER™ CO₂
All Environment Shaker




BT4001	Orbi-Shaker CO ₂ with remote controller & rubber mat platform	\$3265.00	\$2775.25
BT3000-MR	MAGic Clamp universal platform for flasks & tube racks	\$255.50	\$217.18

See page 2 for information on available accessories clamps.



- 3D motion for superior mixing
- Adjustable speed and tilt
- Ideal for staining and general mixing



Optional Stacking Platform

Speed: 5-105rpm
Platform: 12.5 x 13 in.
Tilt Angle: 0 to ±10°


BENCHWAVER™
Undulating 3-D Shaker

B3D5000	BenchWaver Shaker	\$1290.00	\$1096.50
B3D5000-STK	Optional stacking platform	\$130.00	\$110.50
B3D5000-DIMP	Optional dimpled mat	\$86.50	\$73.53



- Temperature control from 0° to 100°C
- Exchangeable blocks, for tube sizes 0.2ml to 50ml
- Precision contoured wells for uniform thermal transfer
- Two models: Heat/Cool or Heat only
- Sequence link function for linking or repeating programs (HC model only)

MULTI-THERM™
COOL - HEAT - SHAKE



H5000-HC	MultiTherm Shaker with heating and cooling	\$2679.00	\$2277.15
H5000-H	MultiTherm Shaker with heating only	\$2039.00	\$1733.15

Please visit www.bsilab.com for a complete listing of all available block options.

Ph: 1-866-674-6294 - Email: info@bsilab.com - Web: www.bsilab.com

Must mention promo code BP1.21 with part number of requested free product.
Valid Through 10.31.21. Pricing, availability, and specifications subject to change without notice.

BenchPress

Vol. 1.1 2021
Use Promo Code: BP1.21



Shakers Incubators Centrifuges Balances Rotators Vortexers DryBaths Rocks

NEW BEAD BLASTER™ 96
BALL-MILL HOMOGENIZER FOR PLATES, TUBES AND JARS



- High-throughput homogenization
2 x 96 well plates 8 x 15/50ml tubes
80 x 2.0ml tubes 48 x 4ml vials
- Powerful mixing action, homogenize bone, teeth and other difficult samples
- For DNA/RNA extractions from plant and animal tissue, QuEChERS, sample preparations for cannabis testing, & capable of wet/dry grinding



Intro Special!
Item: IPD9600 ~~\$10,995~~

Intro Special: **\$9995**

- Includes:
- Microplate holder (2 ea.) Item: IPD9600-96
 - Microtube holder (2ea.) Item: IPD9600-96A2 (Value: \$1700)

Use Promo Code: BP1.21

Wide Range of Accessory Holders

BeadBath Duo™
For
Beads or Water



Item: B2402
2L BeadBath Duo
~~\$689~~ **\$585.65**

No more refills from evaporation

NEW

NEW QüadCount™
Automated Cell Counter

- Fast, accurate and reproducible counting and live/dead determinations
- Intuitive software with auto-focus and touch-screen interface
- Cost-effective QuadSlides™ - disposable slides with four individual chambers

Sample Volume: 20µl Image Detection: 5MP CMOS Sensor
Cell Size Range: 5-60µm Dynamic Range: 1 x 10⁴ - 1 x 10⁷ cells/ml

E7500	QuadCount Cell Counter	\$3,695	\$3,140.75
E7500-S1	QuadSlide, 4 chamber slides	\$149	\$126.65

Replaces manual cell counting...ask for a demo today!

Place your order with BSI Today!

Phone: 866.674.6294 Email: info@bsilab.com Web: www.bsilab.com






SMARTREADER™ 96

- Intuitive, 7 inch, color touch screen control panel
- Filter based, wavelength range from 340-750nm (405, 450, 492, 630nm included)
- Operates as a stand alone system, with USB flash drive for data transfer
- Integrated, multi-speed plate shaker

Your go-to instrument for ELISA and absorbance assays

MR9600	SmartReader 96 Microplate Absorbance Reader	\$5500.00	\$4675.00
MR9600-T	SmartReader 96 Microplate Absorbance Reader with Incubation	\$5885.00	\$4976.75
MR9600-TP	SmartReader Thermal Printer	\$325.00	\$276.25

\$4675




PLATEFUGE™

MICROPLATE MICROCENTRIFUGE



- Quickly spins down droplets in PCR plates
- Unique rotor design prevents sample spillage
- 50% smaller than traditional centrifuges
- Accepts all popular PCR plates (skirted, semi & non) & most microtiter plates



C2000

\$633.25

MORTEXER™

MORE Power, MORE Flexibility... MORE Vortexer!

- Standard cup head PLUS accepts 8 x 2.0ml tubes
- Powerful motor, INSTANT vortex
- Variable speed (touch or continuous)

BV1005

\$246.42



NEW!

ROTO-THERM™ Plus & ROTO-MINI™ Plus

INCUBATED TUBE ROTATOR

Specifications:
 Temp Range: Ambient +5° to 60°C
 Speed Range: 5 to 70rpm (Plus model)
 Max. Capacity: 6 x 50ml, 16 x 15ml, 24 x 0.5ml, 44 x 1.5/2.0ml, 96 x 0.2ml

- Temperature control up to 60°C (H2024 only)
- Includes tube holders for 0.5 to 50ml
- Available with variable speed & mix modes Oscillating ● Pausing ● Rocking

Included Holders
 ● 11 x 1.5-2ml (2/pk.)
 ● 5 x 5-7ml (2/pk.)
 ● 4 x 10-15ml (2/pk.)
 ● 6 x 50ml* (1 ea.)
 *Also accepts 24 x 0.5ml



H2024	Roto-Therm™ Plus w/included tube holders (variable speed)	\$735.00	\$624.75
R2024	Roto-Mini™ Plus, w/included tube holders (variable speed)	\$442.00	\$375.70
R2020-PCR	Optional tube holder for 96 x 0.2 (12 x 8 position PCR strips)	\$72.00	\$61.20

Additional tube holders available, please visit our website for a complete listing.

NEW!

BEADBLASTER™ 24R

REFRIGERATED MICROTUBE HOMOGENIZER

- Completely lyse, grind and homogenize a variety of sample types
- New! **POWERFUL REFRIGERATION** to protect delicate samples
 - True refrigeration - no dry ice or compressed air needed
 - Space saving design
- Process up to 24 samples simultaneously
- Wide variety of grinding media in prefilled 2ml and 5ml tubes

Special Offer

Specifications:
 Capacity: 24 x 2ml tubes
 Speed range: 4.00 - 7.00 m/sec
 Cooling: Compressor
 Temp. Range: -10°C to Ambient
 Temp. Uniformity: ±2°C
 Dimensions: 32 x 51 x 41cm
 Weight: 40 kg
 Warranty: 2 years

D2400-R	BeadBlaster 24R with refrigeration	\$12,990.00	\$10,999.00
D2400-RV5	Optional tube holder for 5ml tubes	\$379.00	\$322.15
D1032-05	Prefilled 2ml tubes, 0.5mm zirconium	\$158.00	\$134.30
D1032-15	Prefilled 2ml tubes, 1.5mm zirconium	\$158.00	\$134.30
D1032-30	Prefilled 2ml tubes, 3mm zirconium	\$158.00	\$134.30
D1033-28	Prefilled 2ml tubes, 2.8mm steel	\$158.00	\$134.30
D1034-30	Prefilled 5ml tubes, 3mm zirconium	\$219.00	\$186.15
D1034-28	Prefilled 5ml tubes, 2.8mm steel	\$211.00	\$179.35

Additional grinding media is available in prefilled tubes and loose bulk. Ask for details.

Also available: D2400 Beadblaster 24 Microtube Homogenizer **\$7,050.75**



INCUSHAKER™

SHAKING INCUBATORS

- Large flask capacity
- **MAGic™** Clamp accessories for instant-magnetic clamp exchange
- Available with CO₂ control or cooling
- **NEW!** Adjustable clamp for 125, 250 and 500ml flasks (4 FREE with purchase)

Limited Time Offer*
FREE!
 4 ea. **MAGic Clamp (Adjustable)**
PATENTED Item: H1000-MR-CMB (Accepts 125ml, 250ml and 500ml Flasks)

Use Promo Code: BP1.21

INCUSHAKER Mini
 Speed Range: 30 to 300rpm
 Temp Range: Amb. +5° to 70°C
 Platform Dims: 9.5 x 11.5 in.
 Universal Platform: H1000-MR
 Flask Capacity: 2 x 1L (H1000-MR-1000), 5 x 500ml (H1000-MR-500), 8 x 250ml (H1000-MR-250), 12 x 125ml (H1000-MR-125)
\$2,600.15
 Item: H1001-M

INCUSHAKER 10L & 10LR
 Speed Range: 30 to 300rpm
 Temp Range: (H1012, Cooling Model) Amb. +5° to 60°C, Amb. -15° to 60°C
 Platform Dims: 17.8 x 17.8 in.
 Universal Platform: H1010-MR
 Flask Capacity: 5 x 2L (H1000-MR-2000), 9 x 1L (H1000-MR-1000), 16 x 500ml (H1000-MR-500), 25 x 250ml (H1000-MR-250), 36 x 125ml (H1000-MR-125)
\$5,300.15 / \$7,497
 Item: H2010 & H2012 (Cooling)

INCUSHAKER Mini CO₂
 Speed Range: 30 to 300rpm
 Temp Range: Amb. +5° to 60°C
 CO₂ Range: 0 to 20%
 Platform Dims: 9.5 x 11.5 in.
 Universal Platform: H1000-MR
 Flask Capacity: 5 x 500ml (H1000-MR-500), 8 x 250ml (H1000-MR-250), 12 x 125ml (H1000-MR-125)
\$7,790.25
 Item: H3501

Z216-M High Speed Microcentrifuges



- Bundles Includes **NEW COMBI-Rotor™** with Quick-Seal™ lid
- Available with refrigeration
- Compact, only 11 inches wide
- High speed, up to **21,631xg**
- Versatile, **10** optional rotors



Z216-CMB	Z216-M Microcentrifuge Bundle w/ COMBI-Rotor	\$4378	\$3700
Z216-K-CMB	Z216-MK Refrigerated Microcentrifuge Bundle w/ COMBI-Rotor	\$7334	\$5800

Max. Speed: 15,000 rpm Max. g-Force: 21,631xg
 Max. Capacity: 44 x 1.5/2.0ml Temp. Range: -20 to 40°C (-MK)
 Dimensions: 11 x 15 x 11" 11 x 21 x 11" (-MK)

Use Promo Code: BP1.21

Place your order with BSI Today!
 Phone: 866.674.6294 Email: info@bsilab.com Web: www.bsilab.com

Place your order with BSI Today!
 Phone: 866.674.6294 Email: info@bsilab.com Web: www.bsilab.com