


Bin Organizer with Lids

Mixed Use Storage Bins

- ◆ Excellent for organizing -20°C and -80°C Freezers & Cold Rooms
- ◆ Can accept a wide variety of tubes, bottles, boxes, etc.
- ◆ Constructed of corrosion resistant stainless steel




Part Number	Description	Dimension L x H x W	Capacity l x h x w	Shipping Weight (lb/Kg)
UFT-S04-5	Bin organizer with lids	421 x 120 x 138 mm 16 ⁹ / ₁₆ x 4 ³ / ₄ x 5 ⁷ / ₁₆ inch	Drawer Area: 16 ¹ / ₈ x 4 ³ / ₄ x 5 ³ / ₁₆ inch 409 x 120 x 132 mm	2.8 / 1.29
UFT-S06-5	Bin organizer with lids	421 x 168 x 138 mm 16 ⁹ / ₁₆ x 6 ⁵ / ₈ x 5 ⁷ / ₁₆ inch	Drawer Area: 16 ¹ / ₈ x 6 ³ / ₈ x 5 ³ / ₁₆ inch 409 x 168 x 132 mm	3.2 / 1.43
UFT-S08-5	Bin organizer with lids	421 x 225 x 138 mm 16 ⁹ / ₁₆ x 8 ⁷ / ₈ x 5 ⁷ / ₁₆ inch	Drawer Area: 16 ¹ / ₈ x 8 ³ / ₄ x 5 ³ / ₁₆ inch 409 x 223 x 132 mm	3.5 / 1.58
UFT-S11-5	Bin organizer with lids	421 x 280 x 138 mm 16 ⁹ / ₁₆ x 11 x 5 ⁷ / ₁₆ inch	Drawer Area: 16 ¹ / ₈ x 10 ¹⁵ / ₁₆ x 5 ³ / ₁₆ inch 409 x 278 x 132 mm	3.8 / 1.73
UFT-S12-5	Bin organizer with lids	421 x 325 x 138 mm 16 ⁹ / ₁₆ x 12 ¹³ / ₁₆ x 5 ⁷ / ₁₆ inch	Drawer Area: 16 ¹ / ₈ x 12 ¹¹ / ₁₆ x 5 ³ / ₁₆ inch 409 x 323 x 132 mm	3.9 / 1.78
UFT-M04-5	Bin organizer with lids	560 x 120 x 138 mm 22 ¹ / ₁₆ x 4 ³ / ₄ x 5 ⁷ / ₁₆ inch	Drawer Area: 21 ⁹ / ₁₆ x 4 ³ / ₄ x 5 ³ / ₁₆ inch 548 x 120 x 132 mm	3.5 / 1.60
UFT-M06-5	Bin organizer with lids	560 x 168 x 138 mm 22 ¹ / ₁₆ x 6 ⁵ / ₈ x 5 ⁷ / ₁₆ inch	Drawer Area: 21 ⁹ / ₁₆ x 6 ³ / ₈ x 5 ³ / ₁₆ inch 548 x 168 x 132 mm	3.9 / 1.76
UFT-M08-5	Bin organizer with lids	560 x 225 x 138 mm 22 ¹ / ₁₆ x 8 ⁷ / ₈ x 5 ⁷ / ₁₆ inch	Drawer Area: 21 ⁹ / ₁₆ x 8 ³ / ₄ x 5 ³ / ₁₆ inch 548 x 223 x 132 mm	4.3 / 1.94
UFT-M11-5	Bin organizer with lids	560 x 280 x 138 mm 22 ¹ / ₁₆ x 11 x 5 ⁷ / ₁₆ inch	Drawer Area: 21 ⁹ / ₁₆ x 10 ¹⁵ / ₁₆ x 5 ³ / ₁₆ inch 548 x 278 x 132 mm	4.7 / 2.11
UFT-M12-5	Bin organizer with lids	560 x 325 x 138 mm 22 ¹ / ₁₆ x 12 ¹³ / ₁₆ x 5 ⁷ / ₁₆ inch	Drawer Area: 21 ⁹ / ₁₆ x 12 ¹¹ / ₁₆ x 5 ³ / ₁₆ inch 548 x 323 x 132 mm	5.0 / 2.26
UFT-L04-5	Bin organizer with lids	680 x 120 x 138 mm 26 ³ / ₄ x 4 ³ / ₄ x 5 ⁷ / ₁₆ inch	Drawer Area: 26 ⁵ / ₁₆ x 4 ³ / ₄ x 5 ³ / ₁₆ inch 668 x 120 x 132 mm	4.2 / 1.89
UFT-L06-5	Bin organizer with lids	680 x 168 x 138 mm 26 ³ / ₄ x 6 ⁵ / ₈ x 5 ⁷ / ₁₆ inch	Drawer Area: 26 ⁵ / ₁₆ x 6 ³ / ₈ x 5 ³ / ₁₆ inch 668 x 168 x 132 mm	4.6 / 2.07
UFT-L08-5	Bin organizer with lids	680 x 225 x 138 mm 26 ³ / ₄ x 8 ⁷ / ₈ x 5 ⁷ / ₁₆ inch	Drawer Area: 26 ⁵ / ₁₆ x 8 ³ / ₄ x 5 ³ / ₁₆ inch 668 x 223 x 132 mm	5.0 / 2.26
UFT-L11-5	Bin organizer with lids	680 x 280 x 138 mm 26 ³ / ₄ x 11 x 5 ⁷ / ₁₆ inch	Drawer Area: 26 ⁵ / ₁₆ x 10 ¹⁵ / ₁₆ x 5 ³ / ₁₆ inch 668 x 278 x 132 mm	5.4 / 2.45
UFT-L12-5	Bin organizer with lids	680 x 325 x 138 mm 26 ³ / ₄ x 12 ¹³ / ₁₆ x 5 ⁷ / ₁₆ inch	Drawer Area: 26 ⁵ / ₁₆ x 12 ¹¹ / ₁₆ x 5 ³ / ₁₆ inch 668 x 323 x 132 mm	5.8 / 2.62
UFT-S04-11	Bin organizer with lids	421 x 120 x 279 mm 16 ⁹ / ₁₆ x 4 ³ / ₄ x 11 inch	Drawer Area: 16 ¹ / ₈ x 4 ⁵ / ₈ x 10 ¹¹ / ₁₆ inch 409 x 118 x 272 mm	5.2 / 2.35
UFT-S06-11	Bin organizer with lids	421 x 168 x 279 mm 16 ⁹ / ₁₆ x 6 ⁵ / ₈ x 11 inch	Drawer Area: 16 ¹ / ₈ x 6 ³ / ₈ x 10 ¹¹ / ₁₆ inch 409 x 166 x 272 mm	5.7 / 2.58
UFT-S08-11	Bin organizer with lids	421 x 225 x 279 mm 16 ⁹ / ₁₆ x 8 ⁷ / ₈ x 11 inch	Drawer Area: 16 ¹ / ₈ x 8 ³ / ₄ x 10 ¹¹ / ₁₆ inch 409 x 223 x 272 mm	6.2 / 2.83
UFT-S11-11	Bin organizer with lids	421 x 280 x 279 mm 16 ⁹ / ₁₆ x 11 x 11 inch	Drawer Area: 16 ¹ / ₈ x 10 ¹⁵ / ₁₆ x 10 ¹¹ / ₁₆ inch 409 x 278 x 272 mm	6.8 / 3.09
UFT-S12-11	Bin organizer with lids	421 x 325 x 279 mm 16 ⁹ / ₁₆ x 12 ¹³ / ₁₆ x 11 inch	Drawer Area: 16 ¹ / ₈ x 12 ¹¹ / ₁₆ x 10 ¹¹ / ₁₆ inch 409 x 323 x 272 mm	7.3 / 3.30
UFT-M04-11	Bin organizer with lids	560 x 120 x 279 mm 22 ¹ / ₁₆ x 4 ³ / ₄ x 11 inch	Drawer Area: 21 ⁹ / ₁₆ x 4 ⁵ / ₈ x 10 ¹¹ / ₁₆ inch 548 x 118 x 272 mm	6.5 / 2.95
UFT-M06-11	Bin organizer with lids	560 x 168 x 279 mm 22 ¹ / ₁₆ x 6 ⁵ / ₈ x 11 inch	Drawer Area: 21 ⁹ / ₁₆ x 6 ³ / ₈ x 10 ¹¹ / ₁₆ inch 548 x 166 x 272 mm	7.1 / 3.20
UFT-M08-11	Bin organizer with lids	560 x 225 x 279 mm 22 ¹ / ₁₆ x 8 ⁷ / ₈ x 11 inch	Drawer Area: 21 ⁹ / ₁₆ x 8 ³ / ₄ x 10 ¹¹ / ₁₆ inch 548 x 223 x 272 mm	7.7 / 3.48
UFT-M11-11	Bin organizer with lids	560 x 280 x 279 mm 22 ¹ / ₁₆ x 11 x 11 inch	Drawer Area: 21 ⁹ / ₁₆ x 10 ¹⁵ / ₁₆ x 10 ¹¹ / ₁₆ inch 548 x 278 x 272 mm	8.3 / 3.76
UFT-M12-11	Bin organizer with lids	560 x 325 x 279 mm 22 ¹ / ₁₆ x 12 ¹³ / ₁₆ x 11 inch	Drawer Area: 21 ⁹ / ₁₆ x 12 ¹¹ / ₁₆ x 10 ¹¹ / ₁₆ inch 548 x 323 x 272 mm	8.8 / 3.99
UFT-L04-11	Bin organizer with lids	680 x 120 x 279 mm 26 ³ / ₄ x 4 ³ / ₄ x 11 inch	Drawer Area: 26 ⁵ / ₁₆ x 4 ⁵ / ₈ x 10 ¹¹ / ₁₆ inch 668 x 118 x 272 mm	7.6 / 3.44
UFT-L06-11	Bin organizer with lids	680 x 168 x 279 mm 26 ³ / ₄ x 6 ⁵ / ₈ x 11 inch	Drawer Area: 26 ⁵ / ₁₆ x 6 ³ / ₈ x 10 ¹¹ / ₁₆ inch 668 x 166 x 272 mm	8.2 / 3.71
UFT-L08-11	Bin organizer with lids	680 x 225 x 279 mm 26 ³ / ₄ x 8 ⁷ / ₈ x 11 inch	Drawer Area: 26 ⁵ / ₁₆ x 8 ³ / ₄ x 10 ¹¹ / ₁₆ inch 668 x 223 x 272 mm	8.8 / 4.00
UFT-L11-11	Bin organizer with lids	680 x 280 x 279 mm 26 ³ / ₄ x 11 x 11 inch	Drawer Area: 26 ⁵ / ₁₆ x 10 ¹⁵ / ₁₆ x 10 ¹¹ / ₁₆ inch 668 x 278 x 272 mm	9.5 / 4.30
UFT-L12-11	Bin organizer with lids	680 x 325 x 279 mm 26 ³ / ₄ x 12 ¹³ / ₁₆ x 11 inch	Drawer Area: 26 ⁵ / ₁₆ x 12 ¹¹ / ₁₆ x 10 ¹¹ / ₁₆ inch 668 x 323 x 272 mm	10.0 / 4.55

Storage Drawers with Lids

Mixed Use Storage Drawers

- ◆ Excellent for organizing -20°C and -80°C Freezers & Cold Rooms
- ◆ Can accept a wide variety of tubes, bottles, boxes, etc.
- ◆ Constructed of corrosion resistant stainless steel



Part Number	Description	Dimension L x H x W	Capacity l x h x w	Shipping Weigh (lb/Kg)
UFTD-SL-5	Drawer style bin organizer	16 ³ / ₄ x 9 ¹⁵ / ₁₆ x 5 ³ / ₈ inch 425 x 252 x 143 mm	Drawer Area: 16 ¹ / ₈ x 4 ³ / ₄ x 5 ³ / ₁₆ inch (Per Drawer) 409 x 120 x 132 mm	8.7 / 3.95
UFTD-SH-5	Drawer style bin organizer	16 ³ / ₄ x 13 ³ / ₄ x 5 ³ / ₈ inch 425 x 349 x 143 mm	Drawer Area: 16 ¹ / ₈ x 6 ⁵ / ₈ x 5 ³ / ₁₆ inch (Per Drawer) 409 x 168 x 132 mm	9.6 / 4.35
UFTD-ML-5	Drawer style bin organizer	22 ¹ / ₄ x 9 ¹⁵ / ₁₆ x 5 ³ / ₈ inch 565 x 252 x 143 mm	Drawer Area: 21 ⁹ / ₁₆ x 4 ³ / ₄ x 5 ³ / ₁₆ inch (Per Drawer) 548 x 120 x 132 mm	10.8 / 4.91
UFTD-MH-5	Drawer style bin organizer	22 ¹ / ₄ x 13 ³ / ₄ x 5 ³ / ₈ inch 565 x 349 x 143 mm	Drawer Area: 21 ⁹ / ₁₆ x 6 ⁵ / ₈ x 5 ³ / ₁₆ inch (Per Drawer) 548 x 168 x 132 mm	11.9 / 5.35
UFTD-LL-5	Drawer style bin organizer	27 x 9 ¹⁵ / ₁₆ x 5 ³ / ₈ inch 685.5 x 252 x 143 mm	Drawer Area: 26 ⁵ / ₁₆ x 4 ³ / ₄ x 5 ³ / ₁₆ inch (Per Drawer) 668 x 120 x 132 mm	12.8 / 5.81
UFTD-LH-5	Drawer style bin organizer	27 x 13 ³ / ₄ x 5 ³ / ₈ inch 685.5 x 349 x 143 mm	Drawer Area: 26 ⁵ / ₁₆ x 6 ⁵ / ₈ x 5 ³ / ₁₆ inch (Per Drawer) 668 x 168 x 132 mm	13.9 / 6.29
UFTD-SL-11	Drawer style bin organizer	16 ³ / ₄ x 9 ¹⁵ / ₁₆ x 11 ³ / ₁₆ inch 425 x 252 x 284 mm	Drawer Area: 16 ¹ / ₈ x 4 ³ / ₄ x 10 ¹¹ / ₁₆ inch (Per Drawer) 409 x 118 x 272 mm	15.7 / 7.14
UFTD-SH-11	Drawer style bin organizer	16 ³ / ₄ x 13 ³ / ₄ x 11 ³ / ₁₆ inch 425 x 349 x 284 mm	Drawer Area: 16 ¹ / ₈ x 6 ⁵ / ₈ x 10 ¹¹ / ₁₆ inch (Per Drawer) 409 x 166 x 272 mm	17.1 / 7.76
UFTD-ML-11	Drawer style bin organizer	22 ¹ / ₄ x 9 ¹⁵ / ₁₆ x 11 ³ / ₁₆ inch 565 x 252 x 284 mm	Drawer Area: 21 ⁹ / ₁₆ x 4 ³ / ₄ x 10 ¹¹ / ₁₆ inch (Per Drawer) 548 x 118 x 272 mm	19.8 / 8.99
UFTD-MH-11	Drawer style bin organizer	22 ¹ / ₄ x 13 ³ / ₄ x 11 ³ / ₁₆ inch 565 x 349 x 284 mm	Drawer Area: 21 ⁹ / ₁₆ x 6 ⁵ / ₈ x 10 ¹¹ / ₁₆ inch (Per Drawer) 548 x 166 x 272 mm	21.3 / 9.65
UFTD-LL-11	Drawer style bin organizer	27 x 9 ¹⁵ / ₁₆ x 11 ³ / ₁₆ inch 685.5 x 252 x 284 mm	Drawer Area: 26 ⁵ / ₁₆ x 4 ³ / ₄ x 10 ¹¹ / ₁₆ inch (Per Drawer) 668 x 118 x 272 mm	23.2 / 10.53
UFTD-LH-11	Drawer style bin organizer	27 x 13 ³ / ₄ x 11 ³ / ₁₆ inch 685.5 x 349 x 284 mm	Drawer Area: 26 ⁵ / ₁₆ x 6 ⁵ / ₈ x 10 ¹¹ / ₁₆ inch (Per Drawer) 668 x 166 x 272 mm	24.8 / 11.23