



Boxes, Cell Dividers, Label Cards

Cardboard Boxes



Part Number	Description	Color	Dimension L x W x H
SB1CC-196	1" water repellent freezer box with 196 cell divider for PCR tubes.	White	5 ¹ / ₄ x 5 ¹ / ₄ x 1 inch 133 x 133 x 26 mm
SB2CC	White or yellow water repellent cardboard boxes for storage of samples or tubes. Supplied with or without cell dividers.		5 ¹ / ₄ x 5 ¹ / ₄ x 2 inch 133 x 133 x 51 mm
SB2CC-XX*			
SB2C	Standard boxes have 4 drain holes in each corner.	Yellow	5 ¹ / ₄ x 5 ¹ / ₄ x 2 inch 133 x 133 x 51 mm
SB2C-XX*			
SB2CC-LP	White water repellent cardboard box with low profile lid.	White	5 ⁵ / ₁₆ x 5 ⁵ / ₁₆ x 1 ¹⁵ / ₁₆ inch 134.8 x 134.8 x 48.8 mm
SB2CC-XX*-LP	Supplied with or without cell dividers.		
SB2CC-H-81	Hold up to repeated freezing & thawing. Water proof fiberboard.	White	5 ¹ / ₄ x 5 ¹ / ₄ x 2 inch 133 x 133 x 51 mm
SB2CC-H-81-P-W	Hinged lid.	Purple	
SB2CC-H-81-G-W	With 81-place cell divider. 3 colors optional.	Green	
SB3C	Large 3" white water repellent freezer boxes.	White	5 ¹ / ₄ x 5 ¹ / ₄ x 3 inch 133 x 133 x 75 mm
SB3C-XX*			
SB3.75C			
SB3.75C-XX*	White water repellent cardboard boxes for storage of samples or tubes. Supplied with or without cell dividers. Standard boxes have 4 drain holes in each corner.	White	5 ¹ / ₄ x 5 ¹ / ₄ x 3 ⁷ / ₈ inch 133 x 133 x 98mm
SB4C	White water repellent cardboard boxes for 15ml or 50ml centrifuge tubes are supplied with or without cell dividers. Boxes have 4 drain holes in each corner.	White	5 ¹³ / ₁₆ x 5 ¹³ / ₁₆ x 4 ¹⁵ / ₁₆ inch 148 x 148 x 125 mm
SB4C-XXX**			
SB10C	Extra large yellow water repellent freezer boxes.	Yellow	10 ¹¹ / ₁₆ x 5 ³ / ₁₆ x 8 ⁵ / ₁₆ inch 271.5 x 132.5 x 211.5 mm
BB3C	Yellow water repellent cardboard boxes for storage of samples or tubes. Supplied with or without cell dividers. Standard boxes have 4 drain holes in each corner.	Yellow	5 ³ / ₄ x 5 ³ / ₄ x 3 inch 146 x 146 x 76 mm
BB3C-100	Large 3" yellow water repellent freezer boxes with 100 cell divider included.		
MB2C	White water repellent mini-cardboard boxes for storage of samples or tubes. Supplied with or without 25-cell dividers. Mini Boxes have drain slots.	White	3 x 3 x 2 inch 75 x 75 x 51 mm
MB2C-25			
MB3C			
MB3C-25			


* To order, replace XX with 100, 81, 64, 49, 36, 25 or 16 for cell dividers.

** To order, replace XXX with 16 or 36 for cell dividers.

Boxes, Cell Dividers, Label Cards

Standard Box Dividers

Standard cardboard cell dividers (CD-XX) fit all 2" and 3" standard cardboard and metal boxes. Dividers for 4" centrifuge tube boxes (CD4-XX) and mini box (MCD-25) are also available.



◆ Sold in packs of 12.


Part Number	Cell Number	Dimension L x D x H	Cell Opening
CD-169	169	4 ⁷ / ₈ x 4 ⁷ / ₈ x 1 ⁵ / ₁₆ inch 123.4 x 123.4 x 33 mm	³ / ₈ inch (8.74mm)
CD-144	144		³ / ₈ inch (9.52mm)
CD-100	100		⁷ / ₁₆ inch (11.5mm)
CD-96	96	7 x 4 ¹¹ / ₁₆ x 1 ¹ / ₂ inch 178 x 118.4 x 38 mm	⁹ / ₁₆ inch (14mm)
CD-64	64	4 ⁷ / ₈ x 4 ⁷ / ₈ x 1 ⁵ / ₁₆ inch 123.4 x 123.4 x 33 mm	⁹ / ₁₆ inch (14.6mm)
CD-49	49		¹¹ / ₁₆ inch (16.7mm)
CD-36	36		¹³ / ₁₆ inch (20mm)
CD-25	25		¹⁵ / ₁₆ inch (24.5mm)
CD-16	16		1 ³ / ₁₆ inch (30.2mm)
CD4-16	16		5 ⁷ / ₁₆ x 5 ⁷ / ₁₆ x 2 ³ / ₄ inch 138 x 138 x 70 mm
CD4-36	36	⁷ / ₈ inch (22.25mm)	
CD4-64	64	³ / ₈ inch (16.45mm)	
CD-49-W	49	4 ⁷ / ₈ x 4 ⁷ / ₈ x 1 ¹ / ₄ inch 124.3 x 124.3 x 31.5 mm	¹¹ / ₁₆ inch (17mm)
CD-81	81	4 ⁷ / ₈ x 4 ⁷ / ₈ x 1 ⁵ / ₁₆ inch 124.5 x 124.5 x 33 mm	¹ / ₂ inch (12.7mm)
CD-100-W	100		⁷ / ₁₆ inch (11.5mm)
MCD-25	25	2 ¹¹ / ₁₆ x 2 ¹¹ / ₁₆ x 1 ⁵ / ₁₆ inch 67.6 x 67.6 x 33 mm	¹ / ₂ inch (12.8mm)

(Continued on next page)

Boxes, Cell Dividers, Label Cards

Polycarbonate Boxes


Polycarbonate 81 or 100 cell 2" boxes hold most popular cryo and cold storage vials and fit into standard freezer racks. Boxes have clear lids with locator grids to facilitate sample identification. Polycarbonate boxes are designed to withstand temperature extremes of -196°C to 121°C. Boxes have drain holes for venting and safety. Keyed lid prevents misalignment of grid.



Part Number	Cell Number	Dimension L x D x H	Cell Opening
PB2-100	100	5 ¹ / ₄ x 5 ¹ / ₄ x 2 ¹ / ₁₆ inch 133 x 133 x 52 mm	1/2 inch (12.5 mm)
PB2-81	81	5 ¹ / ₄ x 5 ¹ / ₄ x 2 ¹ / ₁₆ inch 133 x 133 x 52 mm	1/2 inch (13 mm)
PB4-81	81	5 ¹ / ₄ x 5 ¹ / ₄ x 3 ¹ / ₁₆ inch 133 x 133 x 94 mm	1/2 inch (13 mm)
MPB2-25	25	3 x 3 x 2 ¹ / ₁₆ inch 76 x 76 x 53 mm	3/8 inch (12.8 mm)

Polypropylene Boxes

Polypropylene 81 or 100 cell 2" boxes hold most popular cryo and cold storage vials and fit into standard freezer racks. These boxes offer a convenient and cost effective alternative to polycarbonate boxes for storage of samples. Boxes have clear lids with locator grids to facilitate sample identification. Polypropylene boxes are designed to withstand temperature down to -80°C. Storage of polypropylene boxes in liquid nitrogen is not recommended. Boxes have drain holes for venting and safety. Keyed lid prevents misalignment of grid.



Part Number	Cell Number	Dimension L x D x H	Cell Opening
PPP2-100	100	5 ¹ / ₄ x 5 ¹ / ₄ x 2 inch 133 x 133 x 51 mm	7/16 inch (10.5 mm)
PPP2-81	81	5 ¹ / ₄ x 5 ¹ / ₄ x 2 inch 133 x 133 x 51 mm	1/2 inch (12.5 mm)

100-Cell Hinged Plastic Boxes

100-cell hinged top plastic boxes made of polypropylene offer an alternative storage system for all microtubes up to 2.0 ml and standard 1.2 ml cryovials. The grid is alpha-numerically matrixed with 1-10 and A-J characters hot stamped directly onto the grid for fast, easy sample identification. Boxes fit into standard freezer racks and are available in natural and five different colors.



Part Number	Cell Number	Dimension L x D x H	Cell Opening
PP2-100	100	5 ⁹ / ₁₆ x 6 x 2 ¹ / ₈ inch 141 x 152 x 54mm	1/2 inch (12.5 mm)
PP2-100-A			
PP2-100X*			

* Specify Color: B (Blue), G (Green), O (Orange), P (Pink), or Y (Yellow).


(Continued on next page)

(Continued from last page)

Boxes, Cell Dividers, Label Cards

Metal Boxes - Stainless steel

Stainless steel boxes come with or without cardboard dividers and fit into standard freezer racks.




Part Number	Description	Dimension Wt x Wt x H	Dimension Wb x Wb x Hb
SB2S	Standard 2" Box with Drain Hole	5 1/4 x 5 1/4 x 2 inch 133 x 133 x 51 mm	5 x 5 x 1 15/16 inch 127.7 x 127.7 x 49 mm
SB2S-XX*	Standard 2" Box w/Cell Divider & Drain Hole		
SB3S	Standard 3" Box with Drain Hole	5 1/4 x 5 1/4 x 3 inch 133 x 133 x 75 mm	5 x 5 x 2 15/16 inch 127.7 x 127.7 x 74 mm
SB3S-XX*	Standard 3" Box w/Cell Divider & Drain Hole		

* To order, replace XX with 144,100,81,64,49,36,25 or 16 for cell dividers .

Metal Boxes - Aluminum

Aluminum boxes come with or without cardboard dividers and fit into standard freezer racks.



Part Number	Description	Dimension Wt x Wt x H	Dimension Wb x Wb x Hb
SB2A	Standard 2" Box with Drain Hole	5 1/4 x 5 1/4 x 2 inch 133 x 133 x 51 mm	5 x 5 x 1 15/16 inch 127.7 x 127.7 x 49 mm
SB2A-XX*	Standard 2" Box w/Cell Divider & Drain Hole		
SB3A	Standard 3" Box with Drain Hole	5 1/4 x 5 1/4 x 3 inch 133 x 133 x 75 mm	5 x 5 x 2 15/16 inch 127.7 x 127.7 x 74 mm
SB3A-XX*	Standard 3" Box w/Cell Divider & Drain Hole		


* To order, replace XX with 144,100,81,64,49,36,25 or 16 for cell dividers .

Freezer Rack Labeling and Identification Tags

Specifically identify you racks using our label cards or metal tags. Plastic label cards are used for upright freezer drawer racks and are placed into the card holders. Brass metal tags with either numbers or letters may be used for upright freezer, chest freezer and liquid nitrogen tank racks. Brass tags are available in several colors : natural, blue, red, yellow or green.

Label Cards for Drawer Racks

Pack of 50



Part Number	Description	Dimension
DL-P	Plastic Cards	2 11/16 x 1 3/8 inch 68 x 35 mm

(Continued on next page)