



Upright Freezer Racks

For Standard 2" Boxes

- ◆ Racks are for standard 2" cardboard (white or yellow), plastic (polypropylene or polycarbonate) or metal (stainless or aluminum) boxes.
- ◆ Constructed of corrosion resistant stainless steel.
- ◆ Security lock devices are available for each rack.



Rack Only




Rack w/Cardboard Box



Rack w/Metal Box



Rack w/Polycarbonate Plastic Box




w/DP Lock Device Option

* UF-262 / UF-362 / UF-462 / UF-562 (∅:52.6mm), fitted with 2" cardboard box.



White 2" Cardboard Box





Polycarbonate Plastic Box



Metal Box
S-Stainless Steel
A-Aluminum













[Check page 108 for Boxes]

| Part Number | Description | Dimension L x H x W | Compartment Size l x h x w | Capacity (Box) | Lock Device Option** | Shipping Weight lb(Kg) |
|-------------|-------------|---|--|---|-------------------------|------------------------------|
| UF-222 | Rack Only | 11 1/2 x 4 1/2 x 5 1/2 inch 292 x 114 x 140 mm | 5 3/8 x 2 3/16 x 5 7/16 inch 136x 55 x 138 mm | 2 x 2 =4  | DP-1005 | 2.3(1.04) |
| UF-232 | Rack Only | 11 1/2 x 6 3/8 x 5 1/2 inch 292 x 169 x 140 mm | 5 3/8 x 2 3/16 x 5 7/16 inch 136x 55 x 138 mm | 2 x 3 =6  | DP-1007 | 3.1(1.40) |

(Continued on next page)

Upright Freezer Racks

For Standard 2" Boxes

| Part Number | Description | Dimension L x H x W | Compartment Size l x h x w | Capacity (Box) | Lock Device Option** | Shipping Weight lb(Kg) |
|--------------|---------------------------|---|---|--|-------------------------|--------------------------------|
| UF-242 | Rack Only | 11 ¹ / ₂ x 8 ⁷ / ₈ x 5 ¹ / ₂ inch 292 x 225 x 140 mm | 5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 136x 55 x 138 mm | 2 x 4 =8  | DP-1009 | 3.9(1.76) |
| UF-252 | Rack Only | 11 ¹ / ₂ x 11 x 5 ¹ / ₂ inch 292 x 280 x 140 mm | 5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 136x 55 x 138 mm | 2 x 5 =10  | DP-1011 | 4.5(2.06) |
| UF-262 | Rack Only | 11 ¹ / ₂ x 12 ⁵ / ₈ x 5 ¹ / ₂ inch 292 x 321 x 140 mm | 5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 136 x 52.6 x 138 mm | 2 x 6 =12  | DP-1013 | 5.2(2.35) |
| UF-272 | Rack Only | 11 ¹ / ₂ x 15 ³ / ₈ x 5 ¹ / ₂ inch 292 x 391 x 140 mm | 5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 136x 55 x 138 mm | 2 x 7 =14  | DP-1015 | 6.0(2.72) |
| UF-332 | Rack Only | 16 ⁷ / ₈ x 6 ⁵ / ₈ x 5 ¹ / ₂ inch 429 x 168 x 140 mm | 5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 136x 55 x 138 mm | 3 x 3 =9  | DP-1607 | 4.3(1.97) |
| UF-342 | Rack Only | 16 ⁷ / ₈ x 8 ⁷ / ₈ x 5 ¹ / ₂ inch 429 x 225 x 140 mm | 5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 136x 55 x 138 mm | 3 x 4 =12  | DP-1609 | 5.5(2.49) |
| UF-342-W-XX* | Rack w/2" box and divider | | | | | 5.5(2.49)+12ω ¹ *** |
| UF-352 | Rack Only | 16 ⁷ / ₈ x 11 x 5 ¹ / ₂ inch 429 x 280 x 140 mm | 5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 136x 55 x 138 mm | 3 x 5 =15  | DP-1611 | 6.6(2.97) |
| UF-362 | Rack Only | 16 ⁷ / ₈ x 12 ⁵ / ₈ x 5 ¹ / ₂ inch 429 x 321 x 140 mm | 5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 136 x 52.6 x 138 mm | 3 x 6 =18  | DP-1613 | 7.5(3.38) |
| UF-372 | Rack Only | 16 ⁷ / ₈ x 15 ³ / ₈ x 5 ¹ / ₂ inch 429 x 391 x 140 mm | 5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 136x 55 x 138 mm | 3 x 7 =21  | DP-1615 | 8.5(3.87) |
| UF-432 | Rack Only | 22 ³ / ₁₆ x 6 ⁵ / ₈ x 5 ¹ / ₂ inch 564 x 168 x 140 mm | 5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 136x 55 x 138 mm | 4 x 3 =12  | DP-2107 | 5.8(2.62) |
| UF-442 | Rack Only | 22 ³ / ₁₆ x 8 ⁷ / ₈ x 5 ¹ / ₂ inch 564 x 225 x 140 mm | 5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 136x 55 x 138 mm | 4 x 4 =16  | DP-2109 | 7.1(3.23) |
| UF-442-W-XX* | Rack w/2" box and divider | | | | | 7.1(3.23)+16ω ¹ *** |
| UF-452 | Rack Only | 22 ³ / ₁₆ x 11 x 5 ¹ / ₂ inch 564x 280 x 140 mm | 5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 136x 55 x 138 mm | 4 x 5 =20  | DP-2111 | 8.5(3.86) |
| UF-442-W-XX* | Rack w/2" box and divider | | | | | 8.5(3.86)+20ω ¹ *** |
| UF-462 | Rack Only | 22 ³ / ₁₆ x 12 ⁵ / ₈ x 5 ¹ / ₂ inch 564 x 321 x 140 mm | 5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 136 x 52.6 x 138 mm | 4 x 6 =24  | DP-2113 | 9.7(4.38) |
| UF-462-W-XX* | Rack w/2" box and divider | | | | | 9.7(4.38)+24ω ¹ *** |
| UF-472 | Rack Only | 22 ³ / ₁₆ x 15 ³ / ₈ x 5 ¹ / ₂ inch 564 x 391 x 140 mm | 5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 136x 55 x 138 mm | 4 x 7 =28  | DP-2115 | 11.2(5.06) |
| UF-532 | Rack Only | 27 ⁵ / ₁₆ x 6 ⁵ / ₈ x 5 ¹ / ₂ inch 693 x 168 x 140 mm | 5 ¹ / ₄ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 134 x 55 x 138 mm | 5 x 3 =15  | DP-2707 | 7.0(3.18) |
| UF-542 | Rack Only | 27 ⁵ / ₁₆ x 8 ⁷ / ₈ x 5 ¹ / ₂ inch 693 x 225 x 140 mm | 5 ¹ / ₄ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 134 x 55 x 138 mm | 5 x 4 =20  | DP-2709 | 8.6(3.91) |
| UF-542-W-XX* | Rack w/2" box and divider | | | | | 8.6(3.91)+20ω ¹ *** |
| UF-552 | Rack Only | 27 ⁵ / ₁₆ x 11 x 5 ¹ / ₂ inch 693 x 280 x 140 mm | 5 ¹ / ₄ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 134 x 55 x 138 mm | 5 x 5 =25  | DP-2711 | 10.3(4.66) |
| UF-562 | Rack Only | 27 ⁵ / ₁₆ x 12 ⁵ / ₈ x 5 ¹ / ₂ inch 693 x 321 x 140 mm | 5 ¹ / ₄ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 134 x 52.6 x 138 mm | 5 x 6 =30  | DP-2715 | 11.7(5.32) |
| UF-572 | Rack Only | 27 ⁵ / ₁₆ x 15 ³ / ₈ x 5 ¹ / ₂ inch 693 x 391 x 140 mm | 5 ¹ / ₄ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 134 x 55 x 138 mm | 5 x 7 =35  | DP-2717 | 13.5(6.14) |

* To order, replace XX with either 81 or 100 cell divider.

** Lock device is ordered separately.


*** ω¹=0.22lb(0.1kg) for cardboard box, 0.37lb(0.17kg) for plastic box, 0.84lb(0.38kg) for stainless steel box, 0.53lb(0.24kg) for aluminum box.

**** UF-262, UF-362, UF-462, UF-562, can be fitted with 2" cardboard box only, not applied to plastic box.

Upright Freezer Eco-Racks

For Standard 2" Boxes

- ◆ Designed economically for best cost performance.
- ◆ For standard 2" cardboard (white or yellow), plastic (polypropylene or polycarbonate) or metal (stainless or aluminum) boxes.
- ◆ Constructed of corrosion resistant stainless steel.




Rack Only




Rack w/Cardboard Box



White 2" Cardboard Box



Polycarbonate Plastic Box















Metal Box
S-Stainless Steel
A-Aluminum

[Check page 108 for Boxes]

Upright Freezer Eco-Racks

For Standard 2" Boxes

| Part Number | Description | Dimension L x H x W | Compartment Size l x h x w | Capacity (Box) | Lock Device Option | Shipping Weight lb(Kg) |
|-------------|-------------|--|--|--|-----------------------|------------------------------|
| UFE-332 | Rack Only | 16 ¹³ / ₁₆ x 6 ⁵ / ₈ x 5 ¹ / ₂ inch 427 x 168 x 140 mm | 5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 137 x 55 x 138 mm | 3 x 3 =9  | DP-1607 | 3.1(1.42) |
| UFE-342 | Rack Only | 16 ¹³ / ₁₆ x 8 ⁷ / ₈ x 5 ¹ / ₂ inch 427 x 225 x 140 mm | 5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 137 x 55 x 138 mm | 3 x 4 =12  | DP-1609 | 3.6(1.62) |
| UFE-352 | Rack Only | 16 ¹³ / ₁₆ x 11 x 5 ¹ / ₂ inch 427 x 280 x 140 mm | 5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 137 x 55 x 138 mm | 3 x 5 =15  | DP-1611 | 3.9(1.79) |
| UFE-362 | Rack Only | 16 ¹³ / ₁₆ x 12 ⁵ / ₈ x 5 ¹ / ₂ inch 427 x 321 x 140 mm | 5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 137 x 52.6 x 138 mm | 3 x 6 =18  | DP-1613 | 4.8(2.19) |
| UFE-432 | Rack Only | 22 ¹ / ₈ x 6 ⁵ / ₈ x 5 ¹ / ₂ inch 562 x 168 x 140 mm | 5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 137 x 55 x 138 mm | 4 x 3 =12  | DP-2107 | 4.0(1.81) |
| UFE-442 | Rack Only | 22 ¹ / ₈ x 8 ⁷ / ₈ x 5 ¹ / ₂ inch 562 x 225 x 140 mm | 5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 137 x 55 x 138 mm | 4 x 4 =16  | DP-2109 | 4.5(2.04) |
| UFE-452 | Rack Only | 22 ¹ / ₈ x 11 x 5 ¹ / ₂ inch 562 x 280 x 140 mm | 5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 137 x 55 x 138 mm | 4 x 5 =20  | DP-2111 | 5.0(2.26) |
| UFE-462 | Rack Only | 22 ¹ / ₈ x 12 ⁵ / ₈ x 5 ¹ / ₂ inch 562 x 321 x 140 mm | 5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 137 x 52.6 x 138 mm | 4 x 6 =24  | DP-2113 | 6.1(2.78) |
| UFE-472 | Rack Only | 22 ¹ / ₈ x 15 ³ / ₈ x 5 ¹ / ₂ inch 562 x 391 x 140 mm | 5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 137 x 55 x 138 mm | 4 x 7 =28  | DP-2115 | 6.8(3.10) |
| UFE-532 | Rack Only | 27 ³ / ₁₆ x 6 ⁵ / ₈ x 5 ¹ / ₂ inch 691 x 168 x 140 mm | 5 ⁵ / ₁₆ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 134.8 x 55 x 138 mm | 5 x 3 =15  | DP-2707 | 4.8(2.17) |
| UFE-542 | Rack Only | 27 ³ / ₁₆ x 8 ⁷ / ₈ x 5 ¹ / ₂ inch 691 x 225 x 140 mm | 5 ⁵ / ₁₆ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 134.8 x 55 x 138 mm | 5 x 4 =20  | DP-2709 | 5.4(2.46) |
| UFE-552 | Rack Only | 27 ³ / ₁₆ x 11 x 5 ¹ / ₂ inch 691 x 280 x 140 mm | 5 ⁵ / ₁₆ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 134.8 x 55 x 138 mm | 5 x 5 =25  | DP-2711 | 5.9(2.68) |