



Upright Freezer Racks

For Standard 3" Boxes

- ◆ Constructed of corrosion resistant stainless steel.
- ◆ Security lock devices are available for each rack.




[Check page 108 for Boxes]




Rack Only



Rack w/Cardboard Box



Rack w/Metal Box



w/DP Lock Device Option

* UF-223 / UF-323 / UF-423 / UF-523 (4:82.5mm)

Part Number	Description	Dimension L x H x W	Compartment Size l x h x w	Capacity (Box)	Lock Device Option**	Shipping Weight lb(Kg)
UF-223	Rack Only	11 1/2 x 6 5/8 x 5 1/2 inch 292 x 168 x 140 mm	5 3/8 x 3 1/4 x 5 7/16 inch 136 x 82.5 x 138 mm	2 x 2 =4	DP-1107	2.7(1.23)
UF-233	Rack Only	11 1/2 x 9 3/8 x 5 1/2 inch 292 x 238 x 140 mm	5 3/8 x 3 1/16 x 5 7/16 inch 136 x 78 x 138 mm	2 x 3 =6	DP-1110	3.5(1.60)
UF-243	Rack Only	11 1/2 x 12 2/8 x 5 1/2 inch 292 x 321 x 140 mm	5 3/8 x 3 3/8 x 5 7/16 inch 136 x 79.2 x 138 mm	2 x 4 =8	DP-1113	4.4(2.01)
UF-253	Rack Only	11 1/2 x 15 9/16 x 5 1/2 inch 292 x 395 x 140 mm	5 3/8 x 3 1/16 x 5 7/16 inch 136 x 78 x 138 mm	2 x 5 =10	DP-1116	5.2(2.38)
UF-323	Rack Only	16 7/8 x 6 5/8 x 5 1/2 inch 429 x 168 x 140 mm	5 3/8 x 3 1/4 x 5 7/16 inch 136 x 82.5 x 138 mm	3 x 2 =6	DP-1607	3.9(1.76)
UF-333	Rack Only	16 7/8 x 9 3/8 x 5 1/2 inch 429 x 238 x 140 mm	5 3/8 x 3 1/16 x 5 7/16 inch 136 x 78 x 138 mm	3 x 3 =9	DP-1610	5.0(2.28)
UF-343	Rack Only	16 7/8 x 12 2/8 x 5 1/2 inch 429 x 321 x 140 mm	5 3/8 x 3 3/8 x 5 7/16 inch 136 x 79.2 x 138 mm	3 x 4 =12	DP-1613	6.3(2.84)
UF-353	Rack Only	16 7/8 x 15 9/16 x 5 1/2 inch 429 x 395 x 140 mm	5 3/8 x 3 1/16 x 5 7/16 inch 136 x 78 x 138 mm	3 x 5 =15	DP-1616	7.4(3.37)
UF-423	Rack Only	22 3/16 x 6 5/8 x 5 1/2 inch 564 x 168 x 140 mm	5 3/8 x 3 1/4 x 5 7/16 inch 136 x 82.5 x 138 mm	4 x 2 =8	DP-2207	5.0(2.28)
UF-433	Rack Only	22 3/16 x 9 3/8 x 5 1/2 inch 564 x 238 x 140 mm	5 3/8 x 3 1/16 x 5 7/16 inch 136 x 78 x 138 mm	4 x 3 =12	DP-2210	6.5(2.95)
UF-433-W-XX*	Rack w/ 3" box and divider					6.5(2.95)+12ω ¹ ***
UF-443	Rack Only	22 3/16 x 12 2/8 x 5 1/2 inch 564 x 321 x 140 mm	5 3/8 x 3 3/8 x 5 7/16 inch 136 x 79.2 x 138 mm	4 x 4 =16	DP-2213	8.1(3.66)
UF-453	Rack Only	22 3/16 x 15 9/16 x 5 1/2 inch 564 x 395 x 140 mm	5 3/8 x 3 1/16 x 5 7/16 inch 136 x 78 x 138 mm	4 x 5 =20	DP-2116	9.7(4.40)
UF-523	Rack Only	27 5/16 x 6 5/8 x 5 1/2 inch 693 x 168 x 140 mm	5 1/4 x 3 1/4 x 5 7/16 inch 134 x 82.5 x 138 mm	5 x 2 =10	DP-2707	6.1(2.76)
UF-533	Rack Only	27 5/16 x 9 1/2 x 5 1/2 inch 693 x 241 x 140 mm	5 1/4 x 3 1/16 x 5 7/16 inch 134 x 78 x 138 mm	5 x 3 =15	DP-2710	7.8(3.56)
UF-543	Rack Only	27 5/16 x 12 2/8 x 5 1/2 inch 693 x 321 x 140 mm	5 1/4 x 3 3/8 x 5 7/16 inch 134 x 79.2 x 138 mm	5 x 4 =20	DP-2713	9.8(4.45)


* To order, replace XX with either 81 or 100 cell divider.

** Lock device is ordered separately.

*** ω¹=0.221b(0.1kg) for cardboard box, 0.371b(0.17kg) for plastic box, 0.841b(0.38kg) for stainless steel box, 0.531b(0.24kg) for aluminum box.


Upright Freezer Eco-Racks

For Standard 3" Boxes




[Check page 108 for Boxes]

- ◆ Designed economically for best cost performance.
- ◆ Constructed of corrosion resistant stainless steel.











Rack Only



Rack w/Cardboard Box

** UFE-223 / UFE-423 / UFE-523 (∅.82.5mm)

Part Number	Description	Dimension L x H x W	Compartment Size l x h x w	Capacity (Box)	Shipping Weight lb(Kg)
UFE-323	Rack Only	16 ¹³ / ₁₆ x 6 ⁵ / ₈ x 5 ¹ / ₂ inch 427 x 168 x 140 mm	5 ³ / ₈ x 3 ¹ / ₄ x 5 ⁷ / ₁₆ inch 137 x 82.5 x 138 mm	3 x 2 =6 	3.3(1.51)
UFE-333	Rack Only	16 ¹³ / ₁₆ x 9 ³ / ₈ x 5 ¹ / ₂ inch 427 x 238 x 140 mm	5 ³ / ₈ x 3 ¹ / ₈ x 5 ⁷ / ₁₆ inch 137 x 79 x 138 mm	3 x 3 =9 	3.9(1.78)
UFE-343	Rack Only	16 ¹³ / ₁₆ x 12 ⁵ / ₈ x 5 ¹ / ₂ inch 427 x 321 x 140 mm	5 ³ / ₈ x 3 ¹ / ₈ x 5 ⁷ / ₁₆ inch 137 x 79 x 138 mm	3 x 4 =12 	4.6(2.09)
UFE-423	Rack Only	22 ¹ / ₈ x 6 ³ / ₈ x 5 ¹ / ₂ inch 562 x 168 x 140 mm	5 ³ / ₈ x 3 ¹ / ₄ x 5 ⁷ / ₁₆ inch 137 x 82.5 x 138 mm	4 x 2 =8 	4.3(1.94)
UFE-433	Rack Only	22 ¹ / ₈ x 9 ³ / ₈ x 5 ¹ / ₂ inch 562 x 238 x 140 mm	5 ³ / ₈ x 3 ¹ / ₈ x 5 ⁷ / ₁₆ inch 137 x 79 x 138 mm	4 x 3 =12 	5.0(2.28)
UFE-443	Rack Only	22 ¹ / ₈ x 12 ⁵ / ₈ x 5 ¹ / ₂ inch 562 x 321 x 140 mm	5 ³ / ₈ x 3 ¹ / ₈ x 5 ⁷ / ₁₆ inch 137 x 79 x 138 mm	4 x 4 =16 	5.9(2.69)
UFE-523	Rack Only	27 ³ / ₁₆ x 6 ³ / ₈ x 5 ¹ / ₂ inch 691 x 168 x 140 mm	5 ⁵ / ₁₆ x 3 ¹ / ₄ x 5 ⁷ / ₁₆ inch 134.8 x 82.5 x 138 mm	5 x 2 =10 	5.2(2.36)
UFE-533	Rack Only	27 ³ / ₁₆ x 9 ¹ / ₂ x 5 ¹ / ₂ inch 691 x 241 x 140 mm	5 ⁵ / ₁₆ x 3 ¹ / ₈ x 5 ⁷ / ₁₆ inch 134.8 x 79 x 138 mm	5 x 3 =15 	6.1(2.77)