



# Upright Freezer Racks

## For 100-Cell Hinged Top Plastic Storage Boxes




[ Check page 110 for Boxes ]


- ◆ Made of corrosion resistant stainless steel.
- ◆ Security lock devices are available for each rack.



Rack Only



Rack w/100-Cell Hinged Plastic Storage Box




Rack w/DP Lock Device Option

Part Number	Description	Dimension L x H x W	Compartment Size l x h x w	Capacity (Box)	Lock Device Option*	Shipping Weight lb(Kg)
UFHT-23	Rack Only	13 x 6 <sup>15</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>8</sub> inch 330 x 176 x 149 mm	6 <sup>1</sup> / <sub>8</sub> x 2 <sup>5</sup> / <sub>16</sub> x 5 <sup>13</sup> / <sub>16</sub> inch 156 x 57.6 x 147.5 mm	2 x 3 =6	DP-1307	3.8(1.71)
UFHT-24	Rack Only	13 x 9 <sup>1</sup> / <sub>4</sub> x 5 <sup>7</sup> / <sub>8</sub> inch 330 x 235 x 149 mm	6 <sup>1</sup> / <sub>8</sub> x 2 <sup>5</sup> / <sub>16</sub> x 5 <sup>13</sup> / <sub>16</sub> inch 156 x 57.6 x 147.5 mm	2 x 4 =8	DP-1309	4.7(2.12)
UFHT-25	Rack Only	13 x 11 <sup>9</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>8</sub> inch 330 x 293 x 149 mm	6 <sup>1</sup> / <sub>8</sub> x 2 <sup>5</sup> / <sub>16</sub> x 5 <sup>13</sup> / <sub>16</sub> inch 156 x 57.6 x 147.5 mm	2 x 5 =10	DP-1312	5.6(2.52)
UFHT-26	Rack Only	13 x 13 <sup>13</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>8</sub> inch 330 x 351 x 149 mm	6 <sup>1</sup> / <sub>8</sub> x 2 <sup>5</sup> / <sub>16</sub> x 5 <sup>13</sup> / <sub>16</sub> inch 156 x 57.6 x 147.5 mm	2 x 6 =12	DP-1313	6.5(2.95)
UFHT-33	Rack Only	19 <sup>3</sup> / <sub>16</sub> x 6 <sup>15</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>8</sub> inch 488 x 176 x 149 mm	6 <sup>1</sup> / <sub>8</sub> x 2 <sup>5</sup> / <sub>16</sub> x 5 <sup>13</sup> / <sub>16</sub> inch 156 x 57.6 x 147.5 mm	3 x 3 =9	DP-1907	5.4(2.45)
UFHT-34	Rack Only	19 <sup>3</sup> / <sub>16</sub> x 9 <sup>1</sup> / <sub>4</sub> x 5 <sup>7</sup> / <sub>8</sub> inch 488 x 235 x 149 mm	6 <sup>1</sup> / <sub>8</sub> x 2 <sup>5</sup> / <sub>16</sub> x 5 <sup>13</sup> / <sub>16</sub> inch 156 x 57.6 x 147.5 mm	3 x 4 =12	DP-1909	6.7(3.02)
UFHT-35	Rack Only	19 <sup>3</sup> / <sub>16</sub> x 11 <sup>9</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>8</sub> inch 488 x 293 x 149 mm	6 <sup>1</sup> / <sub>8</sub> x 2 <sup>5</sup> / <sub>16</sub> x 5 <sup>13</sup> / <sub>16</sub> inch 156 x 57.6 x 147.5 mm	3 x 5 =15	DP-1912	8.0(3.63)
UFHT-36	Rack Only	19 <sup>3</sup> / <sub>16</sub> x 13 <sup>13</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>8</sub> inch 488 x 351 x 149 mm	6 <sup>1</sup> / <sub>8</sub> x 2 <sup>5</sup> / <sub>16</sub> x 5 <sup>13</sup> / <sub>16</sub> inch 156 x 57.6 x 147.5 mm	3 x 6 =18	DP-1913	9.3(4.24)
UFHT-43	Rack Only	25 <sup>7</sup> / <sub>16</sub> x 6 <sup>15</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>8</sub> inch 646 x 177 x 149 mm	6 <sup>1</sup> / <sub>8</sub> x 2 <sup>5</sup> / <sub>16</sub> x 5 <sup>13</sup> / <sub>16</sub> inch 156 x 57.6 x 147.5 mm	4 x 3 =12	DP-2407	7.0(3.18)
UFHT-44	Rack Only	25 <sup>7</sup> / <sub>16</sub> x 9 <sup>1</sup> / <sub>4</sub> x 5 <sup>7</sup> / <sub>8</sub> inch 646 x 235 x 149 mm	6 <sup>1</sup> / <sub>8</sub> x 2 <sup>5</sup> / <sub>16</sub> x 5 <sup>13</sup> / <sub>16</sub> inch 156 x 57.6 x 147.5 mm	4 x 4 =16	DP-2409	8.7(3.93)
UFHT-45	Rack Only	25 <sup>7</sup> / <sub>16</sub> x 11 <sup>9</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>8</sub> inch 646 x 293 x 149 mm	6 <sup>1</sup> / <sub>8</sub> x 2 <sup>5</sup> / <sub>16</sub> x 5 <sup>13</sup> / <sub>16</sub> inch 156 x 57.6 x 147.5 mm	4 x 5 =20	DP-2412	10.5(4.75)
UFHT-46	Rack Only	25 <sup>7</sup> / <sub>16</sub> x 13 <sup>13</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>8</sub> inch 646 x 351 x 149 mm	6 <sup>1</sup> / <sub>8</sub> x 2 <sup>5</sup> / <sub>16</sub> x 5 <sup>13</sup> / <sub>16</sub> inch 156 x 57.6 x 147.5 mm	4 x 6 =24	DP-2413	12.2(5.54)

\*Lock device is ordered separately.


# Upright Freezer Drawer Racks

## For 100 Cell Hinged Lid Storage Boxes




[ Check page 110 for Boxes ]

- ◆ Made of corrosion resistant stainless steel.
- ◆ Drawer racks allow easy access to and storage of samples without having to remove the whole rack from the freezer.
- ◆ Security lock devices are available for each rack.



Rack Only



Rack w/100-Cell Hinged Plastic Storage Box



Rack w/DDHT Lock Device Option

Part Number	Description	Dimension L x H x W	Compartment Size l x h x w	Capacity (Box)	Lock Device Option*	Shipping Weight lb(Kg)
UFDHT-23	Rack Only	11 <sup>3</sup> / <sub>4</sub> x 7 <sup>5</sup> / <sub>16</sub> x 6 <sup>7</sup> / <sub>16</sub> inch 299 x 186 x 164 mm	5 <sup>5</sup> / <sub>8</sub> x 2 <sup>3</sup> / <sub>16</sub> x 6 <sup>1</sup> / <sub>8</sub> inch 142.8 x 56 x 156 mm	6	DDHT-06	5.9(2.67)
UFDHT-24	Rack Only	11 <sup>3</sup> / <sub>4</sub> x 9 <sup>3</sup> / <sub>4</sub> x 6 <sup>7</sup> / <sub>16</sub> inch 299 x 247 x 164 mm	5 <sup>5</sup> / <sub>8</sub> x 2 <sup>3</sup> / <sub>16</sub> x 6 <sup>1</sup> / <sub>8</sub> inch 142.8 x 56 x 156 mm	8	DDHT-09	7.6(3.46)
UFDHT-25	Rack Only	11 <sup>3</sup> / <sub>4</sub> x 12 <sup>1</sup> / <sub>8</sub> x 6 <sup>7</sup> / <sub>16</sub> inch 299 x 308 x 164 mm	5 <sup>5</sup> / <sub>8</sub> x 2 <sup>3</sup> / <sub>16</sub> x 6 <sup>1</sup> / <sub>8</sub> inch 142.8 x 56 x 156 mm	10	DDHT-11	9.5(4.29)
UFDHT-26	Rack Only	11 <sup>3</sup> / <sub>4</sub> x 14 <sup>1</sup> / <sub>2</sub> x 6 <sup>7</sup> / <sub>16</sub> inch 299 x 369 x 164 mm	5 <sup>5</sup> / <sub>8</sub> x 2 <sup>3</sup> / <sub>16</sub> x 6 <sup>1</sup> / <sub>8</sub> inch 142.8 x 56 x 156 mm	12	DDHT-14	11.1(5.04)
UFDHT-33	Rack Only	17 <sup>3</sup> / <sub>8</sub> x 7 <sup>5</sup> / <sub>16</sub> x 6 <sup>7</sup> / <sub>16</sub> inch 442 x 186 x 164 mm	5 <sup>5</sup> / <sub>8</sub> x 2 <sup>3</sup> / <sub>16</sub> x 6 <sup>1</sup> / <sub>8</sub> inch 142.8 x 56 x 156 mm	9	DDHT-06	8.0(3.62)
UFDHT-34	Rack Only	17 <sup>3</sup> / <sub>8</sub> x 9 <sup>3</sup> / <sub>4</sub> x 6 <sup>7</sup> / <sub>16</sub> inch 442 x 247 x 164 mm	5 <sup>5</sup> / <sub>8</sub> x 2 <sup>3</sup> / <sub>16</sub> x 6 <sup>1</sup> / <sub>8</sub> inch 142.8 x 56 x 156 mm	12	DDHT-09	10.3(4.66)
UFDHT-35	Rack Only	17 <sup>3</sup> / <sub>8</sub> x 12 <sup>1</sup> / <sub>8</sub> x 6 <sup>7</sup> / <sub>16</sub> inch 442 x 308 x 164 mm	5 <sup>5</sup> / <sub>8</sub> x 2 <sup>3</sup> / <sub>16</sub> x 6 <sup>1</sup> / <sub>8</sub> inch 142.8 x 56 x 156 mm	15	DDHT-11	12.6(5.72)
UFDHT-36	Rack Only	17 <sup>3</sup> / <sub>8</sub> x 14 <sup>1</sup> / <sub>2</sub> x 6 <sup>7</sup> / <sub>16</sub> inch 442 x 369 x 164 mm	5 <sup>5</sup> / <sub>8</sub> x 2 <sup>3</sup> / <sub>16</sub> x 6 <sup>1</sup> / <sub>8</sub> inch 142.8 x 56 x 156 mm	18	DDHT-14	14.9(6.77)
UFDHT-43	Rack Only	23 <sup>1</sup> / <sub>16</sub> x 7 <sup>5</sup> / <sub>16</sub> x 6 <sup>7</sup> / <sub>16</sub> inch 585 x 186 x 164 mm	5 <sup>5</sup> / <sub>8</sub> x 2 <sup>3</sup> / <sub>16</sub> x 6 <sup>1</sup> / <sub>8</sub> inch 142.8 x 56 x 156 mm	12	DDHT-06	10.2(4.62)
UFDHT-44	Rack Only	23 <sup>1</sup> / <sub>16</sub> x 9 <sup>3</sup> / <sub>4</sub> x 6 <sup>7</sup> / <sub>16</sub> inch 585 x 247 x 164 mm	5 <sup>5</sup> / <sub>8</sub> x 2 <sup>3</sup> / <sub>16</sub> x 6 <sup>1</sup> / <sub>8</sub> inch 142.8 x 56 x 156 mm	16	DDHT-09	13.1(5.95)
UFDHT-45	Rack Only	23 <sup>1</sup> / <sub>16</sub> x 12 <sup>1</sup> / <sub>8</sub> x 6 <sup>7</sup> / <sub>16</sub> inch 585 x 308 x 164 mm	5 <sup>5</sup> / <sub>8</sub> x 2 <sup>3</sup> / <sub>16</sub> x 6 <sup>1</sup> / <sub>8</sub> inch 142.8 x 56 x 156 mm	20	DDHT-11	16.1(7.30)
UFDHT-46	Rack Only	23 <sup>1</sup> / <sub>16</sub> x 14 <sup>1</sup> / <sub>2</sub> x 6 <sup>7</sup> / <sub>16</sub> inch 585 x 369 x 164 mm	5 <sup>5</sup> / <sub>8</sub> x 2 <sup>3</sup> / <sub>16</sub> x 6 <sup>1</sup> / <sub>8</sub> inch 142.8 x 56 x 156 mm	24	DDHT-14	18.5(8.64)

\* Lock device is ordered separately.