



Upright Freezer Racks

For 96-Well Microtube Boxes



- ◆ Made of corrosion resistant stainless steel.
- ◆ Locking rods are included.
- ◆ Security lock devices are available for each rack.




Part Number	Dimension L x H x W	Compartment Size l x h x w	Capacity (Box)	Lock Device Option*	Shipping Weight(lb/Kg)
UFTB-33L	11 ⁹ / ₁₆ x 6 ¹³ / ₁₆ x 5 ¹ / ₂ inch 293 x 173 x 140mm	3 ⁹ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	3 x 3 = 9	DP-1107	3.7 / 1.70
UFTB-34L	11 ⁹ / ₁₆ x 9 x 5 ¹ / ₂ inch 293 x 228 x 140mm	3 ⁹ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	3 x 4 = 12	DP-1109	4.7 / 2.12
UFTB-35L	11 ⁹ / ₁₆ x 11 ³ / ₁₆ x 5 ¹ / ₂ inch 293 x 284 x 140mm	3 ⁹ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	3 x 5 = 15	DP-1111	5.6 / 2.53
UFTB-36L	11 ⁹ / ₁₆ x 13 ³ / ₈ x 5 ¹ / ₂ inch 293 x 339 x 140mm	3 ⁹ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	3 x 6 = 18	DP-1113	6.5 / 2.95
UFTB-43L	15 ³ / ₁₆ x 6 ¹³ / ₁₆ x 5 ¹ / ₂ inch 385 x 173 x 140mm	3 ⁹ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	4 x 3 = 12	DP-1407	4.9 / 2.20
UFTB-44L	15 ³ / ₁₆ x 9 x 5 ¹ / ₂ inch 385 x 228 x 140mm	3 ⁹ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	4 x 4 = 16	DP-1409	6.0 / 2.74
UFTB-45L	15 ³ / ₁₆ x 11 ³ / ₁₆ x 5 ¹ / ₂ inch 385 x 284 x 140mm	3 ⁹ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	4 x 5 = 20	DP-1411	7.2 / 3.27
UFTB-46L	15 ³ / ₁₆ x 13 ³ / ₈ x 5 ¹ / ₂ inch 385 x 339 x 140mm	3 ⁹ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	4 x 6 = 24	DP-1413	8.4 / 3.81
UFTB-53L	18 ¹³ / ₁₆ x 6 ¹³ / ₁₆ x 5 ¹ / ₂ inch 478 x 173 x 140mm	3 ⁹ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	5 x 3 = 15	DP-1807	5.9 / 2.69
UFTB-54L	18 ¹³ / ₁₆ x 9 x 5 ¹ / ₂ inch 478 x 228 x 140mm	3 ⁹ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	5 x 4 = 20	DP-1809	7.4 / 3.36
UFTB-55L	18 ¹³ / ₁₆ x 11 ³ / ₁₆ x 5 ¹ / ₂ inch 478 x 284 x 140mm	3 ⁹ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	5 x 5 = 25	DP-1811	8.9 / 4.02
UFTB-56L	18 ¹³ / ₁₆ x 13 ³ / ₈ x 5 ¹ / ₂ inch 478 x 339 x 140mm	3 ⁹ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	5 x 6 = 30	DP-1813	10.3 / 4.68
UFTB-63L	22 ¹ / ₂ x 6 ¹³ / ₁₆ x 5 ¹ / ₂ inch 571 x 173 x 140mm	3 ⁹ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	6 x 3 = 18	DP-2207	7.0 / 3.17
UFTB-64L	22 ¹ / ₂ x 9 x 5 ¹ / ₂ inch 571 x 228 x 140mm	3 ⁹ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	6 x 4 = 24	DP-2209	8.8 / 3.97
UFTB-65L	22 ¹ / ₂ x 11 ³ / ₁₆ x 5 ¹ / ₂ inch 571 x 284 x 140mm	3 ⁹ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	6 x 5 = 30	DP-2211	10.5 / 4.76
UFTB-66L	22 ¹ / ₂ x 13 ³ / ₈ x 5 ¹ / ₂ inch 571 x 339 x 140mm	3 ⁹ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	6 x 6 = 36	DP-2213	12.2 / 5.54
UFTB-73L	26 ¹ / ₈ x 6 ¹³ / ₁₆ x 5 ¹ / ₂ inch 664 x 173 x 140mm	3 ⁹ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	7 x 3 = 21	DP-2607	8.1 / 3.68
UFTB-74L	26 ¹ / ₈ x 9 x 5 ¹ / ₂ inch 664 x 228 x 140mm	3 ⁹ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	7 x 4 = 28	DP-2609	10.1 / 4.59
UFTB-75L	26 ¹ / ₈ x 11 ³ / ₁₆ x 5 ¹ / ₂ inch 664 x 284 x 140mm	3 ⁹ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	7 x 5 = 35	DP-2611	11.9 / 5.42
UFTB-76L	26 ¹ / ₈ x 13 ³ / ₈ x 5 ¹ / ₂ inch 664 x 339 x 140mm	3 ⁹ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₁₆ inch 91 x 55 x 130.7 mm	7 x 6 = 42	DP-2613	14.2 / 6.42


* Lock device is ordered separately.

Upright Freezer Racks

For 50-Cell 0.5ml Storage Boxes



- ◆ Made of corrosion resistant stainless steel.
- ◆ Each rack compartment holds two 50-cell 0.5ml Microtube Storage Boxes.
- ◆ Security lock devices are available for each rack.



w/DP Lock Device Option


*UFM-1605 / UFM-3605 / UFM-4605 (h:52.6mm)

Part Number	Dimension L x H x W	Compartment Size l x h x w	Capacity (Box)	Lock Device Option*	Shipping Weight(lb/Kg)
UFM-1305	6 x 6 ⁵ / ₈ x 5 ¹ / ₂ inch 153 x 169 x 140 mm	5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 137 x 55 x 138 mm	1 x 3 x 2 =6	DP-0507	1.8 / 0.83
UFM-1405	6 x 8 ¹³ / ₁₆ x 5 ¹ / ₂ inch 153 x 224 x 140 mm	5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 137 x 55 x 138 mm	1 x 4 x 2 =8	DP-0509	2.3 / 1.05
UFM-1505	6 x 11 x 5 ¹ / ₂ inch 153 x 280 x 140 mm	5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 137 x 55 x 138 mm	1 x 5 x 2 =10	DP-0511	2.6 / 1.20
UFM-1605	6 x 12 ⁵ / ₈ x 5 ¹ / ₂ inch 153 x 321 x 140 mm	5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 137 x 52.6 x 138 mm	1 x 6 x 2 =12	DP-0512	3.1 / 1.40
UFM-2305	11 ¹ / ₂ x 6 ⁵ / ₈ x 5 ¹ / ₂ inch 292 x 169 x 140 mm	5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 137 x 55 x 138 mm	2 x 3 x 2 =12	DP-1107	3.1 / 1.42
UFM-2405	11 ¹ / ₂ x 8 ¹³ / ₁₆ x 5 ¹ / ₂ inch 292 x 224 x 140 mm	5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 137 x 55 x 138 mm	2 x 4 x 2 =16	DP-1109	4.0 / 1.80
UFM-2505	11 ¹ / ₂ x 11 x 5 ¹ / ₂ inch 292 x 280 x 140 mm	5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 137 x 55 x 138 mm	2 x 5 x 2 =20	DP-1111	4.6 / 2.10
UFM-2605	11 ¹ / ₂ x 12 ⁵ / ₈ x 5 ¹ / ₂ inch 292 x 321 x 140 mm	5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 137 x 52.6 x 138 mm	2 x 6 x 2 =24	DP-1112	5.5 / 2.48
UFM-3305	16 ¹⁵ / ₁₆ x 6 ⁵ / ₈ x 5 ¹ / ₂ inch 430 x 169 x 140 mm	5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 137 x 55 x 138 mm	3 x 3 x 2 =18	DP-1707	4.5 / 2.03
UFM-3405	16 ¹⁵ / ₁₆ x 8 ¹³ / ₁₆ x 5 ¹ / ₂ inch 430 x 224 x 140 mm	5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 137 x 55 x 138 mm	3 x 4 x 2 =24	DP-1709	5.6 / 2.55
UFM-3505	16 ¹⁵ / ₁₆ x 11 x 5 ¹ / ₂ inch 430 x 280 x 140 mm	5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 137 x 55 x 138 mm	3 x 5 x 2 =30	DP-1711	6.6 / 2.98
UFM-3605	16 ¹⁵ / ₁₆ x 12 ⁵ / ₈ x 5 ¹ / ₂ inch 430 x 321 x 140 mm	5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 137 x 52.6 x 138 mm	3 x 6 x 2 =36	DP-1712	7.5 / 3.40
UFM-4305	22 ³ / ₈ x 6 ⁵ / ₈ x 5 ¹ / ₂ inch 569 x 169 x 140 mm	5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 137 x 55 x 138 mm	4 x 3 x 2 =24	DP-2207	5.8 / 2.64
UFM-4405	22 ³ / ₈ x 8 ¹³ / ₁₆ x 5 ¹ / ₂ inch 569 x 224 x 140 mm	5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 137 x 55 x 138 mm	4 x 4 x 2 =32	DP-2209	7.3 / 3.29
UFM-4505	22 ³ / ₈ x 11 x 5 ¹ / ₂ inch 569 x 280 x 140 mm	5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 137 x 55 x 138 mm	4 x 5 x 2 =40	DP-2211	8.6 / 3.88
UFM-4605	22 ³ / ₈ x 12 ⁵ / ₈ x 5 ¹ / ₂ inch 569 x 321 x 140 mm	5 ³ / ₈ x 2 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 137 x 52.6 x 138 mm	4 x 6 x 2 =48	DP-2212	9.9 / 4.50

*Security lock device is ordered separately.


Upright Freezer Racks

For 50-Cell 1.5ml Microtube Storage Boxes



A B C
3 x 5 1/2 x 2 inch
76.2 x 139.0 x 50.8mm

- ◆ Made of corrosion resistant stainless steel.
- ◆ Each rack compartment holds two 50-cell 0.5ml Microtube Storage Boxes.
- ◆ Security lock devices are available for each rack.




w/DP Lock Device Option

Part Number	Dimension L x H x W	Compartment Size l x h x w	Capacity (Box)	Lock Device Option*	Shipping Weight(lb/Kg)
UFM-1315	6 ³ / ₄ x 6 ¹⁵ / ₁₆ x 5 ⁷ / ₈ inch 172 x 176 x 149 mm	6 ¹ / ₈ x 2 ⁵ / ₁₆ x 5 ¹³ / ₁₆ inch 156 x 57.6 x 147 mm	1 x 3 x 2 =6	DP-0607	2.1 / 0.93
UFM-1415	6 ³ / ₄ x 9 ¹ / ₄ x 5 ⁷ / ₈ inch 172 x 235 x 149 mm	6 ¹ / ₈ x 2 ⁵ / ₁₆ x 5 ¹³ / ₁₆ inch 156 x 57.6 x 147 mm	1 x 4 x 2 =8	DP-0609	2.5 / 1.14
UFM-1515	6 ³ / ₄ x 11 ⁹ / ₁₆ x 5 ⁷ / ₈ inch 172 x 293 x 149 mm	6 ¹ / ₈ x 2 ⁵ / ₁₆ x 5 ¹³ / ₁₆ inch 156 x 57.6 x 147 mm	1 x 5 x 2 =10	DP-0612	3.0 / 1.36
UFM-2315	13 x 6 ¹⁵ / ₁₆ x 5 ⁷ / ₈ inch 330 x 176 x 149 mm	6 ¹ / ₈ x 2 ⁵ / ₁₆ x 5 ¹³ / ₁₆ inch 156 x 57.6 x 147 mm	2 x 3 x 2 =12	DP-1307	3.6 / 1.64
UFM-2415	13 x 9 ¹ / ₄ x 5 ⁷ / ₈ inch 330 x 235 x 149 mm	6 ¹ / ₈ x 2 ⁵ / ₁₆ x 5 ¹³ / ₁₆ inch 156 x 57.6 x 147 mm	2 x 4 x 2 =16	DP-1309	4.4 / 2.01
UFM-2515	13 x 11 ⁹ / ₁₆ x 5 ⁷ / ₈ inch 330 x 293 x 149 mm	6 ¹ / ₈ x 2 ⁵ / ₁₆ x 5 ¹³ / ₁₆ inch 156 x 57.6 x 147 mm	2 x 5 x 2 =20	DP-1312	5.3 / 2.39
UFM-2615	13 x 12 ⁵ / ₈ x 5 ⁷ / ₈ inch 330 x 321 x 149 mm	6 ¹ / ₈ x 2 ⁵ / ₁₆ x 5 ¹³ / ₁₆ inch 156 x 52.6 x 147 mm	2 x 6 x 2 =24	DP-1313	5.9 / 2.69
UFM-2715	13 x 14 ³ / ₄ x 5 ⁷ / ₈ inch 330 x 375 x 149 mm	6 ¹ / ₈ x 2 ⁵ / ₁₆ x 5 ¹³ / ₁₆ inch 156 x 52.6 x 147 mm	2 x 7 x 2 =28	DP-1314	6.7 / 3.05
UFM-3315	19 ³ / ₁₆ x 6 ¹⁵ / ₁₆ x 5 ⁷ / ₈ inch 488 x 176 x 149 mm	6 ¹ / ₈ x 2 ⁵ / ₁₆ x 5 ¹³ / ₁₆ inch 156 x 57.6 x 147 mm	3 x 3 x 2 =18	DP-1907	5.2 / 2.34
UFM-3415	19 ³ / ₁₆ x 9 ¹ / ₄ x 5 ⁷ / ₈ inch 488 x 235 x 149 mm	6 ¹ / ₈ x 2 ⁵ / ₁₆ x 5 ¹³ / ₁₆ inch 156 x 57.6 x 147 mm	3 x 4 x 2 =24	DP-1909	6.3 / 2.88
UFM-3515	19 ³ / ₁₆ x 11 ⁹ / ₁₆ x 5 ⁷ / ₈ inch 488 x 293 x 149 mm	6 ¹ / ₈ x 2 ⁵ / ₁₆ x 5 ¹³ / ₁₆ inch 156 x 57.6 x 147 mm	3 x 5 x 2 =30	DP-1912	7.6 / 3.43
UFM-3615	19 ³ / ₁₆ x 12 ⁵ / ₈ x 5 ⁷ / ₈ inch 488 x 321 x 149 mm	6 ¹ / ₈ x 2 ⁵ / ₁₆ x 5 ¹³ / ₁₆ inch 156 x 52.6 x 147 mm	3 x 6 x 2 =36	DP-1913	8.5 / 3.85
UFM-3715	19 ³ / ₁₆ x 14 ³ / ₄ x 5 ⁷ / ₈ inch 488 x 375 x 149 mm	6 ¹ / ₈ x 2 ⁵ / ₁₆ x 5 ¹³ / ₁₆ inch 156 x 52.6 x 147 mm	3 x 7 x 2 =42	DP-1914	9.7 / 4.38


*Security lock device is ordered separately.

Upright Freezer Drawer Racks

For 96-Well Microtube Boxes



- ◆ Made of corrosion resistant stainless steel.
- ◆ Drawer racks allow easy access to fewer microtiter plates at a time thus protecting other samples in the rack from warming up.
- ◆ Security lock devices are available for each rack.



Part Number	Description	Dimension L x H x W	Compartment Size l x h x w	Capacity Microtube Boxes	Lock Device Option*	Shipping Weight (lb/Kg)
UFDMP-33	Rack w/3 drawers (3 compartments per drawer) for Microtube boxes	11 ¹ / ₈ x 7 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 282 x 186 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	9	DD-06	5.4 / 2.45
UFDMP-34	Rack w/4 drawers (3 compartments per drawer) for Microtube boxes	11 ¹ / ₈ x 9 ³ / ₄ x 5 ⁷ / ₁₆ inch 282 x 247 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	12	DD-09	6.9 / 3.13
UFDMP-35	Rack w/5 drawers (3 compartments per drawer) for Microtube boxes	11 ¹ / ₈ x 12 ¹ / ₈ x 5 ⁷ / ₁₆ inch 282 x 308 x 138mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	15	DD-11	8.6 / 3.90
UFDMP-36	Rack w/6 drawers (3 compartments per drawer) for Microtube boxes	11 ¹ / ₈ x 14 ¹ / ₂ x 5 ⁷ / ₁₆ inch 282 x 369 x 138mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	18	DD-14	10.2 / 4.64
UFDMP-43	Rack w/3 drawers (4 compartments per drawer) for Microtube boxes	14 ³ / ₈ x 7 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 372 x 186 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	12	DD-06	6.8 / 3.10
UFDMP-44	Rack w/4 drawers (4 compartments per drawer) for Microtube boxes	14 ³ / ₈ x 9 ³ / ₄ x 5 ⁷ / ₁₆ inch 372 x 247 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	16	DD-09	8.9 / 4.02
UFDMP-45	Rack w/5 drawers (4 compartments per drawer) for Microtube boxes	14 ³ / ₈ x 12 ¹ / ₈ x 5 ⁷ / ₁₆ inch 372 x 308 x 138mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	20	DD-11	10.9 / 4.94
UFDMP-46	Rack w/6 drawers (4 compartments per drawer) for Microtube boxes	14 ³ / ₈ x 14 ¹ / ₂ x 5 ⁷ / ₁₆ inch 372 x 369 x 138mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	24	DD-14	10.7 / 5.86
UFDMP-53	Rack w/3 drawers (5 compartments per drawer) for Microtube boxes	18 ³ / ₁₆ x 7 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 462 x 186 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	15	DD-06	8.2 / 3.70
UFDMP-54	Rack w/4 drawers (5 compartments per drawer) for Microtube boxes	18 ³ / ₁₆ x 9 ³ / ₄ x 5 ⁷ / ₁₆ inch 462 x 247 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	20	DD-09	10.6 / 4.79
UFDMP-55	Rack w/5 drawers (5 compartments per drawer) for Microtube boxes	18 ³ / ₁₆ x 12 ¹ / ₈ x 5 ⁷ / ₁₆ inch 462 x 308 x 138mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	25	DD-11	13.0 / 5.90
UFDMP-56	Rack w/6 drawers (5 compartments per drawer) for Microtube boxes	18 ³ / ₁₆ x 14 ¹ / ₂ x 5 ⁷ / ₁₆ inch 462 x 369 x 138mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	30	DD-14	15.4 / 7.00
UFDMP-63	Rack w/3 drawers (6 compartments per drawer) for Microtube boxes	21 ³ / ₄ x 7 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 552 x 186 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	18	DD-06	9.7 / 4.42
UFDMP-64	Rack w/4 drawers (6 compartments per drawer) for Microtube boxes	21 ³ / ₄ x 9 ³ / ₄ x 5 ⁷ / ₁₆ inch 552 x 247 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	24	DD-09	12.6 / 5.73
UFDMP-65	Rack w/5 drawers (6 compartments per drawer) for Microtube boxes	21 ³ / ₄ x 12 ¹ / ₈ x 5 ⁷ / ₁₆ inch 552 x 308 x 138mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	30	DD-11	15.5 / 7.05
UFDMP-66	Rack w/6 drawers (6 compartments per drawer) for Microtube boxes	21 ³ / ₄ x 14 ¹ / ₂ x 5 ⁷ / ₁₆ inch 552 x 369 x 138mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	36	DD-14	17.9 / 8.14
UFDMP-73	Rack w/3 drawers (7 compartments per drawer) for Microtube boxes	25 ¹ / ₄ x 7 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 642 x 186 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	21	DD-06	11.0 / 4.97
UFDMP-74	Rack w/4 drawers (7 compartments per drawer) for Microtube boxes	25 ¹ / ₄ x 9 ³ / ₄ x 5 ⁷ / ₁₆ inch 642 x 247 x 138 mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	28	DD-09	14.2 / 6.43
UFDMP-75	Rack w/5 drawers (7 compartments per drawer) for Microtube boxes	25 ¹ / ₄ x 12 ¹ / ₈ x 5 ⁷ / ₁₆ inch 642 x 308 x 138mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	35	DD-11	17.4 / 7.91
UFDMP-76	Rack w/6 drawers (7 compartments per drawer) for Microtube boxes	25 ¹ / ₄ x 14 ¹ / ₂ x 5 ⁷ / ₁₆ inch 642 x 369 x 138mm	3 ¹ / ₂ x 2 ³ / ₁₆ x 5 ¹ / ₈ inch 89 x 55.8 x 130 mm	42	DD-14	20.7 / 9.38

*Lock device is ordered separately.