


Upright Freezer Racks

Modifiable Racks for 2" and 3" Standard Boxes


- ◆ Removable clip can be adjusted so that the rack can hold both standard 2" and 3" boxes.
- ◆ Constructed of corrosion resistant stainless steel.



Part Number	Dimension L x H x W	Compartment Size l x h x w(2" boxes)	Compartment Size l x h x w(3" boxes)	Capacity (Boxes)			Shipping Weight(lb/Kg)
				2"	3"	2"and 3"	
UFMC-33-2-3	16 ¹⁵ / ₁₆ x 6 ⁵ / ₈ x 5 ¹ / ₂ inch 430.2 x 168.8 x 140 mm	5 ⁷ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₈ inch 138 x 54.6 x 137 mm	5 ⁷ / ₁₆ x 3 ¹ / ₁₆ x 5 ³ / ₈ inch 138 x 77.7 x 137 mm	9	6	3 + 3	3.8 / 1.73
UFMC-34-2-3	16 ¹⁵ / ₁₆ x 8 ⁷ / ₈ x 5 ¹ / ₂ inch 430.2 x 225 x 140 mm	5 ⁷ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₈ inch 138 x 54.6 x 137 mm	5 ⁷ / ₁₆ x 3 ¹ / ₁₆ x 5 ³ / ₈ inch 138 x 77.7 x 137 mm	12	6	6 + 3	6.7 / 3.04
UFMC-35-2-3	16 ¹⁵ / ₁₆ x 11 x 5 ¹ / ₂ inch 430.2 x 280 x 140 mm	5 ⁷ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₈ inch 138 x 54.6 x 137 mm	5 ⁷ / ₁₆ x 3 ¹ / ₁₆ x 5 ³ / ₈ inch 138 x 77.7 x 137 mm	15	6	9 + 3	5.2 / 2.36
UFMC-43-2-3	22 ³ / ₈ x 6 ⁵ / ₈ x 5 ¹ / ₂ inch 569 x 168.8 x 140 mm	5 ⁷ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₈ inch 138 x 54.6 x 137 mm	5 ⁷ / ₁₆ x 3 ¹ / ₁₆ x 5 ³ / ₈ inch 138 x 77.7 x 137 mm	12	8	4 + 4	4.9 / 2.24
UFMC-44-2-3	22 ³ / ₈ x 8 ⁷ / ₈ x 5 ¹ / ₂ inch 569 x 225 x 140 mm	5 ⁷ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₈ inch 138 x 54.6 x 137 mm	5 ⁷ / ₁₆ x 3 ¹ / ₁₆ x 5 ³ / ₈ inch 138 x 77.7 x 137 mm	16	8	8 + 4	5.8 / 2.63
UFMC-45-2-3	22 ³ / ₈ x 11 x 5 ¹ / ₂ inch 569 x 280 x 140 mm	5 ⁷ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₈ inch 138 x 54.6 x 137 mm	5 ⁷ / ₁₆ x 3 ¹ / ₁₆ x 5 ³ / ₈ inch 138 x 77.7 x 137 mm	20	8	12 + 4	6.7 / 3.02
UFMC-46-2-3	22 ³ / ₈ x 12 ⁵ / ₈ x 5 ¹ / ₂ inch 569 x 321 x 140 mm	5 ⁷ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₈ inch 138 x 54.6 x 137 mm	5 ⁷ / ₁₆ x 3 ¹ / ₁₆ x 5 ³ / ₈ inch 138 x 77.7 x 137 mm	24	12	12 + 12	5.8 / 2.62
UFMC-53-2-3	27 ⁵ / ₁₆ x 6 ⁵ / ₈ x 5 ¹ / ₂ inch 693 x 168.8 x 140 mm	5 ⁵ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₈ inch 135 x 54.6 x 137 mm	5 ⁵ / ₁₆ x 3 ¹ / ₁₆ x 5 ³ / ₈ inch 135 x 77.7 x 137 mm	15	10	5 + 5	5.8 / 2.62
UFMC-54-2-3	27 ⁵ / ₁₆ x 8 ⁷ / ₈ x 5 ¹ / ₂ inch 693 x 225 x 140 mm	5 ⁵ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₈ inch 135 x 54.6 x 137 mm	5 ⁵ / ₁₆ x 3 ¹ / ₁₆ x 5 ³ / ₈ inch 135 x 77.7 x 137 mm	20	10	10 + 5	7.1 / 3.20
UFMC-55-2-3	27 ⁵ / ₁₆ x 11 x 5 ¹ / ₂ inch 693 x 280 x 140 mm	5 ⁵ / ₁₆ x 2 ³ / ₁₆ x 5 ³ / ₈ inch 135 x 54.6 x 137 mm	5 ⁵ / ₁₆ x 3 ¹ / ₁₆ x 5 ³ / ₈ inch 135 x 77.7 x 137 mm	25	10	15 + 5	8.1 / 3.68

For Mixed Storage of 2", 3" and 3.75" Standard Boxes

- ◆ Constructed of corrosion resistant stainless steel.



Part Number	Dimension L x H x W	Compartment Size l x h x w	Capacity (Boxes)			Shipping Weight(lb/Kg)
			2"	3"	2"and 3"	
UF-306	16 ¹³ / ₁₆ x 6 ¹ / ₂ x 5 ¹ / ₂ inch 427 x 165 x 140 mm	5 ³ / ₈ x 6 ⁵ / ₁₆ x 5 ⁷ / ₁₆ inch 136.5 x 160 x 138 mm	3x3 =9	3x2 =6	3x1 + 3x1	3.3 / 1.49
UF-308	16 ¹³ / ₁₆ x 8 ⁹ / ₁₆ x 5 ¹ / ₂ inch 427 x 218 x 140 mm	5 ³ / ₈ x 8 ³ / ₈ x 5 ⁷ / ₁₆ inch 136.5 x 213 x 138 mm	3x4 =12	3x2 =6	3x1 + 3x2	3.8 / 1.71
UF-312L	16 ¹³ / ₁₆ x 12 ³ / ₈ x 5 ⁵ / ₈ inch 427 x 315 x 143 mm	5 ³ / ₈ x 12 ³ / ₁₆ x 5 ⁵ / ₁₆ inch 136.5 x 310 x 135 mm	3x6 =16	3x4 =12	3x3 + 3x2	5.1 / 2.31
UF-406	22 ¹ / ₈ x 6 ¹ / ₂ x 5 ¹ / ₂ inch 562 x 165 x 140 mm	5 ³ / ₈ x 6 ⁵ / ₁₆ x 5 ⁷ / ₁₆ inch 136.5 x 160 x 138 mm	4x3 =12	4x2 =8	4x1 + 4x1	4.2 / 1.92
UF-408	22 ¹ / ₈ x 8 ⁹ / ₁₆ x 5 ¹ / ₂ inch 562 x 218 x 140 mm	5 ³ / ₈ x 8 ³ / ₈ x 5 ⁷ / ₁₆ inch 136.5 x 213 x 138 mm	4x4 =16	4x2 =8	4x1 + 4x2	4.8 / 2.18
UF-412	22 ¹ / ₈ x 9 ³ / ₈ x 5 ¹ / ₂ inch 562 x 238 x 140 mm	5 ³ / ₈ x 9 ³ / ₁₆ x 5 ⁷ / ₁₆ inch 136.5 x 233 x 138 mm	4x4 =16	4x3 =12	4x3 + 4x1	5.1 / 2.3
UF-412L	22 ¹ / ₈ x 12 ³ / ₈ x 5 ⁵ / ₈ inch 562 x 315 x 143 mm	5 ³ / ₈ x 12 ³ / ₁₆ x 5 ⁵ / ₁₆ inch 136.5 x 310 x 135 mm	4x6 =24	4x3 =12	4x3 + 4x2	6.3 / 2.87
UF-506	27 ⁵ / ₁₆ x 6 ¹ / ₂ x 5 ¹ / ₂ inch 691 x 165 x 140 mm	5 ⁵ / ₁₆ x 6 ⁵ / ₁₆ x 5 ⁷ / ₁₆ inch 134.7 x 160 x 138 mm	5x3 =15	5x2 =10	5x1 + 5x1	5.1 / 2.33
UF-508	27 ⁵ / ₁₆ x 8 ⁹ / ₁₆ x 5 ¹ / ₂ inch 691 x 218 x 140 mm	5 ⁵ / ₁₆ x 8 ³ / ₈ x 5 ⁷ / ₁₆ inch 134.7 x 213 x 138 mm	5x4 =20	5x2 =10	5x1 + 5x2	5.9 / 2.66
UF-512L	27 ⁵ / ₁₆ x 12 ³ / ₈ x 5 ⁵ / ₈ inch 691 x 315 x 143 mm	5 ⁵ / ₁₆ x 12 ³ / ₁₆ x 5 ⁵ / ₁₆ inch 134.7 x 310 x 135 mm	5x6 =30	5x3 =15	5x3 + 5x2	7.7 / 3.50