



# Chest Freezer Racks

(Locking Rod Included)



## For 96-Well & 384-Well Microtiter Plates

- ◆ Made of corrosion resistant stainless steel.
- ◆ Locking rods are included.
- ◆ Security lock devices are available for each rack.



Chest Freezer Rack



Rack w/DR Lock Device Option


Part Number	Dimension H x D x W	Compartment Size h x d x w	Capacity (Plate)	Lock Device Option*	Shipping Weight (lb/Kg)
CFMP-5	16 <sup>13</sup> / <sub>16</sub> x 3 <sup>5</sup> / <sub>8</sub> x 5 <sup>1</sup> / <sub>2</sub> inch 427 x 92 x 140 mm	3 <sup>5</sup> / <sub>16</sub> x 3 <sup>7</sup> / <sub>16</sub> x 5 <sup>3</sup> / <sub>16</sub> inch 83.5 x 88 x 131 mm	25 Plates w/lids	DR-15	2.0 / 0.91
			30 Plates w/o lids		
CFMP-6	20 <sup>1</sup> / <sub>8</sub> x 3 <sup>5</sup> / <sub>8</sub> x 5 <sup>1</sup> / <sub>2</sub> inch 511 x 92 x 140 mm	3 <sup>5</sup> / <sub>16</sub> x 3 <sup>7</sup> / <sub>16</sub> x 5 <sup>3</sup> / <sub>16</sub> inch 83.5 x 88 x 131 mm	30 Plates w/lids	DR-19	2.3 / 1.06
			36 Plates w/o lids		
CFMP-7	23 <sup>7</sup> / <sub>16</sub> x 3 <sup>5</sup> / <sub>8</sub> x 5 <sup>1</sup> / <sub>2</sub> inch 595 x 92 x 140 mm	3 <sup>5</sup> / <sub>16</sub> x 3 <sup>7</sup> / <sub>16</sub> x 5 <sup>3</sup> / <sub>16</sub> inch 83.5 x 88 x 131 mm	35 Plates w/lids	DR-22	2.7 / 1.21
			42 Plates w/o lids		
CFMP-8	26 <sup>3</sup> / <sub>4</sub> x 3 <sup>5</sup> / <sub>8</sub> x 5 <sup>1</sup> / <sub>2</sub> inch 679 x 92 x 140 mm	3 <sup>5</sup> / <sub>16</sub> x 3 <sup>7</sup> / <sub>16</sub> x 5 <sup>3</sup> / <sub>16</sub> inch 83.5 x 88 x 131 mm	40 Plates w/lids	DR-25	3.0 / 1.37
			48 Plates w/o lids		
CFMP-9	30 <sup>1</sup> / <sub>16</sub> x 3 <sup>5</sup> / <sub>8</sub> x 5 <sup>1</sup> / <sub>2</sub> inch 763 x 92 x 140 mm	3 <sup>5</sup> / <sub>16</sub> x 3 <sup>7</sup> / <sub>16</sub> x 5 <sup>3</sup> / <sub>16</sub> inch 83.5 x 88 x 131 mm	45 Plates w/lids	DR-29	3.4 / 1.52
			54 Plates w/o lids		

\* Lock device is ordered separately.


# Chest Freezer Racks

(Locking Rod Included)

## For 96 Deep-Well Microtiter Plates



- ◆ Made of corrosion resistant stainless steel.
- ◆ Security lock devices are available for each rack.



Rack w/DR Lock Device Option

Part Number	Dimension H x D x W	Compartment Size h x d x w	Capacity	Lock Device Option*	Shipping Weight(lb/Kg)
CFDP-9	16 <sup>5</sup> / <sub>8</sub> x 3 <sup>5</sup> / <sub>8</sub> x 5 <sup>1</sup> / <sub>2</sub> inch 423 x 92 x 140 mm	1 <sup>13</sup> / <sub>16</sub> x 3 <sup>7</sup> / <sub>16</sub> x 5 <sup>3</sup> / <sub>16</sub> inch 45.5 x 88 x 131.5 mm	9	DR-15	3.1 / 1.39
CFDP-10	18 <sup>7</sup> / <sub>16</sub> x 3 <sup>5</sup> / <sub>8</sub> x 5 <sup>1</sup> / <sub>2</sub> inch 469 x 92 x 140 mm	1 <sup>13</sup> / <sub>16</sub> x 3 <sup>7</sup> / <sub>16</sub> x 5 <sup>3</sup> / <sub>16</sub> inch 45.5 x 88 x 131.5 mm	10	DR-17	3.3 / 1.48
CFDP-11	20 <sup>1</sup> / <sub>4</sub> x 3 <sup>5</sup> / <sub>8</sub> x 5 <sup>1</sup> / <sub>2</sub> inch 515 x 92 x 140 mm	1 <sup>13</sup> / <sub>16</sub> x 3 <sup>7</sup> / <sub>16</sub> x 5 <sup>3</sup> / <sub>16</sub> inch 45.5 x 88 x 131.5 mm	11	DR-19	3.5 / 1.61
CFDP-12	22 <sup>1</sup> / <sub>8</sub> x 3 <sup>5</sup> / <sub>8</sub> x 5 <sup>1</sup> / <sub>2</sub> inch 562 x 92 x 140 mm	1 <sup>13</sup> / <sub>16</sub> x 3 <sup>7</sup> / <sub>16</sub> x 5 <sup>3</sup> / <sub>16</sub> inch 45.5 x 88 x 131.5 mm	12	DR-22	3.8 / 1.73
CFDP-13	23 <sup>15</sup> / <sub>16</sub> x 3 <sup>5</sup> / <sub>8</sub> x 5 <sup>1</sup> / <sub>2</sub> inch 608 x 92 x 140 mm	1 <sup>13</sup> / <sub>16</sub> x 3 <sup>7</sup> / <sub>16</sub> x 5 <sup>3</sup> / <sub>16</sub> inch 45.5 x 88 x 131.5 mm	13	DR-23	4.1 / 1.86
CFDP-14	25 <sup>3</sup> / <sub>4</sub> x 3 <sup>5</sup> / <sub>8</sub> x 5 <sup>1</sup> / <sub>2</sub> inch 654 x 92 x 140 mm	1 <sup>13</sup> / <sub>16</sub> x 3 <sup>7</sup> / <sub>16</sub> x 5 <sup>3</sup> / <sub>16</sub> inch 45.5 x 88 x 131.5 mm	14	DR-25	4.4 / 1.99
CFDP-15	27 <sup>9</sup> / <sub>16</sub> x 3 <sup>5</sup> / <sub>8</sub> x 5 <sup>1</sup> / <sub>2</sub> inch 700 x 92 x 140 mm	1 <sup>13</sup> / <sub>16</sub> x 3 <sup>7</sup> / <sub>16</sub> x 5 <sup>3</sup> / <sub>16</sub> inch 45.5 x 88 x 131.5 mm	15	DR-27	4.7 / 2.12
CFDP-16	29 <sup>7</sup> / <sub>16</sub> x 3 <sup>5</sup> / <sub>8</sub> x 5 <sup>1</sup> / <sub>2</sub> inch 747 x 92 x 140 mm	1 <sup>13</sup> / <sub>16</sub> x 3 <sup>7</sup> / <sub>16</sub> x 5 <sup>3</sup> / <sub>16</sub> inch 45.5 x 88 x 131.5 mm	16	DR-29	5.0 / 2.25

\* Lock device is ordered separately.