

# BenchPress

Valid Through 07.31.24.  
Pricing, availability, and specifications subject to change without notice.

Vol. 1.1 2024  
Promo: BP1.24

Promos and Lab Favorites

Choose one  
**FREE**

With purchase of item below\*

(Free with items: C2000, B5200, H5100-HCT, T5 Series, B2410, MR9600/9610)



Value: \$347

Includes microtube  
& PCR rotors



Value: \$298

Speed Range:  
200-3200 RPM



Value: \$100

Amazon Gift Card

Your choice of:

\*Promo Code BP1.24

**MyFuge™ Mini**

Item: C1008-G

**BenchMixer™**

Item: BV1000

**Gift Card**

Item: AMZ-100

## PLATE FUGE™

MICROPLATE MICROCENTRIFUGE



- Quickly spins down droplets in PCR plates
- Unique rotor design prevents sample spillage
- 50% smaller than traditional centrifuges
- Accepts all popular PCR plates (skirted, semi & non) & most microtiter plates



Before/After Centrifugation During Centrifugation

C2000 PlateFuge™ Microplate Microcentrifuge \$894

## TC 9639

Thermal Cycler



- Intuitive touch screen control
- Gradient Capabilities
- Multi-format block accepts:
  - 0.2ml tubes or strips (96)
  - 0.5ml tubes (39)
  - 96 well plate (1)

Touch screen cycler  
0.2ml tubes or strips (32)



T5000-96 TC 9639 Thermal Cycler, multi-format block \$5,995  
T5000-384 TC 9639 Thermal Cycler, 384 well block \$6,979  
T5005-3205 TC 32 Mini Thermal Cycler \$3,775

## SUREAIR™

PCR WORKSTATION WITH UV AND HEPA FILTRATION



- HEPA filter and UV light reduce cross-contamination during PCR preparation
- Compact, 55 x 40 cm footprint
- Cable port allows easy use of centrifuges, vortexers and other instruments
- Includes pipette holder and shelf for storage of pipettes and tips
- Provides clean air (to ISO 5 standards)

B5200 SureAir™ PCR Workstation \$2,279  
Accessories:  
B5200-FIL SureAir™ Replacement Filter \$307  
B1450-IC UV-C Indicator Cards, 25 ct. \$108  
B1450-ISS UV-C Indicator Stickers, 40 ct. \$96

## BEAD BATH DUO™

THE ONLY BATHS DESIGNED FOR BEADS OR WATER



- Two separate modes for BEADS or WATER
- Quik-Cal™, simple in-lab calibration function
- Offers all of the advantages of a "Bead Bath"
  - Easy (accessory free) securing of tubes and bottles
  - Allows angled incubations
  - Runs constant (no wait for heat-up)



B2410 BeadBath Duo™ 10L \$1,424

## MULTI-THERM™ Touch

COOL - HEAT - SHAKE



- Large color touch screen
- Multiple block options for plates & tubes up to 50ml
- Broad temperature & speed range
- Advanced programming features

Speed range: 200-3,000 rpm (block dependent)  
Shaker Orbit: 3mm  
Temp. range: Ambient -20° to 105°C



Heated lid prevents  
condensation & evaporation

H5100-HCT MultiTherm Touch, Heating/Cooling Vortexer (blocks sold separately) \$4,189  
H5100-HL Optional heated lid \$803  
Available blocks\*  
H5100-05 24 x 0.5ml \$374  
H5100-MP Block, microplate \$374  
H5100-384 Block, 384 well plate \$641  
\*Additional blocks available, see website for details

## SMARTREADER™ 96

Your go-to instrument for  
ELISA and absorbance assays



- Intuitive, 7 inch, color touch screen control panel
- Filter based, wavelength range from 340-750nm (405, 450, 492, 630nm filter included)
- Operates as a stand alone system, with USB flash drive for data transfer
- Integrated, multi-speed plate shaker

Also Available:  
MW9600  
Microplate Washer



MR9600 SmartReader™ 96 Microplate Absorbance Reader \$6,730  
MR9600-T SmartReader™ 96 Microplate Absorbance Reader with Incubation \$7,164  
MR9600-TP SmartReader™ Thermal Printer \$368