



# Freezer Racks & Boxes for Ultra-Low Freezers, Freezers, and Refrigerators



**2025 Cryo Storage Solutions**

*Our Most Popular Premium Racks & Boxes*



# Table of Contents

Standard Side Access Racks for 2 & 3" Boxes . . . . .	4
Adjustable Freezer Racks . . . . .	5
Sliding Tray Freezer Racks for 2 & 3" Boxes . . . . .	6
Premium Drawer Freezer Racks for 2 & 3" Boxes . . . . .	7
Vertical Racks for 2 & 3" Boxes . . . . .	8
Specialty Vertical Racks for 2" Boxes . . . . .	9
Vertical Racks for SBS Plates, Racks/Boxes. . . . .	10
Premium Drawer Freezer Racks for 15mL & 50mL Centrifuge Tubes . . . . .	11
Side Access Racks for 15mL & 50mL Centrifuge Tube Boxes . . . . .	11
Standard Side Access Racks for SBS Plates, 96 Well & 384 Well Plates. . . . .	12
Premium Drawer Racks for 96-Well & 384 Well Plates . . . . .	12
Freezer Boxes . . . . .	13 & 14
Rack Configurator . . . . .	15

Biomedical Solutions, Inc. offers the widest range of Freezer Racks and Boxes to fit most brands/models of Ultra-Low Freezers, Upright Freezers, Chest Freezers and Liquid Nitrogen Tanks. This Mini Catalog highlights our most popular models, please see our website [bsilab.com](http://bsilab.com) or ask your representative for more sizes and configurations.

### All freezer racks:

- Are constructed of high grade, corrosion resistant stainless steel
- Come with an unconditional lifetime warranty on workmanship
- Manufactured using the latest precision technology, robotic welding, CNC machines, and CAD design software
- Feature double welds on all seams
- Standard racks for 2", 3", 4" and 5" Boxes, Microplates, Specialty Tube Racks and 15mL & 50mL Centrifuge Tubes
- Most Racks in stock for quick delivery
- Cryo Storage Solutions for LN2 Tanks available, contact us!



**Don't see what you need? Contact us for easy, economical Custom Racks!**

**We make choosing the right Freezer Racks and Boxes to optimize your storage space easy, check out our website for our Freezer Rack Configurator. <https://www.bsilab.com/freezer-rack-configurator>**

# Freezer Racks & Boxes for Ultra-Low Freezers, Freezers, and Refrigerators

## Standard Side Access Racks for 2" & 3" Boxes



### ◆ Features

- Constructed of high grade corrosion resistant stainless steel.
- Precision craftsmanship with tolerances of +/- 1 mm.
- Compartments accept all cardboard, plastic and metal boxes.
- Handles located on both ends for easy handling.

#### For 2" Freezer Boxes

Item#	Description / Capacity	L x H x W
<b>UF-332</b>	3 x 3 Rack / for 9 Boxes	16-7/8 x 6-5/8 x 5-1/2"
<b>UF-342</b>	3 x 4 Rack / for 12 Boxes	16-7/8 x 8-7/8 x 5-1/2"
<b>UF-352</b>	3 x 5 Rack / for 15 Boxes	16-7/8 x 11 x 5-1/2"
<b>UF-432</b>	4 x 3 Rack / for 12 Boxes	22-3/16 x 6-5/8 x 5-1/2"
<b>UF-442</b>	4 x 4 Rack / for 16 Boxes	22-3/16 x 8-7/8 x 5-1/2"
<b>UF-452</b>	4 x 5 Rack / for 20 Boxes	22-3/16 x 11 x 5-1/2"
<b>UF-462*</b>	4 x 6 Rack / for 24 Boxes	22-3/16 x 12-5/8 x 5-1/2"
<b>UF-472</b>	4 x 7 Rack / for 28 Boxes	22-3/16 x 15-3/8 x 5-1/2"
<b>UF-542</b>	5 x 4 Rack / for 20 Boxes	27-5/16 x 8-7/8 x 5-1/2 "
<b>UF-552</b>	5 x 5 Rack / for 25 Boxes	27-5/16 x 11 x 5-1/2"

#### For 3" Freezer Boxes

Item#	Description / Capacity	L x H x W
<b>UF-323</b>	3 x 2 Rack / for 6 Boxes	16-7/8 x 6-5/8 x 5-1/2"
<b>UF-333</b>	3 x 3 Rack / for 9 Boxes	16-7/8 x 9-3/8 x 5-1/2"
<b>UF-423</b>	4 x 2 Rack / for 8 Boxes	22-3/16 x 6-5/8 x 5-1/2"
<b>UF-433</b>	4 x 3 Rack / for 12 Boxes	22-3/16 x 9-3/18 x 5-1/2"
<b>UF-443</b>	4 x 4 Rack / for 16 Boxes	22-3/16 x 12-5/8 x 5-1/2"
<b>UF-453</b>	4 x 5 Rack / for 20 Boxes	22-3/16 x 15-9/16 x 5-1/2"
<b>UF-523</b>	5 x 2 Rack / for 10 Boxes	27-5/16 x 6-5/8 x 5-1/2"
<b>UF-533</b>	5 x 3 Rack / for 15 Boxes	27-5/16 x 9-1/2 x 5-1/2"

\*UF-462 is designed to hold standard 2" cardboard boxes only, plastic boxes do not fit.

Ask your representative about security lock accessories available for these racks. We have over 1,000+ sizes & styles of racks. Ask your representative if you don't see the size you need.



# Freezer Racks & Boxes for Ultra-Low Freezers, Freezers, and Refrigerators

## Adjustable Freezer Racks for 2" & 3" Boxes



### ◆ Features

- Allows you to store a variety of different boxes in one rack.
- Adjustable & removable stainless steel clips to accommodate a wide variety of box heights.
- Lighter... So it's easier to take in and out of your freezer.

For 2" Freezer Boxes		
Item#	Description / Capacity	L x H x W
UFMC-33-2-3	3 x 3 Adj. Rack / for 9 Boxes	16-15/16 x 6-5/8 x 5-1/2"
UFMC-34-2-3	3 x 4 Adj. Rack / for 12 Boxes	16-15/16 x 8-7/8 x 5-1/2"
UFMC-35-2-3	3 x 5 Adj. Rack / for 15 Boxes	16-15/16 x 11 x 5-1/2"
UFMC-43-2-3	4 x 3 Adj. Rack / for 12 Boxes	22-3/8 x 6-5/8 x 5-1/2"
UFMC-44-2-3	4 x 4 Adj. Rack / for 16 Boxes	22-3/8 x 8-7/8 x 5-1/2"
UFMC-45-2-3	4 x 5 Adj. Rack / for 20 Boxes	22-3/8 x 11 x 5-1/2"
UFMC-46-2-3	4 x 6 Adj. Rack / for 24 Boxes	22-3/8 x 12-5/8 x 5-1/2"
UFMC-54-2-3	5 x 4 Adj. Rack / for 20 Boxes	27-5/16 x 8-7/8 x 5-1/2"
UFMC-55-2-3	5 x 5 Adj. Rack / for 25 Boxes	27-5/16 x 11 x 5-1/2"

For 3" Freezer Boxes		
Item#	Description / Capacity	L x H x W
UFMC-33-2-3	3 x 2 Adj. Rack / for 6 Boxes	16-15/16 x 6-5/8 x 5-1/2"
UFMC-34-2-3	3 x 2 Adj. Rack / for 6 Boxes	16-15/16 x 8-7/8 x 5-1/2"
UFMC-35-2-3	3 x 2 Adj. Rack / for 6 Boxes	16-15/16 x 11 x 5-1/2"
UFMC-43-2-3	4 x 2 Adj. Rack / for 8 Boxes	22-3/8 x 6-5/8 x 5-1/2"
UFMC-44-2-3	4 x 2 Adj. Rack / for 8 Boxes	22-3/8 x 8-7/8 x 5-1/2"
UFMC-45-2-3	4 x 2 Adj. Rack / for 8 Boxes	22-3/8 x 11 x 5-1/2"
UFMC-46-2-3	4 x 2 Adj. Rack / for 12 Boxes	22-3/8 x 12-5/8 x 5-1/2"
UFMC-54-2-3	5 x 2 Adj. Rack / for 10 Boxes	27-5/16 x 8-7/8 x 5-1/2"
UFMC-55-2-3	5 x 2 Adj. Rack / for 10 Boxes	27-5/16 x 11 x 5-1/2"

\*All Above Adjustable Racks Come Complete with Clips.

Additional Clips		
Item#	Description / Capacity	L x H x W
UFMC-L	"L" clip	3-5/8 x 3/8 x 1"
UFMC-T	"T" clip	3-5/8 x 3/8 x 1/2"

Need help with which rack fits your freezer?

Not storing your samples in boxes?  
Ask your representative about our new bin style racks.

# Freezer Racks & Boxes for Ultra-Low Freezers, Freezers, and Refrigerators

## Sliding Tray Freezer Racks for 2" & 3" Boxes



### ◆ Features

- Convenience of accessing each tray individually.
- Premium performance at a value price.
- Space saving design to maximize your storage needs.

For 2" Freezer Boxes		
Item#	Description / Capacity	L x H x W
<b>UFS-332</b>	3 x 3 Slide Rack / for 9 Boxes	16-7/8 x 6-5/8 x 5-1/2"
<b>UFS-342</b>	3 x 4 Slide Rack / for 12 Boxes	16-7/8 x 8-7/8 x 5-1/2"
<b>UFS-352</b>	3 x 5 Slide Rack / for 15 Boxes	16-7/8 x 11 x 5-1/2"
<b>UFS-432</b>	4 x 3 Slide Rack / for 12 Boxes	22-3/8 x 6-5/8 x 5-1/2"
<b>UFS-442</b>	4 x 4 Slide Rack / for 16 Boxes	22-3/8 x 8-7/8 x 5-1/2"
<b>UFS-452</b>	4 x 5 Slide Rack / for 20 Boxes	22-3/8 x 11 x 5-1/2"
<b>UFS-522</b>	5 x 2 Slide Rack / for 10 Boxes	27-1/16 x 4-9/16 x 5-1/2"
<b>UFS-532</b>	5 x 3 Slide Rack / for 15 Boxes	27-1/16 x 6-5/8 x 5-1/2"
<b>UFS-542</b>	5 x 4 Slide Rack / for 20 Boxes	27-1/16 x 8-7/8 x 5-1/2"
<b>UFS-552</b>	5 x 5 Slide Rack / for 25 Boxes	27-1/16 x 11 x 5-1/2"

For 3" Freezer Boxes		
Item#	Description / Capacity	L x H x W
<b>UFS-333</b>	3 x 3 Slide Rack / for 9 Boxes	16-7/8 x 9-1/2 x 5-1/2"
<b>UFS-343</b>	3 x 4 Slide Rack / for 12 Boxes	16-7/8 x 12-11/16 x 5 1/2"
<b>UFS-353</b>	3 x 5 Slide Rack / for 15 Boxes	16-7/8 x 15-13/16 x 5 1/2"
<b>UFS-433</b>	4 x 3 Slide Rack / for 12 Boxes	22-3/8 x 9-1/2 x 5-1/2"
<b>UFS-443</b>	4 x 4 Slide Rack / for 16 Boxes	22-3/8 x 12-11/16 x 5-1/2"
<b>UFS-453</b>	4 x 5 Slide Rack / for 20 Boxes	22 3/8 x 15-13/16 x 5 1/2"
<b>UFS-533</b>	5 x 3 Slide Rack / for 15 Boxes	27-1/16 x 9-1/2 x 5-1/2"
<b>UFS-543</b>	5 x 4 Slide Rack / for 20 Boxes	22-1/16 x 12-11/16 x 5-1/2"

Cardboard boxes are recommended for use with the sliding Tray Racks. Looking for a smaller or larger rack? Ask your representative, as we have additional sizes.

# Freezer Racks & Boxes for Ultra-Low Freezers, Freezers, and Refrigerators

## Premium Drawer Freezer Racks for 2" & 3" Boxes



### ◆ Features

- Fully framed providing even greater sample protection.
- Precision craftsmanship with tolerances of +/- 1 mm.
- Compartment options to accept all cardboard, plastic and metal boxes.
- Label place card holders for easy identification.

For 2" Freezer Boxes		
Item#	Description / Capacity	L x H x W
UFD-332	3 x 3 Drawer Rack / for 9 Boxes	16-3/4 x 7-1/8 x 5-1/2"
UFD-342	3 x 4 Drawer Rack / for 12 Boxes	16-3/4 x 9-7/16 x 5-1/2"
UFD-352	3 x 5 Drawer Rack / for 15 Boxes	16-3/4 x 11-13/16 x 5-1/2"
UFD-432	4 x 3 Drawer Rack / for 12 Boxes	22-5/16 x 7-1/8 x 5-1/2"
UFD-442	4 x 4 Drawer Rack / for 16 Boxes	22-5/16 x 9-7/16 x 5-1/2"
UFD-452	4 x 5 Drawer Rack / for 20 Boxes	22-5/16 x 11-13/16 x 5-1/2"
UFD-462	4 x 6 Drawer Rack / for 24 Boxes	22-5/16 x 14-1/8 x 5-1/2"
UFD-462-S*	4 x 6 Drawer Rack / for 24 Boxes	22-5/16 x 12-5/8 x 5-1/2"
UFD-462-SA**	4 x 6 Drawer Rack / for 24 Boxes	22-5/16 x 13-1/4 x 5-1/2"
UFD-542	5 x 4 Drawer Rack / for 20 Boxes	27 x 9-7/16 x 5-1/2"
UFD-552	5 x 5 Drawer Rack / for 25 Boxes	27 x 11-13/16 x 5-1/2"

For 3" Freezer Boxes		
Item#	Description / Capacity	L x H x W
UFD-323	3 x 2 Drawer Rack / for 6 Boxes	16-3/4 x 6-7/8 x 5-1/2"
UFD-333	3 x 3 Drawer Rack / for 9 Boxes	16-3/4 x 10-1/16 x 5-1/2"
UFD-423	4 x 2 Drawer Rack / for 8 Boxes	22-5/16 x 6-7/16 x 5-1/2"
UFD-433	4 x 3 Drawer Rack / for 12 Boxes	22-5/16 x 10-1/16 x 5-1/2"
UFD-443	4 x 4 Drawer Rack / for 16 Boxes	22-5/16 x 13-1/2 x 5-1/2"
UFD-453	4 x 5 Drawer Rack / for 20 Boxes	22-5/16 x 16-9/16 x 5-1/2"
UFD-523	5 x 2 Drawer Rack / for 20 Boxes	27-1/16 x 6-7/8 x 5-1/2"
UFD-533	5 x 3 Drawer Rack / for 15 Boxes	27-1/16 x 10-1/16 x 5-1/2"

Constructed of high grade corrosion resistant stainless steel with precision craftsmanship. Ask your representative about label cards and security lock accessories for these racks.

\*Top-Loading for Panasonic Freezers for 2" Cardboard Boxes only

\*\*Side-Loading for Panasonic Freezers for 2" Polycarbonate Boxes

# Freezer Racks & Boxes for Ultra-Low Freezers, Freezers, and Refrigerators

## Vertical Racks for 2" & 3" Boxes



### ◆ Features

- Security lock devices are available for each rack.
- Constructed of corrosion resistant stainless steel.
- Locking Rods included for each rack.

### For 2" Freezer Boxes, Stainless Steel Rack

Item#	Description / Capacity	L x H x W
<b>CF-4-2</b>	4 Boxes	9-1/8 x 5-5/8 x 5-1/2"
<b>CF-5-2</b>	5 Boxes	11-5/16 x 5-5/8 x 5-1/2"
<b>CF-6-2</b>	6 Boxes	13-1/2 x 5-5/8 x 5-1/2"
<b>CF-7-2</b>	7 Boxes	15-11/16 x 5-5/8 x 5-1/2"
<b>CF-8-2</b>	8 Boxes	17-7/8 x 5-5/8 x 5-1/2"
<b>CF-9-2</b>	9 Boxes	20-1/16 x 5-5/8 x 5-1/2"
<b>CF-10-2</b>	10 Boxes	22-5/16 x 5-5/8 x 5-1/2"
<b>CF-11-2</b>	11 Boxes	24-1/2 x 5-5/8 x 5-1/2"
<b>CF-12-2</b>	12 Boxes	26-5/8 x 5-5/8 x 5-1/2"
<b>CF-13-2</b>	13 Boxes	28-1/2 x 5-5/8 x 5-1/2"
<b>CF-14-2</b>	14 Boxes	31-1/16 x 5-5/8 x 5-1/2"
<b>CF-15-2</b>	15 Boxes	33-1/4 x 5-5/8 x 5-1/2"

### For 3" Freezer Boxes, Stainless Steel Rack

Item#	Description / Capacity	L x H x W
<b>CF-4-3</b>	4 Boxes	13 x 5-5/8 x 5-1/2"
<b>CF-5-3</b>	5 Boxes	16-3/16 x 5-5/8 x 5-1/2"
<b>CF-6-3</b>	6 Boxes	19-3/8 x 5-5/8 x 5-1/2"
<b>CF-7-3</b>	7 Boxes	22-1/2 x 5-5/8 x 5-1/2"
<b>CF-8-3</b>	8 Boxes	25-11/16 x 5-5/8 x 5-1/2"
<b>CF-9-3</b>	9 Boxes	28-1/2 x 5-5/8 x 5-1/2"
<b>CF-10-3</b>	10 Boxes	32-1/16 x 5-5/8 x 5-1/2"



# Freezer Racks & Boxes for Ultra-Low Freezers, Freezers, and Refrigerators

## Specialty Vertical Racks for 2" Boxes

### ◆ Features

- Made of aluminum.
- Security lock devices are available for each rack.



### For 2" Freezer Boxes, Aluminum Rack

Item#	Description / Capacity	L x H x W
<b>CFA-10-2</b>	10 Boxes	22-5/16 x 5-5/8 x 5-1/2"
<b>CFA-11-2</b>	11 Boxes	24-1/2 x 5-5/8 x 5-1/2"
<b>CFA-12-2</b>	12 Boxes	26-5/8 x 5-5/8 x 5-1/2"
<b>CFA-13-2</b>	13 Boxes	28-9/16 x 5-5/8 x 5-1/2"
<b>CFA-14-2</b>	14 Boxes	31-1/16 x 5-5/8 x 5-1/2"
<b>CFA-15-2</b>	15 Boxes	33-1/4 x 5-5/8 x 5-1/2"

\*Other sizes available

### ◆ Features

- Made of corrosion resistant stainless steel.
- Security lock devices are available for each rack.
- Spring Clips securely hold boxes in rack compartments without using locking rods and allow for easy removal of boxes.



### For 2" Freezer Boxes, Stainless Steel Rack with Spring Clip

Item#	Description / Capacity	L x H x W
<b>CF-9-2-S</b>	9 Boxes	20-1/16 x 5-5/8 x 5-1/2"
<b>CF-10-2-S</b>	10 Boxes	22-5/16 x 5-5/8 x 5-1/2"
<b>CF-11-2-S</b>	11 Boxes	24-1/2 x 5-5/8 x 5-1/2"
<b>CF-12-2-S</b>	12 Boxes	26-5/8 x 5-5/8 x 5-1/2"
<b>CF-13-2-S</b>	13 Boxes	28-1/2 x 5-5/8 x 5-1/2"
<b>CF-14-2-S</b>	14 Boxes	31-1/16 x 5-5/8 x 5-1/2"
<b>CF-15-2-S</b>	15 Boxes	33-1/4 x 5-5/8 x 5-1/2"

### ◆ Features

- Made of corrosion resistant stainless steel.
- Racks are ideal for light load and easy handing of sample storage in freezers or liquid nitrogen tanks.



### For 2" Freezer Boxes, Vertical Half Racks

Item#	Description / Capacity	L x H x W
<b>LR-3-2</b>	3 Boxes	16-1/2 x 2-1/4 x 5-9/16"
<b>LR-4-2</b>	4 Boxes	21-5/16 x 2-1/4 x 5-9/16"
<b>LR-5-2</b>	5 Boxes	26-3/4 x 2-1/4 x 5-9/16"
<b>LR-6-2</b>	6 Boxes	32-13/16 x 2-1/4 x 5-9/16"

# Freezer Racks & Boxes for Ultra-Low Freezers, Freezers, and Refrigerators

## Vertical Racks for SBS Plates, Racks/Boxes



For 15mL & 50mL  
Centrifuge Tube Boxes



For 96-Deep  
Well Plates



For Specialty  
Tube Racks



For 100 Cell Hinged  
Plastic Boxes

### CFLB Series For 15mL & 50mL Centrifuge Boxes

Item#	Description / Capacity	L x H x W
<b>CFLB-4</b>	4 Boxes	20-3/8 x 6 x 6-3/16"
<b>CFLB-5</b>	5 Boxes	25-3/8 x 6 x 6-3/16"
<b>CFLB-6</b>	6 Boxes	30-3/8 x 6 x 6-3/16"

Other sizes and configurations available see our website or ask your representative.

### CFMX Series For Specialty Tube Racks

Item#	Description / Capacity	L x H x W
<b>CFMX-9-2</b>	9 Boxes	22-5/8 x 3-5/8 x 5-1/2"
<b>CFMX-10-2</b>	10 Boxes	25-1/8 x 3-5/8 x 5-1/2"
<b>CFMX-11-2</b>	11 Boxes	27-5/8 x 3-5/8 x 5-1/2"
<b>CFMX-12-2</b>	12 Boxes	30-1/16 x 3-5/8 x 5-1/2"
<b>CFMX-13-2</b>	13 Boxes	32-9/16 x 3-5/8 x 5-1/2"

Other sizes and configurations available see our website or ask your representative.

### CFDP Series For 96-Deep Well Plates

Item#	Description / Capacity	L x H x W
<b>CFDP-9</b>	9 Boxes	16-5/8 x 3-5/8 x 5-1/2"
<b>CFDP-10</b>	10 Boxes	18-7/16 x 3-5/8 x 5-1/2"
<b>CFDP-11</b>	11 Boxes	20-1/4 x 3-5/8 x 5-1/2"
<b>CFDP-12</b>	12 Boxes	22-1/8 x 3-5/8 x 5-1/2"
<b>CFDP-13</b>	13 Boxes	23-15/16 x 3-5/8 x 5-1/2"
<b>CFDP-14</b>	14 Boxes	25-3/4 x 3-5/8 x 5-1/2"
<b>CFDP-15</b>	15 Boxes	27-9/16 x 3-5/8 x 5-1/2"
<b>CFDP-16</b>	16 Boxes	29-7/16 x 3-5/8 x 5-1/2"

Other sizes and configurations available see our website or ask your representative.

### CFHT Series For 100 Cell Hinged Plastic Boxes

Item#	Description / Capacity	L x H x W
<b>CFHT-8</b>	8 Boxes	18-3/4 x 5-7/8 x 6-1/4"
<b>CFHT-9</b>	9 Boxes	21 x 5-7/8 x 6-1/4"
<b>CFHT-10</b>	10 Boxes	23-3/8 x 5-7/8 x 6-1/4"
<b>CFHT-11</b>	11 Boxes	25-5/8 x 5-7/8 x 6-1/4"
<b>CFHT-12</b>	12 Boxes	27-15/16 x 5-7/8 x 6-1/4"
<b>CFHT-13</b>	13 Boxes	30-1/4 x 5-7/8 x 6-1/4"

Other sizes and configurations available see our website or ask your representative.

# Freezer Racks & Boxes for Ultra-Low Freezers, Freezers, and Refrigerators

## Premium Drawer Freezer Racks for 15mL & 50mL Centrifuge Tubes

### ◆ Features

- Available in a wide variety of formats:  
1 or 2 drawer, 15mL, 50mL & 15mL & 50mL.
- Precision craftsmanship with tolerances of +/- 1 mm.
- Easy drawer access to each tube and no need for boxes.



### For 15mL Centrifuge Tubes

Item#	Description / Capacity	L x H x W
<b>UFD-T15-2</b>	120 Tubes (15mL)	16-11/16 x 9-15/16 x 5-1/2"
<b>UFD-LT15-2</b>	160 Tubes (15mL)	22-3/16 x 9-15/16 x 5-1/2"
<b>UFD-XLT15-2</b>	200 Tubes (15mL)	26-15/16 x 9-15/16 x 5-1/2"

### For 50mL Centrifuge Tubes

Item#	Description / Capacity	L x H x W
<b>UFD-T50-2</b>	60 Tubes (50mL)	16-11/16 x 9-15/16 x 5-1/2"
<b>UFD-LT50-2</b>	78 Tubes (50mL)	22-3/16 x 9-15/16 x 5-1/2"
<b>UFD-XLT50-2</b>	102 Tubes (50mL)	26-15/16 x 9-15/16 x 5-1/2"

### For 15mL & 50mL Centrifuge Tubes

Item#	Description / Capacity	L x H x W
<b>UFD-T15/50</b>	60 Tubes (15mL) + 30 Tubes (50mL)	16-11/16 x 9-15/16 x 5-1/2"
<b>UFD-LT15/50</b>	80 Tubes (15mL) + 39 Tubes (50mL)	22-3/16 x 9-15/16 x 5-1/2"
<b>UFD-XLT15/50</b>	100 Tubes (15mL) + 51 Tubes (50mL)	26-15/16 x 9-15/16 x 5-1/2"

Storing a different tube size? We have more racks for those too... Ask your representative.

## Side Access Racks for 15mL & 50mL Centrifuge Tube Boxes



### For 15mL & 50mL Centrifuge Tube Boxes

Item#	Description / Capacity	L x H x W
<b>UFLB-22</b>	2 x 2 Rack / for 4 Boxes	12-3/4 x 10-3/16 x 5-1/2"
<b>UFLB-32</b>	3 x 2 Rack / for 6 Boxes	18-13/16 x 10-3/16 x 5-1/2"
<b>UFLB-42</b>	4 x 2 Rack / for 8 Boxes	24-13/16 x 10-3/16 x 5-1/2"

# Freezer Racks & Boxes for Ultra-Low Freezers, Freezers, and Refrigerators

## Standard Side Access Racks for SBS Plates, 96 Well & 384 Well Plates

### ◆ Features

- Locking rods are included.
- Security lock devices are available for each rack.
- Made of corrosion resistant stainless steel.



Item#	Description / Capacity	L x H x W
<b>UFMP-509L</b>	70 Plates w / lids 80 Plates w / o lids	18-13/16 x 9-3/16 x 5-1/2"
<b>UFMP-511L</b>	85 Plates w / lids 95 Plates w / o lids	18-13/16 x 11 x 5-1/2"
<b>UFMP-609L</b>	84 Plates w / lids 96 Plates w / o lids	22-1/2 x 9-3/16 x 5-1/2"
<b>UFMP-611L</b>	96 Plates w / lids 120 Plates w / o lids	22-1/2 x 11 x 5-1/2"
<b>UFMP-709L</b>	98 Plates w / lids 112 Plates w / o lids	26-1/8 x 9-3/16 x 5-1/2"
<b>UFMP-711L</b>	112 Plates w / lids 140 Plates w / o lids	26-1/8 x 11 x 5-1/2"

Other sizes and configurations available see our website or ask your representative.

## Premium Drawer Racks for 96-Well & 384-Well Plates

### ◆ Features

- Fully framed providing greater sample protection.
- Compartment fits all standard 96-Well & 384-Well Plates
- Label place card holders for easy identification




Item#	Description / Capacity	L x H x W
<b>UFDMP-54</b>	4 drawers 60 Plates w / lids 80 Plates w / o lids	18-3/16 x 9-3/4 x 5-1/2"
<b>UFDMP-55</b>	5 drawers 75 Plates w / lids 100 Plates w / o lids	18-3/16 x 12-1/16 x 5-1/2"
<b>UFDMP-64</b>	4 drawers 72 Plates w / lids 96 Plates w / o lids	21-3/4 x 9-3/4 x 5-1/2"
<b>UFDMP-65</b>	5 drawers 90 Plates w / lids 120 Plates w / o lids	21-3/4 x 12-1/8 x 5-1/2"
<b>UFDMP-66</b>	6 drawers 108 Plates w / lids 144 Plates w / o lids	21-3/4 x 14-1/2 x 5-1/2"
<b>UFDMP-74</b>	4 drawers 84 Plates w / lids 112 Plates w / o lids	25-1/4 x 9-3/4 x 5-1/2"

Other sizes and configurations available see our website or ask your representative.

# Freezer Racks & Boxes for Ultra-Low Freezers, Freezers, and Refrigerators

## Freezer Boxes



### Water Proof Cardboard Boxes

Item#	Description / Capacity	L x W x H	Color
<b>SB2CC-H-81</b>	2" Water Proof Cardboard w/ Hinged Lid w/81 Place Cell Divider inserted	5-1/4 x 5-1/4 x 2"	White
<b>SB2CC-H-81-P-W</b>	2" Water Proof Cardboard w/ Hinged Lid w/81 Place Cell Divider inserted	5-1/4 x 5-1/4 x 2"	Purple
<b>SB2CC-H-81-G-W</b>	2" Water Proof Cardboard w/ Hinged Lid w/81 Place Cell Divider inserted	5-1/4 x 5-1/4 x 2"	Green

### Water Repellent Cardboard Boxes

Item#	Description / Capacity	L x W x H	Color
<b>SB2CC</b>	Standard 2" Cardboard Box	5-1/4 x 5-1/4 x 2"	White
<b>SB2CC-81</b>	Standard 2" Cardboard Box, with 81 Place Cell Divider inserted	5-1/4 x 5-1/4 x 2"	White
<b>SB2CC-100</b>	Standard 2" Cardboard Box, with 100 Place Cell Divider inserted	5-1/4 x 5-1/4 x 2"	White
<b>SB2CC-LP*</b>	Low Profile Water Repellent 2" Cardboard Box	5-5/16 x 5-5/16 x 1-15/16"	White
<b>SB2C*</b>	Standard 2" Cardboard Box	5-1/4 x 5-1/4 x 2"	Yellow
<b>SB3C*</b>	Standard 3" Cardboard Box	5-1/4 x 5-1/4 x 3"	White
<b>SB3.75C*</b>	Standard 4" Cardboard Box	5-1/4 x 5-1/4 x 3-7/8"	White
<b>SB4C*</b>	Standard 5" Cardboard Box	5-13/16 x 5-13/16 x 4-15/16"	White

\*Available with cell dividers

### Standard Box Divider


Item#	Description / Capacity	L x D x W	Cell Opening
<b>CD-64</b>	8 x 8 = 64 Cells	4-7/8 x 4-7/8 x 1-5/16"	9/16 inch
<b>CD-81-W</b>	9 x 9 = 81 Cells	4-7/8 x 4-7/8 x 1-5/16"	1/2 inch
<b>CD-100</b>	10 x 10 = 100 Cells	4-7/8 x 4-7/8 x 1-5/16"	7/16 inch
<b>CD-144</b>	12 x 12 = 144 Cells	4-7/8 x 4-7/8 x 1-5/16"	3/8 inch

Many more sizes available see our website or contact your representative



# Freezer Racks & Boxes for Ultra-Low Freezers, Freezers, and Refrigerators

## Freezer Boxes



### ◆ Features

Stainless steel boxes come with or without cardboard dividers and fit into standard freezer racks.

- 2" and 3" boxes available
- With or without cell dividers
- Fits into standard freezer boxes
- Drain holes



### ◆ Features

Aluminum boxes come with or without cardboard dividers and fit into standard freezer racks.

- Standard 2" Aluminum box with drain holes
- Standard 3" Aluminum box with drain holes
- Cell dividers available for each
- Fits into standard freezer boxes

### PolyCarbonate, Polypropylene and Plastic Boxes

Item#	Description / Capacity	L x H x W
PB2-100*	10 x 10 Polycarbonate Box	5-1/4 x 5-1/4 x 2-1/16"
PB2-81*	9 x 9 Polycarbonate Box	5-1/4 x 5-1/4 x 2-1/16"
PB4-81*	9 x 9 Polycarbonate Box	5-1/4 x 5-1/4 x 3-11/16"
MPB2-25*	5 x 5 Polycarbonate Mini Box	3 x 3 x 2-1/16"
PPP2-100	10 x 10 Polypropylene Box	5-1/4 x 5-1/4 x 2"
PPP2-81	9 x 9 Polypropylene Box	5-1/4 x 5-1/4 x 2"
PP2-100	100-Cell Hinged Top Plastic Box pk1 (Natural)	5-9/16 x 6 x 2-1/8"
PP2-100-A	100-Cell Hinged Top Box pk 5 (5 colors)	5-9/16 x 6 x 2-1/8"
PP2-100X*	100-Cell Hinged Top Box (Specify 5 color: B-Blue, G-Green, O-Orange, P-Pink, Y-Yellow)	5-9/16 x 6 x 2-1/8"

\*Suitable for use to -196°C

### Stainless Steel Boxes

Item#	Description / Capacity	Dimension Wt×Wt×Ht (Lid)	Dimension Wb×Wb×Hb
SB2S	Standard 2" Stainless Steel Box w/Drain Hole	5-1/4×5-1/4×2" 133×133×51mm	5 x 5 x 1-15/16" 127.7×127.7×49mm
SB2S-XX*	Standard 2" Stainless Steel Box w/Divider & Drain Hole	5-1/4×5-1/4×2" 133×133×51mm	5 x 5 x 1-15/16" 127.7×127.7×49mm
SB3S	Standard 3" Stainless Steel Box w/Drain Hole	5-1/4×5-1/4×3" 133×133×75mm	5 x 5 x 2-15/16" 127.7×127.7×74mm
SB3S-XX*	Standard 3" Stainless Steel Box w/Divider & Drain Hole	5-1/4×5-1/4×3" 133×133×75mm	5 x 5 x 2-15/16" 127.7×127.7×74mm

\* To order, replace XX with 144, 100, 81, 64, 49, 36, 25 or 16 for cell dividers.

### Aluminum Boxes

Item#	Description / Capacity	Dimension Wt×Wt×H	Dimension Wb×Wb×Hb
SB2A	Standard 2" Aluminum Box w/Drain Hole	5-1/4×5-1/4×2" 133×133×51mm	5 x 5 x 1-15/16" 127.7×127.7×49mm
SB2A-XX*	Standard 2" Aluminum Box w/Divider & Drain Hole	5-1/4×5-1/4×2" 133×133×51mm	5 x 5 x 1-15/16" 127.7×127.7×49mm
SB3A	Standard 3" Aluminum Box w/Drain Hole	5-1/4×5-1/4×3" 133×133×75mm	5 x 5 x 2-15/16" 127.7×127.7×74mm
SB3A-XX*	Standard 3" Aluminum Box w/Divider & Drain Hole	5-1/4×5-1/4×3" 133×133×75mm	5 x 5 x 2-15/16" 127.7×127.7×74mm

\* To order, replace XX with 144, 100, 81, 64, 49, 36, 25 or 16 for cell dividers.

**Sales & Product Support:**

3727 Greenbriar Dr. Suite 304,  
Stafford, TX. 77477  
866-674-7220 | 281-242-6294  
info@bsilab.com

