

streamlined

SAMPLE PREP

MIX AND STERILIZE

Buy any Lab Line Autoclave
and receive a **FREE Hei-PLATE
Mix 'n' Heat Core 145MM**



Save time and ensure accurate results! Sample preparation is key to ensuring reproducible and accurate results in your laboratory. Streamline your sample preparation with the combination of Heidolph Tuttnauer's premium lab line sterilizers and the reliability of the Hei-PLATE Mix 'n' Heat Core 145mm.

**Industry leading 10 Year Warranty on Hei-PLATE!
3 Year Warranty on all Lab Line Autoclaves!**



A Powerful Combination

Hei-PLATE Mix n Heat Core

The Hei-PLATE Mix 'n' Heat Core is easy to use with two control knobs for speed and temperature. The durable, quick-heating surface fits various samples and accessories for safe and reliable heating and stirring for a variety of applications.

Lab Line Autoclaves

Sterilize everything in your lab including liquids with the versatile lineup of Heidolph Tuttnauer Lab Line Autoclaves. Choose from benchtop, vertical, free standing configurations as you prevent contamination in your laboratory!



HEIDOLPH TUTTNAUER LAB LINE AUTOCLAVES



Liquid



Plastic



Glass



Sealed bottle



Tools



Wrapped Tools



Clothes



Bioreactors

Benchtop Lab Line

Model	Chamber Dimensions (in)	Volume (L)	P/N
2840ELP	11 x 18	28	023210671
3850ELP	15 x 19.7	65	023210681
3870ELP	15 x 27.2	85	023210672
5075ELP	19.7 x 29.5	160	023210675



Vertical Lab Line

Model	Chamber Dimensions (in)	Volume (L)	P/N
3950V	15.3 x 19.8	60	023210691
3970V	15.3 x 27.6	85	023210692
2840ELVP-WR	11 x 18	28	023210744
3870ELVP-WR	15 x 27.2	85	023210745
5075ELVCP	19.7 x 29.5	160	023210990



Free-Standing Lab Line

Model	Chamber Dimensions (in)	Volume (L)	P/N
3870EL-PVG-WS	14 x 27.2	85	023210680
5075EL-PVG-WS	19.7 x 29.5	160	023210685
2840ELC-PVG-P	11 x 18	28	023210683
3870ELC-PVG-P	15 x 27.2	85	023210678
5075ELC-PVG-P	19.7 x 29.5	160	023260002



Hei-PLATE Mix 'n' Heat MAGNETIC STIRRER

Hei-PLATE Mix 'n' Heat Core

Model	Dimensions	Plate Diameter	P/N
Hei-PLATE Mix 'n' Heat Core 145	168 x 101 x 299 mm	145 mm	036110200



Research and Development



Mining



Wastewater Treatment



Pharmaceutical Research



Agar Preparation



Organic Chemistry



heidolph
research made easy

